Leica DM4 B Leica DM4 M Leica DM4 P

Instructions for use
Gebrauchsanweisung
Mode d'emploi

Leica Microsystems CMS GmbH. Instructions for use, 11934109, Revision 1.1, 2017-05-01

Living up to Life

Leica
MICROSYSTEMS

Revision 1.1, published 2017-05-01 by:

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Leica Microsystems CMS GmbH

Ernst Leitz-Straße 17-37

D-35578 Wetzlar (Germany)

<http://www.leica-microsystems.com>

Responsible for contents:

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Leica
MICROSYSTEMS

Copyrights

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The term „Windows“ may appear in the following text without further identification. It is, however, a registered trademark of Microsoft Corporation. The names of companies and products used herein may be trademarks of their respective owners.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information in this manual is subject to modification at any time and without notification.

Contents

1. Important Notes about this Manual.....	7	6.10	Optional accessories	43
2. Intended purpose of the microscope....	10	6.11	Connection to the power supply	44
3. Safety notes.....	12	7. Startup.....	45	
3.1 General safety notes	12	7.1 Functional principle.....	45	
3.2 Electrical safety	13	7.2 Switching on the unit.....	48	
3.3 Laser coupling by the user.....	15	7.3 The display		
3.3.1 Safety hazard at eyepiece	16	(Leica DM4 B/DM4 M/DM4 P).....	49	
3.3.2 Safety hazard in the specimen area.....	16	7.4 The function keys	50	
3.3.3 Safety hazard due to openings		7.5 Köhler illumination.....	51	
at microscope	16	7.5.1 Transmitted light	51	
3.3.4 Notes regarding the safety hazards in		7.5.2 Incident light.....	52	
Section 3.3.1	17	7.6 Checking the phase contrast rings.....	54	
3.3.5 Tubes with integrated		7.7 Setting the motorized polarizers		
interlock switches	17	(DM4 P LED).....	55	
3.3.6. Exception for Class 1 lasers.....	20	7.8 Adjusting the light sources	55	
3.4 Transport and storage	20	8. Operation	61	
3.5 Notes on Handling Light Sources	21	8.1 Switching on the unit.....	61	
3.6 Notes on Handling Immersion Oil.....	21	8.2 Stages and object displacement	61	
3.7 Notes on Handling Acids and Bases.....	21	8.3 Focusing.....	63	
3.8 Disposal.....	21	8.4 Tubes	63	
4. Overview of the instrument.....	22	8.5 Eyepieces.....	65	
5. Unpacking the microscope	26	8.6 Objectives.....	65	
6. Assembling the microscope	29	8.7 Magnification changer	68	
6.1 Specimen stage	29	8.8 HC P 1x/1.6x tube optics with		
6.2 Condenser.....	31	Bertrand lens (encoded)	68	
6.3 Tube and eyepieces	32	8.9 Light sources.....	69	
6.4 Objectives	33	8.10 Aperture diaphragm and		
6.5 Light sources for DM4 B/M and DM4 P 33		field diaphragm	69	
6.6 Light sources for incident light fluores-				
cence axis DM4 B	34			
6.7 Equipping the incident light turret disk .	40			
6.8 Polarizer and analyzer	41			
6.9 DIC prisms.....	42			

Contents

9. Contrast methods for Leica DM4 B/Leica DM4 P.....	70
9.1 Transmitted light.....	70
9.1.1 Bright field (TL).....	70
9.1.2 Phase contrast (TL).....	70
9.1.3 Dark field (TL).....	71
9.1.4 Polarization (TL).....	71
9.1.4.1 Manual method.....	71
9.1.4.2 DM4 P - examinations in polarized transmitted light.....	72
9.1.4.3 Motorized method.....	78
9.1.4.4 Combined methods.....	78
9.1.5 Differential interference contrast (TL) .	78
9.1.5.1 DM4 P.....	78
9.2 Fluorescence.....	80
10. Contrast methods for Leica DM4 M.....	81
10.1 Incident light.....	81
10.1.1 Bright field (IL).....	81
10.1.2 Dark field (IL).....	81
10.1.3 Polarization (IL).....	82
10.1.4 Interference contrast (RL).....	83
10.2 Transmitted light.....	83
10.2.1 Bright field (TL).....	83
10.2.2 Polarization (TL).....	83
11. Troubleshooting.....	84
12. Care of the microscope.....	87
12.1 Dust cover.....	87
12.2 Cleaning.....	87
12.3 Handling acids and bases.....	88
13. Essential wear and spare parts.....	89
14. Abbreviations and pictograms	90
15. Index.....	91
16. EU Declaration of Conformity	93

1. Important Notes about this Manual

Caution!

This Operating Manual is an essential component of the product; it must be read carefully before the product is assembled, put into operation or used, and must be kept for later reference.

This operating manual contains important instructions and information for the operational safety and maintenance of the device.

Text symbols, pictograms and their meanings:

(1.2)

Numbers in parentheses such as “(1.2)”, correspond to illustrations (in the example, Figure 1, Item 2).

→ p.20

Numbers with an arrow, e.g. → p. 20, refer to a specific page of the manual.

WARNING indicates a hazard with a medium degree of risk that, if not avoided, can result in death or serious injury.

CAUTION indicates a hazard with a low degree of risk that, if not avoided, can result in minor or moderate injury.

Caution!

Special safety instructions within this manual are indicated with the triangle symbol shown here, and have a gray background.

Notice! The macroscope and accessories can be damaged when operated incorrectly.

1. Important notes about this manual

Warning of hazardous electrical voltage! Risk of electrical shock!

Danger due to hot surface.

Warning of optical radiation! Never look directly into the light beam! Wear safety goggles!

Warning of electromagnetic field

Warning of UV radiation! Wear safety goggles!

Warning of hazardous laser radiation – take appropriate safety precautions!

Warning of hand injuries

Instructions for disposing of the instrument, its accessories and consumables.

Connection to ground!

Explanatory note.

*

Item not contained in all configurations.

1. Important notes about this manual

IVD manufacturing date, example 11 / 2011 for November 2011.

Device for in vitro diagnostics (IVD).

Manufacturer's address with manufacturing date

China RoHS 50 years EFUP
(Environmentally friendly use period)

2. Intended purpose of the microscope

The Leica microscope DM4 B (11 888 858) is an upright light microscope and intended for use as a general lab microscope for routine examinations of biological specimens.

Directions for use:

Inspection, numbers, classification, identification and monitoring of cell and tissue cultures, examining biological specimens, liquids and sediments that lead to information relating to physiological and/or pathological states, inborn anomalies and the monitoring of therapeutic measures.

The resulting measures are based on the interpretation of the medical professional.

Contraindications:

Not suitable for testing potentially infectious specimens.

IVD

The Leica DM4 B (11 888 858) microscope complies with the EU Directive 98/79/EEC concerning in vitro diagnostics.

The Leica microscope DM4 B (11 888 866) is an upright, inverted light microscope and intended for use as general lab microscopes for routine examinations of biological specimens - not for IVD applications.

Directions for use:

Inspection, counting, classification, identification and monitoring of cell and tissue cultures, examining biological specimens, liquids and sediments.

Contraindications:

Not suitable for testing potentially infectious specimens.

The Leica DM4 M and DM4 P microscopes are intended for materials science, geological or mineralogical examinations and applications in the laboratory area.

The requirements in EU Directives 2006/95/EC and 2014/35/EU regarding the safety of electrical equipment, and 2004/108/EC and 2014/30/EU regarding electromagnetic compatibility have been met for all the instruments mentioned above for use in industrial environments.

Reasonably foreseeable misuse

The following are prohibited:

- Using the microscope for any purpose not in accordance with the Declaration of Conformity (e.g. using it as an in-vitro diagnostic product according to Council Directive 98/79/EC or as a medical product according to Council Directive 93/42/EEC)
- Making use of auxiliary equipment to use the area between the microscope stage and objectives as a clamp or holder (working like a vice).
- Operating the microscope in an inclined position.
- Cleaning the microscope in a way other than specified in the manual.
- Modifying the built-in safety circuits in the instrument.
- Allowing unauthorized personnel to open the instrument.
- Using cables that Leica has not provided or permitted.
- Using combinations with non-Leica components that go beyond the scope of the manual.
- Coupling a laser into not special or unapproved microscopes.

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from using the microscopes for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In such cases the declaration of conformity shall cease to be valid.

Caution!

These (IVD) instruments are not intended for use in the patient environment defined by DIN VDE 0100-710. Nor are they designed to be combined with medical instruments in accordance with EN 60601-1. If a microscope is electrically connected to a medical instrument in accordance with EN 60601-1, the requirements defined in EN 60601-1-1 shall apply.

Not suitable for examining potentially infectious specimens.

Only trained personnel may operate this type of device.

Notes on handling laser devices

In the standard design without additional laser safety precautions, the microscopes are not suitable for coupling laser radiation (such as to the camera ports), since this radiation poses a hazard for the user (eye injury in particular).

For use of the microscope with lasers, Leica Microsystems offers special microscope variants with additional safety features. Laser couplings require corresponding safety devices that have to be inspected and installed by trained personnel.

For further information, please contact your authorized Leica Microsystems representative.

3. Safety notes

3.1 General safety notes

These devices of protection class 1 are built and inspected in accordance with EN 61010-1 / IEC 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use and EN 61010-2-101 / IEC 61010-2-101 Safety requirements for electrical equipment for measurement, control and laboratory use, Part 2 Particular requirements - In vitro diagnostic (IVD) medical equipment.

They also meet the requirements in EN / IEC 61326-1 Electrical equipment for measurement, control and laboratory use – EMC requirements and EN / IEC 61326-2-6, Part 2-6 Electrical equipment for measurement, control and laboratory use – EMC requirements – Particular requirements – In vitro diagnostic (IVD) medical equipment.

The products were also tested in accordance with EN 62471 / IEC 62471 Photobiological safety of lamps and lamp systems, and classified in risk group 1 (low risk).

Caution!

The devices and accessories described in this operating manual have been tested for safety and potential hazards.

The responsible Leica affiliate or the main plant in Wetzlar, Germany, must be consulted whenever the device is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the device or noncompliant use shall void all rights to any warranty claims and void product liability!

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this operating manual.

3.2 Electrical safety

General specifications

Microscopes Leica DM4 B/M/P

For indoor use only.

Supply voltage: 100-240 V AC

Frequency: 50/60 Hz

Power input: max. 80 VA

Fuses: 1.6 A, slow-blowing,
breaking capacity H,
250 VAC, size 5x20 mm

Ambient temperature: 15-35°C

Relative humidity: max. 80% to 30°C
non-condensing

Pollution degree: 2

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the device, or release of the ground wire connection, can cause the device to become hazardous. Intentional ground interruption is not permitted!

Caution!

Through connection to the grounding connection, ancillary equipment with its own and/or extra power supply may be brought to the same ground wire potential. For connections without a ground connector, Leica Service must be consulted.

Never use any fuses as replacements other than those of the types and the current ratings listed here. Bypassing fuse holders is not permitted.

The use of incorrect fuses may result in a fire hazard.

Caution!

The microscope's electrical accessory components are not protected against water. Water can cause electric shock.

3. Safety notes

Caution!

Protect the microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.

Operating temperature: 15-35°C

Caution!

Before exchanging the fuses or lamps, be absolutely certain to switch off the main power switch and remove the power cable.

WARNING

Use only original power cables or alternative cables with a VDE / HAR logo, which at least fulfill the requirement of 3x0.75 mm² and 10A/250V.

WARNING

Per definition, the connection between power supply cable and appliance inlet is the power separator of this device. Thus, the user is responsible for sufficient clearance and accessibility of this area at all times.

Caution!

Do not use the microscope in altitudes exceeding 2000 m ASL/NL.

Caution!

Only insert or unplug data and control circuits when the instrument is switched off; otherwise, the instrument may be damaged.

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

3.3 Laser coupling by the user

Because of a danger for the user due to laser radiation, the Leica DM4 microscopes in the standard version are not suitable without additional laser safety measures for coupling laser radiation¹⁾.

If you are, however, modifying a microscope and couple a laser beam into it, as the user, you assume responsibility as the manufacturer of the complete system.

We do not make provisions for your interference with the device and, therefore, you are fully responsible for the modifications regarding the adherence to basic safety requirements and laser protective measures set forth by the European Union guidelines and by the national legal and technical directives and standards.

If you modify your microscope by coupling a laser beam without authorization, when determining and implementing laser safety measures, also take note of the subsequently described hazards at the microscope, as here the danger exists of serious and irreversible eye damage and skin injuries!

Danger of serious and irreversible eye damage from laser radiation!
Radiation of the eye or skin from direct or indirect laser radiation has to be avoided by all means! The laser light can cause serious eye damage and skin damage!

1) Laser radiation of Class 1 lasers is exempt. Further details regarding this topic in Section 3.3.6

3. Safety notes

3.3.1 Safety hazard at eyepiece

The laser radiation that you coupled into the microscope can disperse from the eyepiece and directly hit the user's eye!

3.3.2 Safety hazard in the specimen area

The laser radiation is accessible without obstruction in the specimen area of the microscope and below the specimen stage!

Specimen area


3.3.3 Safety hazard due to openings at microscope

Make sure, as part of your risk assessment, that all openings at the microscope are tightly sealed and can only be opened with tools to prevent laser radiation from escaping!

These openings include: the camera port on the tube, connection for transmitted light lamp housing, connection of incident light lamp housing/ connection for adapter for Leica EL6000, filter block holder and centering window; seal the filter blocks and centering window light-tight.

Furthermore, the terminals that are not occupied by lenses have to be closed with caps at the objective nosepiece. You also have to cover the front cover for fluorescence and the 3 side openings for various slides.

The material that is used to seal off the openings has to be suitable for the intensity of the coupled laser beam.

- a Transmitted light lamp housing terminal
- b Incident light lamp housing terminal / EL6000 adapter
- c Camera port on the tube
- d Filter block holder
- e Centering window

3.3.4 Notes regarding the safety hazards in Section 3.3.1

The laser radiation coupled into the microscope can escape from the eyepieces.

When using one of the two special tubes (with integrated interlock switches) from Section 3.3.5, the safety hazard on the eyepiece can be eliminated.

The principle of the safety function is that either the laser radiation coupled into the microscope is switched on or the optical axis at the eyepiece is cleared.

Important safety note

In addition, observe the safety hazards described in Section 3.3.2 and 3.3.3! The following tubes eliminate the safety hazard from Section 3.3.1 only.

The following tubes have a built-in interlock switch that you have to connect with a laser beam interrupter or with the laser.

Note that all ports on the tubes must be sealed off by a component (such as a camera) or using a suitable cover.

The safety function is implemented by the tubes as follows:

- When the optical axis between the specimen (object) plane and the eyepiece is released by the selector lever on the tube, the interlock switch opens. The laser radiation coupled into the microscope must then be switched off.
- When the optical axis between the specimen (object) plane and the eyepiece is interrupted by the selector lever on the tube, the interlock switch closes. The laser radiation coupled into the microscope is switched on.

3.3.5 Tubes with integrated interlock switches

Important safety note

The overall responsibility for laser safety of the laser microscope system assembled by you remains with you as the responsible body, even if you use the Leica components listed below in conjunction with the system. Within the context of your system composition and modifications, you have to verify whether the components selected by you from this chapter are effective and suitable for your laser protection concept.

3. Safety notes

When connecting your laser beam interrupter/laser to the interlock connector on the tube, make sure not to exceed the electrical characteristic values of the internal interlock switch, as otherwise the interlock switch will be damaged and no longer function safely and reliably:

Maximum permitted electrical characteristic values of the integrated interlock switches:

$U=30\text{ V DC}$


$I=0.1\text{ A}$

The integrated interlock switch in the microscope and its wiring with a laser beam interrupter/laser ensure that the laser radiation you coupled into the microscope cannot get into the user's eyes through the eyepieces.

Tubes with interlock connector:

Variant 1:

Quasar dual tube output (Leica No. 11888339)

Options for coupling in/out:

- Laser coupling and documentation together on the rear tube port (e.g. Spinning Disk).
A camera can be used on the upper tube port. Otherwise, this port must be sealed off.


or

- Laser coupling into the following openings:
 - Fluorescence axis on the incident light lamp housing terminal
 - Fluorescence axis on the centering window
 - On the transmitted light lamp housing terminal,
connect the camera to the rear tube port.
Otherwise, if no camera is connected, this port must be sealed off.
A camera can be used on the upper tube port. Otherwise, this port must be sealed off.

Variant 2:

SD documentation tube MBDT 25+V100/50/0, mot.


(Leica No. 11888451) in combination with mono attachment (Leica No. 11505161).

Variant 3:

SD documentation tube MBDT 25+V100/50/0, mot.

(Leica No. 11888451) in combination with manual documentation port with warning sign (Leica 11888421)

Options for coupling in/out:

- Laser coupling and documentation together on the tube port (e.g. Spinning Disk)

or

- Laser coupling into the following openings:
 - Fluorescence axis on the incident light lamp housing terminal
 - Fluorescence axis on the centering window
 - On the transmitted light lamp housing terminal, connect the camera to the tube port. Otherwise, if no camera is connected, this port must be sealed off.

Options for coupling in/out:

- Laser coupling and documentation on the rear tube port (Spinning Disk). A camera can be used on the upper tube port. Otherwise, this port must be sealed off.

or

- Laser coupling into the following openings:
 - Fluorescence axis on the incident light lamp housing terminal
 - Fluorescence axis on the centering window
 - On the transmitted light lamp housing terminal, connect the camera to the rear tube port. Otherwise, this port must be sealed off. Otherwise, if no camera is connected, this port must be sealed off. A camera can be used on the upper tube port. Otherwise, this port must be sealed off.

3. Safety notes

3.3.6. Exception for Class 1 lasers

Coupling a laser of Class 1 in accordance with IEC/EN 60825-1 into a microscope does not result in any danger to the user.

The listed safety hazards apply only for lasers of Class 1M, 2, 2M, 3R, 3B and 4.

Note that the rules and regulations of the United States FDA 21CFR1040.10 Class I (CFR – Code of Federal Regulation Title 21, revised as of April 1, 2011) do not distinguish between Class 1 and 1M. Thus the statement that there are no safety hazards on the microscope do not apply for the laser classification in accordance with FDA Class 1.

3.4 Transport and storage

Fully equipped, the microscope weighs more than 18 kg. For transportation, the user has to take care of the corresponding actions. See also p. 21.

Caution!

Transport and storage in a range of $-25^{\circ} - +70^{\circ}\text{C}$ and at a humidity not exceeding 80% (non-condensing).

3.5 Notes on Handling Light Sources

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation). Therefore, lamps have to be operated in closed housings and only after being mounted. Never look directly into the beam path (blinding hazard).

Make sure to follow the safety notes in the operating manual for the light source.

Lamps and lamp housing may be hot! They must be placed at least 10 cm away from the wall and away from flammable substances. Particularly, supply and data lines must not encounter lamp housings!

3.6 Notes on Handling Immersion Oil

Caution!

When using immersion oil, take care to avoid skin contact! Ask the supplier for a safety data sheet!

3.7 Notes on Handling Acids and Bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

Caution!

Be absolutely certain to avoid direct contact with these chemicals.

3.8 Disposal

Once the product has reached the end of its service life, please contact Leica Service or Sales about disposal.

Please observe and comply with the national and federal laws and regulations that are equivalent to EU guidelines such as WEEE.

Note!

Like all electronic devices, the microscope, its accessory components and consumables must never be disposed of with general household waste.

4. Overview of the instrument

Specification	Leica DM4 B	Leica DM4 M Leica DM4 P
Contrast methods	- Transmitted light: BF, DF, PH, Pol - Incident light: Fluorescent	- Transmitted light: DM4 M: BF, DF, PH, ICT, Pol DM4 P: BF, DF, PH, ICT, Pol (conoscopy) - Incident light: DM4 M: BF, DF, ICR, Pol, Fluo DM4 P: BF, ICR, Pol, Fluo
Transmitted light axis	- Automatic illumination manager (mot. aperture diaphragm and field diaphragm, mot. intensity control) - Constant colour temperature with LED - Motorized shutter	
Incident light axis	- Integrated into the stand - Motorized 5x filter turret disk - With FIM (fluorescence intensity manager) for reducing the light intensity in 5 increments - Mechanical booster lens for increasing fluorescence intensity - Motorized shutter	- Integrated into the stand - Motorized 4x filter turret disk - Automatic illumination manager - DM4 M: motorized shutter
Z pinion	Manual	
Objective turret	- Manual, fully encoded 6x/7x with M25 thread	- Manual, fully encoded - DM4 M: 6x with M32 thread - DM4 P: 6x with M25 thread, centerable, encoded - Receptacle for DIC prisms and Pol compensators (for DM4 M: optional)

Specification	Leica DM4 B	Leica DM4 M Leica DM4 P
X/Y stage	• Manual • Replaceable specimen stage • Coaxial drive length: 155 mm	• Manual • DM4 M: • Replaceable specimen stage • Coaxial drive length: 140 mm - DM4 P: • Replaceable Pol stage
Tube	• Manual or motorized (DM4 P: manual) • Optionally with one or two camera outputs • DM4 P: conoscopy module (tube optics HC P1x/1.6x with Bertrand lens, encoded)	
Condenser	• Motorized condenser head • Condenser disk for the light ring, DF-Stop, DIC prisms • Automatic Köhler illumination • Optional polarizer (integrated and motorized)	
Magnification changer (optional)	• Manual • 3x fully encoded • 1x; 1.25x; 1.6x	• Manual • 3x fully encoded • 1x; 1.5x; 2x
Controls	• Operating buttons on the stand for all motorized microscope functions • Additional variable multifunction keys • Focus wheels • Display	
Computer interface	• USB2.0	
Software tools	• Leica Application Suite (LAS or LAS X) for Windows™ 7 64-bit, Windows™ 7 • With plug-ins for: • Microscope and camera configuration • Microscope and camera control • Image acquisition	

Note:

Leica DM4 B/M and Leica DM4 P microscopes are designed for the LH111 LED lamp housing. Connection of lamp housings, other than LED type, are not supported by microscope outlets.

4. Overview of the instrument


Fig. 1 Left side of the stand with Ergonomic Documentation Tube

- | | |
|---------------------------------------|---|
| 1 Eyepiece | 8 Field diaphragm operating buttons |
| 2 Eyepiece tube | 9 Transmitted light / incident light switch |
| 3 Tube | 10 Aperture diaphragm operating buttons |
| 4 Objective turret with objectives | 11 Brightness adjustment buttons |
| 5 Specimen stage with specimen holder | 12 Focus wheel with coarse and fine adjustment |
| 6 Condenser | 13 Variable function keys (preset at the factory) |
| 7 Display | 14 Lamp adjustment window |


Fig. 2 Right side of the stand with Documentation Tube

- 15** Adapter for Leica EL6000 (incident light fluorescence)
- 16** Lamp housing for transmitted light
- 17** Transmitted light filter, optional
- 18** Transmitted light filter, optional
- 19** Variable function keys (preset at the factory)
- 20** x/y coaxial drive, adjustable height
- 21** Handwheel for fine focus
- 22** Motorized filter block exchanger

5. Unpacking the microscope

The microscope is delivered in two packages.

The **stand package** contains the following components:

- Stand with integrated incident light axis and objective turret
- Specimen stage with stage bracket
- Power cable and PC connecting cable
- CD with Leica Application Suite (LAS oder LAS X) software package
- Instructions and list of microscope default settings

The **system package** contains the microscope's accessories:

- Tube
- Eyepieces
- Objectives
- Condenser
- Lamp housings with accessories
- Assembly tools
- Additional microscope accessories such as filter cubes, etc. depending on product configuration

The external ebq 100 supply unit* is delivered in a separate package.

First, carefully remove all components from the transportation and packaging materials.

Note:

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces of optical devices in a short time. See the chapter on „Care of the microscope“ → p. 87 for additional instructions.

Caution!

Do not connect the microscope or peripherals to an AC power source at this time under any circumstances!

Installation location

Work with the microscope should be performed in a dust-free room, which is free of oil vapors and other chemical vapors, as well as extreme humidity. At the workstation, large temperature fluctuations, direct sunlight and vibrations should be avoided. These conditions can distort measurements and micrographic images.

Allowable ambient conditions

Temperature 15-35°C

Relative humidity maximum 80% up to 30°C
non-condensing

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

See the chapter on „Care of the microscope“ → p. 87 for additional instructions.

Caution!

Electrical components must be assembled at least 10 cm away from the wall and from flammable substances.

Caution!

When installing the microscope, make sure the power inlet is freely accessible so that the instrument can be quickly disconnected from the mains if necessary.

Transport

For shipping or transporting the microscope and its accessory components, the original packaging should be used.

As a precaution to prevent damage from vibrations, the following components should be disassembled and packaged separately:

- Unscrew the objectives.
- Remove the condenser.
- Remove the specimen stage.
- Remove the lamp housings.
- Disassemble the burner of 106 z lamp housing.
- Remove all moving or loose parts.

Weight

Fully equipped, the microscope weighs more than 18 kg. For transportation, the user has to take care of the corresponding actions.


Caution!

For transporting it is essential to remove all components listed under “Transport”.

5. Unpacking the microscope

Dimensions (specifications in mm)

Weight, as illustrated: 16 kg

6. Assembling the microscope

The microscope components are logically assembled in this order:

- Specimen stage
- Condenser with condenser head
- Tube
- Eyepieces
- Objectives
- Lamp housings with light sources
- Equipment for the incident light turret disk*

Only a few commonly used screwdrivers and keys are necessary for assembly; these are included in the delivery package.

When using intermediate systems and optical accessories, the sequence may vary. In this case, read chapter, „6.10 Optional accessories“ → p. 43

6.1 Specimen stage


Caution:

Never install objectives before assembling the stage.

- Place the specimen holder on the stage and fasten it with the two screws (3.1).
- Using the condenser height adjuster (3.2), turn the condenser holder completely upwards, i.e. as close to the stage as possible.
- Loosen the stage clamp (3.3) slightly.

Fig. 3 Mechanical object stage

- 1 Locking screws for specimen holder
- 2 Condenser height adjuster
- 3 Stage clamp

6. Assembly

- From above, set the stage clamp onto the dovetail guide (4.2) and push the stage downwards until the upper end of the dovetail guide is tightly fastened to the upper end of the stage clamp.
- Firmly tighten the stage clamp (4.1.).

Note:

For thicker specimens (Leica DM4 M) the stage can be set to a correspondingly lower level.

- Only for DM4 P:

Pol attachable mechanical stage*

Adjust the attachable mechanical stage so that the fastening screw is visible below the holes (4a.1). Set the attachable mechanical stage in the guide holes on the rotating stage and tighten the fastening screw using the hexagonal key.

Attachable mechanical stage*

The attachable mechanical stage can be installed on the left, on the right or on the front (not pictured). The two clamping screws fasten it into place.


Fig. 4 Assembling the stage

- 1 Stage clamp
- 2 Dovetail guide

Fig. 4a Pol rotating stage* and Pol 3 attachable mechanical stage*

- 1 Holes for the fastening screw.
- 2 Lever for the holder for glass slides of various formats, which can be turned inward and outward
- 3 Storage for the centering key
- 4 Locking button pair
- 5 45° click stop
- 6 Clamping system for the stage rotation

6.2 Condenser

- Screw the condenser head into the condenser.
- Using the condenser height adjuster (5.4), turn the condenser holder (5.1) downward as far as it will go.
- Unscrew the clamping screw for the condenser (5.3) far enough so that the condenser can be inserted from the front.
- From the front, insert the condenser into the condenser holder as far as it will go. On the underside of the condenser, there is an orientation pin (6.1) that must be locked into place in the guiding notch (7.1).
- Tighten the condenser's clamping screw (5.3) until the condenser locks into place.
- Connect the condenser over the connection (8.1) with the stand.

Note:

The condenser must be centered before using the microscope.

→ Köhler illumination p. 51.

Fig. 6

Underside of condenser

- 1 Orientation pin


Fig. 7 Condenser holder

- 1 Guiding groove


Fig. 5 Condenser holder

- 1 Condenser holder
2 Condenser centering
3 Clamping screw for the condenser
4 Condenser height adjuster


Fig. 8 Condenser connection

- 1 Condenser cable socket

6. Assembly

6.3 Tube and eyepieces

The tube is mounted on the stand either directly or with the use of intermediate modules. The side clamping screw fastens it into place (9b.1).

- For the MBDT motorized tube only:
Remove the transport anchor (9a.1) on the underside of the tube.
- Partially unscrew the clamping screw (9b.1).
- Insert the tube into the circular receptacle (dovetail ring).
- Retighten the clamping screw (9b.1).
- For the MBDT motorized tube only:
Connect the tube to the stand with the connector bushing (10.1).

- The eyepieces are inserted into the eyepiece tubes on the tube.

- For the BDTP tube only:

The Pol eyepieces are inserted into the eyepiece tubes (using the locking groove).

Fig. 9b Fastening the tube
1 Clamping screw

Fig. 9a Underside of the tube
1 Transport anchor

Fig. 10 Motorized tube connection
1 Connector socket

6.4 Objectives

The receptacles on the objective turrets are numbered (Fig. 11). Based on your equipment, the individual objectives have already been assigned to specific positions at the factory. For details on the exact positions of the objectives, please refer to the enclosed identification sheet.

Caution:

Close vacant threads in the turret with dust protection caps!

6.5 Light sources for DM4 B/M and DM4 P

Caution!

Ensure that the lamp housing has been disconnected from the power supply. Unplug the power plug and the power supply during assembly.

CAUTION

Light sources pose a potential irradiation risk (glare, UV-radiation, IR-radiation). Therefore, lamps have to be operated in closed housings.

CAUTION

Switch on light sources only if they are firmly connected to the microscope. Uncontrolled light output poses a blinding hazard!

Fig. 11

Objective turret with engraved objective receptacles

Fig. 12 Lamp housing LH111 LED
Lifetime: up to 50.000 hrs.

6. Assembly

For DM4 B/M and DM4 P microscopes the lamp housing LH111 LED (Fig. 12) is used for the transmitted light axis. For DM4 M and DM4 P this lamp housing is used for the incident light axis as well. For the incident light fluorescence axis DM4 B see chapter 6.6.

- Place the lamp housing in the transmitted light lamp housing receptacle (13.2) or incident light lamp housing receptacle (13.1) and fasten it with the clamping screw on the side.
- Connect the lamp housing to the power supply for transmitted light (13.3) (symbol: $\triangle$) or incident light (13.4) (symbol: ∇).

Note:

A separate lamp power supply is necessary to connect other lamp housings.

6.6 Light sources for incident light fluorescence axis DM4 B

Caution!

Ensure that the lamp housing has been disconnected from the power supply. Unplug the power plug and the power supply during assembly.

CAUTION

Light sources pose a potential irradiation risk (glare, UV-radiation, IR-radiation). Therefore, lamps have to be operated in closed housings.


CAUTION

Switch on light sources only if they are firmly connected to the microscope. Uncontrolled light output poses a blinding hazard!

Fig. 13 Rear view of the DM4 stand

- 1 Incident light lamp housing receptacle
- 2 Transmitted light lamp housing receptacle
- 3 12V 15 W connection for transmitted LED illumination (symbol: $\triangle$)
- 4 12V 15 W connection for incident LED illumination (symbol: ∇)
- 5 I²C connection
- 6 USB connection, i.e. for PC

**WARNING**

During assembly work on xenon burners, always wear the protective gloves and face protection supplied (Fig. 15) (risk of explosion).

Never touch the glass parts of the burner with bare hands.

Never look directly into the beam path (blinding hazard).

Lamp housing 106/106 z

This lamp housing is suitable for use with a 12 V 100 W halogen lamp or a variety of gas discharge lamps.

**CAUTION**

Make sure to follow the instructions and safety notes of the lamp supplier.

Before changing lamps allow it to cool down for at least 30 min.!


**CAUTION**

The lamp may still be hot.

Fig. 16 106/106 z lamp housing (on the side, open)

- 1 Cover raised
- 2 Collector
- 3 12 V 100 W lamp or gas discharge lamp in mount
- 4 Reflector (mirror)
- 5, 6, 7 Adjusting screw for x-y reflector
- 8 Fastening screw for the lamp mount
- 9 Socket for contact plug

Fig. 15 Protective gloves and mask

6. Assembly

Inserting the 12 V 100W halogen bulb into the 106/106 z lamp housing

- Unscrew the fastening screws of the cover and flip the cover up (16.1).
- Unscrew the fastening screws of the lamp mount (16.8) and pull out the mount (Fig. 17).
- Insert the lamp with the dust cover straight into the socket until it stops.

Caution!

Do not remove the lamp's dust cover until after you have installed the lamp. Avoid fingerprints on the lamp.

- Remove the dust cover.

- Insert the lamp mount, with the burner installed, into the lamp housing and tighten it with the screws (16.8).
- Close the lamp housing and retighten the fastening screws.
- Place the lamp housing in the incident light lamp housing receptacle (18.1) and fasten it with the clamping screw on the side.
- Connect the lamp housing to the power supply for incident light (symbol: ▽) (18.4).
- Now connect the lamp power cable of the microscope (18.5) to the **12V 100W** port of the Leica CTR5000 electronics box (See 6.12 Connection to the CTR5000 electronics box).

Fig. 17 Lamp mount with 12 V 100 W halogen bulb



Fig. 18 Rear view of the DM4 stand
1 Incident light lamp housing receptacle

Inserting gas discharge lamps (Hg and Xe) into the 106/106z lamp housing

Hg and Xe lamps are powered by separate supply units.

Please also read the separate instruction manual provided with these supply units.

The following gas discharge lamps may be used and require different power supplies and lamp mounts (Fig. 19):

Type	Typical bulb life ⁺⁾
50 W high-pressure mercury burner (alternating current)	100 hrs.
100 W high-pressure mercury burner (direct current)	200 hrs.
100 W high-pressure mercury burner, type 103 W/2 (direct current)	300 hrs.
75 W high-pressure xenon burner (direct current)	400 hrs.

⁺⁾ Please observe the data sheets of the lamp manufacturer.

6. Assembly

- To open the 106 z lamp housing, unscrew the fastening screws on the cover.
- Remove the transport anchor (red plastic rod in place of the burner) in the lamp mount. To do so, remove the lower clamp (19.1). Pull up the cooling element (19.3) and turn it to the side. Detach the lower clamp system (19.2) and remove the transport anchor.
- Install the burner in reverse order.

Caution!

Hg 50 burner:




After installation, the labeling must be up-right.

If a glass melt nipple is present (19a.4), position it by turning the burner so that the nipple does not impede the beam path later, but instead is positioned sideways.

Xe 75 burner:

Remove the burner's dust cover (19b.5) after you have installed the burner.

Fig. 19 a- Lamp mounts for gas discharge lamps
1 Upper clamping system, **2** Lower clamping system, **3** Cooling element
4 Melt nipple for the Hg 50 arc lamp, **5** Dust cover for the Xe 75 arc lamp

- Insert the lamp mount, with the burner installed, into the lamp housing and tighten it with the screws (20.8).
- Close the lamp housing and retighten the fastening screws.
- Place the lamp housing in the incident light lamp housing receptacle (21.1) and fasten it with the clamping screw on the side.
- Connect the lamp housing to the external power supply (22.1).

Fig. 20 106/106 z lamp housing (on the side, open)

- 1 Cover raised
- 2 Collector
- 3 12 V 100 W lamp or gas discharge lamp in mount
- 4 Reflector (mirror)
- 5, 6, 7 Adjusting screw for x-y reflector
- 8 Fastening screw for lamp mount
- 9 Socket for contact plug



Fig. 21 Rear view of the stand

- 1 Incident light lamp housing receptacle
- 2 Transmitted light lamp housing receptacle

Fig. 22 Rear panel of the ebq 100 supply unit

- 1 Lamp connection

6. Assembly

6.7 Equipping the incident light turret disk

The positions in the turret disk are numbered. Depending on how they are equipped, the individual filter and/or reflector cubes are set in pre-assigned positions at the factory. For details, check the identification sheet included with your order.

Insert the filter and reflector cubes in the following manner:

- Never fit the incident light turret disk while the microscope is in operation.
- Remove the face plate from the upper part of the microscope (Fig. 25). Press the locking pin (25.2) to turn the turret disk. When the locking pin is released, the turret disk locks into place again.
- With the holder facing you squarely, insert the filter cube or reflector cube into the holder as described in the identification sheet provided. To do so, place the filter or reflector cube on the **right** side and press it toward the **left** into the mounting (Fig. 26).
- Press the locking pin (25.2) and turn the filter turret to the next click stop.
- Make sure that the turret engages (the locking pin springs forward) and insert the next filter and/or reflector cube as described above.
- When all filters and reflector cubes have been inserted, close the front cover plate again.

Fig. 23 Filter cube, front side

Fig. 24 Filter cube, back side

Fig. 25 Removing the front panel

- 1 Filter receptacle
- 2 Locking pin
- 3 Front panel


Fig. 26 Inserting the filter or reflector cubes

- 1 Mounting

6.8 Polarizer and analyzer

Transmitted light polarizer: ICT/P

- Using the left clamping screw, fasten the ICT/P transmitted light polarizer to the underside of the condenser holder (Fig. 27).
- Make sure that the red index point on the front of the polarizer is at 0.
- If necessary, insert the compensators (λ -, $\lambda/4$ plates or elliptical compensators) into the polarizer's receptacle (Fig. 28) or into the compensation slot.

Fig. 27 Assembly of the ICT/P transmitted light polarizer
1 Clamping screw


Fig. 28 Inserting the compensators

Incident light polarizers:

R/P polarizer, rotating polarizer,

L/ICR polarizer, R/ICR polarizer

- Remove the plug cap on the right side of the incident light axis (Fig. 29).
- Insert the polarizer into the receptacle until it locks into place.


Caution:

Push the polarizer into the front receptacle only.

Motorized polarizer

- A motorized polarizer is already installed and ready for operation in the DIC condenser.

Fig. 29 Inserting the polarizer
1 The plug cap is replaced with the polarizer.

6. Assembly

Transmitted light and incident light analyzer

- Remove the plug cap on the left side of the stand.
- Insert the analyzer into the receptacle until it latches in place (Fig. 30).

Motorized analyzer

- Insert the analyzer cube into the corresponding position on the filter turret as described in section „6.7 Equipping the incident light turret disk“ → p. 40. For details on the proper position for the analyzer cube, please refer to the identification sheet provided.

6.9 DIC prisms

- Insert the objective prism slide into the tube slot (31.1). The code letter must match the code letter on the objective.

Fig. 30 Inserting the analyzer

1 The plug cap is replaced with the analyzer.



Fig. 31 Inserting the objective prism slide

1 Objective prism slide

6.10 Optional accessories

ErgoModule

For raising the eye level of the tube opening, the ErgoModule may be used.

It is fastened in place with the side clamping screw.

Mirror housing

- Place the mirror housing directly onto the lamp housing receptacle on the back of the stand and attach it using the side clamping screw.
- Place the lamp housing(s) on the mirror housing and fasten it using the accompanying lateral clamping screw.

Booster lens

- Insert the filter slide into the front receptacle on the right side of the stand (Fig. 32.1).

Fig. 32

1 Inserting the booster lens

6. Assembly

6.11 Connection to the power supply

Leica DM4 B/M and DM4 P only:

After completing the assembly work, connect the stand to the power supply using the power cable (Fig. 34.2).

Fig. 34 Rear side of the Leica DM4 B/M stand

- 1 Power switch
- 2 Power supply

7. Startup

7.1 Functional principle

The microscope's most important functions may be easily accessed using function keys.

- The microscope may be switched between various contrast methods by pressing a button.
- The microscope recognizes the selected objective and associated contrast method. Therefore, the values for intensity (INT), aperture diaphragm (AP) and field diaphragm (FD) are always set correctly.
- The INT, AP and FD values can be adjusted individually. Manual adjustments overwrite the previous settings. The current setting is saved.
- The INT, AP and FD values are always based on the currently activated light axis (transmitted light or incident light).
- In addition to the preset function keys for INT, AP and FD, there are also variable function keys.

Variable function keys:

- At the time of delivery, these function keys are assigned functions suitable to the configuration of your microscope.
- However, the functions can be reprogrammed and/or adapted to your specific requirements.

Note: (reset function)

The microscope can be reset to its factory default programming:

- When the microscope is switched off, press all 3 variable function keys on the left stand section.
- Switch on the power for the stand.
- Hold the buttons until the initialization is complete.
- The standard indicator appears in the display.
- Switch the instrument off and back on. The settings are now saved.

7. Startup

Possible assignments for the function keys

For Leica DM4 B:

Function button	Function
BF	Bright field transmitted light
PH	Phase contrast transmitted light
ICT	Interference contrast, transmitted light
DF	Dark field transmitted light
POL	Polarization transmitted light
CHANGE TL 	Switches through all transmitted light contrast methods
INT ↑	Increases the brightness (transmitted light)
INT ↓	Reduces the brightness (transmitted light)
AP ↑	Opens the aperture diaphragm (transmitted light)
AP ↓	Closes the aperture diaphragm (transmitted light)
FD ↑	Opens the field diaphragm (transmitted light)
FD ↓	Closes the field diaphragm (transmitted light)
SHUTTER TL	Opens/closes the transmitted light shutter
FLUO	Fluorescence (last filter cube)
CUBE 1	Chooses the fluorescence cube on position 1
CHANGE CUBE 	Switches through the fluorescence cube clockwise (1 → 4)
CHANGE CUBE 	Switches through the fluorescence cube counterclockwise (4 → 1)
SHUTTER FLUO	Opens/closes the fluorescence shutter
INT FLUO ↑	Increases the brightness (fluorescence)
INT FLUO ↓	Reduces the brightness (fluorescence)
FD FLUO ↑	Opens the field diaphragm (fluorescence)
FD FLUO ↓	Closes the field diaphragm (fluorescence)
COMBI 	Combination method (PH fluorescence or ICT fluorescence)
CHANGE COMBI 	Switches through all combination methods
CHANGE TUBE	Switches through various types of beamsplitting
100% VIS	100% documentation port
50:50	50% documentation port / 50% camera
100% CAMERA	100% camera

For Leica DM4 M/DM4 P:

Function button	Function
BF	Bright field incident light
ICR	Interference contrast, incident light
DF	Dark field incident light
POL	Polarization incident light
CHANGE RL 	Switches through all incident light contrast methods
INT ↑	Increases the brightness (incident light)
INT ↓	Reduces the brightness (incident light)
AP ↑	Opens the aperture diaphragm (incident light)
AP ↓	Closes the aperture diaphragm (incident light)
FD ↑	Opens the field diaphragm (incident light)
FD ↓	Closes the field diaphragm (incident light)
SHUTTER RL	DM4 M only: Opens/closes the incident light shutter
FLUO	Fluorescence (last filter cube)
CUBE 1	Chooses the fluorescence cube on position 1
CHANGE FLUO	Switches through all filter cubes
CONOS	DM4 P only: The Bertrand lens is in the beam path (conoscopic beam path)
FOCUS FINDER	Finds the smallest incident light field diaphragm and toggles back to the original field diaphragm.
BF TL	Bright field transmitted light
BF-POL	DM4 P only: Polarization transmitted light (conoscopy)
INT ↑	Increases the brightness (transmitted light)
INT ↓	Reduces the brightness (transmitted light)
AP ↑	Opens the aperture diaphragm (transmitted light)
AP ↓	Closes the aperture diaphragm (transmitted light)
FD ↑	Opens the field diaphragm (transmitted light)
FD ↓	Closes the field diaphragm (transmitted light)
COMBI 	Combination methods (BF and BF TL)
CHANGE TUBE	Switches through various types of beamsplitting
100% VIS	100% documentation port
50:50	50% documentation port / 50% camera
100% CAMERA	100% camera

7. Startup

7.2 Switching on the unit

- First, swivel the objective with the lowest magnification into position.
- Switch on the microscope or CTR5000. All motorized microscope components will then run through an initialization phase.

After the initialization is complete, the display on the stand (Fig. 37) shows the current microscope setting

Components such as diaphragms, condensers, light and phase rings have been pre-centered at the factory. It may be necessary to correct the centering after the microscope has been transported and assembled.

Before proceeding with the necessary steps, first familiarize yourself with the stand's display and control panel.

Caution!

After turning on the gas discharge lamp, the burner must be immediately adjusted. Therefore, **do not** turn on the power supply unit yet. First, work in transmitted light in order to familiarize yourself with the microscope's controls.

Note:

The Leica DM4 P microscope is equipped at the factory with an encoded conoscopy module.

- Moving the Bertrand lens activates the conoscopic or orthoscopic beam path. The status is indicated on the display with the methods „Conos“ or „TL Pol“.

Fig. 37 Display after initialization

	TL_BF		
	10x Obj.	1.25x MagCh.	Σ 125x
	INT 100%	AP 33	FD 30
	50%	50%	

7.3 The display (Leica DM4 B/DM4 M/DM4 P)

The screen displays the microscope's current settings. The content of the display depends on the features of the individual microscope. In the first column, corresponding symbols indicate the type of information: contrast method, magnification, light intensity, diaphragms, light splitting for photo tubes.

Please see the abbreviation index for a list of abbreviations used, → p. 45f.

Contrast method

In the first row, you find an indication of the active light axis (incident light or transmitted light) of the current contrast method and, if used, the current filter cube.

The shutter status display for the transmitted light or incident light shutter:

- ↑ Transmitted light shutter is open
- ↑ Transmitted light shutter closed
- ↓ Incident light shutter is open (for fluorescence only)
- ↓ Incident light shutter is closed (for fluorescence only)

Magnification

The current objective magnification (OBJ), sometimes followed by the re-magnification of the magnification changer*, appears along with the total magnification:

$\Sigma = \text{Objective} \times \text{re-magnification} \times \text{eyepiece}$

Light intensity

The current brightness setting is graphically depicted by a beam. In addition, the light intensity is specified in 20 (coarse) or 255 (fine) increments → p. 69.

Diaphragm

The values for the field diaphragm (FD) and the aperture diaphragm (AP) are indicated numerically. The field diaphragm in the incident light may be either round or rectangular. Accordingly, the FD designation is set in parentheses or in brackets: (FD) or [FD].

Note:

When using a digital camera, angular field diaphragms are recommended.

Light distribution

If a motorized tube is present, the light division between ocular (Eye) and photo port (Docu) is indicated in %.

Note:

The display may flash after the initialization phase or even during microscopy. This always occurs when the contrast method selected can not be performed with the microscopic settings. For example, an objective may be rotated in that is not suited to the contrast method chosen. Then check your settings.

7. Startup

7.4 The function keys

There is a row of function keys both on the right and left side of the stand. These can be broken down into fixed and variable buttons. The variable function keys have different functions depending on the features of the individual microscope.

Fixed function keys on the left side of the stand

The **TL/IL** button (38.1) toggles between the incident light and transmitted light axis. The contrast method last used with a given axis is restored when switching.

The **INT** buttons (38.3) adjust the light intensity. The adjustment can be made in coarse or fine steps. Pressing both **INT** buttons at the same time toggles between coarse and fine adjustment. The display indicator changes accordingly.

The **AP** buttons (38.4) for the aperture diaphragm and **FD** (38.2) for the field diaphragm open and close their respective diaphragms.

Note:

Changes to the light intensity as well as aperture and field diaphragm settings are stored for the individual objectives and contrast methods.

Variable function buttons

The variable function buttons are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on button assignments, please refer to the included identification sheet. For information on the abbreviations used, please refer to the list → p. 45f.


Note:

The Leica Application Suite (LAS) module “Set-Up” or LAS X Hardware Configurator is required for changing the key assignments.

Fig. 38 Fixed function keys (left side of stand)

- 1 Toggling between transmitted light / incident light
- 2 Aperture diaphragm
- 3 Light intensity
- 4 Field diaphragm

7.5 Köhler illumination

7.5.1 Transmitted light

Suitable aperture and field diaphragm values have been preset for each objective. The condenser is also already centered at the factory. However, depending on how the condenser is disassembled and reassembled, it may be necessary to re-center the condenser in some cases. Therefore, check the condenser centering.

The following procedure is provided for the transmitted light-bright field illumination.

- Select an objective with moderate magnification (10x-20x).
- Push the **TL/IL** button (38.1) as needed to activate the transmitted light axis. „**TL**” appears in the first line of the display.

Note on DM4 P: Ensure that the Bertrand lens is swiveled outwards.

- Press the **BF** button to activate the bright field contrast method (one of the variable function keys behind the focus wheel). „**TL BF**” appears in the first line of the display.
- Insert the specimen in the stage's specimen holder (39.3).
- Focus on the specimen. The focus wheel on the left side of the stand allows focus adjustment in course and fine increments. On the right side of the stand, there is also a focus wheel for fine focusing.
- Adjust the light intensity with the **INT** buttons (38.3).
- Close the field diaphragm with the **FD** function key (38.2) until the edge of the diaphragm appears in the specimen plane.

Fig. 39 Stage with specimen holder

- 1 Object motion (X direction)
- 2 Object motion (Y direction)
- 3 Specimen holder
- 4 Condenser height adjuster


Fig. 40 Condenser centering

- 1 Centering screws

7. Startup

- Using the condenser height adjuster (39.4), adjust the condenser until the edge of the field diaphragm appears in sharp relief.
- If the image does not appear in the middle of the field of view (41c), the condenser must be moved into the middle of the field of view with the help of the two centering screws (40.1).
- Open the field diaphragm just enough for it to disappear from the field of view (41d).


Caution:

The condenser height adjustment depends on the thickness of the specimen. It may need to be adjusted for each specimen.

Fig. 41 Köhler illumination

- a** Field diaphragm not focused, not centered
- b** Field diaphragm focused, but not centered
- c** Field diaphragm focused and centered
Diameter is too small, however
- d** Field diameter (light) = Field diameter (view)
(Köhler illumination)

7.5.2 Incident light

Suitable aperture and field diaphragm values have been preset for each objective. The incident light module has also been centered at the factory.

However, it may be necessary to readjust the incident light module in some cases after transporting and setting up the stand. Therefore, check the aperture and field diaphragm centering.

The following procedure is provided for the incident light-bright field illumination.

- Select an objective with moderate magnification (10x-20x).
- Activate the incident-light axis with the **TL/IL** button (38.1).
„IL” appears in the first line of the display.
- Select the bright field contrast method by pressing the **IL-BF** button (DM4000 M LED) or select fluorescence by pressing the **FLUO** button (DM4000 B LED, DM5000 B).
These functions can be assigned to the variable function keys on the stand.
IL BF / FLUO appears in the first line of the display.
- Insert the specimen in the stage's specimen holder (39.3).
- Focus on the specimen with the focus wheels.
- Adjust the light intensity using the **INT** buttons (38.3).

Adjusting the field diaphragm

- Close the field diaphragm with the **FD** button (38.4) or manually until the edge of the diaphragm (round or rectangular) appears in the field of view.
- If the border of the field diaphragm does not appear in the middle of the field of view, the field diaphragm must be moved into the middle of the field of view using the two centering screws (42a.1) on the right side of the stand.
- Use the function buttons **FD** (38.4) to open the field diaphragm to the point that they just disappear from the field of view.
- We recommend the use of a rectangular field diaphragm when using a digital camera. Match the size of the diaphragm to the chip size of the camera.

Adjusting the aperture diaphragm (for DM4000 M LED and DM4500 P LED only)

- Remove one eyepiece (e.g. right).
- Close the aperture diaphragm with the **AP** function key (38.2) until the edge of the diaphragm appears in the exit pupil of the objective (aperture diaphragm plane).
- If the image is not in the middle of the field of view of the exit pupil, move the position of the aperture diaphragm to the center of the exit pupil using the two centering screws (42b.2) located on the left side of the stand.
- Open the aperture diaphragm to cover 2/3 of the field of view.

Fig. 42a Adjusting the field diaphragm in the incident light axis

1 Adjusting screws for moving the field diaphragm


Fig. 42b Adjusting the aperture diaphragm in the incident light axis

1 Adjusting screws for moving the aperture diaphragm

7. Startup

7.6 Checking the phase contrast rings

If your microscope is equipped for the use of phase contrast, the light rings that fit the objectives are built into the condenser.

The light rings are already centered in the factory. However, the centering should be rechecked.

Note:

Every objective is assigned its own light ring in the condenser disk. Therefore, a check must be performed for each objective. When rotating in an objective that is suitable for phase contrast, the corresponding light ring is set automatically.

- Press the **BF** (bright field) button (one of the variable function keys, to the left behind the focus wheels).

- In the place of an eyepiece, insert the focusing telescope (Fig. 43) into the observation tube.
- Rotate the phase contrast objective with the lowest magnification into place.
- Focus on the specimen with the focus wheel.
- Focus the ring structure (44a.a) by slightly loosening the clamping ring (43.2) and moving the eyelens (43.1).
- Retighten the clamping ring.
- Press the **PH** (Phase Contrast) button. The ring diaphragm (light ring) in the condenser is turned inward.
- If the light ring and the phase ring are not shown as arranged in Fig. 44a.c, the light ring must be centered.

Fig. 43 Focusing telescope


- 1 Adjustable eyelens
- 2 Clamping ring for fixing the focus position

Fig. 44a Phase contrast centering procedure

PH=phase contrast ring, LR=light ring

- a Condenser in bright field (BF) position
Light ring (LR) not centered
- b Condenser in phase contrast (PH) position
Light ring (LR) not centered
- c Light ring and phase ring centered

- Insert the centering key through the corresponding openings (44b.1) on both sides of the condenser.
- Turn the centering key until the dark ring (phase ring in the objective) is congruent with the slightly narrower bright ring (light ring in condenser) (44a.c).

Note:

When changing the objective, the centering key must no longer be in the corresponding openings.

- Repeat the process for all other phase contrast objectives.
- Always remove the centering keys after the centering procedure.

7.7 Setting the motorized polarizers (DM4 P LED)

- Select the POL method (one of the variable function keys on the stand or on the Leica SmartTouch).
- Insert the centering key into the corresponding openings (44b.2) on the condenser.
- Set the optimum extinction (max. darkness).

7.8 Adjusting the light sources

Transmitted-light axis (TL) with lamp housing 107/2

The lamp housing 107/2 with a 12 V 100 W halogen lamp is fixed. Centering the lamp is not required.

Fig. 44b Light ring centering

- 1 Centering key in the phase contrast centering hole
- 2 Polarization contrast centering hole


Fig. 45 Reflector cube for lamp adjustment

7. Startup

Incident-light axis (IL) with lamp housing 106 z for fluorescence applications

- When a supply unit is used, it is turned on first.
→ Fig. 56.
- Activate the incident-light axis with the **TL/IL** function button. **FLUO** (Leica DM4 B) or **IL** (Leica DM4 M/DM4 P) appears in the display.
- Insert the lamp adjustment reflector (Fig. 45) in the filter turret in place of a filter cube. Make sure to switch off the instrument first.
Make a note of the designation of the replaced filter cube.

Note:

To avoid faulty adjustments, it is a good idea to remove the filter cube located at the left of the reflector cube as well.

- Turn the reflector into the beam path.
The reflector has reached the correct position when the name of the exchanged filter cube is shown in the upper right of the display.

CAUTION

Never look directly into the beam path!
Beware of the glare hazard when switching to reflector BF or Smith!

CAUTION

Light sources pose a potential irradiation risk (glare, UV-radiation, IR-radiation).

In the lamp housing 106z, the direct image of the filament (in halogen lamps) or the arc (in gas discharge lamps) and its reflection are focused separately and adjusted in relation to one another.

On the left side of the microscope, there is an adjustment window (1.14, p. 18) for mapping the light source.

Adjust the lamp as follows while observing the light source in the adjustment window.


CAUTION

The lamp and the lamp housing may still be hot.

Fig. 46 106 z lamp housing

- Lamp height adjustment
- Mirror image height and side adjustment
- Focusing the reflector
- Lamp side adjustment
- Collector (focusing of the lamp image)

Centering the 12 V 100 W halogen bulb

- In the adjustment window, you see the direct filament image and the mirror image, which as a rule are shifted together.
- Focus the direct filament image with the collector (46.6).
- Use the adjusting buttons on the rear side of the lamp housing (46.2, 46.4) to rotate the mirror image of the lamp filament to the side or completely out of the beam path. The lamp filament's focused image remains visible (Fig. 47).
- Adjust the direct filament image using the adjusting knobs (46.1) and (46.5) so that the centering surface is halfway covered (Fig. 48).
- Then rotate the mirror image of the lamp filament using the adjusting knobs (46.2) and (46.4), and focus it using the reflector (46.3).
- Align the mirror image symmetrically to the filament image (Fig. 49). To do so, again use the adjusting knobs (46.2) and (46.4).
- Defocus the image with the collector head (46.6) until the filament image and mirror image are no longer recognizable and the image is uniformly illuminated.
- Replace the lamp adjustment reflector with the original filter cube.

Note:

Make sure to switch off the instrument.

Fig. 47

Direct lamp filament image focused, but not centered
(in reality, the image is less focused)


Fig. 48

Direct lamp filament image in target position
(in reality, the image is less focused)



Fig. 49

Direct lamp filament image and mirror image in target position
(in reality, the image is less focused)

7. Startup

Centering the Hg 50 W Mercury Lamp

- The adjustment window shows the direct image of the arc and its mirror image. These are generally not in alignment with one another.
- Focus the direct image with the collector (46.6).
- Use the adjusting buttons on the rear side of the lamp housing (46.2, 46.4) to rotate the mirror image of the arc to the side or completely out of the beam path. The arc's focused image remains visible (Fig. 50).
- Use the adjusting buttons (46.1) and (46.5) to place the direct arc image right or left on an imaginary center line of the centering plane (Fig. 51).
- Then rotate the mirror image of the arc with the adjusting knobs (46.2) and (46.4), and focus it using the reflector (46.3).
- Orient the mirror image symmetrically to the direct image (Fig. 52). To do so, again use the adjusting knobs (46.2) and (46.4).
- Using the collector, defocus the image with the collector head (46.6) until the arc image and mirror image are no longer recognizable and the image is uniformly illuminated.
- Replace the lamp adjustment reflector with the original filter cube.

Fig. 50

Direct arc image focused but not centered
(in reality, the image is less focused)


Fig. 51

Direct arc image in target position
(in reality, the image is less focused)



Fig. 52

Direct arc image and mirror image in target position
(in reality, the image is less focused)

Centering the Hg 100 W and Xe 75 W mercury lamps

- The adjustment window shows the direct image of the arc and its mirror image. These are generally not in alignment with one another.
- Focus the direct image with the collector (46.6).
- Use the adjusting buttons to pivot the arc's mirror image on the rear side of the lamp housing (46.2, 46.4) to the side or completely out of the beam path. The arc's focused image remains visible (Fig. 53).
- Use the adjusting buttons (46.1) and (46.5) to place the direct arc image in the middle of the centering plane, whereby the bright tip of the arc, the focal spot, should lie slightly outside of center (Fig. 54).
- Then pivot the arc's mirror image with the adjusting knobs (46.2) and (46.4), and focus it using the reflector (46.3).
- Orient the mirror image symmetrically to the direct image (Fig. 55). To do so, again use the adjusting knobs (46.2) and (46.4). The V-shaped irradiation of the direct image and mirror image arcs can be superimposed.

Fig. 53

Direct arc image focused but not centered
(in reality, the image is less focused)


Fig. 54

Direct arc image in target position
(in reality, the image is less focused)



Fig. 55

Direct arc image and mirror image in target position
(in reality, the image is less focused)

7. Startup

The bright tips of the arcs, the focal spot, must never be projected onto each other, as this results in a danger of explosion by overheating.

Caution!

In older lamps, the structure of the arc is no longer clearly recognizable. The image is then more like that of a HG 50 lamp. The image and mirror image can no longer be superimposed exactly. In this case, align both images.

- Using the collector, defocus the image with the knob (46.6) until the arc image and mirror image are no longer recognizable and the image is homogeneously illuminated.
- Replace the lamp adjustment reflector with the original filter cube.

Note:

Make sure to switch off the instrument.

8. Operation

8.1 Switching on the unit

When using a gas discharge lamp, the external supply unit must be turned on separately (56.1). Switch on the microscope at the power switch. All motorized microscope components will then run through an initialization phase. After the initialization is complete, the display on the stand (Fig. 57) shows the current microscope setting.

Fig. 56 Front view of the ebq 100 supply unit

- 1 Power switch
- 2 Lamp status


Fig. 57 Display after initialization

	TL_BF		
	10x Obj. 1.25x MagCh.	Σ 125x	
	INT 100%		
	AP 33	FD 30	
	50%	50%	

8.2 Stages and object displacement

Lengthening the coaxial drive


- For lengthening, pull the lower handle (58.2) downwards. Repeat with the upper handle (58.1).

Setting the movement rate (torque)

The torque for x and y can be individually adjusted using two knurled rings (58.3, 58.4).

Fig. 58 Rotating specimen stage

- 1 Object displacement (Y direction)
- 2 Object displacement (X direction)
- 3 Torque adjuster (Y-direction)
- 4 Torque adjuster (X-direction)
- 5 Focus wheel for fine focusing

8. Operation

Rotating the stage

The swivel range of the rotating stages is 0° - 110° .

- Rotate the stage to loosen the locking screw (59b.1).
- Bring the table into the desired position.
- Retighten the locking screw.

Pol rotating stage (360°)*, Pol attachable mechanical stage *

The specimen can be fastened using either two spring-back stage clips or the Pol multi-format attachable mechanical stage (Fig. 59a). For glass slides of approx. 26 mm (1") in width, the metal plate (59a.2) needs to be pivoted outwards and the specimen needs to be positioned according to Fig. 59a. If commercially available glass slides of 26 mm in width are placed perpendicular to that, then the range of movement of the attachable mechanical stage, which is normally approx. 30 mm x 40 mm, is not fully used. The pair of lock buttons provided allow for locking distances of 0.1, 0.3, 0.5, 1 and 2 mm. Changing them requires forceful, axial removal. When attaching the new lock button, make sure that the catch pins on the inside are properly positioned. For smaller microscopes, the limit stop screw on the underside must be shifted approx. 2 mm inwards toward the stroke limit.

The two verniers allow the angular measurements to be read with a precision of 0.1.


45°click stop:

- Screw in the rotary knob (59a.5) until you can feel it starting to resist.
- Then turn the specimen stage to the next click stop.
- Loosen the rotary knob and look for the starting point of the next click stop (e.g., the object darkness setting).
- Retighten the rotary knob.

Now the rotating stage can be turned using 45° click stop intervals.

Fig. 59a Pol rotating stage* and Pol attachable mechanical stage*

- 1 Holes for the fastening screw.
- 2 Lever for the holder for glass slides of various formats, which can be turned inward and outward
- 3 Storage for the centering key
- 4 Locking button pair
- 5 45° click stop
- 6 Clamping system for the stage rotation

8.3 Focusing

There is a focus wheel on the left side of the stand for coarse and fine focus adjustment (Fig. 59b).

On the right side of the stand, there is also a focus wheel, which is used exclusively for fine focusing (58.4).

The special form of this handwheel makes it possible to simultaneously grasp the coaxial drive with your hand while operating the fine adjustment with one finger.

8.4 Tubes

Note:

Close any unused tube openings, as otherwise stray light can interfere with observation.

Note:

Make sure that the connecting cable is plugged in on the MBDT25+ motorized tube (60.1).

Adjust the interpupillary distance

- Adjust the interpupillary distance of the eyepiece tubes so that a congruent total image is seen (Fig. 60).

Fig. 59b Rotating specimen stage

- 1 Locking screw
- 2 Fine focusing
- 3 Coarse focusing


Fig. 60 Tube setting

↔ Setting the personal interpupillary distance

- 1 Motorized tube connection

8. Operation

Adjusting the viewing angle

- For the AET22 and EDT22 ergo tubes, the viewing angle can be adjusted by tilting the binocular eyepiece in the range of 5° – 32° (Fig. 61).

Adjusting the eyepiece extension to the arm length

- On the AET22 tube, the eyepieces can be extended up to 30 mm (Fig. 61).

Beamsplitting in photo tubes

EDT22 tube:

The light division between the observation and documentation ports has a definite presetting (50:50).

BDT25+ tube:

The light division is set manually by pulling out a control bar.

Control Bar		Observation	Photo
VIS		100%	0%
50/50		50%	50%
PHOTO		0%	100%

MBDT25+ tube:

This tube is similar to the documentation tube BDT25+, but it is motorized.

The control positions are selected using a variable function key on the stand.

HC L 2TU tube:

The light division is set manually by pulling out a control bar.

Control Bar		Observation	Photo
VIS		100%	0%
PHOTO		0%	100%

Fig. 61 Individual settings of the AET22 tube


Fig. 62 BDT25+ tube with digital camera

1 Control bar

8.5 Eyepieces

Note:

The eyepiece's aperture protector must be removed or folded back, during microscopy while wearing glasses.

We recommend removing bifocals and spectacles with progressive-addition lenses when using the microscope.

- For the adjustable tubes with documentation output, choose the 100% VIS position.

Eyepieces with inlaid reticle

- Focus the reticle by adjusting the eyelens in the eyepiece.
- Focus on the object through this eyepiece.
- Then, close that eye and focus on the object by adjusting the second ocular only.

Correction for vision problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.
- Then, with your left eye, view the same position of the specimen and rotate the left eyepiece tube until this position is brought into sharp focus. Do not use the focus wheel.

8.6 Objectives

The objectives are moved into the beam path manually. Be sure that the turret locks into place.

The positions of the objectives in the objective turret have been specified at the factory and must be observed when installing the objectives.

(See Installing objectives → p. 33)

When pivoting an objective inwards, the microscope automatically selects:

- The optimum setting for the field diaphragm
- The optimum setting for the aperture diaphragm and the light intensity for each contrast method.

The objective magnification and the total magnification appear in the display.

- Begin with a small level of magnification. Then switch to the next higher objective.
- For **immersion objectives** use the appropriate immersion medium.
OIL: use optical immersion oil only according to DIN/ISO standards
Cleaning → p. 87
W: Water immersion
IMM: Universal objective for water, glycerol, oil immersion

Caution!

Follow the safety instructions for immersion oil!

8. Operation

For lockable immersion objectives:

- Lock these by pushing the front part upwards until it stops (approx. 2 mm).
- Then, after a gentle turning motion to the right, the objective is locked (Fig. 63).

For objectives with corrective mounts:

- Turn the knurled ring to adjust the objective to the thickness of the cover glass.

Objective centering * (DM4 P)

Caution:

When changing the objective, the centering key must no longer be in the corresponding openings.

When centering the objectives (Fig. 64, 65), use two hexagon wrenches to move the objectives so that the optical axis of the objective (and, therefore, the center of the image) is aligned with the axis of rotation of the objective stage. If the objectives are centered correctly, a programmed specimen position will not drift out of the field of vision when the stage is turned. Therefore, a specimen point located inside the center of the cross-hairs does not change its position when the stage is rotated a full turn. When centering objectives, it is best to use a detailed, high-contrast specimen.

Fig. 63a Immersion objective, released

Fig. 63b Immersion objective, locked

- Switch off the analyzer, the 1.6x tube lens and the Bertrand lens.
- Reduce the aperture diaphragm so that it is very small.
- Insert both objective centering keys above the objectives that need to be centered.
- Focus the specimen.

Two resembling methods can be used for centering objectives:

Method I (Fig. 64)



- Turn the specimen stage and note the position of the specimen that does not move in a circular path. This position of the specimen corresponds to the mechanical axis of rotation of the specimen stage.
- Now move the marked specimen position by shifting the two centering keys to the center of the cross-hairs.
- Turn the specimen stage and refine the centering as needed.

Method II (Fig. 65)

- Move the marked specimen position (65a) into the center M of the cross-hairs.
- Turn the specimen stage until the specimen position is as far from the center of the M cross-hairs as possible (Position A, Fig. 65b). In extreme cases, point A (= maximum deviation of the specimen position) can also be located outside of the field of view.
- Move the image by turning the centering key until position A of the specimen is located in the center (= Pos. B) between pos. A and the center of the M cross-hairs (65c).
- Move position A of the specimen to M and check to see whether A remains in M when the stage is rotated (65d). If necessary, repeat the centering procedure.

The objective centering procedure needs to be repeated for each objective. This ensures that the objectives retain their approximate centering settings when they are removed for cleaning, or other such procedures, and then reinserted into the same holes. If the coarse drive or the height adjustment device is used to change the height of the specimen stage (for example, when viewing thick specimens) the centering precision for all of the objectives may be reduced slightly.

Fig. 64 Centering method I Fig. 65 Centering method II

8. Operation

8.7 Magnification changer

Optionally, a coded magnification changer can be used, which is manually operated.

On the knurled ring, the following magnification factors can be set:

B stand	M stand
1x	1x
1.25x	1.5x
1.6x	2x

The selected factor is indicated in the display and included in the total magnification.

8.8 HC P 1x/1.6x tube optics with Bertrand lens (encoded)

These optics were developed especially for polarizing microscopy; however, they can also be used for all other methods.

- If needed, switch the Bertrand lens off and the **1x** tube factor on.
For the HC P (Pol) tube optics, switching to the 1x tube factor is sufficient.

Setting the tubes and eyepieces → p. 58f.

Features:

- 1x tube factor, can be switched to 1.6x
- Bertrand lens, can be activated, encoded, focused and centered
- Iris diaphragm in the intermediate image for masking out small particles (15 µm for the 100x objective).

Built-in, depolarized quartz plate:

This plate prevents interference colors, which are caused by tube prism polarization effects (pseudopleochroism), when the analyzer is connected and the polarizer is switched on; however, it only prevents this when the **1x** tube factor is being used.

When using the **1.6 fix** tube factor, make sure that the useful magnification for high objective magnifications and apertures (objective aperture x 1000) is not exceeded to the point that overmagnification causes a blurry image impression.

8.9 Light sources

- Adjust the brightness using the function keys (66.5). The **INT** function buttons are always assigned to the currently active transmitted light (TL) or incident light (IL) axis.

- For TL and IL:

The setting can be made in coarse and fine steps. Pressing both **INT** buttons at the same time toggles between coarse and fine adjustment. The light intensity in the display changes accordingly.

0-20

Coarse adjustment: [REDACTED] =====

0-255

Fine adjustment: [REDACTED] -----

- The intensity is individually adjusted and stored for each objective and contrast method.
- For fluorescence:
The brightness is adjusted in 5 defined increments (FIM):
100% / 55% / 35% / 20% / 10%

8.10 Aperture diaphragm and field diaphragm

Both diaphragms have been set to suitable values for the current objective and contrast method at the factory.

- The diaphragms can be adjusted at any time using **AP** (aperture diaphragm) (66.2) and **FD** (field diaphragm) (66.4) function keys.

Caution!

The old values will be overwritten by the current ones!

- The function keys are assigned to the currently active transmitted light (TL) or incident light (IL) axis.

Caution!

When the **PH** or **DF** is being used, the aperture diaphragm is opened completely and cannot be closed to prevent incorrect operations.

Fig. 66 Control panel

- Variable function keys
(also on the right side of the stand)
- Aperture diaphragm
- Transmitted light / incident light
- Field diaphragm
- Light intensity

9. Contrast methods for Leica DM4 B / Leica DM4 P

9.1 Transmitted light

9.1.1 Bright field (TL)

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Select the **BF** (bright field) contrast method. Do so by pressing the **BF** variable key. Alternatively: press the variable button **CHANGE TL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **BF**.
- Insert a transmitted light specimen.
- Rotate an appropriate objective into place.
- Use the focus wheel to bring the image into focus and set the brightness using the **INT** function key.

CAUTION

Do not look into the eyepieces when toggling the contrast method!

During the toggling procedure, the maximum radiant power of the light source may be present on the eyepieces for a short time and can temporarily blind the user!

9.1.2 Phase contrast (TL)

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Select the **PH** (phase contrast) contrast method. Do so by pressing the variable button **PH**. Alternatively: press the variable button **CHANGE TL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **PH**.
- Insert a transmitted light specimen.
- Rotate an appropriate objective into place. Objectives that are suitable for phase contrast are engraved with **PH**.
- Bring the image into focus using the focus wheel and set the brightness using the **INT** function key.

Notes:

The microscope automatically selects the correct light ring in the condenser.

When selecting the phase contrast method, the aperture diaphragm is opened fully and can not be adjusted.

9.1.3 Dark field (TL)

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Select the DF (dark field) contrast method by pressing the variable button **DF**.
Alternatively: press the variable button **CHANGE TL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **DF**.
The dark field ring (dark field stop) is set automatically in the condenser.
- Insert a transmitted light specimen.
- Rotate an appropriate objective into place.
- Bring the image into focus using the focus wheel and set the brightness using the **INT** function key.

Notes:

The maximum objective aperture that can be used for dark field is **0.75**. All objectives with greater apertures are automatically blocked for this method („DF“ flashes in the display).

The microscope automatically selects the correct light ring in the condenser.

When the dark field method is selected, the aperture diaphragm is opened fully and cannot be adjusted.

9.1.4 Polarization (TL)

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Select the **POL** (polarization) contrast method by pressing the variable button **POL**.
Alternatively: press the variable button **CHANGE TL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **POL**.

9.1.4.1 Manual method

- Rotate the polarizer on the underside of the condenser into the beam path (Fig. 67a). Make sure that the red index point on the front of the polarizer is at 0.
- Insert the analyzer up to the stop into the left side of the stand (67b.1).
- Bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Insert a specimen and rotate a suitable objective into place.

9.1.4.2 DM4 P - examinations in polarized transmitted light


One polarizer only

If specimens need to be examined using other transmitted light methods, such as bright field, phase contrast and dark field, instead of with polarizers, then it is usually sufficient to switch off the analyzer or polarizer. If the brightness of the image is not adequate, the polarizer and the analyzer need to be switched off. Stained, birefringent specimens may exhibit brightness or color fluctuations when the specimen stage or the polarizer is rotated (with the analyzer switched off). This is called dichroism or pleochroism, which is an important indication when examining crystals. However, this effect can be simulated by non-polarizing microscopes, because they do not contain depolarized quartz plates, or when an incident light reflector is left in the beam path during the transmitted light method. This also applies to using the 1.6x tube lens in the DM4 P.

Fig. 67 Crossing the polarizers for observation using the Bertrand lens, objective with high aperture, without specimen

- a** crossed precisely
- b** not crossed precisely

If pressure is present in the condenser or in the objective, position „a“ cannot be set at all

Crossed polarizers

According to DIN and ISO ratings the run of the vibration directions of the polarizer is East-West, that of the analyzer with crossed polarizers North-South. However, when polarizers are crossed, these vibrations have the same polarized optical appearances as they do when each of the polarizers are interchanged by 90°.

If the specimen contains a large amount of non-birefringent or opaque particles, then the analyzer is usually turned a few degrees away from the cross position to make these particles somewhat visible (for specimens that remain dark when the polarizers are crossed precisely). Examinations are not usually performed when the polarizers are parallel because this configuration is not sensitive enough to detect birefringent elements.

Changing the brightness when rotating birefringent objectives

When the specimen stage is rotated, the brightness of birefringent (anisotropic) specimens changes periodically. When the specimen is rotated, a total of 4 extinction positions (also called normal position) occur at each 90° interval. Exactly 45° between 2 extinction positions, 4 orientations of maximum intensity occur, the diagonal or 45° positions. In extinction, the specimen vibration directions are parallel to the transmission directions of the polarizers; at maximum intensity, the specimen vibration directions represent the bisectors of the polarizer directions. The crosshairs in the (right) eyepiece on the polarizing microscope can be rotated N – S/E – W; therefore, they can be rotated in the direction of polarization or 45°, or set up according to the specimen vibration directions when in a diagonal position.

Simple overview observations

- Place the transmitted light specimen on the polarizer.
- Turn the condenser head inward.
- Focus through the condenser using a low-power magnifying lens, such as a 5x.

Although this method does not provide good imaging performance; it does make it possible to view rows of specimen very quickly.

$\lambda/4$ - and λ plate, quartz wedge

Depending on their design, $\lambda/4$ and λ plates are installed on the underside of the condenser or, for polarizing microscope, in the 8x condenser disk (the vibration direction γ runs: ↗) or they are inserted in the tube slot. An automatic, spring-mounted dust cover flap closes the tube slot.

For the IC/P analyzer, the λ plate can be activated by turning it so that the mark " λ " is facing upwards.

When the plate is activated, the path difference is increased or reduced according to Fig. 68. The corresponding color changes can be used to determine the vibration direction γ according to the larger refractive index (i.e. refractive index n_1). (The quartz wedge (69.7) allows various color shifts on the polarizing microscope.)


Circular polarization

When the specimen stage is rotated, birefringent specimens show 4 extinction positions.

When a larger number of birefringent specimens are present, a few birefringent specimens always show up randomly in the extinction position. Circular polarization is used for simultaneous interference color observation of all specimens:

- Remove the specimen from the beam path or locate the blank position of the specimen.
- Cross the polarizers precisely. The polarizers must be positioned exactly in the N – S/E – W direction. This means that the analyzer must either be set exactly at 90° or 0° of polarization.

Fig. 68 Interference color with dependence on the path difference or on the thickness and color changes in the addition and subtraction areas of the λ and $\lambda/4$ plates.

9. Contrast methods for Leica DM4 B/DM4 P

- Insert the $\lambda/4$ plate (69.5) into the tube slot.
- Insert the $\lambda/4$ plate (69.1) into the receptacle located on top of the polarizer and rotate it until a specimen-free field of vision appears at its maximum extinction position (make sure to cross the polarizers precisely before doing this).

Compensators for quantitative measurements

Adjustable compensators are used to obtain exact measurements of path differences.

To obtain the usual specimen thickness, d , and the measured path difference, γ (Γ), the birefringence Δn is calculated according to the following formula:

$$\Gamma = d \times \Delta n \text{ [nm]} \text{ or } \Delta n = \frac{\Gamma}{d}$$

When obtaining measurements, the compensator is inserted into the compensator slot and adjusted until the specimen position being measured is at the maximum extinction position. For this, the specimen must be pitched at a diagonal.

For HC P tube optics, use an iris diaphragm to mask out the measuring points.

Detailed information can be found in the instructions for the compensator.

The following compensators are available:

Tilt compensator B with a measurement range of up to 5 orders based on Berek techniques


Compensator (69.7) with an MgF_2 plate for measuring up to approx. 5 orders of path differences in white or monochromatic light. You can read the path difference directly from the provided measurement table by calculating the sum of the two angles that are formed when the compensators are tilted simultaneously.

Tilt compensator K with a measurement range of up to 30 Orders (69.6)

For measuring path differences in white or monochromatic light up to the specified maximum path difference. The compensator plate is made of calcite. The evaluation is created by performing simple calculations using the provided tables and the specified measured constants. Measure in white or monochromatic light.

Fig. 69 Compensators

1, 2 $\lambda/4$ u. λ plate. For polarizing microscopes only: **3** $\lambda/4$ u. λ plate for the 8x turret disk, **4, 5** $\lambda/4$ u. λ plate for the tube slot, **6** Quartz wedge, **7** Tilt compensator

Using conoscopy for crystal structures

Birefringent crystals create interference images, also called axis images or conoscopic images (Fig. 71a,b), in the exit pupil of the objective (e.g. inside the objective). The form of these interference images and the changes that occur in these images when using a compensator make it possible to state how many axes the crystals have (uniaxial or biaxial crystals), how the axes are oriented, and whether or not the birefringence is positive or negative (positive or negative birefringent crystals).

Because the interference images appear in the eyepoint, they are not visible during typical observation (orthoscopy). An improvisational method for observing these images is to remove the eyepiece from the tube and to use a monocular, held a few cm away, for looking into the tube. You can improve the observation by using a focusing telescope for the phase contrast.

However, additional crystals located in the field of view need to be masked out because they disrupt the interference images of a crystal located in the middle of the field of view.

Setting the conoscopy

For conoscopy, the specimen positions that are most suitable are those that have the lowest path differences (table Fig. 68).

Efficient conoscopic observation requires that the objective be centered precisely and that the cross position of the polarizers is exact.

- Rotate an objective with an aperture that is as high as possible (e.g. 40x, 50x or 63x) into the beam path.
- Rotate the condenser head into the beam path.
- Open the aperture diaphragm.
- Place the crystal being examined as close to the center of the field of vision as possible.
- Turn the 1.6x tube lens inward.
- Widen or narrow the iris diaphragm according to the size of crystal, and make the field diaphragm narrower if necessary.

9. Contrast methods for Leica DM4 B/DM4 P

- Insert the Bertrand lens (Fig. 70B) and focus it by rotating the operating button until the interference image or the circular diffuse edge of the eyepoint is in focus.
If needed, center the Bertrand lens: Insert the hexagonal keys into the centering holes in sequence. If necessary, align the right eyepieces so that the cross-hairs correspond approximately to the direction of movement during the centering process.
- Adjust the collector to its optimum setting; use the diffuser if necessary.

Determining the optical character

Uniaxial crystals (Fig. 71a)

When observing uniaxial crystals in the conoscopic (diverging) beam path, a dark cross appears whose center point indicates the optical axis. The cross is surrounded by colored interference bands*. When a variable compensator (quartz wedge or tilt compensator) is used, the rings move toward the center point or outside two opposite quadrants in the cross. The optical character results from the direction of the movement of the rings according to Fig. 71.

Cross sections in which the crystal optical axis is sloped toward the direction of the viewer are suitable for determining the optical character. In addition, an optical character can be determined even if the center point of the cross is located outside of the field of view. Fig. 71 shows that fixed compensators can be used in place of variable compensators for determining the optical character.

Even if only one of the optical axes is within the observer's direction of sight, the optical character can usually be identified. The brightness for specimens oriented in this way changes in the orthoscopic beam path very little or not at all when the objective is rotated. Consequently, only one of the two isogyres are visible in the conoscopic beam path.

Biaxial crystals (Fig. 71b)

The cross sections in which the bisectors of the optical axes run parallel to the direction of sight are especially suitable for determining the optical character (the section is perpendicular to the acute bisectrix).

A dark cross can be identified in the divergent beam path. This cross splits into two hyperbolic lines, also called isogyres, when the specimen stage is turned. The cross and or hyperbolic isogyres are surrounded by colored interference bands. After the compensator has been activated, the optical character can be determined from the direction of movement of these bands according to Fig. 71 or the following rule: The screw axis symmetry of the isogyres must run perpendicular to the γ -direction of the compensator:

Fig. 70 Functions of the HC P Pol tube optics

Control element	1x Orthoscopy	1.6x Orthoscopy	Conoscopy
Tube lens	1x	1.6x	1.6x
Iris diaphragm	user-defined, as it is not in the beam path	adapted for the field of vision	> Specimen
Bertrand lens	on or off	off	on
Polarization	on or off	on or off (not for Dichroism/Pleochroism)	crossed

Biaxial, positive crystals:

When the compensator is activated, the movement of the interference bands proceeds from the convex side to the concave side of the isogyres.

Biaxial, negative crystals:

The interference bands move from the concave side to the convex side.

Possible faults

The polarizers have been damaged (discolored) by strong light sources or they are heavily soiled.

The objectives or the condenser have been deformed by mechanical damage.

There is a beam splitter or filter between the polarizers.

The embedding material for transmitted light specimens is birefringent.

Fig. 71a Determining the optical characters for uniaxial structures

Left: Positive uniaxial crystal, cut perpendicular to the optical axis.

Right: Negative uniaxial crystal, cut perpendicular to the optical axis.

- 1 Illustration of the vibration directions in the specimen and in the compensator
- 2 Interference figure changes when using a $\lambda/4$ plate
- 3 Interference figure changes when using a λ plate

Fig. 71b Table for determining the optical character

Orientation of the condenser plate	Uniaxial		Biaxial			
	+	-	+		-	
* When using $1/4-\lambda$ glimmer plates, black points appear instead of black arcs.						

9. Contrast methods for Leica DM4 B/DM4 P

9.1.4.3 Motorized method

- After the **POL** contrast method has been selected, the position of the polarizer is switched inside the condenser (if the microscope is equipped with these components). The analyzer cube is also automatically brought into the beam path.

9.1.4.4 Combined methods

- For the Leica DM4 B, DM4 P microscopes, it is possible to combine purely mechanical and motorized components.

For DM4 P only:

After the Bertrand lens (links) is turned inwards, **CONOS** appears on the display. The light intensity and aperture and field diaphragms are assigned to the objective currently in use (objective magnification for conoscopy is 40x, 63x or 100x).

9.1.5 Differential interference contrast (TL)

9.1.5.1 DM4 P

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Insert a specimen and rotate a suitable objective into place.
- Select the **DIC** contrast method.
Do so by pressing the **DIC** variable key.
Or by pressing the variable button.
CHANGE TL .
(For details on key assignments, please see the identification sheet.)
ICT appears on the display.
- Turn the polarizer on the underside of the condenser into the light path. Make sure that the red index point on the front of the polarizer is at **0** (Fig. 72).

Fig. 72 Turn the polarizer inward

1 Polarizer


Fig. 73 Inserting the analyzer

1 Analyzer

- Insert the analyzer on the left side of the stand up to the click stop (fig. 73).
- Then insert the objective prism slider into the slot on the objective turret.
- The DIC prism that needs to be used appears on the display.
- For fine adjustment, rotate the knurled screw (74.1) on the objective prism slide.

Fig. 74 Motorized variants
1 Knurled wheel for fine adjustment

9. Contrast methods for Leica DM4 B/DM4 P

9.2 Fluorescence

- Switch to the fluorescence axis / fluorescence cube (FLUO) by pushing the **TL/IL** function key.
- Apply a specimen and rotate a suitable objective into place.
- The current fluorescence filter cube is indicated on the display.
- You may protect your specimen from fading by closing the incident light shutter. Do so by pressing the **SHUTTER** variable button.
(For details on key assignments, please see the identification sheet.)
- Changing the fluorescence filter cube:
By pressing the variable key
Cube ↺ or **Cube** ↻
- Bring the image into focus using the focus wheel and set the brightness using the **INT** function key.
- The fluorescence intensity can be increased using the booster lens on the right side of the stand (Fig. 75). Leica recommends that you insert the booster lens into the front slot.

This symbol appears on the display: ↘

Fig. 75 Inserting the booster lens

10. Contrast methods for Leica DM4 M

10.1 Incident light

10.1.1 Bright field (IL)

- Switch to the incident light axis (IL) by pushing the **TL/IL** function key.
- Select the **BF** (bright field) contrast method
Do so by pressing the **BF** variable key.
Alternatively: press the variable button **CHANGE RL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **BF**.
- Insert a specimen.
- Rotate an appropriate objective into place.
- Bring the image into focus using the focus wheel and set the brightness using the **INT** function key.
The brightness values for each objective are saved.

Do not look into the eyepieces when toggling the contrast method!

During the toggling procedure, the maximum radiant power of the light source may be present on the eyepieces for a short time and can temporarily blind the user!

10.1.2 Dark field (IL)

- Switch to the incident light axis (IL) by pushing the **TL/IL** function key.
- Select the **DF** (dark field) contrast method by pressing the variable button **DF**.
Alternatively: press the variable button **CHANGE RL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **DF**.
The **DF** reflector is rotated into the beam path.
- Insert a specimen.
- Rotate an appropriate objective into place.
- Bring the image into focus using the focus wheel and set the brightness using the **INT** function key.

Notes:

The maximum objective aperture that can be used for dark field is **0.75**. All objectives with greater aperture are automatically blocked for this method („DF“ flashes in the display).

When the dark field method is selected, the aperture diaphragm is opened fully and cannot be adjusted.

10. Contrast methods for Leica DM4 M

10.1.3 Polarization (IL)

- Switch to the incident light axis (IL) by pushing the **TL/IL** function key.
- Select the **POL** (polarization) contrast method by pressing the variable button **POL**.
Alternatively: press the variable button **CHANGE RL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **POL**.

Automatic procedure:

- The ICR filter cube is automatically brought into the beam path.

Manual procedure:

- Manually push the corresponding polarizer (77.3) and the IC/P analyzer (78.1) on the stand into the beam path. Then bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Apply a specimen and rotate an objective with a low magnification level into place.

Fig. 77 Objective prism slide

- 1 Knurled wheel for fine focusing
- 2 Prism slot with inserted prism slide bar
- 3 Insert the polarizer


Fig. 78

- 1 Insert the analyzer

10.1.4 Interference contrast (RL)

- Switch to the incident light axis (IL) by pushing the **TL/IL** function key.
- Insert a specimen and rotate a suitable objective into place.
- Select the **DIC** contrast method
Do so by pressing the **DIC** variable key.
Alternatively: press the variable button **CHANGE RL** .
(For details on key assignments, please see the identification sheet.)
The display indicates **ICR**.
- The ICR filter cube (containing the polarizer and analyzer) is automatically brought into the beam path on the incident light axis. Insert the objective prism slide into the prism slot (Fig. 77.2).

Alternatively:

- Manually push the ICR polarizer (77.3) and the IC/P analyzer (78.1) on the stand into the beam path.
- Then insert the objective prism slider into the slot on the objective turret. (77.2).
- For fine adjustment, rotate the knurled screw (77.1) on the objective prism slide.

10.2 Transmitted light

10.2.1 Bright field (TL)

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Select the **BF** (bright field) contrast method
Do so by pressing the **BF TL** variable key.
(For details on key assignments, please see the identification sheet.)
The display indicates **BF TL**.
- Insert a transmitted light specimen.
- Rotate an appropriate objective into place.
- Bring the image into focus using the focus wheel and set the brightness using the **INT** function key.

10.2.2 Polarization (TL)

- Switch to the transmitted light axis (TL) by pushing the **TL/IL** function key.
- Select the **POL** (polarization) contrast method.
- Push the analyzer on the left side of the stand into the beam path (Fig. 78).
- Turn the polarizer on the underside of the condenser into the light path. Make sure that the red index point on the front of the polarizer is at **0** (Fig. 72).

11. Troubleshooting

Problem	Cause/remedy
Stand	
The microscope does not respond.	‣ Ensure that the AC outlet has power. ‣ Make sure that the microscope is connected to the power supply. ‣ Check the cable connections. ‣ Inform Service and have the supply unit fuse checked.
Illumination	
The image is completely dark.	‣ Open the shutter. ‣ Check the connection to the lamp housing on the microscope. Transmitted light connection: $\triangle$ Incident light (fluorescence) connection: ∇ ‣ Make sure that the lamps are connected to the power supply. ‣ Inform Service and have the ebq 100 supply unit fuse checked.
The image is unevenly or not uniformly illuminated.	‣ Remove all unneeded filters from the light path. ‣ Center the lamp ($\rightarrow$ p. 55ff). ‣ Replace the old lamp ($\rightarrow$ p. 34ff).
The illumination flickers.	‣ Be sure that there is no loose connection at the power supply. ‣ Replace the old lamp ($\rightarrow$ p. 34ff).
The lamp does not illuminate immediately upon being switched on.	‣ The ebq 100 must be switched-on repeatedly. ‣ Hot Hg lamps should cool down before being switching on again.

Problem	Cause/remedy
Bright field	
The specimen can not be brought into focus.	▶ Use the correct immersion medium. ▶ Lay the specimen with the cover glass toward the top. ▶ Make sure that the cover glass thickness is correct and that it meets the specifications on the objective.
Dark field	
No definite DF contrast is possible.	▶ Be sure that a DF objective is being used. ▶ The objective aperture setting is too high (maximum 0.75). If necessary, reduce the objective aperture using the iris diaphragm on the objective. ▶ Check the condenser centering.
The image is unevenly or not uniformly illuminated.	▶ The magnification is too weak. Use a higher magnification.
Undesirable stray light.	▶ Clean the specimen and neighboring lens surfaces (→ p. 87).
Phase contrast	
No phase contrast is possible.	▶ The specimen is too thick. ▶ The cover glass is not placed planar. ▶ Check the centering of the light rings (→ p. 54).

11. Troubleshooting

Problem	Cause/remedy
Polarization	
No polarization contrast is possible.	▶ Bring the polarizer and analyzer into cross position until they reach maximum darkness (without specimen). (→ p. 72). ▶ The embedding material for transmitted light specimen is birefringent.
The image is insufficient.	▶ Polarizers can be damaged by strong light sources. Check to see if the polarizer is damaged (discolored) or heavily soiled.
The image is too dark.	▶ Check to see if a beam splitter or filter is switched on.
No conoscopy image is possible.	▶ Objectives and condensers can be deformed by mechanical damage.

Fluorescence

The image is completely dark (no fluorescence).	▶ Open the shutter. ▶ Select the incident light axis (IL) . ▶ Check the antigen-antibody combination.
The fluorescence is too weak.	▶ Insert the booster (→ p. 43). ▶ Center the lamp (→ p. 55ff). ▶ Insert a new lamp (→ p. 34f).

Display

The display flashes.	▶ Rotate an appropriate objective for the contrast method into place.
FAIL! appears on the display.	▶ Check the installed objectives, filter cubes, etc. ▶ Switch the microscope off and back on.

12. Care of the microscope

Caution!

Unplug the power supply before performing cleaning and maintenance work!
Protect electrical components from moisture!

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

The microscope should be cleaned after each use, and the microscope optics should be kept strictly clean.

12.1 Dust cover

Note:

To protect against dust, cover the microscope and accessories with the dust cover after each use.

Let lamps cool down before covering the stand with a dust cover. The dust cover is not heat-resistant. In addition condensation water may occur.

12.2 Cleaning

Caution:

Residual fiber and dust can create unwanted background fluorescence.

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be removed with little soapy water or ethyl alcohol.

For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Thinners containing acetone, xylene or nitrogen can harm the microscope and thus must not be used.

Test cleaning solutions of unknown composition first on a less visible area of the unit. Be sure that coated or plastic surfaces do not become matted or etched.

Cleaning the stage

Rub the stage with paraffin oil or acid-free Vaseline to remove light spots on the stage.

12. Care of the microscope

Cleaning Glass Surfaces and Objectives

Glass surfaces, and particularly objectives, are always to be cleaned as described in the brochure "Cleaning of Microscope Optics". You can download the information from

<http://www.leica-microsystems.com/products/light-microscopes/life-science-research/upright-microscopes>

Select the type of your microscope and switch to the "Download" page.

You can also contact our Technical Service with any questions.

Removing immersion oil

Caution!

Follow safety instructions for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

12.3 Handling acids and bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

Caution:

Be absolutely certain to prevent the optics and mechanical parts from coming into contact with these chemicals.

13. Essential wear and spare parts

Order No. Material No.	Name	Used for
<u>Replacement lamps</u>		
11 500 974	Halogen bulb 12 V 100 W	107/2 lamp housing
11 500 137	High-pressure mercury burner 50W	106 z lamp housing
11 500 138	High-pressure mercury burner 100W	106 z lamp housing
11 500 321	High-pressure mercury burner 100W (103 W/2)	106 z lamp housing
11 500 139	High-pressure xenon burner 75 W	106 z lamp housing
<u>Screw cap for unused objective receptacles</u>		
020-422.570-000	ScREW cap M 25	Objective turret
<u>Replacement eyecup (antiglare protection) for HC PLAN eyepiece</u>		
021-500.017-005	HC PLAN eyecup	10x/25 eyepiece
021-264.520-018	HC PLAN eyecup	10x/22 eyepiece
021-264.520-018	HC PLAN eyecup	10x/20 eyepiece
<u>Immersion oil conforming to DIN/ISO standards, fluorescence-free</u>		
11 513 859	Type F, ISO 8036, very low autofluorescence, highly recommended for highly recommended for APO objectives, 10 ml	OIL and IMM objectives and oil condenser heads
11 513 860	Type N, ISO 8036, low autofluorescence, 20 ml	
11 513 861	Type N, ISO 8036, low autofluorescence, 250 ml	
<u>Fuses for Leica DM4 B, DM4 M and DM4 P:</u>		
11 362 150 010 160	Fine-wire fuse 5x20, T1.6A H, 250 VAC	

14. Abbreviations and pictograms

	Contrast method
	Magnification
	Light intensity / diaphragms
	Light division in the motorized tube
	Transmitted light shutter is open
	Transmitted light shutter is closed
	Incident light shutter is open
	Incident light shutter is closed
AET	Advanced ergo tube
AP	Aperture diaphragm
BF	Bright field
COMBI	Combination methods
CONOS	Conoscopy
CUBE	Fluorescence cube
DF	Dark field incident / transmitted light
DIC	Differential interference contrast
FD	Field diaphragm
FLUO	Fluorescence axis (incident light)
ICR	Interference contrast (incident light)
ICT	Interference contrast (transmitted light)
IL	Incident light
INT	Intensity
MBDT	Motorized basic documentation tube
PH	Phase contrast
POL	Polarization, incident / transmitted light
RL	Incident light
TL	Transmitted light
90	

15. Index

Symbole

$\lambda/4$ - and λ plate 73

A

Adjusting the light sources 55
Allowable ambient conditions 27
Ambient conditions 27
Analyzer 42, 71, 82
Analyzer cube 78
Aperture diaphragm 24, 49, 69

B

Beam splitting in photo tubes 64
Bertrand lens 68, 78
Booster lens 80
Bright field 70, 81

C

Checking the phase contrast rings 54
Circular polarization 73
Cleaning 87
Coaxial drive 61
Compensators 74
Condenser 23, 24, 31
Condenser centering 51
Condenser connection 31
Condenser height adjuster 29, 31
Condenser holder 31
Conoscopy 75
Conoscopy module 23, 48
Contrast method 22, 70, 81
Correction for vision problems 65

D

Dark field 71, 81
Dark field stop 71
Defined function key 50
DIC prisms 42
Differential interference contrast 78
Dimensions 28
Display 49
Disposal 21

E

ebq 100 supply unit 39, 61
Electrical safety 13
ErgoModule 43
Eyepiece 24

F

Field diaphragm 24, 49, 69
Filter block changer 25
Filter cube 40
Filter cube ICR 82, 83
FIM (fluorescence intensity manager) 22, 69
Fluorescence 80
Fluorescence intensity 80
Focus finder 47
Focusing 51, 63
Focusing telescope 54
Function buttons 45, 50

G

Gas discharge lamps 37, 38

H

Halogen lamp 36, 57
HC P tube optics 68, 76
Hg 50 burner 38
Hg 50 W mercury lamp 58
Hg 100 W and Xe 75 W mercury lamps 59

I

Identification sheet 33, 40, 70, 71, 80, 81, 82, 83
Illumination manager 22
Immersion oil 65, 88
Incident light axis 22
Incident light polarizers 41
Incident light shutter 49
Incident light turret disk 40
Initialization 48
Interference contrast 83
Intermediate systems 29

K

Köhler illumination 31, 51

L

Lamp housing 25
Lamp housing 106 z 35, 56
Lamp housing 107/2 55
Lamp housing receptacle 34, 36, 39
Laser coupling 15
LAS software package 23, 26
Light intensity 49
Light sources 34, 55, 69
Light sources for the transmitted light axis 33
Locking pin 40

M

Magnification changer 23, 68
Mirror housing 43
Motorized analyzer 42
Motorized polarizer 41

O

Objective aperture 71
Objective centering 66
Objective prism slide 83
Objectives 65
Objective turret 22, 24
Operating buttons 24
Optical character 76

P

Phase contrast 54, 70
Polarization 71, 82
Polarizer 71, 72, 82
Polarizer L/ICR 41
Polarizer R/ICR 41
Polarizer R/P 41
Pol attachable mechanical stage 62
Power supply connection 44
Prism slot 83

15. Index

Q

Quartz plate 68
Quartz wedge 73

R

Reflector cube 40
Reflector cube for lamp adjustment 55
Reset function 45
Rotating polarizer 41
Rotating stage 30, 62

S

Shutter 80
Specimen holder 29
Specimen stage 29

T

Transmitted light and incident light analyzer 42
Transmitted light axis 22
Transmitted light filter 25
Transmitted light / incident light toggle 24
Transmitted light polarizer ICT/P 41
Transmitted light shutter 49
Transport 20, 27

V

Variable function buttons 24, 25, 45, 50

W

Weight 27, 28

X

Xe 75 burner 38

16. EU Declaration of Conformity

To download the EU Declaration of Conformity for your product use the link

[http://www.leica-microsystems.com/
products/light-microscopes/life-science-research/upright-microscopes](http://www.leica-microsystems.com/products/light-microscopes/life-science-research/upright-microscopes)

Select the type of your microscope and switch to the Download page.

- Administrative Measures on the Control of Pollution Caused by Electronic Information Products -

部件名称 Name of the part	有毒有害物质或元素 Hazardous substances					
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)
印刷电路板 printed circuit boards	X	0	0	0	0	0
电子元器件 electronic components	X	0	0	0	0	0
机械部件 mechanical parts	X	0	0	X	0	0
光学元器件 optical components	X	0	X	0	0	0
电缆 cables	0	0	0	0	X	X
光源 light sources	0	X	0	0	0	0

- o : 表示该有毒有害物质在该部件中的含量均在SJ/T 11363-2006 标准规定的限量要求以下。
Indicates that the concentration of the hazardous substance in all materials in the parts is below the relevant threshold of the SJ/T 11363-2006 standard.
- x : 表示该有毒有害物质至少在该部件的某一材料中的含量超出SJ/T 11363-2006 标准规定的限量要求。
Indicates that the concentration of the hazardous substance of at least one of all materials in the parts is above the relevant threshold of the SJ/T 11363-2006 standard.

Note: The actual product may or may not include in all the part types listed above

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Copyright © Leica Microsystems CMS GmbH • Ernst-Leitz-Straße • 35578 Wetzlar • Germany 2013 • Tel. (06441)29-0 • Fax (06441)29-2599 LEICA and the Leica logos are registered trademarks of Leica IR GmbH.