From Eye to Insight

LEICA DM8000 M

LEICA DM12000 M

Instructions for use
Gebrauchsanweisung
Mode d'emploi

Leica Microsystems CMS GmbH. Instructions for use, 11933916, Revision 1.2, 2018-05-15

From Eye to Insight

LEICA DM8000 M

LEICA DM12000 M

Instructions for use

Leica Microsystems CMS GmbH. Instructions for use, 11933916, Revision 1.2, 2018-05-15

Copyrights

Copyrights

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. We are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information included in this manual may be changed without prior notice.

Revision 1.2, published 2018-05-15 by:

Leica Microsystems CMS GmbH

Ernst Leitz-Straße 17-37

D-35578 Wetzlar (Germany)

<http://www.leica-microsystems.com>

Responsible for contents:

Marketing CMS

Contents

1. Important Notes about this Manual.....	7	7. Instrument Setup.....	35
1.1 Text symbols, pictograms and their meanings.....	7	7.1 Functional principle.....	35
2. Function of the Microscopes.....	9	7.2 Switching on the microscope	39
3. Safety Notes.....	10	7.3 Front operating panel.....	39
3.1 General Safety Notes.....	10	7.4 Variable function keys on the stand.....	39
3.2 Electrical safety.....	10	7.5 XY joystick.....	40
3.3 Transport und Lagerung	12	7.6 Leica SmartMove	41
3.4 Type plates.....	13	7.7 Leica STP6000	41
3.5 Disposal.....	14	7.8 Adjusting the aperture diaphragm in reflected light.....	42
4. Overview of the Instrument.....	15	7.9 Condenser height adjustment	43
5. Unpacking	21	8. Operation.....	44
6. Assembling the Microscope.....	26	8.1 Switching on the microscope	44
6.1 Inspection stage.....	26	8.2 Stages and object displacement	44
6.2 Condenser.....	27	8.2.1 Manual stage.....	44
6.3 Tube and eyepieces.....	28	8.2.2 Motorized stage	45
6.4 Objectives.....	28	8.3 Toggling between reflected light/transmitted light	46
6.5 Illumination.....	29	8.4 Focusing.....	46
6.6 Equipping the reflected light reflector turret	29	8.5 Tubes	49
6.7 Polarizer and analyzer (reflected light) ...	30	8.6 Eyepieces.....	50
6.8 Filter slide.....	30	8.7 Objectives.....	51
6.9 Objective prism slider.....	30	8.7.1 MACRO overview magnification ...	52
6.10 Optional accessories.....	31	8.8 Light sources.....	53
6.10.1 Camera	31	8.9 Aperture diaphragm.....	54
6.10.2 Footswitch.....	31	8.10 Polarizer and analyzer	55
6.10.3 IR illumination	32	8.11 Selecting the contrast method.....	55
6.11 Wiring.....	33	8.12 Objective prism slider	56
		8.13 Filter slide.....	56

Contents

9. Contrast methods	57	10. Care and Maintenance	64
9.1 Reflected light.....	57	10.1 Dust cover.....	64
9.1.1 Bright field (RL).....	57	10.2 Cleaning	64
9.1.2 Dark field (RL)	58	10.3 Handling acids and bases.....	65
9.1.3 Polarization (RL)	58	10.4 Maintenance interval.....	65
9.1.4 Interference contrast (RL).....	59	11. Essential Wear and Spare Parts	66
9.1.5 Oblique illumination (RL).....	59	12. Abbreviation Glossary.....	67
9.1.6 UV mode	60	13. Index	68
9.2 Transmitted light.....	61	14. EC Declaration of Conformity	70
9.2.1 Bright field (TL)	61		
9.2.2 Polarization (TL)	61		
9.2.3 Infrared* (TL)	62		
9.3 Combination methods.....	63		
9.3.1 Bright field (TL)/ Bright field (RL)...	63		
9.3.2 Bright field (TL)/ Oblique illumination (RL).....	63		

1. Important Notes about this Manual

Caution!

This instruction manual is an essential component of the microscope, and must be read carefully before the microscope is assembled, put into operation or used.

This instruction manual contains important instructions and information for the operational safety and maintenance of the microscope and

1.1 Text symbols, pictograms and their meanings

(1.2)

→ p.20

accessories. Therefore, it must be kept in a safe place.

Separate manuals are provided for operating the Leica STP6000 and the Leica Application Suite (LAS) software.

Numbers in brackets, e.g. (1.2), refer to illustrations, in the example to Fig. 1, item 2.

Numbers with an arrow, e.g. → p.20, refer to a specific page of the manual.

Caution!

Special safety instructions within this manual are indicated with the triangle symbol shown here, and have a gray background.

Warning of hot surface!

Caution! Hot surface. May cause burn. Allow for sufficient time to cool before handling or opening.


1. Important Notes about this Manual

Ground/Earth terminal!

Caution! Optical radiation, UV

Caution! UV Radiation
Switch power off before opening!

Warning of hand injury!

Caution! The microscope and accessories can be damaged when operated incorrectly.

Instructions on disposing of the microscope, accessory components and consumables.

Explanatory note.

*

Item not included in all configurations (optional).

2. Function of the Microscopes

The Leica DM8000 M and DM12000 M microscopes, to which this instruction manual belongs, are intended for industrial routine and research applications in microelectronics and other fields. In addition to examination and process control of wafers or LCD displays, cleanliness analysis of mechanical parts is carried out using filter analysis.

The microscopes specified above conform to the EC machinery directive 2006/42/EC and the Council Directive 2004/108/EC concerning electromagnetic compatibility for use in an industrial environment.

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from, using the microscope for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In such cases, the Declaration of Conformity shall be invalid.

3. Safety Notes

3. Safety Notes

3.1 General Safety Notes

These safety class 1 devices were built and tested in accordance with the safety regulations for electrical measuring, control, regulating and laboratory devices in accordance with EN 61010-1:2002-08,
IEC 61010-1:2001-02,
UL 61010-1:2004,
EN 61326-1:2006,
EN ISO 14121-1:2007-12.

They also comply with EN 62471-2008, Photobiological safety of lamps and lamp systems, and belongs to risk class 1 (low risk).

In addition, the instruments conform to the SEMI standard directives:

S2-0703: Safety Guidelines for Semiconductor Manufacturing Equipment and

S8-1103: Safety Guidelines for Ergonomics Engineering of Semiconductor Manufacturing Equipment

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this instruction manual.

Caution!

The instruments and accessories described in this manual have been safety-tested and checked for possible hazards.

The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the device is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the device or noncompliant use shall void all rights to any warranty claims and void product liability!

3.2 Electrical safety

General specifications

For indoor use only.

Supply voltage: 100-240 V~

Frequency: 50/60 Hz

P max: 100 W

Fuse (power supply): T 3.15 A 250 V

Ambient temperature
without IR illumination: 15-35°C

with optional
IR illumination: 15-28°C

Relative humidity: max. 80% to 30°C

Overvoltage category: II

Pollution degree: 2

**Caution!**

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the instrument, or release of the protective conductor terminal, can cause the instrument to become hazardous. Intentional ground interruption is not permitted!

**Caution!**

By connecting it to the grounding connection (earth screw at the back of the microscope) ancillary equipment with its own and/or extra power supply may be brought to the same ground wire potential. For connections without a ground connector, Leica Service must be consulted.

**Caution!**

Never use any fuses as replacements other than those of the types and the current ratings listed here. Using other fuses or bridging the fuse holder is not permitted. The use of incorrect fuses may result in a fire hazard.

**Caution!**

The microscope's electrical accessory components are not protected against water. Water can cause electric shock.

**Caution!**

Protect the microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.

Operating temperature: 15-35°C without IR illumination, 15-28°C with optional IR illumination.

**Caution!**

Per definition, the connection between power supply cable and appliance inlet is the power separator of this device. Thus, the user is responsible for sufficient clearance and accessibility of this area at all times.

**Caution!**

Do not use the microscope in altitudes exceeding 2000 m ASL/NL.

**Caution!**

Only use original power cord, or alternate VDE- / HAR-certified power cords which meet the minimum requirements of 3x0,75mm² and 10A/250V.

3. Safety Notes

Caution!

Before exchanging the fuses or lamps, be absolutely certain to switch off the main power switch and remove the power cable.

Caution!

For Leica STP6000: Touch the touch screen using your finger only. Never use a pen or other hard, sharp or pointed objects.

3.3 Transport and storage

Caution!

Fully equipped, the microscope Leica DM8000 M weighs approx. 41 kg, the microscope Leica DM12000 M approx. 52 kg. For transportation, the user has to take care of the corresponding actions. See also p. 23.


Caution!

Transport and storage in a range of -25° – $+70^{\circ}\text{C}$ and at a humidity not exceeding 80%.

Fig. 1 Safety label
1 Warning UV radiation

Fig. 2 Safety label
1 Warning UV radiation
2 Caution! Optical radiation, UV

3.4 Type plates

Leica DM8000 M Type plates

Leica Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar

TYPE DM8000 M 11888747 BZ:00
SN 123456
10 / 10 11888XXX BZ:00

U 100V- 240V~ 50/60 Hz
P max 100W
T 3.15A 250V

Made in Singapore

CE

Leica Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar

TYPE DM8000 M 11888748 BZ:00
SN 123456
10 / 10 11888XXX BZ:00

U 100V- 240V~ 50/60 Hz
P max 100W
T 3.15A 250V

Made in Singapore

CE

Leica Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar

TYPE DM8000 M 11888749 BZ:00
SN 123456
10 / 10 11888XXX BZ:00

U 100V- 240V~ 50/60 Hz
P max 100W
T 3.15A 250V

Made in Singapore

CE

Leica Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar

TYPE DM8000 M 11888750 BZ:00
SN 123456
10 / 10 11888XXX BZ:00

U 100V- 240V~ 50/60 Hz
P max 100W
T 3.15A 250V

Made in Singapore

CE

Leica Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar

TYPE DM8000 M 11888761 BZ:00
SN 123456
10 / 10 11888XXX BZ:00

U 100V- 240V~ 50/60 Hz
P max 100W
T 3.15A 250V

Made in Singapore

CE

Leica Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar

TYPE DM8000 M 11888762 BZ:00
SN 123456
10 / 10 11888XXX BZ:00





U 100V- 240V~ 50/60 Hz
P max 100W
T 3.15A 250V

Made in Singapore

CE

3. Safety Notes

Leica DM12000 M Type plates

3.5 Disposal

Once the product has reached the end of its service life, please contact Leica Service or Sales about disposal.

Please observe and comply with the national and federal laws and regulations that are equivalent to EC directives such as WEEE.

Note!

Like all electronic devices, the microscope, its accessory components and consumables must never be disposed of with general household waste.

4. Overview of the Instrument

Specification	Leica DM8000 M/12000 M
Contrast methods	• Reflected light: BF, DF, OBL, POL, DIC (ICR), UV, OUV • Transmitted light: BF, POL, IR • Combination methods: TL BF/RL BF, TL BF/RL OBL
Transmitted light axis	• Integrated LED Illumination • Motorized adjustment of brightness and aperture diaphragm • Condenser focal intercept S10/0.7 • Manual polarizer • IR illumination (optional, retrofittable)
Reflected light axis	• Reflected illuminator VIS or VIS/UV • 2 variants with different height • 2 variants: manual or motorized • 4x reflector turret (motorized for respective variant) • MACRO overview magnification (motorized variant) • 2-position manual filter slide • Holder for manual polarizer slider • Holder for manual analyzer slider • LED status display (motorized variant) • Centerable aperture diaphragm • Fixed field diaphragm • Focus finder
Tube	• Manual operation • With image correction • With variable viewing angle • Optionally with fixed or variable viewing angle • Optionally with 2 camera ports for VIS and UV • Manual change of observation ports
Magnification changer (optional)	• Manual • Magnification levels: 1x; 1.5x; 2x
Objective nosepiece	• Motorized, 6-position for BD objectives with M32 thread • Absolutely encoded • Holder for manual DIC objective prism slider

4. Overview of the Instrument

Specification	Leica DM8000 M/12000 M
X/Y stage	- DM8000 M: 8" stage DM12000 M: 12" stage - Manual with rapid adjustment - Positioning range 8" stage in reflected light and transmitted light: 202 mm x 202 mm or positioning range 12" stage in reflected light: 302 mm x 302 mm in transmitted light: 202 mm x 202 mm - or - Motorized with 2-phase stepper motor - Positioning range 8" stage in reflected light and transmitted light: 200 mm x 200 mm or positioning range 12" stage in reflected light: 300 mm x 300 mm in transmitted light: 202 mm x 202 mm - Maximum travel speed: 240 mm/s with 4 mm spindle - Optional inserts: metal plate, glass plate, wafer holder, vacuum wafer holder, mask holder
Focus	- Manual - Focus handwheel for coarse and fine focusing - Height adjustment - Speed Switch (optional) - Adjustable focus threshold and movement rate - or - Motorized - Focus handwheel - Control of the mot. z-drive via Leica SmartMove or Leica STP6000 - Parfocality adjustment - Focus finder for the reflected light axis - Optional: Laser autofocus

Specification	Leica DM8000 M/12000 M
Controls	• Operating buttons on the stand • Additional variable multifunction keys • Focus wheels • Joystick: Control element for x,y control • Leica SmartMove: Ergonomic control element for x,y,z control with additional variable function keys • Leica STP6000: Ergonomic control element for x,y,z control with 11 additional variable function keys and touch-sensitive LC display • Double footswitch for objective change
Computer Interface	• 2 x USB 2.0
Software	• Leica Application Suite (LAS) for Windows™ 7, XP for: • Microscope and camera configuration • Microscope and camera control • Image acquisition and storage of images • INSPEC IS for Windows™ 7, XP for: • Microscope and camera control • Front-end wafer inspection

4. Overview of the Instrument


Fig. 3 Left side of the stand, Leica DM8000 M

- | | |
|--|--|
| 1 Camera | 7 Swivel lever for transmitted light polarizer |
| 2 Tube port (VIS/UV) | 8 Variable function keys (preset at the factory) |
| 3 Tube | 9 Focus handwheel |
| 4 Centering of aperture diaphragm | 10 Condenser height adjuster |
| 5 6-position objective nosepiece, motorized,
for BD objectives (M32 thread) | 11 Objective prism slider |
| 6 Motorized inspection stage with wafer chuck | 12 Receptacle for analyzer |
| | 13 Joystick for xy movement |


Fig. 4 Right side of the stand, Leica DM8000 M

- | | |
|--|---------------------------------------|
| 1 Eyepiece tube with eyepiece | 10 Coaxial drive |
| 2 LED status display | 11 Locking lever for rapid adjustment |
| 3 Manual inspection stage with wafer chuck | 12 Rapid stage adjustment |
| 4 Aperture diaphragm adjustment | 13 On/off switch |
| 5 Objective changer buttons | 14 Receptacle for polarizer |
| 6 Variable function keys (preset at the factory) | 15 Focus finder |
| 7 Brightness adjustment | 16 Filter slide |
| 8 Focus handwheel | 17 Control slider for beam splitting |
| 9 Variable function keys (preset at the factory) | |

4. Overview of the Instrument

Fig. 5 Leica STP6000

- 1 Touch screen
- 2 Information key
- 3,4 Variable function keys, user-programmable
- 5 Fine focus adjustment
- 6 Coarse focus adjustment
- 7 Movement in y-direction
- 8 Movement in x-direction


Fig. 6 Leica SmartMove

- 1 Movement in x-direction
- 2 Movement in y-direction
- 3 Focus setting
- 4 Individual setting of the knob height position
- 5 Variable function keys (preset at the factory)



Fig. 7 XY joystick

- 1 Movement in xy-direction
- 2 Fastening the x-adjustment in place
- 3 Fastening the y-adjustment in place
- 4 Variable function keys (preset at the factory)

5. Unpacking

Depending on your instrument configuration, the standard delivery contains the following components:

- Stand with integrated reflected light axis and objective turret
- Inspection stage
- Power cable and PC connecting cable
- CD with Leica Application Suite (LAS) software package
- Instructions and list of microscope presets ("Identification Sheet")
- Tube
- Eyepieces
- Objectives
- Condenser head
- Assembly tools
- Additional microscope accessories such as filters, etc. depending on product configuration

First, carefully remove all components from the transportation and packaging materials.

- Lift the stand out of the carton by its carrying handles (8.1).

! Caution!

When lifting the instrument out of the packaging, and when transporting it to its installation site, carry it by its handles only.


Caution!

Note that the stand weighs approx. 29 kg (Leica DM8000 M) or 31 kg (Leica DM12000 M). Two people are required to lift it out of its packaging and carry it to its installation site.

Later, the carrying handles can be unscrewed from the stand. Then, seal the threaded holes using seal plugs.

Fig. 8 Stand with carrying handles

1 Screw-on carrying handles

5. Unpacking

Note:

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces in a short time. See the chapter on "Care of the Microscope" → p. 64 for additional instructions.

Allowable ambient conditions

Temperature (without IR illumination) 15–35°C

Temperature (with IR illumination) 15–28°C

Relative humidity maximum 80% up to 30°C

Microscopes in warm and warm-damp climatic zones require special care in order to prevent fungus contamination.

See the chapter on "Care of the Microscope" → p. 64 for additional instructions.

Caution!

Do not connect the microscope and peripherals to the power supply at this point under any circumstances!

Caution!

Electrical components must be assembled at least 10 cm from the wall and away from flammable substances.

Installation location

Work with the microscope should be performed in a dust-free room, which is free of oil vapors and other chemical vapors, as well as extreme humidity. At the workstation, large temperature fluctuations, direct sunlight and vibrations should be avoided. These conditions can distort measurements and micrographic images.

Caution!

When installing the microscope, make sure the power inlet is freely accessible so that the instrument can be quickly disconnected from the mains if necessary.

Transport

For shipping or transporting the microscope and its accessory components, the original packaging should be used.

As a precaution to prevent damage from vibrations, the following components should be disassembled and packaged separately:

- Remove the tube.
- Unscrew the objectives.
- Remove the condenser head.
- Remove the inspection stage.
- Remove all moving or loose parts.

Weight

The weight of the microscope depends on the particular equipment.

Fully equipped, an approximate value is about 41 kg (Leica DM8000 M) or 52 kg (Leica DM12000 M). The microscope itself weighs approx. 29 kg (Leica DM8000 M) or 31 kg (Leica DM12000 M).

For transportation, the user has to take care of the corresponding actions. Two people are required to carry the microscope.



Caution!

For transporting it is essential to remove all components listed under "Transport"!

5. Unpacking

Dimensions (specifications in mm)

6. Assembly

6. Assembling the Microscope

The microscope components are logically assembled in this order:

- Inspection stage
- Condenser head (optional)
- Stage insert (i.e. metal or glass plate or wafer holder)
- Tube
- Eyepieces
- Objectives

Only a few commonly used screwdrivers and keys are necessary for assembly; these are included in the delivery package.

When using intermediate systems and optical accessories, the sequence may vary.

In this case, read chapter

"6.10 Optional accessories" → p. 31

6.1 Inspection stage

Caution!

Never install objectives before assembling the stage.

- Attach the stage from above and fasten it using the 6 screws (9.1).


Caution!

The inspection stage must be installed by personnel authorized by Leica CMS only.

The reflector turret is already equipped at the factory according to the microscope configuration. For an upgrade of reflector cubes, contact Technical Service.

Fig. 9 Assembling the inspection stage

1 Fastening screws

- Place the insert (metal plate, glass plate, wafer or mask holder) into the stage and fasten it using the 4 screws (10.1).

Caution:

To set the focus position and the lower threshold, see 8.4 Focusing → p. 47.

- For the motorized stage only:
Insert the connections (11.1) and (11.2) into the rear side of the motorized stage and connect the stage cable to the rear side of the stand using the connection (24.5) → p. 34.

Note:

Please observe the provided documents of the scanning stage manufacturer!

6.2 Condenser

- Screw the condenser head into the condenser.

Note:

Before using the microscope, you have to check the condenser height setting, see 7.9 Condenser height adjustment → p. 43.


Fig. 10 Assembly of stage insert

1 Fastening screws

Fig. 11 Motorized stage

1 Connection for x-direction
2 Connection for y-direction

6. Assembly

6.3 Tube and eyepieces

The tube is installed directly on the stand. Optionally, a magnification changer or laser autofocus can be installed between the stand and tube.

- Fasten the tube using the 2 side clamping screws (12.1).
- The eyepieces are inserted into the eyepiece tubes on the tube.

Note:

We recommend teaching in eyepieces that are not included in the scope of delivery using the Leica Application Suite (LAS) software - module: Set-Up.

This ensures that the total magnification specified in the LAS software or on the Leica STP6000 is correct.

6.4 Objectives

The receptacles on the objective nosepiece are numbered (Fig. 13). Based on your equipment, the individual objectives have already been assigned to specific positions at the factory.

For details on the exact positions of the objectives, please refer to the enclosed identification sheet.

Caution:

Close vacant threads in the nosepiece with dust protection caps!

Note:

For instruments with motorized focus drive, we recommend carrying out a parfocality adjustment using the Leica Application Suite (LAS) software - module: Fine Tuning.

Fig. 12 Assembling the tube

1 Fastening screws

Fig. 13

Objective nosepiece with labeled objective receptacles

To prevent contamination on the wafer, you can insert a contamination guard (14.1) above the objective turret. This must be done before you screw the objectives into the objective nosepiece. The transparent Plexiglas contamination guard is pushed onto the bottom of the reflected light axis via the objective nosepiece and screwed on hand-tight (without exerting great force) using the 3 screws provided.

6.5 Illumination

Note:

The microscope has a built-in LED illumination for reflected light (VIS and UV) and transmitted light. The service life of the LED is about 25000 hours. If, despite this, it is necessary to change the LED, this must be done by Leica Microsystems CMS GmbH Technical Service only. See chapter "6.10 Optional accessories" → p. 32 for assembly of the optional IR illumination.

Fig. 14

1 Contamination guard

Caution!

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation). Therefore, lamps have to be operated in closed housings.

6.6 Equipping the reflected light reflector turret

The positions in the reflector turret are numbered. Depending on your equipment, the individual reflector cubes have already been assigned to specific positions at the factory. For details, check the identification sheet included with your order.

Using the stickers provided, label the filter positions on the front side of the stand (15.1).

To insert the reflector cubes, please contact Technical Service.

Note:

Reflector cubes that are not included in the scope of delivery must be taught in using the Leica Application Suite (LAS) software - module: Set-Up.

Fig. 15

1 Filter positions of the reflected light reflector turret, for motorized variant with LED status display

6. Assembly

6.7 Polarizer and analyzer (reflected light)

- Remove the plug cap on the right side of the reflected light axis (Fig. 16.1).
- Insert the polarizer into the receptacle until it locks into place.
- Remove the plug cap on the left side of the stand.
- Insert the analyzer into the receptacle until it latches in place (Fig. 17.1).

6.8 Filter slide

- For activating and deactivating filters, the two sliders (16.2) on the right side of the stand are provided.

Various filter inserts (e.g. yellow filters) with a diameter of 25 mm are available.

6.9 Objective prism slider

- The objective prism slider, used for Differential Interference Contrast (DIC), is inserted into the holder (16.3).
- For fine adjustment use the adjusting screw (17.2) above the objective nosepiece.

Fig. 16

- 1 Polarizer mount
- 2 Filter slide
- 3 Slot for objective prism slider

Fig. 17

- 1 Analyzer slot
- 2 Fine adjustment for objective prism slider

6.10 Optional accessories

6.10.1 Camera

Up to two cameras can be adapted via a tube attachment with two camera ports.

- Attach the tube attachment to the top port of the tube and fasten it tightly using the side clamping screw (18.1).
- Screw on the camera(s).

Note:

The use of a tube attachment should be taught in using the Leica Application Suite (LAS) software - module: Set-Up.

- When using a Type DC digital camera, the camera is connected to the PCI card of the PC.

Note:

Observe the separate instruction manual for the digital camera!

Note:

When upgrading to the tube with dual attachment, both ports should be centered by Technical Service.

6.10.2 Footswitch

To swing in the objectives, a double footswitch (Fig. 19) can be connected to the rear side of the stand (24.3).

Fig. 18

1 Fastening screw

Fig. 19

Double footswitch

6. Assembly

6.10.3 IR illumination

When silicon wafers are inspected under IR light they appear transparent. This effect is used, for example, for aligning single layers on top of each other. As IR light cannot be seen with the naked eye, a special camera is used in combination with LAS software and a monitor. This monochromatic screen image enables different planes in the silicon wafer to be focused.

If your system configuration already contains the IR illumination unit on delivery, it has already been assembled at the factory and only has to be connected up to the external power supply (23.3).

The IR illumination is installed before the stage is assembled and the objectives are screwed in. To retrofit an IR illumination unit therefore, the objectives and the stage have to be disassembled first.

The IR illumination comprises

- the IR lamphousing with the lamp cord
- the external power supply
- the TL/IR wafer holder
- the UV/IR filter kit for the DFC camera

Note:

The UV/IR filter kit may only be retrofitted to the DFC camera by the Technical Service of Leica Microsystems CMS GmbH.

- Insert the IR lamphousing into the stage from above (20.1).

Fig. 20 Assembly of IR illumination

1 IR lamp housing receptacle



Fig. 21 Assembly of IR illumination

1 IR lamp housing

- Slacken the three screws (22.1) on the base plate of the stand. Do not remove the screws completely, and only loosen the base plate, do not remove!
- Guide the lamp cord between the bottom of the microscope and the base plate toward the front.
- Plug the lamp cord into the external power supply (23.3).

Fig. 22 Microscope stand from underneath
1 Fixing screws of the base plate



Fig. 23 External power supply
for IR illumination

- 1 On/off switch
- 2 Connection power supply
- 3 Connection lamp cable
- 4 Earth screw

6. Assembly

6.11 Wiring

Depending on the configuration, various operating elements have to be connected to the rear side of the stand. To do so, use the cables provided.


- To connect a computer, connect the USB interface of the PC to the USB interface (24.1) on the rear side of the stand.
- The second USB port (24.2) is provided for the Leica STP6000.
- The footswitch is connected to the port (24.3).
- When using a motorized stage, the stage cable is connected to the port (24.5).
- Connect the XY joystick to the port (24.7).
- The RS232 interface (24.8) is intended for connecting the Leica SmartMove.
- After completing the assembly work, connect the microscope to the power supply using the power cable supplied (24.6).

For the optional IR illumination:

- Connect the external lamp power supply to the power supply using the power cable supplied (23.2).

Fig. 24 Rear view of the stand

- 1 USB port, e.g. for PC
- 2 USB port, e.g. for Leica STP6000
- 3 Connection for footswitch
- 4 Earth screw
- 5 Motorized stage connection
- 6 Power supply connection
- 7 Joystick for XY movement
- 8 RS232 port for Leica SmartMove
- 9 3 x I²C port
- 10 Maintenance indicator light (p. 52)

7. Instrument Setup

7.1 Functional principle

Based on its intelligent automation, the Leica DM8000 M or DM12000 M can be operated using various controls.

1. Intelligent automation

- Switch between different contrast methods at the touch of a button.
- The microscope recognizes the selected objective and associated contrast method. Therefore, the values for intensity (INT) and aperture diaphragm (AP) are always set properly.
- The INT and AP values are always based on the currently activated light axis (transmitted light or reflected light).
- The INT and AP values can be adjusted individually. Manual adjustments overwrite the previous settings. The current setting is stored and remains in effect even after the microscope is switched off.

2. Controls

- Footswitch for objective change
- XY joystick for stage control
- Leica SmartMove for stage and focus control
- Leica STP6000: For stage and focus control and control of the microscope functions via a touchscreen
- Setting wheels on the stand for INT and AP
- Variable function keys on the stand, XY joystick, Leica SmartMove, Leica STP6000
These function buttons have functions suitable to the configuration of your microscope assigned to them at the factory. The functions can be reprogrammed and/or adapted to your specific requirements, however.
Switching between transmitted light and reflected light axis possible by pressing a button
- Complete control of microscope and camera via software
(Leica Application Suite, Leica application software, INSPEC IS inspection software)

7. Instrument Setup

Note: (Reset function)

The table on the following page provides an overview of which microscope components can be controlled via the respective controls.

The microscope can be reset to its factory default programming:

- Press the two objective changer buttons (25.1) and the variable function key (25.2) on the front operating panel simultaneously and switch off the microscope.
- Switch on the power for the microscope.
- Hold the buttons until the initialization is complete.

Fig. 25 Front operating panel

- 1 Objective changer buttons
- 2 Variable function key

7. Instrument Setup

Function	Touchscreen Leica STP6000	Setting wheels or fixed function keys Stativ	Variable function keys stand Leica joystick, Leica STP6000	Rotary knobs Leica SmartMove or Leica STP6000	Software LAS
Switching contrast methods	+	-	+	-	+
Switching between transmitted light/reflected light axis	+	-	+	-	+
Positioning objectives 1) Teaching in parfocality	+ -	+ -	+ -	- -	+ +
Illumination manager	+	+	+	-	+
Magnification changer*	manually only				
Switch observation port	manually only				
Focusing 2) Setting the thresholds Setting the thresholds Changing step size (Coarse/Fine)	- + + +	- - - -	- + + +	+ - - -	+ + + +
XY stage positioning ³⁾ (motorized stage only) Changing step size Stage positions (store/position)	- + +	- - -	- + -	+ - -	+ + +

+ always possible

- not possible

¹⁾ Alternative for objective change via footswitch

²⁾ Alternative for focusing via focus handwheels

³⁾ Alternative for stage positioning via joystick

7. Instrument Setup

Possible assignments of the variable function keys

For Leica DM8000/12000 M:

Function button	Function
BF (RL)	Reflected light bright field
ICR (RL)	Interference contrast reflected light
DF (RL)	Dark field reflected light
POL (RL)	Polarization reflected light
UV (RL)	UV reflected light
OBL (RL)	Oblique reflected light illumination
MACRO (RL)	MACRO UV reflected light
CHANGE RL ●	Switches through all reflected light contrast methods
INT ↑	Increase brightness (intensity)
INT ↓	Decrease brightness (intensity)
AP ↑	Open the aperture diaphragm (reflected light)
AP ↓	Closes the aperture diaphragm (reflected light)
BF	Bright field transmitted light
RL/TL	Switch between transmitted light and reflected light
BF-BF	Bright field reflected light - bright field transmitted light combination methods
OBJn	Select objective at position n
CHANGE OBJ CW	Switch through objectives clockwise
CHANGE OBJ CCW	Switch through objectives counterclockwise
Z FINE	Enable fine focus
Z COARSE	Enable coarse focus
MEMn	Move to memory position n
CAM	Camera button - capture image
-	Function key not assigned

7.2 Switching on the microscope

- Switch on the microscope at the on/off switch (27.1). All motorized microscope components will then run through an initialization phase.

Note:

If you have connected a PC, switch on the microscope first, then the computer.

After you switch the microscope on, the reflected light axis is active first. The LED of the **TL/RL** variable function key is illuminated. To toggle to the transmitted light axis, press the **TL/RL** key. The LED is not illuminated.

The aperture diaphragm is already precentered at the factory. It may be necessary to correct the centering after the microscope has been transported and assembled. Before performing the required steps, please familiarize yourself with the front operating panel and the controls.

7.3 Front operating panel

The following assignments are made for all stand variants: the left setting wheel for the aperture diaphragm (26.1), the right setting wheel for the light intensity (26.2) and the objective changer buttons (26.3).

Other functions, such as certain contrast methods, can also be assigned to the variable function keys (26.4).

7.4 Variable function keys on the stand

In addition to the variable function keys on the front operating panel, different versions also have 2 other freely assignable function keys on the left and right side of the stand, each behind the focus handwheel. The variable function keys are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on key assignments, please refer to the included identification sheet.

For information on the abbreviations used, please refer to the list → p. 38.

Fig. 26 Front operating panel

- Aperture diaphragm adjustment
- Brightness adjustment
- Objective changer buttons
- Variable function keys with integrated LED (factory preassigned)

Fig. 27

- On/off switch
- Variable function keys (preset at the factory)

7. Instrument Setup

Note:

Other variable function keys can be defined for the XY joystick, the Leica SmartMove and the Leica STP6000.

Note:

The Leica Application Suite (LAS) module "Set-Up" is required for changing the key assignments.

7.5 XY joystick


To move the motorized inspection stage, you can use the XY joystick (Fig. 28).

In addition, variable function keys (28.4) can be assigned. The variable function keys are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on key assignments, please refer to the included identification sheet.

For information on the abbreviations used, please refer to the list → p. 38.

Fig. 28 XY joystick

- 1 Movement in XY direction
- 2 Fastening the x-adjustment in place
- 3 Fastening the y-adjustment in place
- 4 Variable function keys (preset at the factory)

7.6 Leica SmartMove

The rotary knobs (29.1, 29.2) move the motorized inspection stage in the x and y-direction. The image is focused using the rotary knob (29.3).

You can adjust the height of the rotary knobs by turning the wheel (29.4) to a working height that is convenient for you.

In addition, variable function keys (29.5) can be assigned.

The variable function keys are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on key assignments, please refer to the included identification sheet. For information on the abbreviations used, please refer to the list → p. 38.

Note:


The "Set-Up" module of the Leica Application Suite (LAS) software is required for changing the key assignments.

7.7 Leica STP6000

When using the Leica STP6000 external control panel, observe the separate instruction manual.

Fig. 29 Leica SmartMove

- 1 Movement in x-direction
- 2 Movement in y-direction
- 3 Focus setting
- 4 Individual setting of the knob height position
- 5 Variable function keys (preset at the factory)

7. Instrument Setup

7.8 Adjusting the aperture diaphragm in reflected light

The aperture diaphragm is already centered at the factory.

As a result of transport and setup of the stand, however, in some cases centering may be required again. Therefore check the aperture diaphragm centering.

The following procedure is provided for the reflected light-bright field illumination.

- Select an objective with moderate magnification (10x-20x).
- If necessary, activate the reflected light axis by pressing the **TL/RL** button (one of the variable function keys).
- Press the **RL-BF** button to activate the bright field contrast method (one of the variable function buttons).
- Now place a specimen on the stage.
- Focus on the specimen using the Leica Smart-Move, Leica STP6000 or the focus handwheels.
- Adjust the light intensity using the setting wheel (26.2).
- Remove one eyepiece (e.g. right). To make adjusting the aperture diaphragm easier, you can use the focusing telescope (Fig. 31) instead.
- Close the aperture diaphragm with the setting wheel (26.1) until the edge of the diaphragm appears in the exit pupil of the objective (aperture diaphragm plane).
- If the image is not in the middle of the field of view of the exit pupil, move the position of the aperture diaphragm to the center of the exit pupil using the two centering screws (30.1) located on the left side of the stand.
- Open the aperture diaphragm to cover 2/3 of the field of view.

Fig. 30 Adjusting the aperture diaphragm in the reflected light axis

1 Adjusting screws for moving the aperture diaphragm

Fig. 31 Focusing telescope (order number 11505070)

- 1 Adjustable eyelens
- 2 Clamping ring for fixing the focus position

7.9 Condenser height adjustment

Before using the microscope, you have to check the condenser height setting.

- If necessary, activate the transmitted-light axis by pressing the **TL/RL** button (one of the variable function keys).
- Remove one eyepiece.
- Swing in the objective 50x.
- Set the light intensity to a moderate value using the setting wheel (26.2) to avoid the risk of glare.
- Use the 3 mm Allen key provided to adjust the screw (32.1).
Turn the screw until the illuminated area covers the pupil completely. At this point the diameter of the illuminated area is maximized.
- Insert the eyepiece again.

Fig. 32 Condenser height setting

1 Adjusting screw

8. Operation

8. Operation

8.1 Switching on the microscope

- Switch on the microscope at the on/off switch (33.1). All motorized microscope components will then run through an initialization phase.

Note:

If you have connected a PC, switch on the **microscope first**, and **then the computer**.

After you switch the microscope on, the reflected light axis is active first. The LED of the **TL/RL** variable function key is illuminated. To toggle to the transmitted light axis, press the **TL/RL** key. The LED is not illuminated.

In addition, all of the user's most recent illumination settings are restored.

Note: (Reset function)

The microscope can be reset to its factory default programming:

- Press the two objective changer buttons (33.3) and the variable function key (33.2) on the front operating panel simultaneously and switch of the microscope.
- Switch on the power for the microscope.
- Hold the buttons until the initialization is complete.

8.2 Stages and object displacement

Caution!

When moving the stage, particularly the motorized stage, there is a crushing hazard between the stage and stand!

Fig. 33

- 1 On/off switch
- 2 Variable function key (preset at the factory)
- 3 Objective changer buttons

8.2.1 Manual stage

The stage is moved via the rapid adjustment (34.5) or, for fine positioning, via the coaxial drive (34.1, 34.3).

- Before using the coaxial drive, lock the rapid adjustment using the lever (34.4).

8.2.2 Motorized stage

Object displacement

The motorized stage is moved via the joystick (35.1), the rotary knobs (36.1, 36.2) on the Leica SmartMove or Leica STP6000 (37.7, 37.8)

Fig. 34 Inspection stage

- 1 Movement in x-direction
- 2 Torque adjustment in x-direction
- 3 Movement in y-direction
- 4 Locking mechanism for rapid adjustment
- 5 Rapid adjustment

Fig. 36 Leica SmartMove

- 1 Movement in x-direction
- 2 Movement in y-direction
- 3 Focus setting
- 4 Individual setting of the knob height position
- 5 Variable function keys (preset at the factory)


Fig. 35 XY joystick

- 1 Movement in xy-direction
- 2 Fastening the x-adjustment in place
- 3 Fastening the y-adjustment in place
- 4 Variable function keys (preset at the factory)



Fig. 37 Leica STP6000

- 1 Touch screen
- 2 Information key
- 3,4 Variable function keys, user-programmable
- 5 Fine focus adjustment
- 6 Coarse focus adjustment
- 7 Movement in y-direction
- 8 Movement in x-direction

8. Operation

Programming the step sizes

You can change the travel speed of the stage by switching between the step sizes.

If the value **Coarse** is selected, the travel speed is the same for all objectives.

The value **Fine** is adapted to the respective objective.

Simultaneously, the focus speed is also adapted.

► Switching between Coarse and Fine using

- Variable function keys on the stand, Leica SmartMove and joystick
- Leica STP6000
- Leica Application Suite software - module: Acquire

► Storing positions temporarily and move to them using

- Leica STP6000
- Leica Application Suite software - module: Acquire

8.3 Toggling between reflected light/transmitted light

To toggle between the reflected light axis (**RL**) and the transmitted light axis (**TL**), you normally use one of the freely configurable function keys on the front operating panel.

If your microscope is equipped with the TL axis, this key is already predefined at the factory.

Saving and moving to stage positions

Up to 14 stage positions can be stored using the Leica Application Suite software - module: Acquire or using the Leica STP6000. The xy-position is stored, but not the current z-position.

In addition, a load position (Load) and a stage center position can be determined. As well as the xy position, the corresponding z position is then also permanently stored. The load position is the user-defined point at which the wafer is transferred from the wafer loader to the microscope. When the "Load" key is pressed, the stage travels to the pre-defined x, y and z positions. The Stage Center position is the center position of the stage in xy direction, also freely determined beforehand. The stage can be moved to these positions using the LAS software.

► Toggling between the illumination axes using:

- Variable function keys on the stand, SmartMove and joystick
- Leica STP6000
- Leica Application Suite software - module: Acquire

8.4 Focusing

Note:

Leica objectives are parfocal with each other in wide ranges. Small deviations in parfocality can result from screwing in the objectives during assembly. For this purpose, stands with motorized focus offer the ability to teach in the parfocality. We recommend checking the parfocality before setting the thresholds and, if necessary, teaching them in again using the Leica Application Suite software - module: Fine Tuning or in the Acquire module.

Manual focusing

Coarse and Fine Focusing

Focus knobs for coarse and fine focusing are located on both sides of the stand (Fig. 38, 39).

Height adjustment of the focus knobs

! Caution!

During the height adjustment of the focus knobs the stage is moving - the sample is moving towards the objective. So a collision between sample and objective may be possible. Remove the sample before adjustment and select an objective with lower magnification (5x or 10x).

- Defocus the image by moving the stage down with a full turn of the coarse focus handwheel (38.2, 39.2).

- Grasp the right-hand and left-hand focus knobs at the same time and press the knobs gently upward or downward into the desired position.

- Refocus the image.

Speed switch (optional)

Beside the coarse focusing two further speeds are available for fine focusing. We distinguish between the fine and super-fine switching stages. To switch speeds, press the left focus knob to the right or vice versa.

Adjusting the focus threshold

The current position can be set as the focus threshold by locking the knurled wheel (39.3) at the right-hand focus knob. It will then no longer be possible to travel past that position during coarse focusing.

To do so, rotate the knurled wheel clockwise. Rotate it in the opposite direction to release the wheel once again.

Fig. 38 Focus knob

- 1 Fine focusing
- 2 Coarse focusing
- 3 Torque adjuster

Fig. 39 Focus knob

- 1 Fine focusing
- 2 Coarse focusing
- 3 Setting the focus threshold

8. Operation

Setting the resistance

The torque of the focus drive can be adjusted using the knurled wheel (38.3) at the left focus knob.

! Caution!

Ensure that the action is not too light. Otherwise, the stage can drift downward unintentionally.

Motorized focus

Focusing the image

Focus the image using the rotary knob (36.3) on the Leica SmartMove or the rotary knobs (37.5, 37.6) on the Leica STP6000.

When the microscope is switched on, before stage initialization, the z-drive moves to a position that is 4 mm below the last programmed z-position.

Setting thresholds

You can use the Leica Application Suite software to set the lower focus threshold.

In addition, a **focus position** that cannot be exceeded can be set. This position should be defined for the dry objective with the highest magnification level. For all other objectives, it is set automatically, taking into account the parfocality adjustment and the working distance.

► Setting the thresholds using:

- Variable function keys on the Leica SmartMove, Leica STP6000 or joystick
- Leica Application Suite software - module: Fine Tuning

When the thresholds have been set, the LED (40.1) is illuminated. For the Leica STP6000, this is likewise displayed:

► Moving to the thresholds using:

- Variable function keys on the stand, Leica SmartMove and joystick
- Leica STP6000
- Leica Application Suite software - module: Fine Tuning or Acquire module

Note:

When moving to the threshold using the **Focus position** key, you have to press and hold the corresponding key until the position is reached.

Fig. 40 LED status display

- 1 LED: Focus set
- 2 Focus finder

Programming the step sizes

You can toggle between two step sizes, **Fine** and **Coarse**.

The value **Fine** is adapted to the respective objective. The values are predefined logically. You can change the allocation using the Leica Application Suite software - module: Fine Tuning.

If the value **Coarse** is selected, the travel speed is the same for all objectives. **Coarse** corresponds to the maximum speed.

Note:

The allocation of a certain step size to an objective applies not only for the z-drive, but also for defining the stage step size.

► Toggling between Fine and Coarse via

- Variable function keys on the stand, Smart-Move and joystick
- Leica STP6000
- Leica Application Suite software - module: Acquire

Focus finder for reflected light

For surfaces that have only few structures on which focusing is possible or none at all, you can use the focus finder.

- Swing the focus finder into position using the rotary knob (40.2).

- Then, in the field diaphragm plane, focus on the edge of the dot mask that is now visible.
- Swing the focus finder back out.

8.5 Tubes


Note:

Close any unused tube openings, as otherwise stray light can interfere with observation.

Adjusting the interpupillary distance

- Adjust the interpupillary distance of the eyepiece tubes so that a congruent total image is seen (Fig. 41).

Fig. 41 Tube setting
↔ Setting the personal interpupillary distance

8. Operation

Adjusting the viewing angle

- For the Ergo tubes, the viewing angle can be adjusted by tilting the binocular eyepiece (Fig. 42).

Beam splitting

The light division is set manually by pulling out a control bar (43.1).

Available tubes are:

- 11565024 Ergo-Tube with upright image for visible light only
- 11565025 Ergo-Tube with upright image for VIS and UV light

Control Bar		Observation	Camera
VIS		100 %	0 %
CAM		0 %	100 %

- 11565050 Ergo-Tube with upright image for VIS and UV light

Control Bar		Observation	Camera
VIS		50 %	50 %
CAM		0 %	100 %

Fig. 42 Individual settings on the ErgoTube

Note:

When using a dual attachment on the camera port, the light division between both cameras is fixed (50:50).

8.6 Eyepieces

Note:

The eyepiece's aperture protector must be removed during microscopy while wearing eyeglasses, or at least folded back.

We recommend removing bifocals and spectacles with progressive-addition lenses when using the microscope.

- For the adjustable tubes with documentation port, choose the 100% VIS position.

Fig. 43 Beam splitting on the tube
1 Control bar

Eyepieces with inlaid graticule

- Focus the graticule by adjusting the eyelens in the eyepiece.
- Focus on the object through this eyepiece.
- Then, close that eye and focus on the object by adjusting the second ocular only.

Correction for vision problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.
- Then, with your left eye, view the same position of the specimen and rotate the left eyepiece tube until this position is brought into sharp focus. While doing so, do not change the adjustment of the z-position!

Note:

We recommend teaching in eyepieces that are retrofitted or are not included in the scope of delivery using the Leica Application Suite (LAS) software - module: Set-Up. This ensures that the total magnification specified in the LAS software or on the Leica STP6000 is correct.

8.7 Objectives

Changing objectives

The objective change takes place by motorized action via fixed function keys on the front operating panel of the stand (44.3) or the variable function keys on the joystick, Leica SmartMove or Leica STP6000. Alternatively, objectives can also be swung into the beam path by rotating the objective nosepiece manually. During the manual objective change, be sure that the nosepiece locks into place.

Optionally, to change the objectives, you can also connect a double footswitch to the stand.

The positions of the objectives in the objective turret have been specified at the factory and must be observed when installing the objectives. (See Installing objectives → p. 28).

When positioning the objective (manually or by motorized action), the microscope configures the following parameters automatically:

- The optimal setting for the aperture diaphragm
- The light intensity in the respective contrast method

The objective magnification and the total magnification appear in the LAS software or on the Leica STP6000.

8. Operation

► Positioning the objectives using:

- Fixed function keys on the stand
- Variable function keys on the joystick and Leica SmartMove
- Leica STP6000
- Double footswitch
- Leica Application Suite software - module: Acquire
- Manual operation possible

Note:

If objectives are retrofitted, they have to be taught in using the Leica Application Suite software - module: Set-Up. Afterwards, for stands with motorized focus, the parfocality should also be taught in again.

8.7.1 MACRO overview magnification

Only for motorized stands and bright field reflected light

The MACRO overview magnification enables you to illuminate a field with a diameter of approx. 40 mm. This is done using a special objective (MACRO objective) and the MACRO filter cube.

- If necessary, switch to the reflected light axis (RL) by pushing the **TL/RL** function key.
- Swing in the MACRO objective.
Press the button that is preassigned for **MACRO**.
(For details on key assignments, please see the identification sheet.)
Alternative
Press the objective changer buttons.
- Insert a specimen.
- The MACRO objective is focused by trained Leica personnel when the instrument is installed. Due to the great field depth of the MACRO objective, the image of the specimen is pre-focused. Refocusing with the focus wheel is not necessary.

Note:

The MACRO objective is parfocal with the Leica objectives.

- Adjust the brightness using the setting wheel (44.2) on the front operating panel.

Fig. 44 Front operating panel

- 1 Aperture diaphragm adjustment
- 2 Brightness adjustment
- 3 Objective changer buttons
- 4 Variable function keys (preset at the factory)

8.8 Light sources

Transmitted light and reflected light

- The brightness is adjusted using the setting wheel (44.2) on the front operating panel. The setting is always assigned to the currently active axis for transmitted light (TL) or reflected light (RL).

Note:

If the reflected light axis is activated, the LED of the **TL/RL** button is illuminated. To toggle to the transmitted light axis, press the **TL/RL** key. The LED is not illuminated.

- The intensity is individually adjusted and stored for each objective and contrast method.

► Adjusting light via

- setting wheel on the stand
- Variable function keys on the stand, Leica SmartMove and joystick
- Leica STP6000
- Leica Application Suite software - module: Acquire

UV light

- If necessary, switch to the reflected light axis by pushing the **TL/RL** function key.
- For the manual stand:
Set the BF/UV position on the nosepiece wheel on the right side of the stand.
- Press the function key to which you have assigned the **UV** function.
- The aperture is set via the dial (44.1) on the front control panel. This setting should generally be about 2/3 of the total value.
- The brightness is adjusted using the setting wheel (44.2) on the front operating panel.

! Caution!

As soon as the **UV** key on the stand is pressed, the motorized objective nosepiece moves the UV objective 150x into the light path. Remember to check that the working distance of the objective to the specimen is adequate before switching, as otherwise the objective and/or specimen may be damaged.

Refer also to **9.1.6 UV mode** on p. 60.

8. Operation

IR illumination

- Switch on the external power supply with the On/Off switch (45.1).
- Set the required intensity with the rotary switch (46.1).

Caution!

Lamp housing and stage plate may be hot!

See also **9.2.3 Infrared* (TL)**.

Fig. 45 External power supply for the IR-illumination

1 On/off switch



Fig. 46 External power supply for the IR-illumination

1 Intensity adjustment

8.9 Aperture diaphragm

Transmitted light and reflected light

The aperture diaphragm has been set to suitable values for the current objective and contrast method at the factory.

- You can change the diaphragm at any time using the setting wheel (44.1) on the front operating panel.

The function buttons are always assigned to the currently active axis for transmitted light (TL) or reflected light (RL, UV) .

Note:

If the reflected light axis is activated, the LED of the **TL/RL** button is illuminated. To toggle to the transmitted light axis, press the **TL/RL** key. The LED is not illuminated.

Caution:

The old values will be overwritten by the current ones!

Caution:

When using **DF**, the aperture diaphragm is open all the way and should not be closed.

► Adjusting diaphragm via

- setting wheel on the stand
- Variable function keys on the stand, Leica SmartMove and joystick
- Leica STP6000
- Leica Application Suite software - module: Acquire

8.10 Polarizer and analyzer

Only manual operation of the polarizer and analyzer is possible.

Transmitted light

- Swing the polarizer in or out using the lever (50.1).

Reflected light

- Where applicable, remove the plug cap on the right side of the stand.
- Insert the polarizer into the receptacle until it latches in place (51.1).
- Where applicable, remove the plug cap on the left side of the stand.
- Insert the analyzer into the receptacle until it latches in place (52.1).

Fig. 50

- 1 Swing-in lever for TL polarizer

8.11 Selecting the contrast method

Manual stand

- The front side of the stand has a label specifying which position (1 - 4) of the reflector disk must be set to carry out a certain method. The individual filter positions are numbered accordingly with 1 - 4.
Turn the knurled ring on the right side of the stand to set the desired position.

Fig. 51

- 1 Polarizer mount
2 Filter slide
3 Slot for objective prism slider

Fig. 52

- 1 Analyzer slot
2 Fine adjustment for objective prism slider

8. Operation

Motorized stand

- To select a contrast method or the MACRO overview magnification (one of the freely assignable function keys), press the corresponding function key.
The necessary filter position and components are positioned automatically.
- Alternatively, this can also be controlled using the Leica STP6000 or the LAS software.

► Changing the contrast method - motorized stand

- Variable function keys on the stand, Leica SmartMove and joystick
- Leica STP6000
- Leica Application Suite software - module: Acquire

8.12 Objective prism slider

- Insert the objective prism slider for DIC into the fixture (51.3). The identifying letter, e.g. D, must match the identifying letter of the objective engraving (top row). The number after the code letter only specifies one variant, e.g. D1 = also applies for pupil position D. For fine adjustment use the adjusting screw (52.2) above the objective prism slider.

8.13 Filter slide

- For using filters, the two levers (51.2) on the right side of the stand are provided.
Insert the desired filter into the receptacle until it locks into place.

9. Contrast methods

For the motorized stand, all contrast methods can be selected and operated via the variable function keys on the stand, Leica SmartMove and joystick, as well as using the Leica STP6000 and the Leica Application Suite software. The only exceptions are methods involving components that have to be operated manually (e.g. those with manual objective prism slider). The following describes operation via the function keys on the stand. Refer to the separate instructions for operation via the Leica STP6000 and software. For the manual stand, the required filter position of the reflector disk is adjusted manually by turning the knurled ring on the right side of the stand.

9.1 Reflected light

9.1.1 Bright field (RL)

- If necessary, switch to the reflected light axis (RL) by pushing the **TL/RL** function key.
- Select the **BF** (bright field) contrast method
Motorized stand:
 To do so, press the button preassigned for bright field (**BF**).
 Alternatively: Press the button to which the **CHANGE RL**  function is assigned.
 (For details on key assignments, please see the identification sheet.)
Manual stand
 Rotate the reflector disk manually until the position labeled on the front operating panel for **BF** is reached.
- Insert a specimen.
- Position an appropriate objective.
- Focus the image using the rotary knob on the Leica SmartMove, Leica STP6000 or the focus handwheel and adjust the brightness using the setting wheel or the **INT** function keys.

9. Contrast Methods

9.1.2 Dark field (RL)

- If necessary, switch to the reflected light axis (RL) by pushing the **TL/RL** function key.
- Select the **DF** (dark field) contrast method
Motorized stand:
To do so, press the button preassigned for dark field (**DF**).
Alternatively: Press the button to which the **CHANGE RL**  function is assigned.
(For details on key assignments, please see the identification sheet.)
The DF reflector is rotated into the beam path.
Manual stand
Swing the DF reflector into the beam path. To do so, rotate the reflector disk manually until the position labeled on the front operating panel for **DF** is reached.
- Insert a specimen.
- Position an appropriate objective.
- Focus the image using the rotary knob on the Leica SmartMove, Leica STP6000 or the focus handwheel and adjust the brightness using the setting wheel or the **INT** function keys.

Note:

When the dark field method is selected, the aperture diaphragm is opened fully and should not be adjusted.

9.1.3 Polarization (RL)

- If necessary, switch to the reflected light axis (RL) by pushing the **TL/RL** function key.
- Insert a suitable specimen and rotate an appropriate objective into place.
- Select the **POL** (polarization) contrast method
Motorized stand:
To do so, press the button preassigned for polarization (**POL**).
Alternatively: Press the button to which the **CHANGE RL**  function is assigned.
(For details on key assignments, please see the identification sheet.)
The corresponding reflector cube (containing the polarizer and analyzer) is automatically brought into the beam path on the reflected light axis.
Manual stand
Rotate the reflector disk manually until the position labeled on the front operating panel for **POL** is reached.

Alternative procedure:

- Configure the bright-field position **BF**.
(Using the corresponding function key for the motorized stand or by turning the reflector turret for the manual stand.)
- Manually push the corresponding polarizer and the analyzer on the stand into the beam path. Then bring the polarizer and analyzer into cross position until they reach maximum darkness.

9.1.4 Interference contrast (RL)

- If necessary, switch to the reflected light axis (RL) by pushing the **TL/RL** function key.
- Place a specimen on the stage and select a suitable objective.

- Select the **DIC** (Differential Interference Contrast) contrast method.

Motorized stand:

- To do so, press the button preassigned for interference contrast (DIC).
- Insert the objective prism slider into the slot on the objective nosepiece (51.3).
- For fine adjustment, rotate the knurled screw (52.3) on the objective prism slider.

Alternatively: Press the button to which the **CHANGE RL**  function is assigned.

(For details on key assignments, please see the identification sheet.)

Manual stand

- Rotate the reflector disk manually until the position labeled on the front operating panel for **DIC** is reached.
- Insert the objective prism slider into the slot on the objective nosepiece (51.3).
- For fine adjustment, rotate the knurled screw (52.3) on the objective prism slider.

Alternative procedure:

- Adjust the bright-field position **BF**.
(Using the corresponding function key for the motorized stand or by turning the reflector turret for the manual stand.)
- Manually push the corresponding polarizer and the analyzer on the stand into the beam path. Then bring the polarizer and analyzer into cross position until they reach maximum darkness.

9.1.5 Oblique illumination (RL)

- If necessary, switch to the reflected light axis (RL) by pushing the **TL/RL** function key.

- Place a specimen on the stage and select a suitable objective.

- Select one of these methods: **BF**, **POL**, **DIC** or **UV**.

Motorized stand:

To do so, press the corresponding function key or press the key to which the **CHANGE RL**  function is assigned.

(For details on key assignments, please see the identification sheet.)

Manual stand

To do so, rotate the reflector disk manually until the position labeled on the front operating panel for this method is reached.

- Focus the image using the rotary knob on the Leica SmartMove, Leica STP6000 or the focus handwheel and adjust the brightness using the setting wheel or the **INT** function keys.

- Press the **OBL** (Oblique) function key.

- The fine adjustment of the angle of light incidence is carried out using the aperture setting wheel. If necessary, adapt the intensity adjustment.

- Pressing the **OBL** function key again switches the oblique illumination off again, so that by pressing this key multiple times, you can compare the images in the selected method with and without oblique illumination.

**Note:**

When changing objectives, the oblique illumination is switched off automatically.

9. Contrast Methods

9.1.6 UV mode

The following components are necessary for microscopy in UV mode:

- The VIS/UV reflected illuminator (365nm i-line)
- 150x UV objective
- BF/UV reflector cube (also suitable for bright field (RL))
- C-mount, 1x without glass
- Camera with UV filter kit
- PC with monitor

- If necessary, switch to the reflected light axis (I) by pushing the **TL/RL** function key.

- Configure the filter position BF/UV.

Motorized stand:

Press the button that is preassigned for **BF/UV**. Alternatively: Press the button to which the **CHANGE RL**  function is assigned.

(For details on key assignments, please see the identification sheet.)

Manual stand

Rotate the reflector disk manually until the position labeled on the front operating panel for **BF/UV** is reached.

- Press the button that is preassigned for **UV**. The visual illumination is switched off and the UV illumination is switched on. The 150x UV objective is positioned.

Caution!

As soon as the **UV** key on the stand is pressed, the motorized objective nosepiece moves the UV objective 150x into the light path. Remember to check that the working distance of the objective to the specimen is adequate before switching, as otherwise the objective and/or specimen may be damaged.

- If necessary, pull out the control bar (43.1) on the right side of the tube so that 100% of the light is allocated to the camera port.

Note:

In UV mode, no UV light reaches the eyepieces. Therefore, the image can be viewed on the monitor of a PC only.

- Insert a specimen and adjust the intensity, aperture and focus as usual.
- Pressing the **UV** function button switches the UV illumination back off. The visual illumination is switched back on. The objective 150x UV remains positioned and can be used for the bright field (BF), dark field (DF) and interference contrast (DIC) methods.

Note:

When changing objectives, the UV illumination is switched off automatically.

Note:

If you press the **OBL** button in UV mode, this provides the **OUV** method (oblique illumination in UV mode).

9.2 Transmitted light

9.2.1 Bright field (TL)

- If necessary, switch to the transmitted light axis (TL) by pushing the **TL/RL** function key.
- Select the **BF** (bright field) contrast method.
Motorized stand:
 To do so, press the button preassigned for bright field (**BF**).
 (For details on key assignments, please see the identification sheet.)
Manual stand
 To do so, rotate the reflector disk manually until the position labeled on the front operating panel for the **empty position** is reached.
 Alternatively, the **BF** filter position can also be configured.
- Insert a transmitted light specimen.
- Position an appropriate objective.
- Focus the image using the rotary knob on the Leica SmartMove, Leica STP6000 or the focus handwheel and adjust the brightness using the setting wheel or the **INT** function keys.

9.2.2 Polarization (TL)

- If necessary, switch to the transmitted light axis (TL) by pushing the **TL/RL** function key.
- Select an empty position in the filter turret.
Motorized stand:
 To do so, press the button preassigned for the **empty position**.
 (For details on key assignments, please see the identification sheet.)
Manual stand
 To do so, rotate the reflector disk manually until the position labeled on the front operating panel for the **empty position** is reached.
- Rotate the polarizer into the beam path using the lever (50.1) on the left side of the stand.
- Insert the analyzer on the left side of the stand up to the click stop (52.1).
- For adjustable analyzer: Bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Insert a suitable specimen and rotate an appropriate objective into place.

9. Contrast Methods

9.2.3 Infrared* (TL)

When silicon wafers are inspected under IR light they appear transparent. This effect is used, for example, for aligning single layers on top of each other. As IR light cannot be seen with the naked eye, a special camera is used in combination with LAS software and a monitor. This monochromatic screen image enables different planes in the silicon wafer to be focused.

- Switch on the external power supply with the On/Off switch (53.1).
- Set the required intensity with the rotary switch (54.1).

Caution!

Remember to reduce the intensity before removing a wafer!

- Focus the required layer plane.

Caution!

The lamphousing and the stage plate may be hot!

Note:

The generated image is only visible via the camera on the monitor, not in the eyepieces.

Note:

The contrast technique brightfield (RL) (**BF/RL**) is possible at the same time.

Note:

When the Infrared*(TL) technique is used, no other transmitted light contrast method is possible. Therefore, this technique is only available for reflected light microscopes.

Fig. 53 External power supply for the IR-illumination

1 On/off switch



Fig. 54 External power supply for the IR-illumination

1 Intensity adjustment

9.3 Combination methods

9.3.1 Bright field (TL)/ Bright field (RL)

- If necessary, assign the **BF-BF** function to one of the variable function keys using the LAS software.
- Press the **BF-BF** button. The LED of this button is illuminated. The light is turned on for both axis (TL and RL).
- To toggle between the transmitted light axis and the reflected light axis press the **TL/RL** button.
If the reflected light axis is activated, the LED of the **TL/RL** button is illuminated. The LED is not illuminated, if the transmitted light axis is activated.
- If the LED is not illuminated (TL axis active), the settings of intensity and aperture diaphragm via the setting wheels (44.1) and (44.2) are assigned to the TL axis.
- If the LED is illuminated (RL axis active), the settings of the aperture diaphragm and intensity via the setting wheels (44.1) and (44.2) are assigned to the RL axis.

9.3.2 Bright field (TL)/Oblique illumination (RL)

- If necessary, assign the **BF-BF** function and the **OBL** function to variable function keys using the LAS software.
- Press the **BF-BF** button. The LED of this button is illuminated. The light is turned on for both axis (TL and RL).
- If necessary, press the **TL/RL** button to activate the reflected light axis. The LED of the **TL/RL** button is illuminated.
- Press the **OBL** button. The LED of the **OBL** button is illuminated.
- Now the fine adjustment of the angle of light incidence is carried out using the aperture setting wheel (44.1).
- Press again the **TL/RL** button to switch to the transmitted light axis. The settings of intensity and aperture diaphragm via the setting wheels (44.1) and (44.2) are now assigned to the TL axis.
- To toggle between Bright field (TL) and Oblique illumination (RL) press the **TL/RL** button successively.

10. Care and Maintenance

Caution!

Unplug the power supply before performing cleaning and maintenance work!
Protect electrical components from moisture!

Microscopes in warm and warm-damp climatic zones require special care in order to prevent fungus contamination.

The microscope should be cleaned after each use, and the microscope optics should be kept painstakingly clean.

10.1 Dust cover

Note:

To protect against dust, cover the microscope and accessories with the dust cover after each use.

Caution!

Allow the microscope to cool down first. The dust cover is not heat-resistant. In addition, condensation may occur.

10.2 Cleaning

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with a low-concentrated soap solution, petroleum ether or ethyl alcohol.

For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Thinners containing acetone, xylen or nitrogen can harm the microscope and thus must not be used.

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure that coated or plastic surfaces do not become matted or etched.

Cleaning glass surfaces

Remove dust on glass surfaces with a fine, dry and grease-free brush made from hair, by blowing with a squeeze blower, or vacuum suction.

Carefully remove stubborn dirt on glass surfaces with a clean cloth moistened with distilled water. If the dirt still can not be removed, use ethyl alcohol or petroleum ether instead of water.

Cleaning objectives

Caution!

The objective may not be unscrewed during cleaning. If damage appears on inner surfaces, the objectives must be sent to your Leica subsidiary for repair. We also advise against cleaning the inside surfaces of the eyepieces.

The front lenses of objectives are cleaned as described under "Cleaning glass surfaces". The upper lens is cleaned by being blown off with a pneumatic pump.

Note:

For cleaning glass surfaces, particularly objectives, download the brochure "Cleaning of Microscope Optics" at

<http://www.leica-microsystems.com/products/light-microscopes/industrial-materials/upright-microscopes/details/product/leica-dm6000-m/downloads/>

In case of questions, please contact our Technical Service.

10.3 Handling acids and bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

Caution:

Never allow the optics and mechanical parts to come into direct contact with these chemicals.

10.4 Maintenance interval


The intervals between the recommended service inspections depend on how heavily the instrument is used. The rear side of the stand has a maintenance indicator light (53.1) that indicates when maintenance is recommended.

At first, the LED illuminates in green. After about 1 year (depending on how heavily the instrument is used), the color changes to orange. After another 8 months or so, the LED is illuminated in red. Then, a service inspection should be carried out.

In case of questions, please contact Technical Service.

Fig. 53

1 Maintenance indicator light

11. Essential Wear and Spare Parts

11. Essential Wear and Spare Parts

Order No. Material No.	Name	Used for
<u>Screw cap for unused objective receptacles</u>		
11-020-422.557-000	Screw cap M 32	Objective nosepiece
<u>Replacement eyecup (antiglare protection) for HC PLAN eyepiece</u>		
11-021-500.017-005	HC PLAN eyecup	10x/25 eyepiece
11-021-264.520-018	HC PLAN eyecup	10x/22 eyepiece
11-021-264.520-018	HC PLAN eyecup	10x/20 eyepiece
<u>Fuses</u>		
11-302-150.001-203	T 3.15 A 250 V	Fuse (power supply)

12. Abbreviation Glossary

AP	Aperture diaphragm
BF	Bright field
CAM	Camera button - capture image
COARSE	"Fast" stage step size and coarse focus
COMB	Combination methods
DF	Dark field reflected light
DIC (ICR)	Differential interference contrast (DIC)
FINE	"Precise" stage step size and fine focus
ICR	Interference contrast, reflected light
ICT	Interference contrast, transmitted light
INT	Intensity
IR	Infrared
MACRO	MACRO overview magnification
MEM	Memory position
OBL	Oblique illumination
OBJ	Objective
OUV	Oblique illumination in UV mode
POL	Polarization, reflected / transmitted light
RL	Reflected light
TL	Transmitted light
UV	UV illumination/UV mode

13. Index

A

Abbreviation Glossary 67
Allowable ambient conditions 22
Ambient conditions 22
Analyzer 30, 55, 58, 61
Aperture diaphragm 39, 42, 54

B

Beam splitting 50
Bright field (RL) 57
Bright field (TL) 61

C

Cables 34
Camera 31
Care 64
Changing objectives 51
Cleaning 64
Coaxial drive 44
Combination methods 15, 63
Condenser 27
Condenser head 27
Condenser height adjustment 43
Connecting the computer 34
Connection for motorized stage 27, 34
Contamination guard 29
Contrast methods 15, 55, 57
Controls 17
Correction for vision problems 51

D

Dark field (RL) 58
Declaration of conformity 70
Dimensions 24
Disposal 14
Dual attachment 50
Dust protector 64

E

External power supply 33, 54, 62
Eyepieces 28, 50

F

Filter slide 30, 56
Focus 48
Focus finder 49
Focusing 46
Focusing telescope 42
Focus position 48
Focus threshold 47, 48
Footswitch 31, 34
Front operating panel 39
Function keys 38, 39

G

Graticule 51

H

Handling acids and bases 65
Height adjustment
of the focus knobs 47

I

Illumination 29
Infrared* (TL) 62
Inspection stage 26
Installation location 22
Intensity Setting 39, 53
Interface 17
Interference contrast (RL) 59
Interpupillary distance 49
IR illumination 32, 54

L

LED illumination 29
Leica SmartMove 20, 34, 41
Leica STP6000 20, 41
Light sources 53
Load position 46

M

MACRO overview magnification 52
Magnification changer 28
Maintenance indicator light 34, 65
Maintenance interval 65

O

Object displacement 44, 45
Objective changer buttons 39
Objective nosepiece 16, 28
Objective prism slider 30, 56
Objectives 28
Oblique illumination (RL) 59
Operating temperature 11

P

Polarization (RL) 58
Polarization (TL) 61
Polarizer 30, 55, 58, 59, 61
Power supply 34

R

Rapid adjustment 44, 45
Reflected light 57
Reflected light axis 15
Reflected light reflector turret 29
Reflector cube 29
Reflector turret 29
Reset function 36, 44
Resistance (focus knob) 48
RS232 interface 34

S

Safety Notes 10
Software 17
Spare parts 66
Specifications 10
Speed switch 47
Stage Center Position 46
Stage insert 27
Stage, manual 44
Stage, motorized 45
Stage positions 46
Step sizes 46, 49
Switching on 39
Symbols 7

T

Technical data 10
Toggling between reflected light/
transmitted light 46
Transmitted light 61
Transmitted light axis 15
Transport 12, 23
Tube 28, 49, 50
Tube attachment 31
Type plate 13, 14

U

USB port 34
UV light 53
UV mode 60

V

Variable function keys 39
Viewing angle 50

W

Warnings 12
Weight 21, 23

X

XY joystick 20, 34, 40

14. EC Declaration of Conformity

14. EC Declaration of Conformity

To download the EU Declaration of Conformity for your product use the link

[http://www.leica-microsystems.com/
products/light-microscopes/industrial-materials/microelectronics](http://www.leica-microsystems.com/products/light-microscopes/industrial-materials/microelectronics)

Select the type of your microscope and switch to the Download page.

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Copyright © Leica Microsystems CMS GmbH · Ernst-Leitz-Straße · 35578 Wetzlar · Germany 2018
Tel. (0 6441) 29-0 · Fax (06441) 29-2599 LEICA and the Leica logos are registered trademarks of Leica IR GmbH.