Leica FS M

Operating Manual

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Mode d'emploi

Leica Microsystems CMS GmbH. Manual, 11933989, Revision 1.1. from December 15, 2013

Living up to Life

Leica
MICROSYSTEMS

Revision 1.1, published December 15, 2013 by:

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Leica
MICROSYSTEMS

Copyrights

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Nevertheless, no liability of any kind shall be accepted for the accuracy of this manual's contents. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information included in this manual may be changed without prior notice.

Contents

1	Important Notes about this Manual.....	6	9.	Accessories	25
2.	Intended Purpose of the Device.....	8	9.1	Motorized work table.....	25
3.	Safety Notes.....	10	9.2	Stackable stages	26
3.1	General safety notes.....	10	9.3	Slide holder.....	28
3.2	Electrical safety	10	9.4	Universal rotating mechanism	35
3.3	Disposal.....	12	9.5	Light panels	35
4.	Installation Location	13	9.6	Filter for color contrast.....	36
5.	Unpacking the Macroscope	14	9.7	Additional illumination options.....	36
6.	Overview of the Instrument.....	16	9.8	Stage scales and calibration standard.....	38
6.1	Instrument identification	17	10.	Technical Data.....	39
7.	Assembly	18	11.	Care of the Instrument.....	42
7.1	Comparison bridge.....	18	11.1	Dust cover.....	42
7.2	Illumination	18	11.2	Cleaning	42
7.3	Objectives	19	11.3	Handling acids and alkalines	43
7.4	Tube	19	12.	Essential Wear and Spare Parts.....	44
7.5	Eyepieces.....	19	13.	EC Declaration of Conformity	45
7.6	Connection to the power supply.....	20			
8.	Operation	20			
8.1	Switching on the microscope	20			
8.2	Stages.....	21			
8.3	Tube	21			
8.4	Eyepieces.....	22			
8.6	Focusing.....	22			
8.7	Objectives	22			
8.8	Magnification changer*.....	23			
8.9	Adjusting the height of the stand.....	23			
8.10	Diaphragm setting	24			
8.11	Magnification calibration.....	24			

1 Important Notes about this Manual

Important!

This operating manual is an essential component of the microscope; it must be read carefully before the product is assembled, put into operation or used, and must be kept for later reference.

This operating manual contains important instructions and information for the operational safety and maintenance of the comparison microscope and accessories. Therefore, it must be kept in a safe place.

Text symbols and their meanings:

(1.2)

Numbers in parentheses such as "(1.2)", correspond to illustrations (in the example, Figure 1, Item 2).

→ p.20

Numbers with an arrow, e.g. → p. 20, refer to a specific page of the manual.

WARNING indicates a hazard with a medium degree of risk that, if not avoided, can result in death or serious injury.

CAUTION indicates a hazard with a low degree of risk that, if not avoided, can result in minor or moderate injury.

Notice!
Special safety instructions within this manual are indicated with the triangle symbol shown here, and have a gray background.

Notice! The macroscope and accessories can be damaged when operated incorrectly.

Warning of hazardous electrical voltage! Risk of electrical shock!

Danger due to hot surface.

Warning of hazardous laser radiation – take appropriate safety precautions!

Warning of optical radiation! Never look directly into the light beam! Wear safety goggles!

Warning of electromagnetic field

Warning of UV radiation! Wear safety goggles!

Warning of hand injuries

Instructions for disposing of the instrument, its accessories and consumables.

Connection to ground!

Explanatory note.

*

Item not contained in all configurations.

Manufacturer's address with manufacturing date

China RoHS 50 years EFUP
(Environmentally friendly use period)

2. Intended Purpose of the Device

Comparisons of small surfaces or profiles can be performed with sufficient precision only if the specimens of comparison can be observed simultaneously.

If only a macroscope with just one imaging beam path is available for such comparisons, then, depending on the interruption of observation by a change of objective, high demands are made on the memory of the observer. Therefore, it is impossible to rule out incorrect interpretation in the case of specimens with very minor differences in the surface or profile.

The Leica FS M macroscope with two image beam inputs is designed for comparing structures of items such as those listed below in split images and

- Evidence on fired ammunition pieces
- Toolmarks
- Documents

composite images in forensic laboratories, document printing facilities, mints and banks, and documenting the results on photographs and videos: A wide variety of illumination equipment is available for various applications.

The device meets the Council Directives 2006/95/EC concerning electrical apparatus and 2004/108/EC concerning electromagnetic compatibility for use as a laboratory instrument.

Reasonably foreseeable misuse

The following are prohibited:

- Using the comparison macroscope for any purpose not in accordance with the Declaration of Conformity (e.g. using it as an in-vitro diagnostic product according to Council Directive 98/79/EC or as a medical product according to Council Directive 93/42/EEC)
- Making use of auxiliary equipment to use the area between the microscope stage and objectives as a clamp or holder (working like a vice).
- Altering the proper adaptation of the original stage on the system (tilting, rotating, lifting, etc.)
- Using the stage surface as storage location for objects (utensils, fluid containers, etc.)
- Operating the comparison macroscope in an inclined position.
- Cleaning the macroscope in a way other than specified in the manual.
- Modifying the built-in safety circuits in the instrument.
- Allowing unauthorized personnel to open the instrument.
- Using cables that Leica has not provided or permitted.
- Using combinations with non-Leica components that go beyond the scope of the manual.

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from, using the microscope for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In such cases, the Declaration of Conformity shall be invalid.

WARNING

Notes on handling laser devices

The comparison macroscope is not suitable for coupling laser radiation (such as to the camera ports), since this radiation poses a hazard for the user (eye injury in particular).

Laser couplings require corresponding safety devices that have to be inspected and installed by trained personnel.

For further information, please contact your authorized Leica Microsystems representative.

3. Safety Notes

3.1 General safety notes

This safety class 1 device is constructed and tested in accordance with EN 61010-1/IEC 61010-1, Safety regulations for electrical measuring, control, and laboratory devices. It also meets the requirements for EN 61326-1, Electrical equipment for measurement, control and laboratory use – EMC requirements.

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this operating manual.

Caution!

The instruments and accessories described in this manual have been safety-tested and checked for possible hazards.
The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the device is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the instrument or noncompliant use shall void all rights to any warranty claims and product liability!

3.2 Electrical safety

General specifications

For indoor use only.	
Supply voltage:	100–240 VAC
Frequency:	50-60 Hz
Power input:	max. 120 W
Fuses:	T2.5A H 250 VAC
Ambient temperature:	15–35°C
Relative humidity:	90% to 30 °C non-condensing
Overvoltage category:	II
Pollution degree:	2

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the instrument, or release of the protective conductor terminal, can cause the instrument to become hazardous. Intentional ground interruption is not permitted!

WARNING

Use only original power cables or alternative cables with a VDE / HAR logo, which at least fulfill the requirement of 3x0.75 mm² and 10A/250V.

WARNING

By definition, the main circuit breaker of this instrument is the connection between the power cable and device port. The user must ensure unobstructed access to the main circuit breaker at all times.

Caution!

Ancillary devices connected to the comparison microscope that have their own and/or a separate power supply can be brought to the same ground wire potential by being connected to the ground. For connections without a ground connector, Leica Service must be consulted.

Caution!

Never use any fuses as replacements other than those of the types and the current ratings listed here. Bypassing the fuse holder is not permitted.

Caution!

The comparison microscope's electrical accessory components are not protected from water. Water can cause electric shock.

Do not assemble the comparison microscope and its accessory components in the immediate vicinity of a water connection or other areas where the possibility of water exposure exists.

Caution!

Protect the comparison microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.

Caution!

Before exchanging the fuses or lamps, be absolutely certain to switch off the main power switch and remove the power cable.

3. Safety Notes

Caution!

This microscope may not be used at elevations over 2000 m above sea level.

Caution!

Transport and storage at -20 °C – +85 °C and max. 90% humidity (non-condensing).

CAUTION

At a full load, the microscope weighs more than 18 kg. The user must take corresponding precautions for transporting the instrument.

Caution!

This instrument has been developed and tested in accordance with CISPR 11 Class A. It can cause radio interference in a household environment. In this case, you may have to take measures to eliminate the interference.

3.3 Disposal

Once the product has reached the end of its service life, please contact Leica Service or Sales about disposal.

Please observe and comply with the national and federal laws and regulations that are equivalent to EC directives such as WEEE.

Note!

Like all electronic devices, this instrument, its accessory components and consumables must never be disposed of with general household waste.

4. Installation Location

Caution!

Ensure that the workplace area is free of oil and chemical fumes. Vibrations, direct sunlight and strong temperature fluctuations interfere with measurements as well as microphotography and image documentation.

The basic prerequisite is a stable instrument table of optimum height (70–80 cm). Combined with an ergonomically correct, adjustable chair, these are the basic requirements for fatigue-free microscopy.

Permitted ambient conditions

Temperature	15–35 °C
Relative humidity	maximum 90% up to 30 °C

Macroscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

For additional instructions, see the chapter entitled "Care of the Instrument" Page 42

CAUTION

Electrical components must be placed at least 10 cm away from the wall and away from flammable substances.

CAUTION

The macroscope has to be placed so that the power socket is not obstructed in order to quickly disconnect the instrument from the power if necessary.

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

Caution!

This instrument has been built for use in a managed electromagnetic environment. This means that the use of radio transmitting devices such as cell phones can cause distortion when evaluating pictures if used in the immediate vicinity of the macroscope.

5. Unpacking the Macroscope

The stand is delivered in a solid box with packaging inserts.

The delivery includes the following components:

- Stand with base as well as the stages, which are already installed. The stand is equipped with two carrying handles.
- Comparison bridge
- Tube
- One or two cold light sources
- Accessory case with slide holders
- Optical accessories
- Oblique incident illuminator
- Optional accessories

When unpacking the equipment, follow these instructions:

- After opening the box, first remove the comparison macroscope and then the cardboard boxes, which hold the accessories.
- Lift the stand out of the carton by its carrying handles (1.1).

**CAUTION**

When lifting the instrument out of the packaging, and when transporting it to its installation site, carry it by its handles only.

**CAUTION**

Note that the stand weighs approx. 32 kg. Two people are required to lift it out if its packaging and carry it to its installation site.

Fig. 1 Stand with carrying handles
1 Removable carrying handle

Later, the carrying handles can be unscrewed and removed from the upper part of the macroscope. Then, seal the threaded holes using seal plugs.

- Carefully remove all components from the transportation and packaging materials.

**Note:**

Keep the packaging materials for storage purposes and transportation of the individual components.

**Caution!**

If at all possible, avoid touching the lens surfaces of the objectives. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces of optical surfaces in a short time. For additional instructions, see the chapter entitled "Care of the Instrument" → Page 42

- Check the delivery carefully against the packing list, delivery note or the invoice. We highly recommend keeping a copy of these documents with the manual in order to have delivery date and delivery package information available when (for example) placing repeat orders or making service calls. Please make sure that no small parts are lost in the packaging materials.

**Caution!**

Do not yet connect the comparison macroscope and peripherals to the power supply at this point!

6. Overview of the Instrument


Fig. 2 Leica FS M

- 1 Comparison bridge
- 2 Leica DFC digital camera
- 3 Tube with eyepieces
- 4 Bridge mode setting dial, right
- 5 Elements for positioning the dividing line
- 6 Adjusting screw for the dividing line
- 7 Objective turret with objectives
- 8 LED spot illumination
- 9 Coarse and fine focusing, right
- 10 Synchronous movement of both cross-stages in the X direction
- 11 Rocker switch for automatic height adjustment of the comparison bridge
- 12 Coarse and fine focusing, left
- 13 Specimen stage with coaxial drive
- 14 Receptacle for filter slides
- 15 Bridge mode setting dial, left

6.1 Instrument identification

Leica FS M Type label

Dimensions and weight

Leica FS M:

Width x Height X Depth: Approx. 70 x 84.5 x

53 cm

Weight: approx. 45 kg

7. Assembly

Only a few commonly used screwdrivers and keys are necessary for assembly; these are included in the delivery package.
A holder is provided on the rear side of the stand for storing the adjusting key for setting the dividing line.

! Caution!

Unscrew the transport anchor screws (3.1) on the bottom of the comparison bridge and screw them into the threading provided for storage (3.2) (rest position).

7.1 Comparison bridge

- From above, set the comparison bridge onto the dovetail guide and push the bridge downwards until the upper end of the dovetail guide is tightly fastened.
- Fasten the comparison bridge using the side clamping screw.

7.2 Illumination

- If present:
- Install the rotating mechanism with joint holder (Fig. 4) for holding the light on the stage.
- For more options for illumination equipment, see the chapter entitled "9. Accessories" → Page 25

Fig. 3 Transport anchor for the comparison bridge

- 1 Transport anchor screws
- 2 Rest position for the transport anchor screws


Fig. 4 Illumination equipment

- 1 LED spot illumination
- 2 Joint holder
- 3 Rotating mechanism

CAUTION

Make sure to follow the safety notes in the operating manual for the light source.

7.3 Objectives

The position of the individual objectives is preset at the factory. The objective turret collars are labeled (5.2).

! Caution!

The objectives for the left and right objective turrets are packaged separately. Observe the side assignment when screwing in the objectives.

- Screw the objective into the objective turret.

Position	Objective
50	0.4x
100	1x
200	2x
500	4x
1000	8x

Note:

For the 0.4x objective, the SL 0.3:1 auxiliary lens (order No. 11 532 301), which is attached to the objective, is available. 1:1 imaging is achieved using the auxiliary lens.

7.4 Tube

- Unscrew the clamping screw on the tube holder
- Insert the tube into the tube holder.
- Retighten the clamping screw.

Note:

When using the optional magnification changer, it is used between the tube holder and the tube and fixed in place using the clamping screw.

7.5 Eyepieces

- The eyepieces are inserted into the eyepiece tubes on the tube.

Fig. 5 Objective turret with objectives

- 1 Filter receptacle
- 2 Objective positions

8. Operation

7.6 Connection to the power supply

- Insert the power cable at the back of the comparison microscope (6.1).
- Then connect the comparison microscope to the power supply.

The Leica FS M comparison microscope offers the following functions:

- Simultaneous (in the X direction) or separate movement of the left and right stage
- Left and right-side focusing
- Left and right-side objective turret
- Height-adjustable stand

8.1 Switching on the microscope

The stand is switched on using the switch (7.1). The illumination is switched on separately.

If applicable, only switch on the light source if the light guide is firmly connected to the microscope. Uncontrolled light output from the light guide poses a blinding hazard! Follow the safety notes in the operating manual for the light source.

Fig. 6 Rear side of the Leica FS M
1 Power supply connection

8.2 Stages

The stages on the left and right side can be moved either independently of each other with the respective coaxial drive (7.2), or simultaneously in the X direction (7.3).

8.3 Tube

Note:

Close any unused tube openings, as otherwise stray light can interfere with observation.

Adjusting the interpupillary distance

- Adjust the viewing distance of the eyepiece tubes so that a congruent total image is seen.

For the ErgoTube:

Adjusting the viewing angle

- The viewing angle can be adjusted from 5° to 32° by tilting the binocular eyepiece.

Fig. 7 Right stand side

- 1 On/off switch
- 2 Coaxial drive
- 3 Synchronous movement of both stages in the X direction

8.4 Eyepieces

Note:

The eyepiece's aperture protector must be removed or folded back during microscopy while wearing eyeglasses.

Eyeglasses with multifocal lenses (bifocals and smooth view glasses) must be removed while operating the microscope.

Correction for vision problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.

If both eyepieces have an adjustable eyelens, the eyelens of the right eyepiece can be brought to the position for normal vision, which is indicated by an encircling line.

When using a graticule in the right eyepiece, the eyelens must first be focused on the graticule.

- Then, view the same specimen using the left eye and rotate the eyelens of the left eyepiece until the specimen appears in sharp focus. While doing so, do not use the focus dial.

8.5 Brightness adjustment

To set the brightness, follow the "LED Controller Coolcontrol" operating instructions provided separately.

Note:

The Leica FS M comparison microscope gives you a variety of illumination options, from LED illumination to cold light illumination with light guides all the way to light panels for diffuse illumination and special illumination equipment for UV illumination.

Refer also to "9. Accessories" → Page 25 and the respective special manual for each.

8.6 Focusing

Coarse and fine focusing is done separately for the left and right side using the respective focus knobs (8.1 and 8.3).

8.7 Objectives

The desired objectives are swiveled into place manually.

All objectives have an iris diaphragm. The iris diaphragm allows you to continuously adjust the numerical aperture. Five ratchets make it easy to adjust to identical settings on both sides.

Fig. 8 Front view of the Leica FS M

- 1 Coarse and fine focusing, left
- 2 Rocker switch for adjusting the height of the stand
- 3 Coarse and fine focusing, right

8.8 Magnification changer*

The magnification changer is switched on and off manually.

Magnification factors: 1x, 1.5x, 2x

When the magnification changer is switched on, the current magnification of both images is increased by the selected factor.

The magnification changer controls the visual and photo outputs simultaneously.

8.9 Adjusting the height of the stand


The rocker switch (8.2) adjusts the height of the stand individually. This is provided for aligning the height of specimens or specimen holders which vary in height.

There is a risk of pinching and crushing injuries near the base when the lifting column of the Leica FS M is moving downwards. Ensure that there are not objects or body parts in this area!

Fig. 9 Bridge mode controls

- 1 Diaphragm setting, left
- 2 Diaphragm setting, right
- 3 Adjusting screw for the dividing line (parallelism)
- 4 Positioning the dividing line
- 5 Center position for accessing the adjustment screws
- 6 Adjusting screw for the dividing line (width)

8.10 Diaphragm setting

The diaphragms allow transition from a superimposed image to the separate view of the right and left image, as well as setting of the dividing line.

Four diaphragm positions come predefined and can be selected using the rotary knobs (9.1 and 9.2):

- Mix The left and right images are superimposed (composite image method). Both knobs at position ○
- LR **The left and right image appear side-by-side.** A line in the middle of the field of vision separates the two images (split-image method).
Left knob at position ●
Right knob at position ●
- L **Only the left image is displayed.**
Left knob at position ○
Right knob at position ●
- R **Only the right image is displayed.**
Left knob at position ●
Right knob at position ○

The dividing line is moved to the left or to the right by rotating the knob (9.4) and has an impact on the size of the fields in split image mode. The center position (9.5) puts the dividing line back in the middle.

The adjusting screw (9.6) is used to change the width of the dividing line. The preferential setting should be a thin, black hairline to highlight the distinction between the two fields. If the edges of the dividing line do not run parallel to each other, you can adjust the parallelism using the adjusting screw (9.3).

In order to access the adjusting screws, the pointer on the knob (9.4) must point at the center position (9.5). The width of the dividing line is then set using an Allen key located in a holder on the back of the stand.

Caution!

Make sure that the adjusting key is returned to its storage location after adjusting the dividing line. In particular, the adjusting key must not be in the adjustment opening when turning the knob (9.4) for positioning the dividing line.

8.11 Magnification calibration

You do not have to calibrate the magnification on the left and right, such as for zoom-based optics, or calibrate the optics using stage micrometers since the maximum deviation is 0.4%.

9. Accessories

A specially developed work table and correctly sized work chair allow ergonomically optimized work.

In addition to various optical accessories, such as filters and stage scales, a number of different slide holders can be used.

9.1 Motorized work table

The work table with motorized height adjustment is designed for optimum ergonomic conditions when working at the comparison macroscope. Its table surface measures 120 cm x 56 cm and it is mounted symmetrically on the height-adjustable column.

Caution!

When setting up the comparison macroscope, make sure that the stand is positioned exactly over the adjustable-height column.

Fig. 30 Motorized work table
 1 Buttons for height adjustment
 2 Adjustable-height feet

After the power supply is connected, the work table can be continuously adjusted within a range of approx. 62 to 92 cm over floor height. To do so, use the buttons (30.1).

The lifting speed is approx. 12 mm/sec. The lifting capacity is approx. 2000 N.

If the floor on which the work table is installed is not completely level, differences in level can be compensated for by screwing the feet (30.2) in or out.

9.2 Stackable stages

The macroscopes are equipped with manual cross-stages.

Stage surface: 140 mm x 202 mm
Travel range of the stage: 50 mm x 50 mm


On the stage surface, there are threaded drill holes as receptacles for the adjustable holder (universal holder), Ø 118 mm rotating stages, large sample stages and rotating device for tool-mark examinations.

Ø 118 mm rotating stages

Rotating stages, which can be stacked on the basic equipment cross-stages (Fig. 32), are fastened using two Allen screws. They are suitable for incident light and transmitted light examinations. The rotating stages each have one opening in the center, which is sealed with a removable glass plate (32.3). The threaded bore holes on the stage surfaces are provided for screwing on the adjustable holder. Additionally, there are also holes for inserting the stage clips. The rotating stage surface can be locked into place using a fastening screw (32.2).

In the bottom parts of the rotating stages, there are receptacles for inserting the polarizers (32.1) with λ -plates (red = OK). These can be rotated 90°.

Fig. 32 Ø 118 mm rotating stage
1 Filter slot for polarizer
2 Fastening screw
3 Glass plate

Large sample stages

The stackable large sample stages can be screwed directly onto the cross-stages. They are primarily used for holding documents. Documents can be securely held in place and stretched using the magnetic clamps provided (33.1).

Of course, the large sample stages can also be used to hold other large-surface carrier material, such as materials bearing toolmarks.

Ø 75 mm inclining rotating stage

The inclining rotating stages are also screwed directly onto the cross-stages. After releasing the lever (34.1), you can rotate them horizontally and incline them up to 45° using the built-in ball-and-socket joint (34.2). After positioning the comparison specimens, be sure to close the clamping system using the lever.

The surface of the stages (which are 75 mm in diameter) has a ridged surface (34.3) for fastening specimens using plasticine. Four drill holes are provided, into which stage clips (34.4) can be inserted.

In the center of the stage plate, there is a conical drill hole (34.5), into which all specimen receptacles can be inserted (in a similar fashion as the adjustable holder).

Fig. 33 Large sample stage (assembly)

1 Magnetic clamps


Fig. 34 Ø 75 mm inclining rotating stage

1 Clamping lever, 2 Ball-and-socket joint, 3 Ridged surface

4 Drill hole for stage clips, 5 Conical drill hole

9.3 Slide holder

Adjustable holder with integrated rotating stage

Fasten the adjustable holders (universal holders) (35.1) to the cross-stage using the two screws (35.4).

After the clamping screw (35.3) is released, the universal holder can be moved sideways. This is especially useful for long, horizontally clamped pieces of evidence.

You can adjust the adapter sleeve (35.6) in the segment guide to an angle of up to 90° after releasing the clamping screw (35.2). The adapter sleeve itself has a rotation device with gradation (35.7) and can be fixed into place using the clamping screw. The receptacles for pieces of evidence listed below can be inserted into the conical bore of the adapter sleeve (35.6).

After the clamping screw (35.5) is unscrewed, the entire unit can be rotated around its own axis.

Fig. 35 Adjustable holder

- 1 Adjustable holder
- 2 Locking device: angle adjustment of the adapter sleeve
- 3 Locking device: lateral adjustment of the holder
- 4 Fastening onto stage
- 5 Locking device: rotation of the holder
- 6 Adapter sleeve
- 7 Rotation device with gradation and clamping screw

Bullet Holder

To assemble the bullet holder (36.1), slide the adapter sleeve in the segment guide downwards to the "0" marking. Guide the bullet holder vertically to the adapter sleeve and fasten it with the knurled screw (36.2). Adjust the adapter sleeve horizontally again (slide upwards to the "90" marking) and press the rubber-coated pressure plate (37.3) into the conical bore of the adapter sleeve. After unscrewing the clamping screw (37.2), you can move the centering insert (37.3) horizontally in the bullet holder, thus adjusting it to the length of the projectiles currently being compared.

Caution!

Before removing the unit, be sure there is sufficient distance between it and the objective.

Fig. 36 Adjustable holder with bullet holder (assembly)

- 1 Bullet holder
- 2 Fastening screw



Fig. 37 Adjustable holder with bullet holder

- 1 Rubber-coated pressure plate
- 2 Clamping screw
- 3 Centering insert with spring system
- 4 Bullet

Paste-on cylinder for bullets

Using adhesive – or using piceine – you can stick (for example) deformed bullets that fit neither into the bullet receptacles nor into the receptacles for deformed pieces of ammunition onto the paste-on cylinder (38.1). There are paste-on cylinders with diameters of 5, 6, 8 and 10 mm.

Receptacles for deformed bullets

Receptacles for deformed bullets (39.1) are available for holding deformed bullets, which cannot be held by the receptacles described above. Their cavities sufficiently hold such specimens in place.


Fig. 38 Adjustable holder with paste-on cylinder for bullets

- 1 Paste-on cylinder
- 2 Bullet

Fig. 39 Adjustable holder with receptacle for deformed bullets

- 1 Receptacle for deformed bullets

Plug receptacles for bullet casings

For fastening undamaged bullet casings of various calibers, plug receptacles (40.1) are provided. Their plastic holding devices are perfectly suited for holding bullet casings (40.2). The various calibers are marked on the rear sides of the plug receptacles. Worn-out holding devices can be replaced. Corresponding replacement parts are included in the delivery.


The brush holder set (40a) is available for deformed shell casings.

Fig. 40a Brush holder set for shell casings
(including deformed casings)

Fig. 40 Adjustable holder with plug receptacle
for shell casings

- 1 Plug receptacle
- 2 Shell casing


Plug receptacles for air rifle ammunition

The stainless steel pyramids (41.1) are intended for attaching air rifle ammunition.

Fig. 41 Adjustable holder with plug receptacle for air rifle
ammunition

- 1 Plug receptacle


Shotgun shell casing holder

The shotgun shell casing holder (42.1) has three spring rods rounded at the edges (42.2), onto which shotgun shell casings (43.1) of various calibers can be attached. The conical base (43.2) is provided for centering the shell casings.

Fig. 42 Adjustable holder with shotgun shell casing holder
1 Shotgun shell casing holder
2 Spring rods

Fig. 43 Adjustable holder with shotgun shell casing holder
1 Shotgun shell casing
2 Conical base

Wire holder

In the wire holder (44.1), round parts between 0.6 mm and 8 mm in diameter can be clamped in order to examine features such as cutting marks.

Receptacles for locking cylinders

Using an underlay of plasticine, these receptacles (45.1) can be used to securely fasten lock cylinders of safety locks bearing toolmarks.

Fig. 44 Adjustable holder with wire holder
1 Wire holder

Fig. 45 Adjustable holder with receptacle for locking cylinders
1 Receptacle for locking cylinder

Ø 20 mm stages

You can fasten pieces of evidence of various dimensions to the ridged surface of the stages with Ø 20 mm (46.1). As the adhesive, you can use plasticine or wax.

Ø 60 mm stages

Specimens with a circumference of up to 55 mm and a thickness of up to 15 mm can be fastened with clamps (47.2) to the Ø 60 mm stages (47.1). To do so, release the screws (47.3), raise the clamps, attach the specimens and press the clamps back in before tightening the screws. The additional threading (47.4) on the border of the Ø 60 mm stages is for arranging the clamps asymmetrically. If you remove these, the Ø 60 mm stages can also be used as large paste-on stages. You can fasten pieces of evidence of various dimensions to the ridged surface of the stages using plasticine.


Fig. 46 Adjustable holder with stage Ø 20 mm

- 1 Ø 20 mm stage
- 2 Joint holder

Fig. 47 Adjustable holder with Ø 60 mm stage

- 1 Ø 60 mm stage
- 2 Clamp
- 3 Fastening screw
- 4 Thread

Ø 25 mm stages with ball-and-socket joint

The Ø 25 mm stages with ball-and-socket joint (48.1) are inserted into the adapter sleeve of the adjustable holder.

After releasing the clamping screw (48.2) on the ball-and-socket joint, the stage surface, which has a ridged surface, can be inclined arbitrarily. After aligning the specimen, retighten the clamping screw (48.2).

Joint holder

If, for example, you want to compare the casing sleeves of shell casings, the profiles of the strike pin impressions on the percussion caps or the impressions of the extractor one after the other in one clamping, it is necessary to use joint holders (49.1). Its joint is fastened with a knurled screw (49.2).

The joint holders are provided for use as connecting pieces for the above-mentioned receptacles, with the exception of the bullet holder. In addition to the adjustment options provided by the adjustable holders, they offer an additional option for rotating and inclining the specimens. The incline of the specimens can be locked in place using the knurled screw (49.2).


Fig. 48 Adjustable holder with Ø 25 mm stage with ball-and-socket joint

- 1 Ø 25 mm stage with ball-and-socket joint
- 2 Clamping screw

Fig. 49 Adjustable holder with joint holder

- 1 Joint holder
- 2 Fastening screw

9.4 Universal rotating mechanism

The rotating mechanism (order No. 581 088) (Fig. 52) consists of a stage plate, which can be mounted directly onto the cross-stage, and an articulated illumination arm, which allows you to rotate the illumination around the specimen.

Either the inclining rotating stage $\varnothing$ 75 mm or the universal slide holder can be installed on the rotating plate.

9.5 Light panels

For homogeneous illumination of larger specimens, you can use the light panel (Fig. 53a, 53b). The light panel (Fig. 53a) is delivered with an articulated arm that can be attached to the universal specimen holder or onto the stage bracket.

Alternatively, the light panel (Fig. 53b) can be used in conjunction with a cold light source and a light guide.

Fig. 52 Universal rotating mechanism

Fig. 53a Light panel

1 On/off switch


Fig. 53b Light panel

9.6 Filter for color contrast

For the differentiation of color deviations of two specimens being compared, the two filters CRA (red) and CGA (green) in sliders (Fig. 54) are provided. They are inserted into the filter receptacles in the macroscope holders (55.1) above the objectives.

For adjusting the composite image with color contrast, start by setting composite imaging to **MIX**.

If the specimens of comparison overlap in their essential structures, these only appear in their original color where matching details are present.

Depending on the assignments of the filters to the two beam paths, the deviations appear in either red or green.

Fig. 54 Color filter

9.7 Additional illumination options

A wide assortment of other illumination options are available in addition to LED spot illumination (→ Page 18) and light panels (→ Page 35).

LED illumination

- LED power spotlight with one LED (3.5 W)
- LED spotlight with 19 LEDs
- LED ring illuminator (Fig. 57) with partial activation of the LED segments for near vertical illumination

To set the brightness, follow the "LED Controller Coolcontrol" operating instructions provided separately.

Fig. 55 Objective turret with objectives

1 Filter receptacle

UV illumination

The UV ring illuminator (Fig. 58) is available for UV illumination.

Cold light illumination

The cold light source, such as the Leica CLS250 LED A (Fig. 56), can be placed next to the stand or onto a shelf screwed to the base of the comparison microscope.

A double-armed, flexible light guide is used to feed in the light.

The cold light source is controlled using the rotary knobs (56.2).

Fig. 56 LED illuminator, e.g. Leica CLS250 LED A
1 Connection for light guide

Fig. 57 LED ring illuminator

Fig. 58 UV ring illuminator

9.8 Stage scales and calibration standard

During assembly of the comparison microscope, the macro-objectives have been chosen so that they match with regard to their optical properties and permissible tolerances. For identical objectives, the image scales generated thus fall within the permissible tolerances.

The stage scale (stage micrometer) 11 519 963 or the calibration standard 11 581 080 can be used to verify this. To do so, adjust the section image for identical specimen magnifications and check that the image scales on both sides match exactly.

Fig. 59 Stage micrometer

10. Technical Data

Stand

- Stable and torsion-free cast-iron stand with motorized height adjustment for the carrier of the comparison bridge (stroke 255 mm)

Focus

- Two focusing drives
- Stroke: 25 mm

Stages

- Cross-stages
- Stage surface: 140 mm x 202 mm
- Travel range: 50 mm x 50 mm
- Ergonomic operating knobs for cross and synchronized movement
- Bore hole receptacles for the adjustable holder, the rotating stages, the large sample stages or the bullet holder

Ø 118 mm stackable rotating stages

- Ø 50 mm openings
- Removable glass stage plates
- Clamping device for the stage rotation
- Stage carrier with receptacle for rotating polarizer

Ø 75 mm inclining stages

- Can be locked into place
- Inclines in all directions up to 45°
- Surface with ridged surface

Large sample stages

- Stage surface: 210 mm x 300 mm
- Metal plates for thin specimens such as documents
- Can be stacked on the cross-stage

Comparison bridge with ergonomic tube

- For superimposed image or split image comparisons with variable dividing line
- Variable width of dividing line
- Combination of split image and superimposed image possible
- Color variation of deviation in traces possible in superimposed image observations
- Maximum deviation of the magnification on the left and right: 0.4%
- Receptacles for macro-objectives 0.4x, 1x, 2x, 4x, 8x with adjustable iris diaphragms
- Rotating fasteners for the holders for oblique incident illumination
- Slots for inserting filter slides
- Distance between the optical axes: 450 mm

Tube factor

1x, 1.5x, 2x with magnification changer

Field number

FOV 22

Image orientation

Portrait or landscape

10. Technical Data

Objectives

Macro-objectives 0.4x, 1x, 2x, 4x, 8x

Eyepieces

HC PLAN S 10x/22
HC PLAN S 10x/20
HC PLAN S 10x/18

Magnifications and object fields

Magnification changer in position	Objectives	Available working distance	FOV with 10x eyepiece	Total magnification with 10x eyepiece
	0.33:1	100 mm	166.0 mm	1.32x
	0.4x	60 mm	55.0 mm	4x
1.5x	0.4x	60 mm	36.6 mm	6x
2x	0.4x	60 mm	27.5 mm	8x
	1x	60 mm	22.0 mm	10x
1.5x	1x	60 mm	14.6 mm	15x
	2x	60 mm	11.0 mm	20x
1.5x	2x	60 mm	7.3 mm	30x
	4x	60 mm	5.5 mm	40x
1.5x	4x	60 mm	3.6 mm	60x
	8x	48 mm	2.75 mm	80x
1.5x	8x	48 mm	1.8 mm	120x
2x	8x	48 mm	1.37 mm	160x

Motorized work table

- Height adjustment: 619 mm + 300 mm stroke
- Lifting capacity: 2000 N
- Lifting speed: approx. 12 mm/s
- Load capacity: 200 kg
- Control system: up and down switches and foot switch
- Stage surface: 1200 mm x 560 mm
- Universal power supply

Dimensions without camera

Height: 845 mm
(maximum; Z column extended)

Width: 700 mm
(max. extended cross-stages)

Depth: 530 mm
(incl. front control knobs)

Observation height for ballistic specimens in the universal holder: 590 mm - 510 mm

Weights

- Basic equipment: 45 kg
 Comparison bridge: 15 kg
 Macroscope stand with table: 32 kg
 Motorized work table: 50 kg

General specifications

- For indoor use only.
 Supply voltage: 100–240 VAC
 Frequency: 50–60 Hz
 Power input: max. 120 W
 Fuses: T2.5A H 250 VAC
 Ambient temperature: 15–35°C
 Relative humidity: 90% to 30 °C
 non-condensing
 Overvoltage category: II
 Pollution degree: 2

11. Care of the Instrument

Caution

Unplug the power supply before performing cleaning and maintenance work!
Protect electrical components from moisture!

Macroscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus.

The instrument should be cleaned after each use, and the optics must be kept extremely clean.

11.1 Dust cover

Note:

To protect against dust, cover the instrument and accessories with the dust cover after each use.

11.2 Cleaning

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with a low-concentrated soap solution, petroleum ether or ethyl alcohol.

For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Acetone, xylene or nitro-containing thinners can harm the instrument and thus must not be used.

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure that coated or plastic surfaces do not become matted or etched.

Cleaning the specimen stage

Remove light-colored spots on the specimen stage by rubbing with paraffin oil or acid-free Vaseline.

Cleaning glass surfaces

Remove dust on glass surfaces with a fine, dry and grease-free brush made from hair, by blowing with a squeeze blower, or vacuum suction.

Carefully remove stubborn dirt on glass surfaces with a clean cloth moistened with distilled water. If the dirt still can not be removed, use pure ethyl alcohol, chloroform or petroleum ether instead of water.

Cleaning objectives**Caution**

The objective may not be unscrewed during cleaning. If damage appears on inner surfaces, the objectives must be sent to your Leica subsidiary for repair. We also advise against cleaning the inside surfaces of the eyepieces.

The front lenses of objectives are cleaned as described under "Cleaning glass surfaces". The upper lens is cleaned by being blown off with a pneumatic pump.

Removing immersion oil**Caution**

Follow safety notes for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

11.3 Handling acids and alkalines

For examinations using acids or other aggressive chemicals, particular caution must be taken.

**Caution:**

Never allow the optics and mechanical parts to come into direct contact with these chemicals.

12. Essential Wear and Spare Parts

Order No. Material No.	Name	Used for
<u>Fuses:</u> 11 362 150 010 250	T2.5A H 250 VAC	2 x fuse (power supply)

13. EC Declaration of Conformity

To download the EC Declaration of Conformity, use this link

<http://www.leica-microsystems.com/products/light-microscopes/forensic/>

Then select the corresponding instrument and go to the "Download" page.

- Administrative Measures on the Control of Pollution Caused by Electronic Information Products -

部件名称 Name of the part	有毒有害物质或元素 Hazardous substances					
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)
印刷电路板 printed circuit boards	X	0	0	0	0	0
电子元件 electronic components	X	0	0	0	0	0
机械部件 mechanical parts	X	0	0	X	0	0
光学元件 optical components	X	0	X	0	0	0
电缆 cables	0	0	0	0	X	X
光源 light sources	0	X	0	0	0	0

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Indicates that the concentration of the hazardous substance in all materials in the parts is below the relevant threshold of the SJ/T 11363-2006 standard.

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Indicates that the concentration of the hazardous substance of at least one of all materials in the parts is above the relevant threshold of the SJ/T 11363-2006 standard.

Note: The actual product may or may not include in all the part types listed above

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