From Eye to Insight

K5

Instructions for use

Leica Microsystems CMS GmbH. Instructions for use, 11934174, Revision 1.0 from 2019-06-29

Copyrights

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other methods (including electronic systems) is not allowed without written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Still, we are always grateful for comments and suggestions regarding potential improvements within this documentation.

The term „Windows“ may appear in the following text without further identification. It is, however, a registered trademark of Microsoft Corporation. The names of companies and products used herein may be trademarks of their respective owners.

The information in this manual is subject to modification at any time and without notification.

Contents

- 1. Introduction.....3**
- 2. Intended Use of the Device5**
- 3. Safety Notes.....6**
 - 3.1 EU-Declaration of Conformity6
 - 3.2 Technical Specifications7
 - 3.3 Tolerable Ambient Conditions.....10
 - 3.4 Notes on Cleaning.....11
 - 3.5 Disposal11
- 4. Notes for the User12**
- 5. Components13**
 - 5.1 Standard Components.....13
 - 5.2 Connectors14
 - 5.3 Computer Requirements15
- 6. Installation16**
 - 6.1 Installing Leica LAS X Software16
 - 6.2 Installing the Camera Framegrabber Card 17
 - 6.3 Important Installation Information for USB 3.0 Camera18
 - 6.4 Installation of Leica Camera Hardware.....19
- 7. Troubleshooting21**
- 8. Service23**

Revision 1.0, published 2019-06-27 by:

Leica Microsystems CMS GmbH
Ernst-Leitz-Straße 17-37
D-35578 Wetzlar (Germany)
<http://www.leica-microsystems.com>

Responsible for contents:
Marketing CMS

1. Introduction

Attention!

This User Manual is an essential component of the instrument and must be read carefully before assembly, commissioning or use.

This instructions for use contains important instructions and information for the operational safety and maintenance of the instrument and accessories. Therefore, it must be kept in a safe place.

Advantages of the K5

The K5, is a high-value sCMOS camera that offers 80%QE, high 4.2megapixel resolution at 6.5µm/pixel, fast 40fps imaging, 2.1e- read noise, stable passive cooling, and a full-well capacity of 45000e-. The K5 can be triggered with additional experimental factors to minimize phototoxicity, and improve physiologically relevant observation. These features come together to provide an excellent camera for routine fluorescent imaging applications.

Signs and Symbols

(1.2) Numbers in parentheses correspond to illustrations (in the example, Figure 1, Item 2)

→ p.20 Numbers with pointer arrows point to a certain page of this manual

Caution!
Special safety instructions within this Instructions for use are indicated with the triangle symbol shown here, and have a gray background.

Notice! The camera and accessories can be damaged when operated incorrectly

Instructions on disposing of the device, accessory components and consumables.

Explanatory note

Item not included in all configurations (optional)

Warning of electromagnetic field

MM/YYYY

Manufacturer's address with manufacturing date

China RoHS 50 years EFUP
(Environmentally friendly use period)

Warning of hazardous electrical voltage! Risk of electrical shock!

WARNING indicates a hazard with a medium degree of risk that, if not avoided, can result in death or serious injury.

CAUTION indicates a hazard with a low degree of risk that, if not avoided, can result in minor or moderate injury.

2. Intended Use of the Device

This camera system is designed for use by technicians, engineers and scientists. It is a scientific measuring instrument, which provides images. The camera may only be used according to the instructions of this manual. The disclosures and operating conditions in these operating instructions installation must be respected. Unauthorized modifications and changes of the device are forbidden for safety reasons.

Caution!

This instructions for use is an essential component of the camera, and must be read carefully before the camera is assembled, put into operation or used, and must be kept for later references

Directives of European Community (EC-Directives)

The camera K5 complies with the requirements in EU Directive 2014/30/EU regarding electromagnetic compatibility for use in industrial environments.

Reasonably foreseeable misuse

- ▶ Cleaning the camera in a way other than specified in the manual
- ▶ Allowing unauthorized personnel to open the camera
- ▶ Installing the camera in a way other that Leica Microsystems CMS GmbH has specified in the manual
- ▶ Using cables that Leica Microsystems CMS GmbH has not provided or permitted

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from using the device for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In the event of such, the declaration of conformity shall become null and void.

This device is not intended for use in the patient environment defined by DIN VDE 0100-710. Neither is it intended for combining with medical devices according to EN 60601-1. If this device is electrically connected to a medical device according to EN 60601-1, the requirements defined in EN 60601-1-1 shall apply.

3. Safety Notes

This camera is constructed and tested in accordance with EN 61326-1/ IEC 61326-1, Electrical equipment for measurement, control and laboratory use – EMC requirements.

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this operating manual.

The instruments and accessories described in this manual have been safety-tested and checked for possible hazards.

The responsible Leica affiliate or the main plant in Wetzlar must be consulted whenever the instrument is altered, modified or used in conjunction with non-Leica components that are outside of the scope of this manual.

Unauthorized alterations to the device or noncompliant use shall void all rights to any warranty claims and product liability!

3.1 EU-Declaration of Conformity

To download the EU Declaration of Conformity use the link

<http://www.leica-microsystems.com/products>
Select your product and switch to the Download page.

3.2 Technical Specifications

Power to the camera via an external power supply.


For indoor use only.

Image Sensor	
type of sensor	customized sCMOS
resolution (h x v)	2048 x 2048 pixel
pixel size (h x v)	6.5 μm x 6.5 μm
sensor format / diagonal	13.3 x 13.3 mm / 18 mm
shutter mode	Rolling Shutter
fullwell capacity (typ.)	45000 e ⁻
readout noise (typ.)	2.1 e ⁻ med / 2.3 e ⁻ rms
dynamic range	up to 80%
spectral range	370 nm ... 1100 nm
dark current (typ.)	15 e ⁻ /pixel/s @ 21°C
DSNU	0.5 e ⁻ rms @ 21 °C
PRNU	0,60 %
Camera	
frame rate	40 fps @ full resolution
exposure / shutter time	10 μs ... 5 s
dynamic range A/D	16 bit
A/D conversion factor	e ⁻ /count 0.65
region of interest	64x16 minimum ROI
binning (h x v)	1x1 ... 4x4
non linearity	< 0.6 %
cooling method	passive cooling
trigger input signals	frame trigger, acquire (SMA connectors)
trigger output signals	exposure, busy (SMA connectors)
data interface	USB 3.1 Gen1
timestamp	in image (1 μs resolution)

3. Safety Notes

General	
power delivery	via USB Type-C
power consumption	< 4.5 W
weight	420 g
operating temperature	+10°C ... +40°C
operating humidity range	10% ... 90% (non-condensing)
storing temperature range	-10°C ... +60°C
optical interface	C-mount, F-mount (optional)
CE / FCC certifie	yes

3.2.1 Quantum efficiency curv


3.2.2 Mechanical dimensions

**Note:**

C-Mount flange focal distance: 17,52 mm

All dimensions given in millimeter.

3. Safety Notes

3.3 Tolerable Ambient Conditions

Temperature: +10°C ... 40°C
Relative humidity: 10% ... 90%
non-condensing

Caution!

Do not use the device in altitudes exceeding 2000 m ASL/NL.

Transport and storage in a range of -10° – +60°C and at a humidity not exceeding 70% (non-condensing).

Make sure that a free air circulation is provided

Caution!

The electrical components are not protected against water. Water can cause electric shock.

Electrical components must be assembled at least 10 cm from the wall and away from flammable substances

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

3.4 Notes on Cleaning

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Exterior Body Cleaning

1. Unplug the power supply before performing any cleaning.
2. Protect the window of the camera from damage or liquid splashes.
3. The exterior surfaces can be cleaned by gently wiping down with a lint-free wipe moistened with water.

Window Cleaning

1. Blow off any loose contamination with clean filtered air under a low flow rate.
2. Add a drop of methanol to a Cleanroom grade swab (eg Texwipe TX761), taking care not to touch the swab tip with fingers. Shake to remove any excess solvent.
3. With minimal pressure, gently clean the window from the centre outwards in short quick strokes taking care not to touch the glue at the perimeter of the window. If contact is made with the glue, replace the swab immediately.
4. Repeat the process as often as required, replacing the swabs after approximately ten seconds or ten wipes to prevent smearing marks.

Caution!

Acetone, xylene or nitro-containing thinner can harm the device and thus may not be used.

Test clean solutions of unknown composition first on a less visible area of the unit. Be sure that coated or plastic surfaces do not become matted or etched.

3.5 Disposal

To dispose of the product at the end of its service life, please contact Leica Service or Sales.

Please observe national laws and regulations, such as those implementing and enforcing the WEEE EU directive.

Note!

Like all electronic instruments, the device, its components and expendables may not be disposed of as general household waste!

4. Notes for the User

Caution!

This Instructions for use is an essential component of the product; it must be read carefully before the product is assembled, put into operation or used, and must be kept for later reference.

Further pay attention to the following points:

- ▶ Follow the directions for safety and care.
- ▶ Read the user manual for your microscope or macroscope.
- ▶ Ensure that the K5 Camera is used only by persons authorized to do so.
- ▶ Ensure that the K5 Camera user manual is always available.
- ▶ Check on a regular basis that the users are complying with safety requirements.
- ▶ If a defect that could cause personal injury is detected, inform your local Leica agency.
- ▶ Ensure that electrical connections are plugged in tightly during use.

Caution!

Only insert or unplug data and control circuits when the instrument is switched off; otherwise, the instrument may be damaged.

Note:

Illustrations, descriptions and technical data are not binding and may be changed without notice. You will find the latest release of this user manual and other information on our home page:

www.leica-microsystems.com

5. Components

5.1 Standard Components

Description		Pcs
	K15 Camera • C-mount optical connection • F-mount adapter optional • For standard C-mount / F-mount lenses and adapters	1

Description		Pcs.
	USB 3.1 Interface Card • 2x USB Type A socket • PCI Express x1 V2.0	1
	Camera to PC Interface cable USB 3.0	1 x 3 m (USB)

5. Components

5.2 Connectors

(at the rear side of the camera head)

1 Status LEDs

(LED1: green / red; LED 2: orange)

Color	Description
Red / green continuous	Camera is booting
Green	Camera is ready for operation
Orange continues	Arm camera / rec state off
Orange blinking	Recording on
Red continues	Error

2 Interface Connect

USB 3.1 plug Type-C screwable; Power delivery for the camera over USB cable

3 SMA Input

Input (1=Exposure Trigger; 2=Acquire Enable)	
Type	Digital
Level	3.3 V LVTTTL (5 V tolerant)
Coupling	DC
Impedance	1 kΩ
Slew rate	> 1 V/ms

4 SMA Output

Output (3=Status Busy; 4=Status Expos)	
Type	Digital
Level	3.3 V LVTTTL
Coupling	DC
Load current	Max. 100 mA

5.3 Computer Requirements

Processor	Intel® Core™ i7
Operating System	Windows 7/10 (64 bit)
Ports	USB 3.1 Gen 1
Random access memory	Min. 8 GB DDR3
Media	Hard disk to install the Leica LAS X software.
Monitor	Full-HD resolution display

6. Installation

Note:

Before installing any software, please ensure that your computer meets the minimum system requirements for the K5 Camera ("5.3 Computer Requirements" p. 15).

Install service packs to update your operating system, add memory, and, if necessary, clean up your hard drive.

Install the software and hardware only on a stable and well organized computer system. To prevent any data loss, always make backup copies of all your important data before you start installing.

Note:

Before you can use the K5 camera, check that the camera is connected properly to a port on your computer.

Caution!

The K5 only functions properly with (USB 3.0 PCI card and USB 3.0 cable) specified by Leica.

For Windows® 7, you are required to have Administrative or Power-User privileges to be able to install this software. If in doubt consult your System Administrator.

6.1 Installing Leica LAS X Software

Download the latest version of microscope imaging software from Leica Microsystems:

Link: <https://webshare.leica-microsystems.com/latest/>

User name: releases

Password: Leica1997!

The installation wizard guides you through the installation process.

During this process you are requested to install the camera drivers. Select the K5 camera.

After complete installation restart the computer.

Note:

The camera driver can be installed separately at any time later.

6.2 Installing the Camera Framegrabber Card


Caution!

Ensure your PC is powered off. Ensure anti-static precautions using anti-static strap provided. Keep all screws safe. The electronic circuits on the PCI card and the circuits of your computer are sensitive to static electricity and surges.

1. Unplug all cables from the rear of the computer.
2. Remove the PC enclosure so as to gain access to the expansion slots.
3. Locate a PCIe x16, x8 or x4 slot.
Framegrabber card operates in a x16, x8 and x4 slot ONLY.
If you are unsure which PC slot is correct, please consult the PC user manual.
4. Remove the filler bracket corresponding to the slot you intend to use.
5. Making sure that the card's mounting bracket is flush with any other mounting brackets or filler brackets to either side of it, secure the controller card in place.
6. Replace the cover of the computer and secure it with the mounting screws if applicable.
7. Reconnect any accessories you were using previously.

1 Fasten the Card to the slot with a screw

6.3 Important Installation Information for USB 3.0 Camera

- ▶ For USB 3.0 cameras, the drivers for the PCIe card must be installed manually by the user as they are not part of the software installation. Refer to the instructions on the driver CD supplied with the USB 3.0 Card (Install Renesas driver).
- ▶ Within the provided installation files USB_HBA, open the folder U3X4-PCIE4XE101, U3X4-PCIE1XE101, U3-PCIE1XG202.
- ▶ Open the subfolder Driver and run RENESAS-USB3-Host-Driver-30230-setup.exe.
- ▶ If your current OS is Win7/8 and the User Account Control is enabled, a dialog prompts whether you wish to start the setup: accept with Yes.

6.4 Installation of Leica Camera Hardware

6.4.1 Attaching the Camera to the Microscope

Use a Suitable Video Objective

The K5 Camera can be attached to any microscope with a C-mount interface.

Suitable C-mount is 1x

- ▶ Screw the appropriate C-mount (or video objective for stereo microscopes) on the camera
- ▶ Plug-in all cables on the backside of the camera (see p. 20) (USB3.0 and if necessary the optional trigger cable).

Caution!

Only insert or unplug data and control circuits when the instrument is switched off; otherwise, the instrument may be damaged.

- ▶ Fasten the C-mount with the attached camera at the documentation port (e.g. phototubes or sideports).
- ▶ Position the camera so that the cables point backwards and do not touch heat emitting devices like light sources etc.

Note:

Via the Leica software LAS X the image can be flipped in x and y direction if necessary.

Positioning the K5 Camera

To view the image upright and the right way round on the monitor, you must first position the camera correctly.

The K5 Camera fits onto the side of the HD videotube or phototube and onto the top of all other videotubes or phototubes (HV, HU and A) and trinocular tube.

6.4.2 Connecting the camera and switching on

Caution!

Ensure that the PC is switched OFF.

- ▶ Securely connect the Camera to PC cable between the appropriate connector(s) on the rear of the camera and the corresponding port on the PCIe card. (See Figure below).

Notes:

- ▶ For the K5 camera, use only the supplied USB 3.0 cable to ensure performance.
- ▶ If required connect the trigger cable to the 15 way I/O connector.
- ▶ Turn on the PC.
- ▶ Start your software.
- ▶ Start acquiring data.

- 1 Status LED
- 2 USB 3.1 plug Type-C screw-able
- 3 SMA Input
- 4 SMA Output

7. Troubleshooting

This chapter is intended to provide you with potential solutions in case the imaging system does not react the way you expect. Please contact your local Leica supplier if you need further support.

The live image window is black in LAS X

- ▶ Check that all cables are correctly connected.
- ▶ Ensure that the port to which the camera is connected and the port you have selected in the Camera Configuration of LAS X match.
- ▶ Check the light conditions and the microscope settings. Ensure that light reaches the camera.
- ▶ Inspect the shutter settings, if necessary increase exposure time.

Image is too dark or too bright

- ▶ Check the exposure time, if necessary increase or reduce the exposure time (LAS X software).
- ▶ Check the microscope settings.
- ▶ Check the lighting settings.

Poor image quality

- ▶ Check the microscope settings.
- ▶ Check the exposure settings to avoid under-exposure or over-exposure (LAS X software).
- ▶ Check the image histogram settings. The image histogram should spread across the middle of the field (LAS X software).

Poor image quality - noise

- ▶ Avoid stray light through the eyepieces. For inverted microscopes cover the sample if necessary.
- ▶ Adjust the illumination correctly.
- ▶ Check the Gain setting. If you increase the Gain setting, you should expect an increase in noise. If you reduce the Gain setting, you should increase the exposure time (LAS X software).

The Camera is not displayed as selectable camera

- ▶ Check that all cables are connected properly
- ▶ Ensure that the camera is assigned to the appropriate port in the Hardware Configurator of the Leica LAS X
- ▶ Check if the camera is displayed in the hardware manager of Windows®. If not, you might install the USB3.0 driver again and install the Leica LAS X software again.

7. Troubleshooting

Camera buzzer does not sound on start-up

- ▶ The camera buzzer should be audible momentarily when the camera is switched on.
- ▶ If this does not occur, ensure that power is connected to the camera and the On/Off switch is set to On

Camera is not recognized by PC

- ▶ Ensure the camera is switched on
- ▶ For USB 3.0 models, check that USB 3.0 cable is connected between the camera and the USB 3.0 slot on the PC.
- ▶ Check that the frame grabber card and drivers have been correctly installed (see p. 18)

Buzzer sounds continuously

This indicates that an “over temperature” condition has occurred within the camera. Follow the instructions below to rectify this situation:

- ▶ Power the camera off and allow it to cool down
- ▶ Check the camera is operating within the specified environmental conditions (see p. 7)

Camera does not cool to the required temperature

- ▶ Check that the operating ambient temperature is within allowable limits (see p. 7) when cooling the sensor to 0°C.
- ▶ Check that the camera vents are not blocked and have sufficient clearance to allow air flow.

8. Service

Service request

If technical problems arise or you need help using your K5 Camera, contact the local Leica representative.

To answer your questions effectively, please provide as much detailed information as possible. The following information is particularly important:

- Computer: type, model, speed, memory
- Operating System: version
- Software: version, language
- Screenshot of the error message
- Description of exactly when and where the problem arises
- Log file (LAS X) if available

Caution!

Mechanical, electronic and/or optical maintenance may only be carried out by Leica-trained maintenance technicians.

Only original Leica spare parts may be used.

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○：表示该有害物质在该部件所有均质材料中的含量均在 GB/T 26572 规定的限量要求以下。
Indicates that said hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement of GB/T 26572.

×：表示该有害物质至少在该部件的某一均质材料中的含量超出 GB/T 26572 规定的限量要求。
Indicates that said hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement of GB/T 26572

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The “Hazardous Substance Marking Table” covers the here listed products.

显微镜	控制	光源	光学和照相机	电源和服务模块
Microscopes	Controls	Light Sources	Optics and Cameras	Power Supply and Service Modules

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