
From Eye to Insight

LEICA LAS X

Imaging and Analysis Software for Widefield Microscopes

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
From Eye to Insight

LEICA LAS X

Imaging and Analysis Software for Widefield Microscopes
Instructions for use

Leica Microsystems CMS GmbH. Instructions for use, 11934003, Revision 1.2, 2017-01-05

COPYRIGHTS

The instructions contained in the following document reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Nevertheless, no liability of any kind may be assumed for the accuracy of this manual's contents. If you have any comments on these operating instructions or on any of our other documentation we would be pleased to receive them. The information in this manual is subject to modification at any time and without notification.

All rights to this documentation are held by Leica Microsystems CMS GmbH. Adaptation, translation and reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

Programs such as LAS and LAS X are protected by copyright laws. All rights reserved. Reproduction, adaptation or translation of these programs is prohibited without prior written permission from Leica Microsystems CMS GmbH.

The term "Windows" may be used in the following text without further identification. It is a registered trademark of the Microsoft Corporation. Otherwise, no inference with regard to the free usability of product names may be drawn from the use of those names. All other brand names and product names in this document are brands, service marks, trademarks or registered trademarks of the respective manufacturers.

Rev 1.2, published 2017-01-05 by:

Leica Microsystems CMS GmbH
Ernst-Leitz-Straße 17-37
D-35578 Wetzlar (Germany)
www.leica-microsystems.com

Responsible for contents:
Marketing CMS

CONTENTS

1. Important Notes about this Manual.....	4	6. First Experiment.....	27
2. Intended Use.....	6	6.1 Configuration.....	27
3. The Leica LAS X.....	7	6.2 Switch on the system	27
3.1 Description	7	and start the software	27
3.2 Dongle	7	6.3 Selecting and focusing specimen section	28
4. Starting up Leica LAS X.....	8	6.4 Camera resolution	30
4.1 Open Leica LAS X	8	6.5 Acquisition mode.....	31
4.2 Calling up the online help.....	8	6.5.1 EL6000 variant	31
4.3 Full-text search	9	6.6 Image optimization	34
with logically connected search terms ..	9	6.7 Setting color (LUT).....	35
4.4 General structure of the user interface	10	6.8 Settings for additional fluorochromes ..	36
5. Parameters for Multidimensional Image Acquisition.....	12	6.9 Capturing 3D image stack.....	37
5.1 Light Path Settings	13	7. Index.....	40
5.2 Image menu.....	17		
5.2.1 Camera	17		
5.2.2 EM CCD cameras.....	19		
5.3 z-menu	19		
5.3.1 Setting the image stack with Wide Focus.....	20		
5.3.2 Setting the image stack with Fine Focus*	20		
5.4 t-menu	20		
5.5 Mark and Find	23		
5.6 Tile Scan.....	24		
5.6.1 Camera adjustment for Tile Scan.....	24		
5.6.2 Tile Scan Acquisition	25		

1. Important Notes about this manual

1. IMPORTANT NOTES ABOUT THIS MANUAL

Attention!

These operating instructions are an essential component of the Leica LAS X and must be read carefully before the system is put into operation or used. Read and comply also with all operating instructions accompanying the individual microscope, components and peripheral devices. The emphasis is on the safety instructions, which must be strictly observed during all work with the Leica LAS X .

Keep these instructions in a convenient, safe place.

These operating instructions present the first steps for activating the system and describe the basic functions of the Leica Application Suite X (LAS X) software.

For further descriptions of the software functions, please refer to the online help of the LAS X, which presents the most up-to-date explanations of the corresponding software functions.

Please read Chapter 4 LAS X to familiarize yourself with the structure and use of the online help.

Text symbols, pictograms and their meanings:

(1.2)

Numbers in parentheses, such as “(1.2)”, correspond to illustrations (in the example, Figure 1, Item 2).

p.20

Numbers with an arrow, e.g. → p.20, refer to a specific page of the manual.

Attention!

Special safety instructions within this manual are indicated with the triangle symbol shown here, and have a gray background.

Attention! The microscope and accessories can be damaged when operated incorrectly.

Explanatory note.

Item not contained in all configurations.

2. Intended Use

2. INTENDED USE

Application

This software has been developed for observing and documenting magnified depictions of samples acquired using Leica microscopes in a variety of routine work and research applications. The product is not considered a medical product or in-vitro diagnostic product.

Do not use this product for any use other than the intended application.

Attention!

The manufacturer assumes no liability for damage caused by, or any risks arising from, using the software for purposes other than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

Quality directives

It is expected that a laboratory that uses LAS X software will work in accordance with the directives encompassing the following (and other) areas: User training, assay development and validation and external audits of laboratory test performance.

Application methods must be reviewed before implementation and documented. Standard operating procedures (SOPs) for analysis methods should be introduced.

3. THE LEICA LAS X

3.1 Description

The Leica LAS X (Leica Application Suite X) offers solutions for different application levels. Each solution is upgradable and allows a wide range of possible combinations with hardware and software modules from Leica Microsystems CMS GmbH.

The complete system unifies the latest technology for image acquisition, perfectly coordinated – for fast results in highest quality.

**Note:**

The manufacturer assumes no responsibility or liability for any use outside of the intended use or use outside of the specifications from Leica Microsystems CMS GmbH or any risks resulting from such use. In such cases, the EU Declaration of Conformity shall be invalid.

**Attention!**

Be sure to follow the included operating instructions for the Leica microscope.

**Attention!**

Do not install any other hardware or software on the workstation, as otherwise serious damage to the system or loss of data can result. If such work is required, contact the Technical Service of Leica Microsystems CMS GmbH.

**Attention!**

To prevent loss of data, make regular backups of your data and images on a suitable data carrier.

Also, be sure to back up your data before any service or repair work is performed! Leica Microsystems CMS GmbH shall not be liable for any loss of data.

3.2 Dongle

To enable the function of the software and the system, a USB dongle that unlocks the software is included with the system. Remove the USB dongle only if you want to prevent anyone from starting the software.

4. Starting up Leica LAS X

4. STARTING UP LEICA LAS X

4.1 Open Leica LAS X

Attention!

Do not start the software until the microscope and all other components have been switched on!

- Click the LAS X icon on the desktop to start the software:

Fig. 1: LAS X icon on the desktop

Now the Leica LAS X system is ready for operation.

4.2 Calling up the online help

The online help can be called up in three ways:

- In the respective context (context-sensitive)
- Via the Help menu
- With the key combination CTRL + F1

In the respective context (context-sensitive)

Click on the small question mark located in the top right corner of every dialog window. Online help opens directly to the description for the corresponding function.

Via the Help menu

Click on the Help menu on the menu bar. The menu drops down and reveals, among others, the following search-related options:

Contents This dialog field contains the table of contents in form of a directory tree that can be expanded or collapsed.

Double-click an entry of the table of contents to display the corresponding information.

Index Enter the term to be located. The online help displays the keyword that is the closest match to the specified term.

Select a keyword. View the corresponding content pages by double-clicking the key word or selecting it and then clicking the Display button.

Search Enter the term or expression to be located and click on LIST TOPICS. A hierarchically structured list of topics is displayed.

About Opens the **User Configuration** dialog box, where you can select the language and other items, in which the online help is shown.

4.3 Full-text search with logically connected search terms

Click on the triangle to the right of the input field on the Search tab to view the available logical operators.


- Select the desired operator.
- Enter the second search term you would like to associate with the first search term behind the operator:

Examples	Result
Pinhole AND sections	This phrase finds help topics containing both the word "pinhole" and the word "sections".
Pinhole OR sections	This phrase finds help topics containing either the word "pinhole" or the word "sections" or both.
Pinhole NEAR sections	This phrase finds help topics containing the word "pinhole" and the word "sections" if they are located within a specific search radius. This method also looks for words that are similar in spelling to the words specified in the phrase.
Pinhole NOT sections	This phrase finds help topics containing the word "pinhole" and not containing the word "sections".

4. Starting up Leica LAS X

4.4 General structure of the user interface

The user interface of the LAS X is divided in five areas:

(1) Menu row and the area for selection of the **Operating Mode/Wizard**:

The various menus for calling up functions are available here.

- File
- Help

and the **Wizard** with:

- Main
- Live Data Mode
- FRET SE
- HCS A (MatrixScreener) for High Content Screening

(2) Arrow symbols: Operating step with the individual functions. These operating steps mirror the typical sequence of image acquisition and subsequent image processing. The functions are arranged appropriately in these operating steps.

- Configuration
- Acquire
- Process
- Quantify

(3) Tab area: Different tabs belong to every operating step (arrow symbol) in which the settings for the experiment can be made or image data can be selected.

Acquire Experiments: Directory tree of the files stored in memory (RAM)

Setup: Hardware settings for the current experiment

Acquisition: Parameter settings for the image acquisition

Process Experiments: Directory tree of the files stored in memory (RAM)

Tools: Directory tree with all the functions available in the respective operating step

Quantify Experiments: Directory tree of the files stored in memory (RAM)

Tools: Tab with the functions available in this operating step

Graphs: Graphical display of values measured in regions of interest (ROI)

Statistics: Display of statistical values that were determined in the drawn in regions of interest (ROI)

(4) Working area: This area under the Acquire operating step provides the Light Path Settings dialog window in which the control elements for setting the recording parameters are located. For image editing, previews are created and displayed in this area.

(5) Viewer window: Displays the recorded images. In the standard setting, the Viewer window consists of the image window in the center and the buttons for image viewing and selection of regions of interest (5a) and channel display (5b).

5. Parameters for Multidimensional Image Acquisition

5. PARAMETERS FOR MULTIDIMENSIONAL IMAGE ACQUISITION

This chapter describes the software menus used for defining the acquisition parameters of the Leica LAS X. The following **acquisition mode** is available when you click the **Acquire** arrow and **Acquisition** in the tab area:


Fig. 2: Acquisition Mode

Buttons **x** and **y** are always active. Both buttons refer to the **Image** menu, which is used for setting camera parameters.

Button **z** activates the **z-menu**, which defines three-dimensional image stacks.

Button **t** activates the **t-menu**, which defines the time series.

The following two buttons, **Tile Scan** and **Mark and Find** refer to the motorized specimen stage:

This button activates the **Tile Scan menu**, which is required for generating composite overview images:

This button activates the **Mark and Find menu**, which is required for defining various stage positions for multiposition experiments:

Note:

A black button means that the function is disabled and is therefore not used for that experiment. The button can be enabled by clicking it with the left mouse button. A red button means that the function is enabled.

5.1 Light Path Settings

The Light Path Settings menu is displayed on the left side of the Leica LAS X user interface. The light channels for the various experiments can be defined in this menu. This menu is important and is therefore always displayed:

Fig. 3: Light Path Settings

5. Parameters for Multidimensional Image Acquisition

Note:

When defining an experiment always start with setting the image channels in the **Light Path Settings** menu.

To define a channel the associated button must be enabled first. Start with channel 1:


Select the LUT (look-up table) wanted for this channel from the color palette, which can be opened with the following button:

The desired contrast method is selected in the **Contrast Method** window, and the filter cube to be used for this channel is selected in the **Filter Cube** field.

Fig. 4 External adapter coupling with two short filter sliders

1 Cable connection number

Using Leica external filter wheels* :

The Leica LAS X can control up to four external filter wheels, which are operated simultaneously. (see fig 4).

The filter positions are defined in the **Light Path Settings** menu.

If the system is equipped with the external **FFIM**, the transmission is always controlled via the intensity slider in the **Image** menu (in five steps). The internal filter wheel of the microscope is set to 100% transmission in this case.

Fig. 5: Adjusting transmission for FFIM

Steps for FFIM:

100%, 50%, 25%, 12.5%, 6.25%

* Requires high-speed control box


5. Parameters for Multidimensional Image Acquisition

Selecting the filter cube in the LAS X Light Path Settings menu activates the corresponding filter wheel.

Attention: The suitable filter wheel must be positioned in the light path before!

Now select the appropriate emission or excitation filter in the **Emission Filter** field or **Excitation Filter** field. The filter wheels move to a blank position (100% transmission), if no filter wheel is activated.

Fig. 6: Select excitation filter

5. Parameters for Multidimensional Image Acquisition

The external filter wheels can be operated in **high speed mode**. In this case the option **Use sequencer board** is marked (Fig 15). If the whole experiment is controlled by the high speed board, the square in front of **Sequencer can control the whole sequence** turns from red to green.

Another channel can be created by clicking the "+" button. Channels can be deleted with the "-" button. If you wish to change the **sequence of channels**, the position of the channel can be changed by clicking the relevant **FCr button** with the right mouse button (in this example FCr 1).

Note:

You can define a maximum of eight channels.

The channel can be saved by clicking the **FCr button** with the right mouse button. Saved channels can be called in the **list field** (to the right of the **FCr button**).

The combination of all channels including their **acquisition parameters** can be saved in the **Load/Save Settings** field and also opened from there.

Fig. 7: Example for a three fluorescence channel experiment

5.2 Image menu

The **Image** menu is used to set the camera parameters:

Fig. 8: Image menu

The most important parameters are **Binning** and **Exposure**.

Binning Camera pixels can be read out together to improve the signal/noise ratio. For example, in **Binning 2x2** 4 pixels (2 pixels in x-direction and 2 pixels in y-direction) are combined. This reduces the image resolution by a factor of four. **Binning 2** or **None** are typical selections.

Exposure The **Exposure** slide is used to set the integration time separately for each image channel. To be able to carry out image acquisition in fast realtime mode (**high speed**), the integration time and also **Gain** and **Intensity** should be set to the same value for all channels.

Gain Allows the electronic camera gain to be set.

Intensity Allows the light intensity in the illumination beam path to be set. The internal filter wheel of the Leica microscope is used for light reduction.

Press the  button to activate a **live histogram**. The live histogram is used with objectives with a compensation ring as an aid for correcting image errors.

5.2.1 Camera

If the **high-speed control box** is used in the **Leica LAS X**, the camera can be operated in **high speed mode** (trigger-controlled with the aid of the sequencer).

Prerequisite: The camera must be supported by the high-speed control box!

However, note that for all cameras used in high speed mode, gain and binning need to be the same for all channels. The intensity setting may be different when using the Leica external filter wheel FFIM in high speed mode. The option **Use sequencer board** is marked (Fig 9).

5. Parameters for Multidimensional Image Acquisition

If the whole experiment is controlled by the high speed board, the square in front of **Sequencer can control the whole sequence** in the **Light Path Settings** menu turns from red to green.


Fig. 9: Experiment Settings

The **Experiment Settings** dialog box is opened with the **Acquire** arrow icon in the **Acquisition** tab.

This is relevant when working with multiple fluorochromes (FCrs). The **Same Exposure for all probes** option is available to prevent accidental use of different exposure times. A change of the exposure time for one fluorochrome is automatically applied to all fluorochromes:

Fig. 10: Same Exposure for all probes

Furthermore **Single image mode** can be activated in the **Experiment Settings** dialog box. In this case the **Single image** button at the bottom of the working area appears. It is used to capture one image of the current channel, whereas the **Capture Image** button takes one image of each defined channel.

For shutter operation mark the appropriate option at **Shutter control of sequencer***.

! **Attention!**

Please follow the included operating instructions for the individual components.

* Requires high-speed control box

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To adjust the EM CCD camera use the **Gain** and/or the **EM Gain** slider in the **Image** menu.


Fig. 11: Image menu for EM CCD cameras

When using the high speed mode (= trigger-controlled with the aid of the sequencer), note that for these cameras gain, EM gain, binning and intensity must be the same for all images of the sequence. Intensity can be different if the Leica external filter wheel FFIM is used for attenuation.

5.3 z-menu

The **z-menu** is enabled with the **Z** button and is used for **defining three-dimensional image stacks**:


Fig. 12: Definition of the z-Stack

The **z-Position** can be changed with the microscope z-drive.

For a faster and more precise fine focus the Leica LAS X uses either the **Super Z Galvo Focus*** (movement of the specimen in the z-direction independently of the objective) or the **Piezo Focus*** (movement of the objective in the z-direction).

Wide Focus (microscope z-drive) or, if available, **Fine Focus** (SuperZ Galvo Focus or Piezo Focus) can be selected in the dropdown window of the z-menu.

* Requires high-speed control box

5. Parameters for Multidimensional Image Acquisition

5.3.1 Setting the image stack with Wide Focus

z-Position is moved by the **z-drive** on the microscope or the **z-drive** of the **SmartMove**. As soon as the lowest image plane appears in the live camera image this **bottom z-plane** can be marked with the **Begin** button. If the bottom z-plane is saved, the **Begin** button is shown in red. The **top z-plane** of the image stack is marked with the **End** button. The **Begin** and **End** buttons can be used for both the top z-plane and the bottom z-plane. The **Focus** button can be used to save a **middle image plane**.

Note:

The image planes are saved as soon as the corresponding button (**Begin**, **End** and **Focus**) are shown in red. The **current z-plane** is shown by the ocher-colored area in the z-menu (see fig.12).

The **overall height of the image stack** is displayed as the **z-Volume** in the bottom half of the z-menu. The preset **System Optimized** button suggests the best **Step in z** in accordance with the scanning theorem depending on the objective in use. The user can also overwrite the **step** or **number of z-planes** by enabling the **no. of steps** button instead of **System Optimized**.

5.3.2 Setting the image stack with Fine Focus*

The image stack is set with **Fine Focus** in the same way as with **Wide Focus** (see chapter 6.3.1). Instead of the microscope z-drive and the SmartMove the **Fine Focus** is adjusted with the **mouse pointer** only. The left mouse button is used to drag the ocher-colored area in the z-menu to the desired position.

The **Go to** button can be used to approach the previously defined focus planes (bottom, top and center).

5.4 t-menu

The **t-menu** is enabled with the button and is used for **defining time series**:


Fig. 13: t-menu

A fixed number of time points can be defined with the **Cycle** button. If **Duration** is selected, the total duration of the experiment can be specified. If **Acquire until Stopped** is selected there is no specific end time for the experiment. The user can stop the experiment at any time.

* optional accessories. Requires either Super Z Galvo or Piezo Focus

5. Parameters for Multidimensional Image Acquisition

The time difference between the time points is defined in the **Time Interval** input field. If the **Minimize** button is activated, the minimum time that is required for acquisition of a time point is selected as the time interval.

If the **High Speed*** button is enabled, the complete time series is run in realtime with the aid of the high-speed control box (fastest possible mode without pauses between the images). The following requirements must be specified:

- All image channels must have the same values for **Gain** and **Intensity**. Intensity can be different if the Leica external filter wheel FFIM is used for attenuation.
- All image channels must use the same filter cube or the same contrast method.
- **Fine Focus**** is enabled in the **z-menu** (if a **z-image stack** is to be acquired).
- **Minimize** must be enabled in the **t-menu**.

If these requirements are met, the **High Speed** field can be enabled. In the **Light Path Settings** menu a green field shows that the **high-speed control box** controls the complete experiment:


Fig. 14: Display field
"Sequencer can control the whole experiment"

After acquisition of every image the associated image data must be transferred from the camera to the PC memory. This may take 1 to 150 ms (depending on the camera type, binning, resolution, ROI1 size and ROI1 position).

The high-speed control box uses this transfer time to move all actuators to the position required for the next image. In most acquisition situations the camera transfer time is sufficient for reaching the new position(s), in all other cases the next image acquisition is delayed accordingly.

* only available on systems equipped with a high-speed control box.

** optional accessories. Requires either Super Z Galvo or Piezo Focus

5. Parameters for Multidimensional Image Acquisition

The following components can be controlled by the high-speed control box (the typically required setting times are shown in parentheses):

- Leica external filter wheels
(switching time for adjacent filter positions 27 ms)
- SuperZ FineFocus
typical stepsize $1\text{ }\mu\text{m} = 40\text{ ms}$
- Piezo FineFocus
(typical stepsize $1\text{ }\mu\text{m} = 3\text{ ms}$)
- Leica AM TIRF MC
- Leica EL6000 high speed shutter (6 ms)

To enable high speed acquisition, only components that can be controlled by the high-speed control box may be changed in a sequence. Therefore, the following components must remain constant in sequence-controlled series:

- Objective
- Camera type/port
- Post-magnification (MagChanger)
- Contrast method
- Fluorescence cube
- Light intensity
(except if FFIM with an external filter wheel is used)
- Binning and camera gain
- Internal microscope shutter and filter wheels
- Stage position

If different values are defined within the experiment for one of these components, the system automatically switches to the **Normal** mode and thus allows acquisition of any desired sequences.

The green display field (**Light Path Settings** menu) indicates that the defined experiment can be completely controlled by the high speed board (Fig. 14).

Notes:

If the field has a red background, you can click the left mouse button on the red display field to query which settings must be changed to set up an experiment that can be completely controlled by the high-speed control box.

In complex experiments the system automatically divides it into partial series that can be controlled by the high-speed control box.

Example: If an experiment is defined from multiple z-series using different fluorescence cubes, the status field has a red background but the individual z-series are controlled by the high-speed control box, i.e. acquired time-optimized.

5.5 Mark and Find

Mark and Find is enabled with the  button and is used to define different stage positions for **multiposition experiments**.

When a stage position is approached, its **x**, **y** and **z-coordinates** can be saved with the **Mark** button:

If multiple positions are saved, they can be approached again using the dropdown menu. When the experiment is started, a time series of all positions is approached and acquired at every time point.


Fig. 15: Mark and Find

5. Parameters for Multidimensional Image Acquisition

5.6 Tile Scan

5.6.1 Camera adjustment for Tile Scan

Note:

A prerequisite for Tile Scan is a correctly calibrated camera!

Depending on the type of camera and the port used, the camera must be mounted in the correct orientation (see following table):

Camera type	Camera orientation on left, right and bottom port eyepiece image = camera image
Andor iXon 897	cable to the back
Hamamatsu 9100-13	black plastic plate in front
Hamamatsu 9100-02	black plastic plate in front
DFC365 FX	cable to the back
DFC450 (C)	cable to the back
DFC310 FX	cable to the back

5. Parameters for Multidimensional Image Acquisition

Aligning the camera

If you find that the edges of the live image do not run parallel to the edges of the image window, you can now fine-adjust the camera. (See chapter 7: "First experiment" for displaying a live image).

Proceed as follows:

- Switch to the objective with the highest magnification.
- Focus the microscope image.
- Select a prominent specimen detail and move it to the top right corner of the image window.
- Move the stage in y direction so that the specimen detail in the image window moves to the bottom.
- Check that the distance to the edge of the image window remains the same.
If this is not the case, the camera must be rotated accordingly.

5.6.2 Tile Scan Acquisition

The **Tile Scan** menu is enabled with the  button and is used to define an area for generation of composite overview images.

The Tile Scan area is specified by two oppositely placed positions which uniquely define the Tile Scan area. This will typically be the top right corner and the bottom left corner of the area that is to be highlighted.

As in the **Mark and Find** menu, the positions are defined with the **Mark** button:

If the Tile Scan area is set, the **Scan Field** shows the number of images that will be required to cover the previously defined scan area. The **Merge Image** button should always be enabled to merge the acquired images to a Tile Scan image:


Fig. 16: Tile Scan

5. Parameters for Multidimensional Image Acquisition

If the single images are not merged correctly, it may be necessary to swap the images horizontally or vertically.

To do this click the **Configurations** button.

This opens the **Stage Configuration** menu (fig. 17).

In **Set merge settings** to select the option **Custom** and mark the appropriate option **Swap images horizontal** or **Swap images vertical**.


Fig. 17: Stage configuration

6. FIRST EXPERIMENT

The "First Experiment" chapter gives an overview of the power of the Leica LAS X and shows you how fast you can acquire the first image. In the following example the specimen "Convallaria" is used. It is stained with Safranin and Fast Green. We will go through the simplest settings and run a fast xyzλ scan. In the first experiment an image stack is acquired at one stage position for both fluorochromes.

The FineFocus (Piezo or SuperZ)* is used to set the z-Positions. The general procedure for such an experiment is shown below. For a detailed description of the dialog boxes and the software functions see the LAS X online help.

Note:

Note that all settings are reference values. Depending on the specimen the values and settings may vary.

6.1 Configuration

Use a 40x or 63x objective for the experiment. The LiveCell HCX PL APO 63x/1.30 GLYC CORR 37°C objective is the optimum solution. Use the BGR Cube 11600197 filter system B/G/R; HP; size "K" as the filter cube.

6.2 Switch on the system and start the software

- Switch on the system.
Follow the specified power-up sequence (see Chapter 5 Starting up).
- Start the LAS X software by clicking the icon on the desktop or by selecting it in the start menu.

- At the following prompt select **Yes** to initialize the xy-stage:


Fig. 18: LAS X stage initialization

* optional accessories. Requires high-speed control box

6. First Experiment

- When you see the following dialog box the software has been started and is ready for use:


Fig. 19: LAS X user interface

6.3 Selecting and focusing specimen section

- Select TL-BF (bright field) as the contrast method.


Fig. 20: Contrast method

- Click the **Setup** tab.
- If available the objective **HCX PL APO 63x/1.30 GLYC CORR 37°C** should be selected, otherwise select a 40x or 63x objective:

Fig. 21: Objective selection

- Place the specimen on the xy-stage (use an immersion medium if necessary).
- Switch the light to the eyepieces and adjust the image directly at the microscope.

Note:

To switch the light to the camera, so that an image is directly shown on the monitor, the **Live** button must be clicked.

6. First Experiment

6.4 Camera resolution

- Click on **Configuration** in the **Light Path Settings** dialog box:


Fig. 22: Camera configuration

- Under **Set Digitization to** (23.1) select the value **12 bits** (when using the cascade camera select 16 bits). Close the window by clicking on the **X**. The settings are automatically applied.


Fig. 23: Camera configuration

6.5 Acquisition mode

- Click the **Acquire** arrow symbol:

Fig. 24: Acquire arrow symbol

- Click the **Acquisition** tab:


Fig. 25: Acquisition tab

- In **Acquisition Mode** select in addition to **x**, **y** the **z** button. The **x** and **y** acquisition parameters are always enabled.

Selected and therefore active fields are underlined in red in this example. If you click the same button again, it is given a black underline. Black icons are not selected and therefore are not enabled (in this example the **t** button is not enabled).


Fig. 26: Acquisition mode xyz

6.5.1 EL6000 variant

When using the EL6000 the wavelength is set with a filter cube, the filter wheel integrated into the microscope (IFW) or with external filter wheels (EFW).

- In the **Light Path Settings** menu in **Contrast Method** select **FLUO** from the list field (27.1).

Filter selection for EFW:

- As **Filter Cube** select **TRI**.
- Selecting the filter cube in the **Light Path Settings** menu (27.4) activates the corresponding filter wheel.
Attention: The suitable filter wheel must be positioned in the light path before!
- Now select the **G** excitation filter in the **Excitation Filter** field (29).

6. First Experiment


Fig. 27: Light Path Settings menu

Filter selection for IFW:

- Press the **IFW** button to be able to define the filter wheel position:


Fig. 28: Image dialog

- As **Filter Cube** select **BGR** and select intensity position **G**:

Fig. 29: Image dialog

If the system is equipped with the external **FFIM**, the transmission is always controlled via the intensity slider in the **Image menu** (in five steps). The internal filter wheel of the microscope is set to 100% transmission in this case.


Fig. 30: Adjusting transmission for FFIM

Steps for FFIM:

100%, 50%, 25%, 12.5%, 6.25%

Attention!

Set the light intensity of the EL6000 to the lowest level to protect the specimen.

Set the values in the **Image dialog**:

Binning:

Exposure:

Gain:

Fig. 31: Image dialog

- Click **Live** to be able to view the image directly on the monitor:


Fig. 32: Live or Stop button

6. First Experiment

- Now you see the live image on your screen in the **Viewer** window (see fig. 33).

Fig. 33: Acquisition view in the viewer

Note:

Move the specimen to the desired position with SmartMove*. For subsequent image processing it is advantageous if you use relatively small cells with many starch granules (bright dots inside the cells). If you cannot find any starch granules, focus a different plane of the specimen and try again.

6.6 Image optimization

- Optimize the image quality by adjusting the exposure time by moving the **Exposure** slider in the **Image** menu slowly left and right. Experience shows that a new Convallaria specimen has exposure values of 100 ms to 200 ms. Observe the changes in the image quality:


Fig. 34: Image menu

Attention!

Reduce the exposure time to protect your specimen.

- When you have reached the desired image quality click **Stop**:

Fig. 35: Live or Stop button

* optional accessories

6.7 Setting color (LUT)

Note:

The **Convallaria** specimen is stained with two stains (**Safranin/Fast Green**). For this reason two fluorochrome settings (FCr1 and FCr2) are explained in this example.

- In this example in **FCr1** the stain **Green** is selected. You can click the button for the color palette (36.1) to select the desired color (36.2). Confirm your selection with double click on left mouse button.


Fig. 36: Color palette

- Click in the text field and assign a name for the **FCr1** (in this example **GFP**):

Fig. 37: Text field

- To save all settings for the **FCr1** click the **right mouse button** on the **FCr1** field (38.1) and then on **Save** (38.2):

Fig. 38: Saving settings

6. First Experiment

- You are prompted to assign a name again. Use the same name as in fig. 37 and confirm it with **OK**:


Fig. 39: Confirming name

6.8 Settings for additional fluorochromes

- To specify settings for an additional fluorochrome, click the **+** button. This creates the FCr2:


Fig. 40: Adding fluorochromes

Proceed with FCr2 as with FCr1:

- **FLUO** remains selected as the **Contrast Method**.
- If EFW is used, the **Filter Cube TRI** remains selected.
- If IFW is used, **BGR** remains selected as the **Filter Cube**.
- In the **Image** menu select intensity setting **R** (for IFW) or select the filterwheel position **R** (EFW) in the **Light Path Settings** menu.
- Click **Live** to be able to view the image directly on the monitor.
- If you wish to optimize the image quality, proceed as described in the section on 6.6 Image Optimization.
- Select the color for display (red in this example, 41.1) and assign a name for **FCr2** (41.2):


Fig 41: Setting color and name

- Save the settings by clicking the **right mouse button** on **FCr2** (42.1) and then clicking **Save** (42.2):


Fig. 42: Saving settings

- Enter a name and confirm with **OK**.
- To start image acquisition click **Capture Image**.

This defines FCr1 and FCr2 and saves the settings for the two fluorochromes. Now you can start the image acquisition.

6.9 Capturing 3D image stack

In the list field of the z-menu you can select either **Wide Focus** (microscope z-drive) or, if available, **Fine Focus** (Super z-Galvo Focus or Piezo Focus)*:


Fig. 43: z-Stack

Example for Fine Focus:

- In the list field select **Fine Focus** (43.6).
- Click **Live** to be able to view the image directly on the monitor.
- Make sure that the image is in focus.
- Click the **Focus** button (43.1) to save the current **z-Position**.
- To define the upper position of your z-stack move the mouse pointer to the ocher-colored area (43.2) and drag it upwards while holding down the left mouse button. The change of the z-position can be followed in the live image.

* optional accessories. Requires high-speed control box

6. First Experiment

- As soon as the desired top image plane appears in the live camera image highlight it by clicking the **Begin** button (43.5).
- Now drag the ocher-colored area (43.2) with the left mouse button held down to the desired bottom position and highlight it by clicking the **End** button (43.3).

Example for z-Wide:

- In the list field select **z-Wide**.
- Click **Live** to be able to view the image directly on the monitor.
- Make sure that the image is in focus.
- Click the **Focus** button (43.1) to save the current **z-Position**.
- To define the upper position of your z-stack move the microscope z-drive to the desired upper limit. The change of the z-position can be followed in the live image.
- As soon as the desired top image plane appears in the live camera image highlight it by clicking the **Begin** button (43.5).
- Now move the microscope z-drive to the desired lower limit and highlight it by clicking the **End** button (43.3).

Note:

When using the 63x objective we recommend 150 to 200 image planes (sections) for definition of the image stack. We recommend 100 to 150 image planes when using a 40x objective.

Fig. 44: z-Stack

- Select the **System Optimized** option (44.1).
- Click the **Start** button to start image acquisition.

Note:

The saved image series can be found under **Acquire** → **Experiments** (fig. 45). The image data created by image processing is also saved here.


Fig. 45: Experiments

Note:

The acquisition is started by clicking **Start** and ends automatically (if that is set). The acquisition can be stopped at any time by clicking **Stop**.

Note:

Save the captured images directly after acquisition before starting a new acquisition, processing the image data or before shutting down the LAS X. In addition to saving the data on the hard disk, the data should also be backed up regularly on a suitable data medium.

Note:

This means that the image data will not be lost if the software crashes.

7. INDEX

C

Copyrights 2

F

First experiment 27

Acquisition mode 31

Camera resolution 30

Capturing 3D image stack 37

Configuration 27

Image optimization 34

Selecting and focusing specimen section 28

Setting color (LUT) 35

Settings for additional fluorochromes 36

Starting the software 27

Switching on the system 27

I

Intended use 6

L

LAS X 8

Calling up the online help 8

Full-text search 9

Structure of the user interface 10

Leica LAS X 7

M


Multidimensional image acquisition 12

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Copyright © Leica Microsystems CMS GmbH • Ernst-Leitz-Straße • 35578 Wetzlar • Germany 2016 • Tel. (06441)29-0 • Fax (06441)29-2599
LEICA and the Leica logos are registered trademarks of Leica IR GmbH.