From Eye to Insight

LIBS

User Manual

Leica Microsystems CMS GmbH • Manual 11934161, Revision 1.1, 2018-01-05

From Eye to Insight

LIBS

User Manual

Leica Microsystems CMS GmbH • Manual 11934161, Revision 1.0, 2017-10-17

Copyright

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information included in this manual may be changed without prior notice.

Revision 1.1, published on 2018-01-05 by:



Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17–37
D-35578 Wetzlar (Germany)
<http://www.leica-microsystems.com>
Responsible for contents:
Marketing CMS

Contents

1. Important Notes about this Manual	5
1.1 Text symbols and their meanings.....	5
2. Function	7
3. Safety Notes	9
3.1 General safety notes	9
3.2 Electrical safety	9
3.3 Laser safety	10
3.4 Safety features.....	10
3.4.1 Safety notes for handling the laser	12
3.4.2 Safety labels	12
3.5 Transport and Storage	13
3.6 Ambient conditions	13
3.7 Disposal.....	13
4. Installation Location.....	14
5. Unpacking	14
6. Description and Assembly of the Components.....	15
6.1 Stand.....	15
6.2 Objective 20x/0.4 LIBS NUV	15
6.2.1 Changing objectives.....	16
6.3 LIBS module	16
6.4 LIBS filter cubes	17
6.5 Cables.....	18
6.6 Software Components	20
7. Operation	21
7.1 Switching on.....	21
7.2 Starting the software	21
7.3 Prepare measurement.....	22
7.4 Start measurement.....	24
7.5 Measurement result.....	25
8. Care.....	28
8.1 Cleaning	28
8.2 Handling acids and bases	28
9. Essential Consumable and Spare Parts.....	29
10. EC Declaration of Conformity.....	29

The LIBS (Laser-Induced Breakdown Spectroscopy) module in combination with the light microscope is designed for visual inspection and chemical analysis in one step. The following combinations are possible.

Matrix of combination options

1. Important Notes about this Manual

Caution!

This User Manual is an essential component of the product; it must be read carefully before the product is assembled, put into operation or used, and must be kept for future reference.

This User Manual includes important instructions and information related to the operating safety and maintenance of the instrument.

In addition to these instructions, depending on the equipment configuration, other instructions are included for the individual components and must be observed:

- User manual for microscope
- User manual for software
- User manual for accessories

1.1 Text symbols and their meanings

(1.2)

Numbers in parentheses refer to illustrations; for example, (1.2) refers to Fig. 1, item 2.

→ page 20

Numbers with a pointer arrow (for example → page 20), point to a certain page of this manual.

WARNING indicates a potentially hazardous situation with a medium degree of risk that, if not avoided, could result in death or serious injury.

CAUTION indicates a potentially hazardous situation with a low degree of risk that, if not avoided, may result in minor or moderate injury.

Caution!

In this manual, additional safety notes are indicated with the triangle symbol shown here, and have a gray background.

Caution! The instrument and accessories can be damaged when operated incorrectly.

1. IMPORTANT NOTES ABOUT THIS MANUAL

Warning of hazardous electrical voltage! Danger of electric shock!

Warning of optical radiation! Never look directly into the light beam! Wear safety goggles!

Warning of electromagnetic field

Warning of UV radiation! Wear safety goggles!

Warning of hazardous laser radiation – take appropriate safety precautions!

Warning of hand injuries

Instructions for disposing of the instrument, its accessories and consumables.

Connection to ground

Explanatory note

*

Item not contained in all equipment configurations

Manufacturing date, e.g. 11 / 2011 for November 2011

China RoHS 50 years EFUP
(Environmentally friendly use period)

2. Function

The Leica DM6/6000 M system with the LIBS module, for which this operating manual has been written, is intended for industrial routine and research applications.

The Leica DM6/6000 M system with the LIBS module is used for inspection and analysis during materials science examinations. It can be used universally. All contrast methods, such as bright field, dark field, polarization, DIC or fluorescence, are integral components of the microscope and can be adapted or exchanged quickly and easily.

Directives of the European Community (EU Directives)

The microscope described here fulfills the requirements of EU Directive 2014/35/EU concerning the safety of electrical equipment and 2014/30/EU concerning electromagnetic compatibility and has been declared accordingly.

The complete system does not conform to the requirements of an IVD device according to 98/79/EC.

According to the VO1194/2012 regulation regarding the implementation of the Ecodesign Directive 2009/125/EC, microscopes are considered special products. The light sources in this instrument are special lamps which are only suitable for microscopes. They are not suitable for general lighting purposes or room lighting.

Reasonably foreseeable misuse

The following are prohibited:

- To use the system for any purpose not in accordance with the Declaration of Conformity (e.g. the use of medical products in accordance with *EU Directive 93/42/EEC*)
- Making use of auxiliary equipment to use the area between the microscope stage and objectives as a clamp or holder (working like a vice)
- Operating the system in an inclined position
- Cleaning the system in a way other than specified in the manual
- Modifying the built-in safety circuits in the instrument
- Allowing unauthorized personnel to open the instrument
- Using cables that Leica Microsystems CMS GmbH has not provided or permitted
- Using the touch screen of the microscope as a transport handle
- Using combinations with non-Leica components that go beyond the scope of the manual
- Using power supplies not approved by Leica.

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from, using the microscope for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH. In such cases, the Declaration of Conformity shall be invalid.

Notes on handling laser devices

In the standard design without additional laser safety precautions, the microscopes are not suitable for coupling laser radiation (such as to the camera ports), since this radiation poses a hazard for the user (eye injury in particular).

For use of the microscope with lasers, Leica Microsystems CMS GmbH offers special microscope variants with additional safety features. Laser couplings require corresponding safety devices that have to be inspected and installed by trained personnel. For further information, please contact your authorized Leica Microsystems CMS GmbH representative.

Due to the special sources of danger of the laser system, only those persons may use the system who have received special training in working with the system, observing the specific hazards of a laser system.

3. Safety Notes

3.1 General safety notes

This safety class 1 device is constructed and tested in accordance with EN 61010-1: 2010, Safety requirements for electrical equipment for measurement, control and laboratory use.

It also meets the requirements for EN 61326-1, Electrical equipment for measurement, control and laboratory use – EMC requirements.

The product has also been tested in accordance with EN 60825-1 / IEC 60825-1 Safety of laser products
Part 1: Equipment classification and requirements.

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this user manual.

The instruments and accessories described in this user manual have been tested for safety and potential hazards.

The responsible Leica Microsystems CMS GmbH affiliate or the main plant in Wetzlar must be consulted whenever the instrument is altered, modified or used in conjunction with non-Leica components that go beyond the scope of this manual.

Unauthorized alterations to the instrument or noncompliant use shall void all rights to any warranty claims!

In addition to this manual, the safety notes of the other manuals provided must be observed!

The system may be set up, installed and put into operation only by authorized Leica personnel!

Caution!

Only personnel who have been instructed by authorized Leica personnel may operate the microscope!

Caution!

The microscope's electrical accessory components are not protected against water. Water can cause electric shock.

! Caution!

Do not plug in or unplug data lines and control circuits unless the instrument is switched off; otherwise, the instrument may be damaged.

3.2 Electrical safety

General specifications of the LIBS module

For indoor use only.

Supply voltage:	24 V DC
Power consumption:	max. 50 W
Fuses:	5x20, 3.15 A, slow-blowing, Breaking capacity H, 250 V DC
Ambient temperature:	15° to 30°C
Relative humidity:	90% up to 30°C, non-condensing

Use only a permitted power supply unit. Other power supplies must not be used.

If the original power supply fails or is damaged, it must be replaced. Repair is not permitted.

Original power supplies are available from your Leica branch office or Leica dealer.


The specifications of the microscope and CTR box as well as other components can be found in the corresponding user manuals.

Specifications of the power supply unit for the LIBS module

Power supply type:	SINPRO SPU80-108
Supply voltage:	100 – 240 VAC
Frequency:	47 – 63 Hz
Power:	max. 80 W
Output voltage:	24V DC
Output current:	3.8 A max.
Ambient temperature:	0° to 40°C
Relative humidity:	10% to 95% R.H., non-condensing

3.3 Laser safety

The LIBS module is a component of the Leica DM6/6000 M system with LIBS module. The LIBS module can be coupled with a microscope with the configuration options described on page 4. In all configuration options shown on page 4, the system fulfills Laser Class 1 in accordance with the requirements of IEC60825-1.

The LIBS module includes a Class 3B laser. Due to the laser safety measures, only radiation of laser class 1 is found in the specimen area of the microscope.

As long as the LIBS module and the laser protection hood are correctly installed on the incident light axis, you do not have to take any special safety precautions.

Specifications of the LTB MNL100 nitrogen laser used in the LIBS module:

Wavelength:	337 nm
Spectral bandwidth:	0.1 nm
Pulse half-width – FWHM:	typically 3 ns
Guaranteed pulse energy (90% after 60 million):	≥ 80 µJ
Typical pulse energies (typically 70% after 100 million):	≥ 90 µJ
Pulse power:	typically 30 kW
Laser Class:	3B

Power supply:	24 V DC
Current draw:	max. 2.4 A
Power consumption:	max. 60 W
Operating temperature range:	+15 to +38 °C
Storage temperature range:	-10 to +60 °C
Max. rel. humidity (non-condensing)	85% R.H.

3.4 Safety features

The laser may be operated only when the complete system is fully installed.

The laser beam is completely shielded until it exits the objective.

Accordingly, only the objectives are potential laser exit points.

Only persons authorized by Leica CMS are permitted to open the LIBS module.

3. SAFETY NOTES

To ensure laser safety, the laser protection hood (Fig. 1 and Fig. 2) must be installed on the incident light axis. Make sure that the laser protection hood has clicked into place.


Fig. 1 Front view of the microscope
1 Laser protection hood


Fig. 2 Laser protection hood
1 Safety switch

Safety switch (interlock)

The laser protection hood is connected to a safety switch (Fig. 2). Removing the laser protection hood automatically shuts off the laser radiation.


Fig. 3 Left side of the stand
1 Power supply connection
2 Fuse
3 USB port
4 Interlock connector for laser protection hood

Protective filter

Additionally, a protective filter is permanently installed in the stand cover (Fig. 4).

Fig. 4 Stand cover
1 Protective Filter

Note:

The laser is activated exclusively through the LAS X software.

3.4.1 Safety notes for handling the laser

WARNING

The user is not permitted to open the instrument!
The laser system is encapsulated. It may be repaired or replaced by personnel authorized by Leica CMS only. Incorrect handling of the laser can cause injuries!

WARNING

If other operation or adjustment devices than those listed here are used or other procedures are carried out, dangerous beam exposure may result.
Risk of injury!

WARNING

To be classified as a Class 1 laser product according to IEC/EN60825-1 and fulfill the electrical safety regulations, all safety devices, safety interlocks and safety systems of the laser product must be in operational condition.

Deactivating or damaging these safety devices or any intervention in any of these safety devices may lead to serious eye injuries, physical injuries or property damage. In these cases, Leica Microsystems CMS shall not assume any liability.

According to IEC/EN60825-1: Safety of laser products, Part 1: Equipment classification, requirements and user's guide," and applicable country-specific regulations, the owner/operator is required to designate a Laser Safety Officer or a Laser Protection Advisor.

WARNING

Never deactivate the laser protection devices.

WARNING

Never look directly into a laser beam or a reflection of the laser beam. Avoid all contact with the laser beam.

WARNING

Do not introduce any reflective objects into the laser beam path.

WARNING

All unoccupied positions in the objective nosepiece must be closed using the supplied caps.
The system must be switched off when the turret is being replaced or fit with an objective.

3.4.2 Safety labels

The safety labels are attached at various places (Fig. 5).

Fig. 5 Rear side of the LIBS module

3.5 Transport and Storage

Transport

Caution!

Transport and storage at -20 °C – +85 °C and max. 90% humidity.

For shipping or transporting the system and its accessory components, the original packaging should be used.

Weight

At a full load, the system weighs more than 18 kg. The user must take corresponding precautions for transporting the instrument.

3.6 Ambient conditions

Caution!

This microscope may not be used at elevations over 2000 m above sea level.

Caution!

Protect the microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.
Operating temperature: 15–30 °C.

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend operating the instrument in ESD-protected areas as well assessing the electromagnetic environment and giving corresponding instructions before operating this instrument.

This instrument has been developed and tested in accordance with CISPR 11 Class A. It can cause radio interference in a household environment. In this case, you may have to take measures to eliminate the interference.

3.7 Disposal

Once the product has reached the end of its service life, please contact Leica Microsystems GmbH Service or Sales about disposal.

Please observe and comply with the national and federal laws and regulations that are equivalent to EC directives such as WEEE.

Note!

Like all electronic devices, the microscope, its accessory components and consumables must never be disposed of with general household waste.

4. Installation Location

Work with the microscope should be performed in a dust-free room, which is free of oil vapors, other chemical vapors and extreme humidity. At the workstation, large temperature fluctuations, direct sunlight and vibrations should be avoided. These conditions can distort measurements and longer microscopic scans.

Permitted ambient conditions:

Temperature	15–27 °C
Relative humidity	maximum 90% up to 30°C non-condensing

Microscopes in warm and warm-damp climatic zones require special care in order to prevent fungus contamination. Also make sure to follow the installation instructions in the microscope's manual!

Caution!

The microscope must be set up so that the power voltage input is not obstructed and the instrument can quickly be removed from the power supply if necessary.

Caution!

Electrical components must be placed at least 10 cm away from the wall and from flammable substances.

5. Unpacking

Please compare the contents of the package you receive carefully with the packing slip, the delivery note or the invoice. We highly recommend keeping a copy of these documents with the manual in order to have delivery date and delivery package information available when (for example) placing repeat orders or making service calls. Please make sure that no small parts are lost in the packaging materials. Some of our packaging materials have symbols to help in environmentally responsible recycling.

Note:

Keep the packaging materials for storage purposes and transportation of the microscope and its accessory components.

Caution!

If at all possible, avoid touching the optics. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces of optical surfaces in a short time.

Caution!

Do not connect the laser and peripherals to the power socket before assembly is complete!

6. Description and Assembly of the Components

Caution!

The system may be set up, installed and put into operation only by authorized Leica personnel!

Do not connect the electronics box, laser or peripherals to the power socket before assembly is complete!

6.1 Stand

LIB spectroscopy is possible in combination with the Leica DM6/6000 M microscope.

- Also observe the User Manual for operating the microscope.

Fig. 6 Configuration example with Leica DM6 M

6.2 Objective 20x/0.4 LIBS NUV

The special 20x/0.4 LIBS NUV objective (Fig. 7) was developed for LIB spectroscopy.

The receptacles on the objective turret are numbered. Depending on your equipment, a specific position has already been assigned to the objective at the factory. For details on the exact positions of the objectives, please refer to the enclosed identification sheet.

CAUTION

Close vacant threads in the turret with dust protection caps!

Note:

We recommend carrying out a parfocality adjustment via the Leica Microsystems LAS X software.

Fig. 7 Objective 20x/0.4 LIBS NUV

Caution!

If the objective receptacle provided for the LIBS objective is fitted with a different objective, proper function of the device cannot be guaranteed. Use of an incorrect objective can result in damage to the objective and the LIBS module.

**Note:**

To protect the objective when the laser is in use, the objective protection glass set (Order No. 11518150) is available.

- Insert the protective glass into the holder.
- Pull the holder over the objective and screw the holder into place on the objective.

6.2.1 Changing objectives

The laser protection hood has to be removed to change the objective. When this happens, the laser is automatically shut off by the laser protection equipment.

- Replace the objective.
- Reinstall the laser protection hood.

6.3 LIBS module

The LIBS module (Order No. 11 560 200) consists of the laser module as well as the laser protection hood and a stand cover with a built-in protective filter (Fig. 8).

**Caution!**

Only persons authorized by Leica Microsystems are permitted to install the LIBS module.

- When the LIBS module is upgraded, the existing stand cover is replaced with the stand cover with a built-in protective filter.
- The laser module is installed on the incident light axis.
- Push the laser protection hood into the guide of the incident light axis.
- Ensure that the laser protection hood clicks into place.


Fig. 8 Front view of the microscope

- 1 Laser protection hood
- 2 Stand cover with protective filter
- 3 Laser module

**Note:**

The magnification changer and conoscopy module cannot be combined with the LIBS module.

The user is not permitted to open the instrument!
The laser system is encapsulated. It may be repaired or replaced by personnel authorized by Leica CMS only. Incorrect handling of the laser can cause injuries!

If other operation or adjustment devices than those listed here are used or other procedures are carried out, dangerous beam exposure may result.
Risk of injury!

Never deactivate the laser protection devices.

Never look directly into a laser beam or a reflection of the laser beam. Avoid all contact with the laser beam.

6.4 LIBS filter cubes

The positions on the turret disk are numbered. Depending on your equipment, specific positions have already been assigned to the individual filter cubes at the factory. For details, check the identification sheet included with your order.

Caution!

Use only the provided filter cubes. In combination with LIBS, only the cubes with order number 11561105 (Fig. 9) are used.

To insert the filter cubes, follow the instructions in the microscope's manual.

Fig. 9 LIBS filter cubes

Caution!

If the filter holder of the turret disk provided for the LIBS filter cube is fitted with a different filter cube, proper function of the device cannot be guaranteed. Use of an incorrect filter cube can result in damage to the filter cube and the LIBS module.

6.5 Cables

After completing the assembly work, connect the LIBS module to the power supply using the provided power supply unit.

- Connect the power supply unit to the port (Fig. 11.1) on the rear side of the LIBS module.
- Connect the power supply unit to the power supply.

WARNING

Use only a permitted power supply unit. Other power supplies must not be used.
If the original power supply fails or is damaged, it must be replaced. Repair is not permitted.
Original power supplies are available from your Leica branch office or Leica dealer.

Caution!

By definition, the main circuit breaker of this instrument is the connection between the power cable and device port. The user must ensure unobstructed access to the main circuit breaker at all times.


Fig. 10 Connection of the LIBS module to the computer

- 1 Power supply connection
- 2 Fuse
- 3 USB port
- 4 Interlock connector for laser protection hood


Fig. 11 Connection of the LIBS module to the power supply
1 Power supply connection

Specifications of the power supply unit for the LIBS module

Power supply type:	SINPRO SPU80-108
Supply voltage:	100 – 240 VAC
Frequency:	47 – 63 Hz
Power:	max. 80 W
Output voltage:	24V DC
Output current:	3.8 A max.
Ambient temperature:	0° to 40°C
Relative humidity:	10% to 95% R.H., non-condensing

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation. We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

6. DESCRIPTION AND ASSEMBLY OF THE COMPONENTS

Caution!

This instrument has been developed and tested in accordance with CISPR 11 Class A.

It can cause radio interference in a household environment. In this case, you may have to take measures to eliminate the interference.

Fig. 12 Power supply unit of the LIBS module

6.6 Software Components

The system is shipped fully configured as standard. All hardware components, such as the laser, are preconfigured, and the Leica LAS X software required for LIB spectroscopy is already installed.


The LIBS software is license-free and is activated when the LIBS module is activated in the **Hardware Setup** module of the **LAS X Hardware Configurator**.

The installation of the software is described in the LAS X Operating Manual, which is included on the installation CD.

The microscope is also configured using the LAS X software.

- Make sure the microscope and all components are switched on.
- Start the **LAS X Hardware Configurator** and select the **Hardware Setup** module.
- Select **Leica DM6/6000** as the connected microscope, **USB** as the connection type and **Leica** for the stage.
- Select the **LIBS** application and click **Apply**.

Fig. 13 LAS X Hardware Configurator – Hardware Setup module

7. Operation

7.1 Switching on

The following switch-on sequence must be adhered to:

- Switch on the computer.
- Switch on the microscope using the CTR box.
- Connect the LIBS module to the power supply.

7.2 Starting the software


- Start the LAS X program by clicking the corresponding icon on the desktop or by selecting it in the start menu.

The Leica Microsystems LAS X program starts with the **Acquire** step.

If the LIBS application was activated in the Hardware Setup module, the LIBS button appears in the bar at the left edge of the window.

Fig. 14 LIBS application activated

Note:

Before you start the measurement, be sure to check the parfocality and parcentricity and, if necessary, perform a calibration.

7.3 Prepare measurement

- Adjust the sample using the controls of the microscope.
- Open the **Laser Spectroscopy** control panel field by clicking the button



Fig. 15 LAS X working window - After clicking the LIBS button and "Toggle reticule"

- Click **Toggle reticule**. A crosshair with a circle appears in the center of the image window (Fig. 15). The circle indicates the place where the laser strikes.
- Position the sampling site to be analyzed under the circle.
- The preinstalled database is specified automatically.
- Click **Prepare Measurement**. The LIBS objective is swung in. Since the LIBS objective is corrected specifically for transmission of the laser wavelength 337.2 nm, the optical contrast when the LIBS objective is swung in is less than for the "normal" objectives used for the visual inspection. (Fig. 16).

The **Prepare Measurement** button becomes **Start Measurement**.

- Check the image setting again.


Note:

A new database can be created using the external "SpectralTools" software. Then the active database can be selected from the database list. The data can be applied after


the  button is clicked.

Fig. 16 LAS X working window - After clicking "Prepare measurement"


7.4 Start measurement

- Click **Start Measurement**.

- The LIBS filter cube is swung in.
- The laser is activated. During the first measurement of a work step, the laser starts with a short delay (warm-up phase of the laser).
- After the shot, the filter cube is swung out again and the spectrum is displayed (Fig. 17).
- The **Start Measurement** button remains active. Another analysis can be repeated immediately after the results window is closed.


Fig. 17 LAS X working window - After the measurement

7.5 Measurement result


- The result is depicted graphically. The five highest matches of the measured signal with the spectra stored in the database are specified.
- The **Quality of Match** value indicates the agreement of a measured spectrum with the spectra in the database. The maximum value of 1000 corresponds to a 100% agreement, which can be achieved only in theory due to the normal measurement noise. Values less than 650 are generally not to be used for unambiguous identification.
- The different colors of the displayed spectra mean:

Red	=	Raw spectrum
Blue	=	Corrected spectrum
Green	=	Database spectrum with the highest agreement
- To export or save the result, enter a name for **Analysis name**.
Otherwise the following error message appears:

- **Signal to noise ratio** is a measure for the technical quality of the measured signal that is superimposed by a noise signal (substrate).


Fig. 18 Results window


7. OPERATION

- **Export raw data**

To export the raw data, enter the name of the target folder:

The created and selected folder contains five files (in this example, for the analysis name "Test"):

Analysis name_data_classifications.csv = Agreement results table

Analysis name_data_spectrum.csv = Data points of the measured spectrum (wavelength and amplitude)

Analysis name_info.txt = Analysis name, date, S/N

Analysis name_snapshot.jpg = Snapshot after the analysis

Analysis name_spectrum.jpg = Snapshot of the spectrum

- **Save for Report**

To enable the results to be displayed later in a Word report, click **Save for Report**. If an acquisition folder has been defined in LAS X, the analyses are automatically saved in this. If an acquisition folder has not yet been defined, the created analysis files are stored in the "unsaved data" folder.

- **Generating a report**


Right-click the analysis file. A display window opens with the menu item **Export** → **AS LIBS Report** (Fig. 19).

To create a report with multiple analyses, select these accordingly.

- **Save to database**

If the spectrum of a measured substance is not yet in the database as a reference spectrum, the obtained spectrum can be taught into the database to serve as a reference spectrum for new measurements. To do so, click **Save to database**. The spectrum is entered into the database being used.

Fig. 19 Generating a report

8. Care

Caution!

Unplug the power supply before performing cleaning and maintenance work!

Protect electrical components from moisture!

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus. The microscope should be cleaned after each use, and the microscope optics should be kept painstakingly clean.

8.1 Cleaning

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with a low-concentrated soap solution, petroleum ether or ethyl alcohol. For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

! Caution:

Thinners containing acetone, xylene or nitrogen can harm the microscope and thus must not be used.

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure that coated or plastic surfaces do not become matted or etched.

8.2 Handling acids and bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

! Caution:

Never allow the optics and mechanical parts to come into direct contact with these chemicals.

Removing immersion oil

Caution!

Follow safety notes for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

9. Essential Consumable and Spare Parts

Order No. Part number	Designation	Use for
<u>Screw cap for unused objective receptacles</u>		
020-422-570-000	Screw cap M 25	Objective turret
11 020-422-580-028	Screw cap M 32	Objective turret
<u>Dust and light stop for analyzer opening</u>		
11 020-530-100-013	Analyzer opening cover	Microscope Stand
11 020-437-101-013	Polarizer opening cover	Microscope Stand
<u>Dust and light stop for camera port openings</u>		
11 020-387-556-009	Port opening cover	Microscope Stand
<u>Replacement eyecup (antiglare protection) for HC PLAN eyepiece</u>		
021-500-017-005	HC PLAN eyecup	10x/25 eyepiece
021-264-520-018	HC PLAN eyecup	10x/22 eyepiece
021-264-520-018	HC PLAN eyecup	10x/20 eyepiece
<u>LIBS objective</u>		
11 518 149	20x/0.4 LIBS NUV	Objective turret
<u>Objective protection glass set</u>		
11 518 150	10 protective glasses with holder	LIBS objective
<u>LIBS filter cubes</u>		
11 561 105		Incident light turret disk

Additional spare parts such as cover caps, stickers, screws, dust covers etc., on request.

10. EC Declaration of Conformity

To download the EC Declaration of Conformity, use this link

<http://www.leica-microsystems.com/products/light-microscopes/industrial-materials/upright-microscopes>

Select the microscope type and go to the "Download" page.

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