Leica LMD6 Leica LMD7

Laser Microdissection Systems
Lasermikrodissektionssysteme
Systèmes de microdissection laser

Instructions for use
Gebrauchsanweisung
Mode d'emploi

Leica Microsystems CMS GmbH. Instructions for use, 11934117, Revision 1.1, 2016-01-15

Leica
MICROSYSTEMS

Revision 1.1, published 2016-01-15 by:

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Leica Microsystems CMS GmbH

Ernst-Leitz-Straße 17-37

35578 Wetzlar (Germany)

<http://www.leica-microsystems.com>

Leica LMD6 Leica LMD7

Laser Microdissection Systems

Instructions for use

Leica Microsystems CMS GmbH. Instructions for use, 11934117, Revision 1.1, 2016-01-15

Leica
MICROSYSTEMS

Copyrights

All rights to this documentation are owned by Leica Microsystems CMS GmbH. Copying of text and illustrations – in full or in part – by printing, photostat, microfilm or other techniques, including electronic systems, is only permitted subject to the express written consent of Leica Microsystems CMS GmbH.

The information contained in the following documentation represents the latest stage of technology. We have composed the texts and illustrations with great care. Nevertheless, we are always grateful to be notified of any errors.

The information in this manual may be altered without prior notice.

Contents

1. Important notes on this manual	6	7. Getting started.....	33
2. Intended Use.....	10	7.1 Installing microdissection software	33
3. Safety notes	12	7.2 Configuring the Microscope	34
3.1 General safety notes	12	7.3 Switching on the Leica LMD6	43
3.2 Electrical safety.....	12	7.4 Switching on the Leica LMD7	44
3.3 Laser safety.....	15	8. Care and maintenance.....	45
3.4 Protection filter for possible laser emission points.....	17	9. Main spare parts and consumables.....	46
3.5 Technical Data of the built-in laser system	17	10. Index	48
3.6 Transport and storage	18	11. EU Declaration of Conformity.....	49
3.7 Disposal	18		
4. Installation site	19		
5. Unpacking.....	19		
6. Description and assembly of the components	20		
6.1 Stage	21		
6.2 Condenser	23		
6.3 Tube and Eyepieces.....	24		
6.4 Camera.....	24		
6.5 Objectives.....	25		
6.6 Lamphousing.....	26		
6.7 Leica EL6000 compact light source.....	26		
6.8 Equipping the Incident Light Turret Disc	27		
6.9 Inserting the UV stray light shield	28		
6.10 Leica CTR6/CTR6 LED electronics box ...	29		
6.11 The Smartmove control.....	29		
6.12 Cabling	30		

1. Important notes on this manual

Attention!

This manual is an integral part of the Leica LMD6/7 and must be read carefully before switching on and using the system!

It contains important instructions and information for safe operation and maintenance of the product and must therefore be kept in a safe place.

This manual refers to both system variants Leica LMD6 and LMD7. Both systems differ in the built-in laser system.


The Leica LMD6 incorporates a laser of the type FTSS355-50, the Leica LMD7 a more powerful laser of the type Explorer 349.

The expression Leica LMD6/70 is used for the designation of both system variants.

Fig. 1 Full view of Leica LMD7 with Leica DM6 B (PC not illustrated)

1. Important notes on this manual

1. Important notes on this manual

Text symbols, pictograms and their meanings:

(1.2)

Numbers in parentheses such as “(1.2)”, correspond to illustrations (in the example, Figure 1, Item 2).

→ p.20

Numbers with an arrow, e.g. → p. 20, refer to a specific page of the manual.

WARNING indicates a hazard with a medium degree of risk that, if not avoided, can result in death or serious injury.

CAUTION indicates a hazard with a low degree of risk that, if not avoided, can result in minor or moderate injury.

Caution!

Special safety instructions within this manual are indicated with the triangle symbol shown here, and have a gray background.

Notice! The macroscope and accessories can be damaged when operated incorrectly.

Warning of hazardous electrical voltage! Risk of electrical shock!

Danger due to hot surface.

Warning of optical radiation! Never look directly into the light beam! Wear safety goggles!

Warning of electromagnetic field

1. Important notes on this manual

Warning of UV radiation! Wear safety goggles!

Warning of hazardous laser radiation – take appropriate safety precautions!

Warning of hand injuries

Instructions for disposing of the instrument, its accessories and consumables.

Connection to ground!

Explanatory note.

*

Item not contained in all configurations.

IVD manufacturing date, example 11 / 2011 for November 2011.

Device for in vitro diagnostics (IVD).

Manufacturer's address with manufacturing date

China RoHS 50 years EFUP
(Environmentally friendly use period)

Warning of hand injuries

2. Intended Use

Laser microdissection is an important, preparatory intermediate step in the area of biomedical research and molecular diagnostics for diseases like cancer that is also used to improve molecular examination methodologies.

This includes examining specimens taken from the human body for the purpose of gaining information about physiological or pathological conditions or inborn anomalies, or testing for safety and compatibility for potential recipients, or for monitoring therapeutic measures.

Contraindications:

Not suitable for testing potentially infectious specimens.

Examination methodologies and possible applications:

The examination methodologies include the PCR methodology, whereby small amounts of specific genetic material (DNA, RNA) can be amplified. Proteins can also be analyzed (proteomics).

Laser microdissection is usually used for histological sectioned materials and is tasked with isolating specific portions of tissue (such as tumorous material) from single cells to cell groups with precision and without any contamination in accordance with morphological criteria. Afterwards, this material is immediately ready for further analysis.

LCC (live cell cutting) provides the option of dissecting living cells from cell cultures without any risk of contamination.

The optional AVC (auto vision control) module provides an automated microdissection solution where single cells are detected, sectioned and collected automatically.

11888869 (LMD6 LED)

11888871 (LMD6 halogen)

11888877 (LMD7 LED)

11888876 (LMD7 halogen)

IVD

The four systems listed above satisfy EU Directive 98/79/EC regarding in vitro diagnostics.

The laser microdissection systems:

11888868 (LMD6 LED)

11888870 (LMD6 halogen)

11888875 (LMD7 LED)

11888878 (LMD7 halogen)

are not approved for IVD applications:

Examining specimens taken from the human body for the purpose of gaining information about physiological or pathological conditions or inborn anomalies, or testing for safety and compatibility for potential recipients, or for monitoring therapeutic measures is prohibited.

These systems are intended for the routine examination of biological samples.

Contraindications:

Not suitable for testing potentially infectious specimens.

The requirements in EU Directives 2006/95/EC and 2014/35/EU regarding the safety of electrical equipment, and 2004/108/EC and 2014/30/EU regarding electromagnetic compatibility have been met for all the systems mentioned above for use in industrial environments.

**Attention!**

The manufacturer assumes no liability for damage caused by, or any risks arising from using the microscopes for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In such cases the conformity declaration shall cease to be valid.

**Attention!**

These (IVD) devices are not intended for use in the patient environment defined by DIN VDE 0100-710. Neither are they intended for combining with medical devices according to EN 60601-1. If a microscope is electrically connected to a medical device according to EN 60601-1, the requirements defined in EN 60601-1-1 shall apply.

Owing to the special risks of the laser system, operators of the Leica LMD6/7 must have received specific training on the Leica LMD6/7 under consideration of the specific dangers of a laser system.

Reasonably foreseeable misuse

The following are prohibited:

- To use the microscope for any purpose not in accordance with the Declaration of Conformity
- Making use of auxiliary equipment to use the area between the microscope stage and objectives as a clamp or holder (working like a vice)
- Operating the microscope in an inclined position
- Cleaning the microscope in a way other than specified in the manual
- Modifying the built-in safety circuits in the instrument
- Allowing unauthorized personnel to open the instrument
- Using cables that Leica Microsystems CMS GmbH has not provided or permitted
- Using combinations with non-Leica components that go beyond the scope of the manual
- Operating light sources if they are not connected.
- Using the Touch Screen of the microscope as a transport handle
- Do not couple the laser to housing openings/ports not intended or not approved for such use.

3. Safety notes

3.1 General safety notes

These devices of protection class 1 are built and inspected in accordance with EN 61010-1 / IEC 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use and EN 61010-2-101 / IEC 61010-2-101 Safety requirements for electrical equipment for measurement, control and laboratory use, Part 2 Particular requirements - In vitro diagnostic (IVD) medical equipment.

They also meet the requirements in EN / IEC 61326-1 Electrical equipment for measurement, control and laboratory use – EMC requirements and EN / IEC 61326-2-6, Part 2-6 Electrical equipment for measurement, control and laboratory use – EMC requirements – Particular requirements – In vitro diagnostic (IVD) medical equipment.

LMD6/LMD7 have also been tested in accordance with IEC/EN 60825-1, Safety of laser products, and classified as laser class 1.

Attention!

In order to maintain this condition and to ensure safe operation, the user must note and adhere to the following instructions and warnings.

Attention!

The instruments and accessories described in this manual have been safety-tested and checked for possible hazards.

Before modifying the instrument in any way or combining it with non-Leica products not dealt with in this manual, it is essential to consult the Leica agency for your area or the main factory in Wetzlar!

Any unauthorized interference with the instrument or use of the instrument for applications for which it is not designed will automatically cancel any warranty claim and product liability!

The safety information in the Leica DM6 B manual must also be observed!

3.2 Electrical safety

Ambient temperature Leica LMD6 and LMD7 system: 18°-35°C

General specifications

Microscope

For indoor use only.

Supplied by Leica CTR electronics box.

Power input: See CTR6/CTR6 LED

Fuses: See CTR6/CTR6 LED

Ambient temperature: 15-35°C

Relative humidity: max. 80% up to 30°C non-condensing

Pollution degree: 2

Leica CTR6 electronics box

For indoor use only.

Supply voltage: 100-240 V AC
 Frequency: 50/60 Hz
 Power input: max. 290 VA
 Fuses: 5x20, T6,3A H 250V AC
 (IEC 60127-2/3)
 Ambient temperature: 15-35°C
 Relative humidity: max. 80% to 30°C
 non-condensing
 Overvoltage category: II
 Pollution degree: 2

Leica CTR6 LED electronics box

For indoor use only.

Supply voltage: 100-240 V AC
 Frequency: 50/60 Hz
 Power input: max. 150 VA
 Fuses: 5x20, T3,15A H 250V AC
 Ambient temperature: 15-35°C
 Relative humidity: max. 80% to 30°C
 non-condensing
 Overvoltage category: II
 Pollution degree: 2

Leica EL6000

For indoor use only.

Supply voltage: 100-240 V AC ($\pm 10\%$)
 Frequency: 50/60 Hz
 Power consumption: max. 210 W
 Fuses: 5x20, T2,5A H 250V AC
 → EL6000 manual
 Ambient temperature: 0°-40°C
 Relative humidity: 10-90%
 non-condensing
 Overvoltage category: II
 Contamination class: 2
 (see enclosed manual)

Laser

For indoor use only.

FTSS355-50 LASER SYSTEM

Supply voltage: 100-240 V AC
 Frequency: 50-60 Hz
 Power consumption: < 100 W
 Ambient temperature: 18°-38°C
 Relative humidity: max. 80%
 non-condensing
 External power supply
 Type: Cincon TRG70A240

Attention!

Use the original power supply only. Original power supplies are available from your Leica branch office or Leica dealer.

E-Box LMD7000 for

Explorer 349 LASER SYSTEM

Supply voltage: 100-240 V AC
 Frequency: 50/60 Hz
 Power consumption: 140 VA
 Ambient temperature: 18°-35°C
 Relative humidity: max. 80%
 non-condensing
 Pollution degree: 2

3. Safety notes

WARNING

The mains plug may only be inserted into a grounded socket.

Any extension cable used must be similarly grounded. Any break in the ground conductor inside or outside the instrument or loosening of the ground conductor connection can render the unit dangerous. Intentional severance is forbidden.

WARNING

Use only original power cables or alternative cables with a VDE / HAR logo, which at least fulfill the requirement of $3 \times 0.75 \text{ mm}^2$ and 10A/250V.

WARNING

By definition, the main circuit breaker for this system is the connection between the power cable and CTR box device connection and the external power supply (for the LMD6) or the CTR box and the LMD7000 E-box (for LMD7). The user must ensure unobstructed access to the main circuit breaker at all times.

Attention!

Using the ground connection, any accessories connected to the microscope which have their own and/or a different power supply can be given the same ground conductor potential. Please consult our servicing personnel if you intend to connect units without a ground conductor.

Attention!

Only fuses of the specified type and rating may be used. Never bridge the fuse holders.

Attention!

The electric accessories of the microscope are not waterproof. If water gets inside them, it may cause electrical shock.

Do not put the microscope and its accessories too near a water supply or anywhere else where water may get inside them.

Before changing fuses or lamps, always turn the mains switch off and disconnect the mains cable.

Attention!

Protect the microscope from major temperature fluctuations. These may lead to condensation which can damage the electric and optical components.

Achtung!

Do not use the microscope in altitudes exceeding 2000 m ASL/NL. 2.

Caution!

Transport and storage in a range of $-20^{\circ} - +85^{\circ}\text{C}$ and at a humidity not exceeding 90%.

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

3.3 Laser safety

The Leica LMD6 incorporates a FTSS355-50 LASER SYSTEM, the Leica LMD7 an Explorer 349 LASER SYSTEM.

According to IEC/EN 60825-1:2007, the built-in laser belongs to Laser Class 3b and the system as a whole to Laser Class 1.

Warnings in connection with laser radiation for both system variants are indicated with labels in the places described below:

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Fig. 2a Left side of Leica LMD6


3.4 Protection filter for possible laser emission points

The laser in the Leica LMD6/7 is only used when the system is fully assembled. Therefore, laser light may be emitted at the following exits: eye-pieces (2b.1), specimen area covered with the UV hood (2b.2) and the bottom of the condenser (2b.3).

All observation ports are equipped with UV shields.

The system as a whole complies with the specifications for Laser Class 1.

3.5 Technical Data of the built-in laser system

FTSS355-50 LASER SYSTEM

Wavelength:	355 nm
Repetition rate:	80 Hz
Pulse width:	~1 ns
Pulse energy:	> 70 µJ

Explorer 349 LASER SYSTEM

Wavelength:	349 nm
Repetition rate:	< 5 kHz
Pulse width:	< 4 ns
Pulse energy:	> 120 µJ

The user is not entitled to open the laser unit!
The laser system is encapsulated and may only be repaired or exchanged by Leica CMS authorised personnel. Improper use of the laser may lead to injury!

Fig. 2b Possible Laser emission points

3. Safety notes

If using other operating or adjusting equipment or other methods than the ones described here, you may expose yourself to dangerous radiation.
Risk of injury!

3.6 Transport and storage

Attention!

Transport and storage in a range of $-25^{\circ} - +70^{\circ}\text{C}$ and at a humidity not exceeding 80%.

CAUTION

Fully equipped, the microscope weighs more than 18 kg. For transportation, the user has to take care of the corresponding actions.

3.7 Disposal

After the end of the product's life, please contact Leica Service or Leica Sales on how to dispose of it.

Please observe the national laws and ordinances which, for example, implement and ensure compliance with EU directive WEEE.

Note:

Like all electronic instruments, the microscope, its components and expendables may not be disposed of as general household waste!

4. Installation site 5. Unpacking

The Leica LMD6/7 should be used in a room that is free from dust, oil, chemical fumes and extreme humidity. Vibrations, direct sunlight and major temperature fluctuations should also be avoided at the workplace.

Ambient conditions:

Temperature 20°-35°C

Relative humidity 20-80% up to 35°C

Microscopes in warm and humid climates need special care to prevent fungus growth.

Please also read the instructions on the installation site in the microscope manual!

Attention!

When installing the system, make sure the power inlet is freely accessible so that the instrument can be quickly disconnected from the mains if necessary.

Attention!

Electrical components including lamphousings and lasers must be installed at least 10 cm away from walls and other inflammable objects.

Attention!

The Leica LMD6/7 must not be placed in close proximity to an air-conditioning system and it must not be exposed to any direct air draft.

Please compare the delivery carefully with the packing note, delivery note or invoice. We strongly recommend you keep a copy of these documents with the manual, so that you have information on the time and scope of the delivery later when ordering more equipment or when the microscope is serviced. Make sure that no small parts are left in the packing material. Some of our packing material has symbols indicating environment-friendly recycling.

Note:

Keep the packing material for storage and transport of the Leica LMD6/7 and its accessories.

Attention!

Avoid touching optical surfaces if possible. Any fingerprints on glass surfaces should be removed with a soft leather or linen cloth. Even slight traces of finger perspiration can soon attack the surfaces of optical instruments.

On no account should you connect the laser and peripherals to the mains before the whole system has been assembled!

6. Description and assembly of the components

First of all, take all the components out of the transport and packing material.

On no account should you connect the electronics box, the laser or the peripherals to the mains before completing the assembly of the whole system!

The LMD module with the laser is already pre-assembled on the Leica DM6 B microscope. The remaining microscope components are logically assembled in this order:

- Stage with collector device
- Condenser with condenser head
- Tube
- Eyepieces
- Objectives
- Lamp housings with light sources
- Equipping the incident light turret disc

Only a few commonly used screwdrivers and keys are necessary for assembly, which are included in the delivery package.

Note:

Please also read the assembly and operation instructions in the Leica DM6 B microscope manual!

Attention!

Only insert or unplug data and control circuits when the instrument is switched off; otherwise, the instrument may be damaged.

6. Description and assembly of the components

6.1 Stage

! Attention!

Before assembling the stage, make sure no objectives are installed!

There are two different stages available:
a motor stage (Fig. 6) and a scanning stage (Fig. 5).

- Using the condenser height adjuster (4.2), turn the condenser holder completely upwards, i.e. as close to the stage as possible.
- Loosen the stage clamp (3.1) slightly. A ball-headed screwdriver is supplied for this purpose in the case of the scanning stage.

- From above, set the stage clamp onto the dovetail guide (3.2) and push the stage downwards until the stage securing screw (4.1) seats on the upper end of the dovetail.
- Firmly tighten the stage clamp (3.1).
Again, use the ball-headed screwdriver for the scanning stage.

In spite of the low speed, there is a risk of pinching and crushing injuries when moving the scanning stage.
Never put objects or body parts in the travel path!

Fig. 3 Assembling the stage

- 1 Stage clamp
- 2 Dovetail guide

Fig. 4

Microscope stand with stage bracket and motorized condenser

- 1 Stage securing screw
- 2 Condenser height adjuster
- 3 Connection for motorized condenser

6. Description and assembly of the components

- Slide in the collection device as far as the stop (Fig. 5, 6, 7).

In the case of the scanning stage, the stage must be moved to its foremost position first.

- Place the specimen holder (Fig. 8) onto the stage.

Fig. 5

Motor stage with motorized collection device

1 4-fold tube holder with PCR tubes

Fig. 6

Scanning stage with motorized collection device

1 4-fold tube holder with PCR tubes


Fig. 7



Inserting the motorized collection device into the scanning stage. The stage must be moved to its foremost position first.

Fig. 8 Removable specimen holder

a 3-fold holder for scanning stage

b holder for motor stage

6. Description and assembly of the components

6.2 Condenser


You can use either a motorized or a manual condenser. Both are equipped with condenser head S28.

! Attention!

Assembly is only permitted to Leica CMS authorised personnel.

- Screw the condenser head into the condenser.
- Using the condenser height adjuster (10.4), turn the condenser holder (10.1) completely downwards.
- Unscrew the clamping screw for the condenser (10.3) far enough so that the condenser can be inserted from the front.
- From the front, insert the condenser into the condenser holder as far as it will go. On the underside of the condenser, there is an orientation pin, which must be located in the guiding notch (9.1).

Fig. 9 Condenser holder
1 Guiding notch

- Pull the condenser's clamping screw (10.3) so that the condenser is locked in place.

! Attention:

- It is necessary to turn the condenser holder (10.1) up again using the condenser height adjuster (10.4). The microscope must not be switched on before.
- Connect the motorized condenser by connecting the condenser cable (4.3) with the stand.


Note:

The condenser must be centered before using the microscope.

→ DM6 B manual
(Köhler illumination).

Fig. 10 Condenser holder
1 Condenser holder
2 Condenser centering
3 Clamping screw for condenser
4 Condenser height adjuster

6. Description and assembly of the components

6.3 Tube and Eyepieces

The tube is mounted to the stand either directly or with the use of intermediate modules. It is fastened in place with the side clamping screw.

- Only for the motorized tube MBDT:
Remove the transportation lock (11.1) from the bottom side of the tube.
Connect the tube to the stand with the connector socket (12.1).
- The eyepieces are inserted into the eyepiece tubes on the tube.

Note:

For eyepieces that are not included in shipment, we recommend to learn them in with the software Leica Application Suite - module: Acquire. This ensures that the information about total magnification on the Leica SmartTouch is correct.

6.4 Camera

You can choose the camera LMD CC7000 or one of the digital cameras from the Leica DFC range (e.g. Leica DFC7000 T).

LMD CC7000

Beside the camera the delivery comprises an external power supply, the camera connection cable and a PC card (Gigabit-Ethernet).

- Normally the PC card has been built into the PC at the factory. Otherwise, i.e. for retrofitting the LMD CC7000 camera, mount the PC card into your computer.
- Screw the TV adapter into the c-mount thread of the camera and mount the whole unit onto the tube (13a, 13b.1).

The Leica logo must point **to the front**, that means towards the eyepieces.

Fig. 11 Bottom side of the tube

1 Transportation lock

Fig. 12 Motorized tube connection

1 Connector socket

6. Description and assembly of the components

Leica DFC

- Screw the TV adapter into the c-mount thread of the camera and mount the whole unit onto the tube (13b.1).
- Affix the TV adapter with the clamp screw on the side.

The camera is not adjusted until the system is started up.

Fig. 13a Camera LMD CC7000
1 C-mount 0.55x



Fig. 13b

- 1 Camera connection
- 2 Laser module
- 3 Microscope adapter for Leica EL6000 compact light source
- 4 Transmitted light illumination

6.5 Objectives

The receptacles on the objective turrets are numbered (Fig. 14). The individual objectives have already pre-assigned positions at the factory according to their configuration. A list of the exact objective positions is provided in shipment ("Identification Sheet").

! Attention!

Cover unoccupied threads on the turret with dust protector caps!

Note:

We recommend to perform a parfocality adjustment with the software Leica Application Suite, module: Fine Tuning.

Fig. 14
tubes

Objective turret with labeled objective recep-

6. Description and assembly of the components

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation). Therefore, lamps have to be operated in closed housings and only after being mounted.

6.6 Lamphousing

- Mount the lamphousing (13.4) to the back of the stand and secure it with the clamp screw on the side.
- Instructions on inserting the lamp are given in the Leica DM6 B microscope manual.

6.7 Leica EL6000 compact light source

- Instructions on inserting the lamp are given in the Leica EL6000 compact light source manual.
- The microscope adapter for connecting the light guide is attached to the back of the stand (13.3).

Attention!

Connect the light guide to the microscope first to prevent exposing the user to the high-energy light output of the Leica EL6000 compact light source.

- Insert the light guide (15.1) into the microscope adapter (13.3, 15.2) and secure it with the clamping screw. Insert the opposite input of the light guide into the port (16.1) on the rear panel of the compact light source. The light guide must click into place.

! Attention!

Take care not to kink or otherwise damage the light guide when connecting it to the light source or microscope adapter. Do not overtighten the clamping screw.

Use only light guides with Storz long connectors to prevent damage to the unit and danger to the user (blinding hazard).

Fig. 15 Light guide with adapter

- 1 Light guide
- 2 Adapter for Leica microscopes


Fig. 16 Rear panel with connectors

- 1 Light port
- 2 Remote control port
- 3 AC input

6. Description and assembly of the components

Connect the light guide **at both ends** (light source/adaptor) **before** opening the shutter or attenuator.

The emitted light may result in eye or skin injuries or damage to materials.

Never look directly into the light emitted by the light guide.

Do not switch on the light source until the light guide is connected to the microscope.

Uncontrolled light output from the light guide poses a blinding hazard!

There is also a burn hazard at the light guide output!

Further information → Leica EL6000 manual.

Fig. 17 Filter cube front side

Fig. 18 Filter cube back side

Fig. 19 Removing the front cover

- 1 Filter receptacle
- 2 Retention pin
- 3 Front cover

6.8 Equipping the Incident Light Turret Disc

The receptacles on the turret are numbered. According to your equipment, the individual filter and/or reflector cubes have already pre-assigned positions. A list is provided along with your shipment ("Identification sheet").

Fig. 20 Inserting the filter or reflector cubes

- 1 Mounting

6. Description and assembly of the components

Insert the filter and reflector cubes in the following manner:

- Equip the incident light turret only when the microscope is switched off.
- Remove the front cover from the upper part of the microscope (Fig. 19). Push the retention pin (19.2) to move the turret.
Releasing the retention pin locks the turret.
- Insert the filter or reflector cube into the mounting in front of you according to the identification sheet provided.
To do so, place the filter or reflector cube on the **right** side and press it to the **left** into the mounting (Fig. 20).

Note:

Filter cubes that are not included in shipment have to be learned in with the Software LMD Hardware Configurator.

6.9 Inserting the UV stray light shield

- Slide the UV stray light shield into its groove from the front.
Make sure it clicks into position. Otherwise it will not be possible to operate the laser.

Note:

At 5-fold filter turret, the numbers are right below the mounting.

- Push the retention pin (19.2) and continue to turn the filter turret until you reach the next locking position.
- Again make sure that the turret engages (retention pin unlocks) and insert the next filter and/or reflector cube as described above.
- When all filters and reflector cubes have been inserted, close the front cover plate again.

6.10 Leica CTR6/CTR6 LED electronics box

The Leica CTR6/CTR6 LED electronics box (Fig. 22) is used to control the LMD system and the Leica DM6 B microscope.

The most convenient position for the electronics box is to the left of the microscope.

6.11 The Smartmove control

The Smartmove control is used for controlling the specimen stage (23.1, 23.2) and the image focus (23.3) by using the turning knobs. The variable function keys (23.5) and (23.6) are assigned magnification change and fine/coarse focusing before delivery. The keys are labeled accordingly.

The motorized stage can be moved to a special loading position for easily changing the specimens. Therefore the objective change keys (23.5) must be pressed simultaneously. The stage will be lowered and the slide will be moved to a front/right position. Pressing both keys again moves the stage to the former position.

For a comfortable operation the height position of the above knobs can be set to an individual level by turning the wheel (23.4).

Fig. 22


Leica CTR6 electronics box

- 1 Mains switch
- 2 Pilot indicator

Fig. 23 Smartmove control

- 1 Movement in x direction
- 2 Movement in y direction
- 3 Setting the focus
- 4 Individual knob height position
- 5 Predefined with magnification change up/down
- 6 Predefined with fine and coarse focusing

6. Description and assembly of the components

6.12 Cabling

Connecting up to the electronics box

For information about connecting the electronics box see also the enclosed manual of the Leica CTR Electronics box.

- Connect the port (24.1) at the microscope to the "Microscope" socket at the electronics box using the cable with the 50-pin plug.
- Connect the microdissection stage to the "X/Y stage" port at the electronics box.
- The motorized collection device is connected to the "Well Handling" port at the electronics box.
- The Smartmove control is connected to "SmartMove" at the electronics box.
- For the laser, port "EXT" at the electronics box is connected to port (25.2) on the right of the microscope, while port "Trigger 1" at the electronics box is connected to port (25.1) for the laser trigger.
- Connect the lamp power cable of the microscope (24.8) to the socket "12V 100W" at the electronics box (CTR6 only).
- To connect the PC, connect the USB port of the PC to the "USB" port at the microscope.

Note:

The laser can be activated via the software. To do this, the serial interface of the PC has to be connected to the RS232 port of the laser box (33.6, p.39 or 27.6).

Fig. 24 Rear side of stand

- 1 Connection to the Leica CTR6/CTR6 LED electronics box
- 2 Connections EXT and USB
- 3 FAN port (LMD7 only)
- 4 LASER HEAD 1 port
- 5 LASER HEAD 2 port
- 6 Connection incident light
- 7 Connection for transmitted light
- 8 Lamp power cable of the microscope (CTR6 only)


Fig. 25 Right side of the microscope

- 1 Laser trigger port
- 2 Laser port

Other ports on the microscope stand

- Connect the cable (24.4) and (24.5) at the microscope to the ports "LASER HEAD" (26.2) and "LASER HEAD CURRENT" (26.3) for the Leica LMD6 or to the ports "LASER HEAD 2" (27.4) and "LASER HEAD 1" (27.5) for Leica LMD7 at the back of the laser box.
- Connect the lamphousing 107/2 to the connection for transmitted light (24.7) and for the LED lamphousing to connection for transmitted light of the CTR6 LED electronics box.


Fig. 26 Rear side of the laser box for LMD6

- 1 Power supply
- 2 LASER HEAD port
- 3 LASER HEAD CURRENT port

Fig. 27 Rear side of the laser box for LMD7

- 1 Power supply
- 2 Connection external laser trigger
- 3 On-/off switch
- 4 LASER HEAD 2 port
- 5 LASER HEAD 1 port
- 6 Connection PC (Com 3)

6. Description and assembly of the components

Connecting the camera

LMD CC7000

- Connect the provided external camera power supply to the port (28.1).

Attention!

Use the original power supply only. Other power supplies must not be used.

In case of damage or malfunction of the original power supply, an exchange is mandatory. You must not repair the power supply. Original power supplies are available at your Leica local branch or dealer.

- Connect the port (28.2) to the PC card (Gigabit-Ethernet).

Leica DFC

- Connect the camera to the USB port of your PC.

Connection to the Power Supply

- Connect the electronics box to the power supply and the laser box to the power supply using the power cable.
- For the fluorescence version only:
Connect the EL6000 compact light source to the power supply (29.3), too.
- For the LMD CC7000 camera only:
Connect the external camera power supply to the mains socket.

Fig. 28 Backplane of the camera LMD CC7000

- 1 Power supply
- 2 Connection for PC card (Gigabit-Ethernet)
- N.C. not connected


Fig. 29 Rear panel of the EL6000 compact light source

- 1 Light port
- 2 Remote control port
- 3 AC input

7. Getting started

The Leica LMD6/7 is normally delivered in its full configuration. All the hardware components, such as the laser and the camera, are preset and the microdissection software is ready installed.

On changing the configuration, installing updates, using a different computer or camera, etc. it will be necessary to reconfigure and set individual components or reinstall software.

You will find all the software components needed for this on the supplied CD.

The following is a description of the individual steps for installing the software and setting the hardware components.

7.1 Installing microdissection software

First make sure to close all active programs.

- Insert the installation CD. The CD starts automatically, giving installation instructions. Close the text window.
- Use the windows explorer to switch to the LMDxxxx directory of the installation CD.
- Start the **SETUP.EXE** program.
- Follow the instructions of the installation program.
- The microdissection software is recorded as **Leica Laser Microdissection** in the **Programs** menu and put on the desktop as a program icon.

n.b.:

The differentiation between the Leica LMD6500 system and the Leica LMD7000 system takes place in the **Options/Settings/Misc** dialog at **Type of connected laser**.

n.b.:

The **Leica Hardware Configuration** program is installed together with the microdissection software. This program is used to configure your microscope.

n.b.:

If you use the LMD CC7000 camera and a firewall is activated on your computer, the microdissection software must be added to the list of exceptions to admit the communication between the program and the LMD CC7000 camera.


n.b.:

Optional software is also installed together with the microdissection software. To use the software insert the USB dongle in port of the PC.

7. Getting started


7.2 Configuring the Microscope

Make sure that the microscope and all the components are switched on. Start the **LMD Mic Configurator** and select the **Hardware Setup** module.


Choose the **Leica DM6** as the connected microscope, the type of connection (**USB**) and the used stage controller (**Leica**) and click on **Apply**.

If the firmware of certain hardware modules is not up to date, the Firmware Update Wizard gives the possibility to update the firmware now. Therefore click on **Yes**.

The following window is opened:


Select the option **for Out of date microscope modules** and click on **Start**. The firmware of these hardware modules is updated now.

7. Getting started

In the next step set the current configuration by choosing **Configure** in the **Setup** box and then the particular microscope component in the **Components** box. On the right side of the window the current configuration is shown in the upper part of the window and the available options and settings below in the database field.

DM6B


The microscope values of your stand are listed.

Ports


Select **Ports** in the **Components** box.

Tube

Select **MANUAL-TUBE** in the list box on the top right and choose your tube in the database list below.

Select Simulator for the CC7000 camera or a DFC camera according to its tube port.

7. Getting started

Eyepiece

Select **VISUAL** in the list box and choose the mounted eyepiece.

The screenshot shows the Nikon software interface. On the left, there are tabs for 'File', 'Cameras', and 'Help'. Below these are icons for 'Hardware Setup', 'Firmware Update', 'Configure', and 'Software Updater'. The 'Components' section includes icons for 'DM6000B', 'Cameras', 'Ports', 'IL Turret', 'Dic Turret', 'Nosepiece', 'Zdrive', 'Stage', 'Condenser', and 'Function Keys'. The 'Subcomponents' section has icons for 'MANUAL-TUBE', 'VISUAL', and 'Doc Ports'. The main area is titled 'VISUAL' and 'Setup Component Data'. It contains a table with the following data:


Position	ArticleNo	Magnification	Name
1	11507802	10	VISUAL

Below the table, there is a search filter: 'Search filter: ACTIVE, entries found (12/12)'. At the bottom, there is a detailed table of eyepieces:

ArticleNo	Name	Magnification
1 11507808	Eyepiece HC PLAN s 10x/25 Br. M	10
2 11557804	Eyepiece HC PLAN s 10x/25 Br.M POL	10
3 11507807	Eyepiece HC PLAN s 10x/22 Br. M	10
4 11557803	Eyepiece HC PLAN s 10x/22 Br.M POL	10
5 11507801	Eyepiece HC PLAN 10x/20 BR.	10
6 11507802	Eyepiece HC PLAN 10x/20 BR.M	10
7 11506515	Eyepiece HC PLAN 12.5x/16 BR. M	12.5
8 10445301	Eyepiece 16x/14B, adjustable	16
9 10445302	Eyepiece 25x/9.5B, adjustable	25
10 13613530	Eyepiece 10x/20 Nonfocusing with Eyeguard	10
11 13613531	Eyepiece 10x/20 Nonfocusing with Eyeguard & Pointer	10
12 13613532	Eyepiece 10x/20 Focusing with Eyeguard	10

C-Mount

Select **DOC-TOP** in the list box and choose the recommended c-mount adapter for the used camera.

7. Getting started

IL Turret

Filter cube(s)

Select **IL Turret** in the **Components** box to assign to your LMD filter the position in the incident light turret disc in which it has been inserted.

The screenshot displays the IL Turret software interface. On the left, the 'Components' panel lists various hardware components: DM6000B, Cameras, Ports, IL Turret (selected), Dic Turret, Nosepiece, Zdrive, Stage, Condenser, and Function Keys. The main area shows the 'IL Turret IL TURRET(8-POS)' setup. A table lists the 8 positions of the turret disc, with Position 1 currently selected and labeled 'Empty BF EMP'. To the right of this table is a circular diagram of the turret disc with 8 numbered positions. Below the turret setup, a search filter is set to 'ACTIVE', showing 142/142 entries. A table lists the available filter cubes with their names, excitation ranges, and excitation filters.

Filtercube	McFilterCubeName	ExcitationRange	ExcitationFilter	Dichromatic
1 A	A	UV	BP 340 - 380	400
2 A4	A4	UV	BP 360/40	400
3 D	D	UV/violet	BP 355 - 425	455
4 E4	E4	violet/blue	BP 436/7	455
5 H3	H3	violet/blue	BP 420 - 490	510
6 B	B	blue	BP 450 - 490	510
7 K3	K3	blue	BP 470 - 490	510
8 L5	L5	blue	BP 480/40	505
9 M2	M2	green	BP 546/14	580
10 N2.1	N21	green	BP 515 - 560	580
11 N3	N3	green	BP 546/12	565
12 G/R	G/R	blue: green	BP 490/20; BP 575/30	505; 600

Stage

Select **Stage** in the **Components** box and **STAGE** in the list box on the top right. Choose the micro-dissection stage (scanningstage or motorstage) out of the database list below.


The screenshot shows the Asten software interface. On the left, the 'Components' sidebar lists various hardware components, with 'Stage' highlighted. The main window displays the 'STAGE' configuration for article 11501420. It includes tabs for 'Setup', 'Component Data', and 'Fine Tuning'. Below these tabs is a search filter set to 'ACTIVE' with 36 entries found. A table lists available stages with columns for ArticleNo, Name, and Description.

	ArticleNo	Name	Description
1	11522057	Motorstage DMi 127x83	-
2	11522021	Motorstage DMi 35x35	-
3	11501222	Motorstage DMi 40x40	-
4	11522601	Motorstage DMi 100x40	-
5	11501240	Motorstage DMi 75x50	-
6	11522068	Motorstage DMi 127x83	-
7	11522069	Motorstage DMi 35x35	-
8	11501419	Motorstage LMD 6000	-
9	11522603	Motorstage DMi 100x40	-
10	11600182	Scanningstage DM 100x80	-
11	11532302	Scanningstage DM 100x100	-
12	12724177	Scanningstage DM 215x80	-
13	11532411	Scanningstage DM 150x150	-

7. Getting started

Condenser

Select **Condenser** in the **Components** box and **CONDENSER** in the list box on the top right. Choose the **LMD Condenser** out of the database list below.

The screenshot shows the Leica Hardware Configurator software. On the left, the 'Components' panel has 'Condenser' selected. The main panel shows the 'CONDENSER' configuration with article number 11023124014000 and S28. Below this is a table of condenser options.

ArticleNo	Name	MichName	Description
1 11522004	Smart Condenser S1/0.90 dry with mot. flip	S1-0.90	-
2 11522005	Manual Condenser S1/0.90 dry with man. flip	S1-0.90	-
3 11522006	Manual Condenser S1/0.90 dry w/o flip	S1-0.90	-
4 11522007	Smart Condenser S1/0.90 dry w/o flip	S1-0.90	-
5 11522004	Smart Condenser S1/1.40 oil with mot. flip	S1-1.40	-
6 11522005	Manual Condenser S1/1.40 oil with man. flip	S1-1.40	-
7 11522006	Manual Condenser S1/1.40 oil w/o flip	S1-1.40	-
8 11522007	Smart Condenser S1/1.40 oil w/o flip	S1-1.40	-
9 11522004	Smart Condenser S23/0.53 with mot. flip	S23	-
10 11522005	Manual Condenser S23/0.53 with man. flip	S23	-
11 11522006	Manual Condenser S23/0.53 w/o flip	S23	-
12 11522007	Smart Condenser S23/0.53 w/o flip	S23	-
13 11522004	Smart Condenser S28/0.55 with mot. flip	S28	-

n.b.:

The Leica Hardware Configurator (version 2010.4.0.3373 or higher) offers the facility to learn the operating range of the stage in z direction in order to avoid the collision risk between the objective and the sample. The maximum range is defined for each particular objective. See **Online Help** for further information.

n.b.:

For further information about camera alignment and settings see also **Online Help**, as well as the separate manual of the digital camera.

7.3 Switching on the Leica LMD6

- Switch on the Leica CTR6/CTR6 LED Electronics Box at the on/off switch (32.1). When in operation, the pilot lamp will light up green (32.2). All motorized microscope components first undergo an initialization phase. After initialization is completed, the Leica SmartTouch shows the current microscope setting.
- Switch on the PC.

Caution!

The service plug (33.5) must be plugged in.

- Switch on the laser:
After this, the laser is switched on with the key switch (33.1), i.e. the switch (33.1) must be set at 1.
The laser now goes through a warming up phase (approx. 5-10 min). The left-hand status LED (33.2) lights up orange, the right-hand status LED (33.4) is off.
- At the end of the warm-up phase, both status LEDs light up green. The laser is now ready for operation and can be activated with the push button (33.3). After activation, the right-hand status LED (33.4) lights up red.
- Make sure to deactivate the laser when not in use by pushing button (33.3). The right-hand status LED (33.4) then lights up green. If you reactivate the laser (by pressing button (33.3) again), it is immediately ready for use without another warming up phase. Otherwise its life span will be shortened!


Fig. 32 Leica CTR6 electronics box

- 1 On/off switch
- 2 Pilot lamp

Fig. 33 Front side of the laser box for LMD6

- 1 Key-operated switch
- 2 Status LED: Operating status of laser
- 3 Push-button: Laser on/off
- 4 Status LED: Laser active/inactive
- 5 Service plug
- 6 RS232 port

7. Getting started

Alternatively:

The laser can be activated/deactivated via the software. To do this, the serial interface of the PC has to be connected to the RS232 port on the front of the laser box (33.6).

Additionally, the option "Laser connected to computer" in the microdissection software has to be selected in the operating step "Configure".

The laser is then automatically activated for cutting and deactivated afterwards if it is not used. The laser timeout can be set by the user in the operating step "Configure".

Fig. 34 Rear side of the laser box for LMD6

1 Power supply

Fig. 35 Rear side of the laser box for LMD7

1 Power supply
2 Power switch

7.4 Switching on the Leica LMD7

- Switch on the Leica CTR6 Electronics Box at the on/off switch (32.1). When in operation, the pilot lamp will light up green (32.2). All motorized microscope components first undergo an initialization phase.

After initialization is completed, the Leica SmartTouch shows the current microscope setting.

- Switch on the PC.
- Switch on the laser:
The power switch (35.2) is switched on. After this, the laser is switched on with the key switch, i.e. the switch must be set at 1.
The laser now goes through a warming up phase (approx. 5-10 min).

For using fluorescence:

Switch on the power switch for the compact light source Leica EL6000 (36.1).

You can now start the microdissection program and make a basic settings.

Before you start, read the "Microdissection software" manual.

Fig. 36 Leica EL6000 compact light source

1 Power switch
2 Intensity control

8. Care and maintenance

Follow the instructions in the Leica DM6 B manual!

Specimen holders and collection devices can be sterilized by autoclaving or baking.

9. Main spare parts and consumables

Order no.	Description
SPARE LAMPS	
11-500.974	Halogen lamp 12 V 100 W Lamphousing 107/2
11-504.120	HXP 120 L spare lamp for EL6000
SCREW CAP FOR VACANT OBJECTIVE THREADS	
11-020-422.570-000	Screw cap M 25 for objective nosepiece
SPARE EYECUPS (ANTI-GLARE PROTECTION) FOR HC PLAN EYEPIECE	
11-507.808	Eyecup HC PLAN Eyepiece 10x/25
11-507.807	Eyecup HC PLAN Eyepiece 10x/22
11-507.802	Eyecup HC PLAN Eyepiece 10x/20
FLUORESCENCE ACCESSORY	
for simultaneous cutting and observation during Laser Microdissection	
11-513.911	Triple band filter RGB cube
11-532.715	LMD-GFP/LP, size "K"
11-532.651	GFP for LMD, ET, size "K"
11-505.260	CY3 filter set for LMD, ET, size "K"
11-505.261	Dapi filter set for LMD, ET, size "K"
11-505.263	Fluorescence filter LMD-CFP
11-505.265	Filter system LMD-GFP/CY3, ET size "K"
11-505.266	Filter system LMD-YFP, ET size "K"
11-505.269	Filter system LMD-CY5, ET size "K"
PARTS FOR MOTOR STAGE 11888826	
<u>Slide holders</u>	
11-532.732	Holder for 1 slide 25 x 76 mm
11-505.214	Holder for big slide 50 x 76 mm
11-505.257	Holder for Petri dish
<u>Collection devices</u>	
Accessories for PCR tubes collection device:	
11-090-615-026-000	Lever for PCR tubes, 0.2 ml, 1 piece
11-090-615-025-000	Lever for PCR tubes, 0.5 ml, 1 piece
11-505.131	Levers for PCR tubes, 0.2 ml, set of 4 pieces
11-505.130	Levers for PCR tubes, 0.5 ml, set of 4 pieces
11-505.258	Lever for 1x 8-well strips (accessory for Live Cell Cutting)

Order no.	Description
-----------	-------------

PARTS FOR SCANNING STAGE 11888827

Slide holders

11-505.226	Holder for 3 slides 25 x 76 mm
11-505.227	Holder for Petri dish and big slide 50 x 76 mm
11-505.255	Holder for 18-well double stack ibidi slide

Collection devices

11-505.229	Collector for PCR tubes, 4 x 0.2 ml
11-505.228	Collector for PCR tubes, 4 x 0.5 ml
11-505.276	Universal collector for 2x 8-well strips/2x8-well strip caps/ LOC (Lab on a chip)/multi-well slides

CONSUMABLES

Membrane slides

11-505.158	Glass, PEN-membrane (2 µm) 50 slides/box
11-505.189	Glass, PEN-membrane (2 µm) 50 slides/box, RNase-free
11-600.288	Glass, PEN-membrane (4 µm) 50 slides/box
11-505.273	Glass, PPS-membrane, (1.2 µm), 50 slides/box
11-505.151	Steel frames, PET-membrane (1.4 µm) 50 slides/box
11-505.190	Steel frames, PET-membrane (1.4 µm) 50 slides/box, RNase-free
11-505.188	Steel frames, POL-membrane (0.9 µm) 50 slides/box
11-505.191	Steel frames, POL-membrane (0.9 µm) 50 slides/box, RNase-free
11-600.289	Steel frames, PEN-membrane (4 µm) 50 slides/box
11-505.268	Steel frames, PPS-membrane, (1.2 µm), 50 slides/box
11-505.215	Large slides for neuroscience, Glass PEN-membrane (2 µm) 100 slides/box
11-600.264	Demo Slides LMD (set)

Consumables for Live Cell Cutting (LCC)

11-505.172	Petri dish (Ø 5 cm), PEN membrane, sterile, 20 dishes/box
11-505.240	µ-Slide with PEN membrane, 18 well, ibidi, sterile, 15 slides/box
11-532.638	Stackable MembraneRing, Petri dishes (Ø 5 cm), PEN-Membrane, sterile, 4 pieces/box

Slide accessories

11-532.325	Plexiglas frame support for steel-frame slides to facilitate mounting of tissue section on the membrane
------------	--

Fuses

11 362 150 010 160	T1,6 A, 250 VAC, Breaking capacity H
--------------------	--------------------------------------

10. Index

4-fold tube holder 22
12 V 100 W 30

A

Activation of the laser 43, 44
Ambient conditions 19

C

Cabling 30
Camera 24, 32
Collection device 22, 30
Compact light source Leica EL6000
26, 32, 44
Condenser 23
Condenser height adjuster 21, 23
Condenser holder 23

D

Disposal 18

E

Electronics box Leica CTR6/CTR6
LED 29, 30, 43, 44
Explorer 349 LASER SYSTEM
13, 15, 17
Eyepieces 24

F

Filter cube 28
FTSS355-50 LASER SYSTEM 15, 17

I

Identification sheet 25, 27
Incident Light Turret Disc 27
Initialization 43, 44
Installing microdissection software
33

K

Key switch 43, 44
Köhler illumination 23

L

Lamphousing 107/2 26, 31
Lamp power cable 30
Laser connected to PC 44
Laser emission points 17
Laser safety 15
Laser timeout 44
Laser trigger 30
Leica CTR6/CTR6 LED 29, 43, 44
Leica EL6000 13, 26, 32, 44
Leica LMD6 6, 43
Leica LMD7 6, 44
Leica SmartTouch 24, 43, 44
Light guide 26

M

Microdissection stage 30
Motor stage 22

O

Objectives 25

P

Power supply 32
Protection filter 17

R

Reflector cube 28
Retention pin 28
RS232C 44

S

Safety labels 15
Scanning stage 22
Service plug 43
Smartmove 29
Smartmove control 30
Specimen holder 22
Stage 21
Stage securing screw 21
Storage 18

T

Technical data 12, 17
Transport 18
Tube 24
Tube MBDT 24

U

UV stray light shield 28

11. EU Declaration of Conformity

To download the EU Declaration of Conformity for your product use the link

<http://www.leica-microsystems.com/products/light-microscopes/life-science-research/laser-micro-dissection/details/product/>

Select the type of your system and switch to the Download page.

If you have any further questions or require technical support, please directly contact your country's Leica branch office or your local contact person. The appropriate contacts can be found on the Internet under:

<http://www.leica-microsystems.com>

- Administrative Measures on the Control of Pollution Caused by Electronic Information Products -

部件名称 Name of the part	有毒有害物质或元素 Hazardous substances					
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)
印刷电路板 printed circuit boards	X	O	O	O	O	O
电子元件 electronic components	X	O	O	O	O	O
机械部件 mechanical parts	X	O	O	X	O	O
光学元件 optical components	X	O	X	O	O	O
电缆 cables	O	O	O	O	X	X
光源 light sources	X	X	X	O	O	O

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Indicates that the concentration of the hazardous substance in all materials in the parts is below the relevant threshold of the SJ/T 11363-2006 standard.

x : 表示该有毒有害物质至少在该部件的某一材料中的含量超出SJ/T 11363-2006 标准规定的限量要求。
Indicates that the concentration of the hazardous substance of at least one of all materials in the parts is above the relevant threshold of the SJ/T 11363-2006 standard.

Note: The actual product may or may not include in all the part types listed above

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	ArticleNo	Name	Magnification	
1	11507808	Eyepiece HC PLAN s 10x/25 Br. M	10	
2	11557804	Eyepiece HC PLAN s 10x/25 Br.M POL	10	
3	11507807	Eyepiece HC PLAN s 10x/22 Br. M	10	
4	11557803	Eyepiece HC PLAN s 10x/22 Br.M POL	10	
5	11507801	Eyepiece HC PLAN 10x/20 BR.	10	
6	11507802	Eyepiece HC PLAN 10x/20 BR.M	10	
7	11506515	Eyepiece HC PLAN 12.5x/16 BR. M	12.5	
8	10445301	Eyepiece 16x/148, adjustable	16	
9	10445302	Eyepiece 25x/9.58, adjustable	25	
10	13613530	Eyepiece 10x/20 Nonfocusing with Eyeguard	10	
11	13613531	Eyepiece 10x/20 Nonfocusing with Eyeguard & Pointer	10	
12	13613532	Eyepiece 10x/20 Focusing with Eyeguard	10	

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The screenshot shows the IL Turret software interface. On the left, the **Components** panel lists various hardware components. The **IL Turret** component is highlighted. The main panel displays the **IL Turret IL TURRET(8-POS)** setup, showing a circular turret with 8 positions. The **Position** and **Name** columns are visible. Below the turret, a search filter is set to **ACTIVE**, showing 142 entries found. A table lists the filtercube configurations.

	Filtercube	McFilterCubeName	ExcitationRange	ExcitationFilter	Dichromatic
1	A	A	UV	BP 340 - 380	400
2	A4	A4	UV	BP 360/40	400
3	D	D	UV/violet	BP 355 - 425	455
4	E4	E4	violet/blue	BP 436/7	455
5	H3	H3	violet/blue	BP 420 - 490	510
6	B3	B3	blue	BP 450 - 490	510
7	K3	K3	blue	BP 470 - 490	510
8	L5	L5	blue	BP 480/40	505
9	M2	M2	green	BP 546/14	580
10	N2.1	N21	green	BP 515 - 560	580
11	N3	N3	green	BP 546/12	565
12	G/R	G/R	blue: green	BP 490/20; BP 575/30	505; 600

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The screenshot shows the Stage software interface. On the left, the 'Components' menu is open, displaying various hardware components like DM6000B, Cameras, Ports, IL Turret, Dic Turret, Nosepiece, Zdrive, Stage, Condenser, and Function Keys. The 'Stage' component is highlighted. On the right, the 'STAGE' configuration window is shown, displaying a table of stage configurations. The table has columns for ArticleNo, Name, and Description. The search filter is set to 'ACTIVE', and 36 entries are found.

	ArticleNo	Name	Description
1	11522057	Motorstage DMI 127x83	-
2	11522021	Motorstage DMI 35x35	-
3	11501222	Motorstage DMI 40x40	-
4	11522601	Motorstage DMI 100x40	-
5	11501240	Motorstage DM 75x50	-
6	11522068	Motorstage DMI 127x83	-
7	11522069	Motorstage DMI 35x35	-
8	11501419	Motorstage LMD 6000	-
9	11522603	Motorstage DMI 100x40	-
10	11600182	Scanningstage DM 100x80	-
11	11532302	Scanningstage DM 100x100	-
12	12724177	Scanningstage DM 215x80	-
13	11532411	Scanningstage DM 150x150	-

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The screenshot shows the Leica Hardware Configurator software interface. On the left, the 'Components' panel is active, displaying various hardware options. The 'Condenser' option is selected, and the 'Subcomponents' list below it shows 'CONDENSER' and 'CondenserTu (REL COND. TURRET) (7-POS)'. The main window displays the 'CONDENSER' configuration page for article number 11023124014000 S28. It includes a table with component data and a search filter set to 'ACTIVE', showing 37 entries found.

Position	ArticleNo	Name
1	11023124014000	11023124014000 S28

ArticleNo	Name	McName	Description
1 11522004	Smart Condenser S1/0.90 dry with mot. flip	S1-0.90	-
2 11522005	Manual Condenser S1/0.90 dry with man. flip	S1-0.90	-
3 11522006	Manual Condenser S1/0.90 dry w/o flip	S1-0.90	-
4 11522007	Smart Condenser S1/0.90 dry w/o flip	S1-0.90	-
5 11522004	Smart Condenser S1/1.40 oil with mot. flip	S1-1.40	-
6 11522005	Manual Condenser S1/1.40 oil with man. flip	S1-1.40	-
7 11522006	Manual Condenser S1/1.40 oil w/o flip	S1-1.40	-
8 11522007	Smart Condenser S1/1.40 oil w/o flip	S1-1.40	-
9 11522004	Smart Condenser S23/0.53 with mot. flip	S23	-
10 11522005	Manual Condenser S23/0.53 with man. flip	S23	-
11 11522006	Manual Condenser S23/0.53 w/o flip	S23	-
12 11522007	Smart Condenser S23/0.53 w/o flip	S23	-
13 11522004	Smart Condenser S28/0.55 with mot. flip	S28	-

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The screenshot shows the MTA software interface. The top menu bar includes 'File', 'Cameras', and 'Help'. The left sidebar contains icons for 'Hardware Setup', 'Firmware Update', 'Configure', and 'Software Updater'. The main panel is titled 'VISUAL' and contains a 'Setup' section with a table of components. The 'VISUAL' button in the sidebar is highlighted.


Position	ArticleNo	Magnification	Name
1	11507802	10	VISUAL

Search filter: ACTIVE, entries found (12/12)

ArticleNo	Name	Magnification
11507808	Eyepiece HC PLAN s 10x/25 Br. M	10
11557804	Eyepiece HC PLAN s 10x/25 Br.M POL	10
11507807	Eyepiece HC PLAN s 10x/22 Br. M	10
11557803	Eyepiece HC PLAN s 10x/22 Br.M POL	10
11507801	Eyepiece HC PLAN 10x/20 BR.	10
11507802	Eyepiece HC PLAN 10x/20 BR.M	10
11506515	Eyepiece HC PLAN 12.5x/16 BR. M	12.5
10445301	Eyepiece 16x/148, adjustable	16
10445302	Eyepiece 25x/9.58, adjustable	25
13613530	Eyepiece 10x/20 Nonfocusing with Eyeguard	10
13613531	Eyepiece 10x/20 Nonfocusing with Eyeguard & Pointer	10
13613532	Eyepiece 10x/20 Focusing with Eyeguard	10

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The screenshot shows the JVC software interface with the 'Doc. Ports' window open. The window title is 'Doc. Ports' and the subtitle is 'HC C-Mount 0.55x DOC-TOP Simulator'. The 'Setup Component Data' section shows a table with the following data:

Position	ArticleNo	Magnification	Name
1	11541544	0.55	HC C-Mount 0.55x DOC-TOP Simulator
2	11543669	0.63	C-mount adapter 0.63x DOC-REAR DFC295-317552710

Below the table, there is a search filter: 'ACTIVE, entries found (30/30)'. The 'Subcomponents' section shows a list of components with the following data:

ArticleNo	Name	Magnification
1 10445772	C-mount adapter 0.5x	0.5
2 11541510	C-mount adapter 1x HC f. 1"	1
3 11541537	C-mount adapter 0.63x HC f. 2/3"	0.63
4 11541511	C-mount adapter 0.5x HC f. 1/2"	0.5
5 11541512	C-mount adapter 0.35x HC f. 1/3"	0.35
6 11541538	TV-Optics HC 0.5x for C/B/F-mount head	0.5
7 11543706	C-mount head 1x f. TV optics 0.5x HC	1
8 11543702	B-mount head 1x f. TV optics 0.5x HC	1
9 11541539	B-mount-Head 1.25x f. TV optics 0.5x HC	1.25
10 11541540	F-mount head 1x f. TV optics 0.5x HC	1
11 11541541	F-mount head 1.25x f. TV optics 0.5x HC	1.25
12 11541006	C-mount adapter 1x, delta	1
13 11541016	C-mount adapter 0.5x, delta	0.5

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The screenshot shows the IL Turret software interface. On the left, the 'Components' panel lists various hardware components, with 'IL Turret' selected. The main window displays the 'Empty BF EMP' setup for the 'IL Turret IL TURRET(8-POS)'. A table lists the filter configurations for each of the 8 positions. A circular diagram on the right shows the turret layout with positions 1 through 8. Below the table, a search filter is set to 'ACTIVE', showing 142/142 entries found.

	Filtercube	MicFilterCubeName	ExcitationRange	ExcitationFilter	Dichromatic
1	A	A	UV	BP 340 - 380	400
2	A4	A4	UV	BP 360/40	400
3	D	D	UV/violet	BP 355 - 425	455
4	E4	E4	violet/blue	BP 436/7	455
5	H3	H3	violet/blue	BP 420 - 490	510
6	B	B	blue	BP 450 - 490	510
7	K3	K3	blue	BP 470 - 490	510
8	L5	L5	blue	BP 480/40	505
9	M2	M2	green	BP 546/14	580
10	N2.1	N21	green	BP 515 - 560	580
11	N3	N3	green	BP 546/12	565
12	G/R	G/R	blue: green	BP 490/20; BP 575/30	505; 600

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The screenshot shows the MIP software interface. On the left, the 'Components' menu is open, and 'Stage' is selected. The main window displays the 'STAGE' configuration for '11501420 STAGE'. Below this, a search filter is set to 'ACTIVE', showing 36/36 entries. A table lists various stage components with their ArticleNo, Name, and Description.

	ArticleNo	Name	Description
1	11522057	Motorstage DMI 127x83	-
2	11522021	Motorstage DMI 35x35	-
3	11501222	Motorstage DMI 40x40	-
4	11522601	Motorstage DMI 100x40	-
5	11501240	Motorstage DM 75x50	-
6	11522068	Motorstage DMI 127x83	-
7	11522069	Motorstage DMI 35x35	-
8	11501419	Motorstage LMD 6000	-
9	11522603	Motorstage DMI 100x40	-
10	11600182	Scanningstage DM 100x80	-
11	11532302	Scanningstage DM 100x100	-
12	12724177	Scanningstage DM 215x80	-
13	11532411	Scanningstage DM 150x150	-

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The screenshot shows the Leica software interface for configuring the Condenser. The sidebar on the left has a 'Components' section with icons for DM6000B, Cameras, Ports, IL Turret, Dic Turret, Nosepiece, Zdrive, Stage, Condenser, and Function Keys. Below this is a 'Subcomponents' section with icons for CONDENSER and CondenserTu (REL COND. TURRET(7-POS)). The main area is titled 'CONDENSER' and shows the setup for component 11023124014000 S28. A table lists the components:

ArticleNo	Name	MicName	Description
11522004	Smart Condenser S1/0.90 dry with mot. flip	S1-0.90	-
11522005	Manual Condenser S1/0.90 dry with man. flip	S1-0.90	-
11522006	Manual Condenser S1/0.90 dry w/o flip	S1-0.90	-
11522007	Smart Condenser S1/0.90 dry w/o flip	S1-0.90	-
11522004	Smart Condenser S1/1.40 oil with mot. flip	S1-1.40	-
11522005	Manual Condenser S1/1.40 oil with man. flip	S1-1.40	-
11522006	Manual Condenser S1/1.40 oil w/o flip	S1-1.40	-
11522007	Smart Condenser S1/1.40 oil w/o flip	S1-1.40	-
11522004	Smart Condenser S23/0.53 with mot. flip	S23	-
11522005	Manual Condenser S23/0.53 with man. flip	S23	-
11522006	Manual Condenser S23/0.53 w/o flip	S23	-
11522007	Smart Condenser S23/0.53 w/o flip	S23	-
11522004	Smart Condenser S28/0.55 with mot. flip	S28	-

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Copyright © Leica Microsystems CMS GmbH • Ernst-Leitz-Straße • 35578 Wetzlar • Germany 2016 • Tel. (06441)29-0 • Fax (06441)29-2539 LEICA and the Leica logos are registered trademarks of Leica IR GmbH.