
OPERATING INSTRUCTIONS

FOR

CLS 100 LED

power supply can change depending on availability
(can be either desktop OR wall plug power supply)

- 1 .. LED green (POWER)
- 2 .. LED red (TEMP.)
- 3 .. Brightness control
- 4 .. Lightguide connection
- 5 .. ESD connection
- 6 .. 24DC input
- 7 .. Type plate
- 8 .. Socket for mains cable
- 9 .. 24VDC plug connector
- 10 .. On/Off switch
- 11 .. Remote (optional)

OPERATING INSTRUCTIONS FOR CLS 100 LED COLD LIGHT SOURCE

Please read all the operating instructions completely and carefully before operating this device!

TECHNICAL DATA

Lamps	High-power LED
Range control	min. 30:1
Color temperature typ.	5,500 K / daylight
Colour reproduction index	80
Power input typ.	22W
Power supply	24V DC +/- 5%, 1000mA Socket f. hollow pin 5.5 x 2.1 mm
On/Off switch	✓
Brightness control	Infinitely variable using rotary potentiometer, with Index marks 1-6
Indicator lamps	LED RED (Temperature status) LED GREEN (Functional status)
ESD socket	for 4mm banana plug for potential equalisation
Cooling	without fan by convection
Housing	vents on top and side
Remote	only with "R" models
Lightguide connection	Photonic standard dia 15mm (galvanically isolated)
Operating environment for normal operation	10-30 °C and 30-70 % relative humidity
Storage/transportation	0-40 °C and 30-70 % relative humidity
Dimension height x width x depth	96 x 110 x 175 mm without controls
Weight	2500g
Type of protection (Protection class)	III
Conformity	CE (EMC Directive 89/336/EEC), RoHS

(*) LEDs are suitable for use at temperatures up to 40 °C. However, increasing temperature and humidity results in a reduction in service life. For optimal use, the ambient temperature should be 25 °C (or lower) and the relative humidity 50 % (or lower).

Description

This new type of LED illumination device (Light Emitting Diode) has been developed for stereo-microscopy illumination purposes in industrial and laboratory settings. In purchasing this product, you have acquired a high-quality LED illumination unit for all applications requiring intensive light with minimal heat generation. LEDs are employed as the illuminant.

Only use spareparts and accessories which are approved for use in medical applications. Contact your dealer for more information.

The following features are characteristic of this type of illumination:

- ▼ Color temperature 5,500 K (daylight!)
- ▼ Ripple and flicker free
- ▼ Long LED service life
- ▼ Quiet vibration-free operation
- ▼ Extremely compact construction

This light source has a light guide connection (4) for fibre optic products such as ring lights, goose necks, transmitted light illuminators and flexible light guides.

The modular design of the LED illumination device provides in combination with fibre-optic

- ▼ ringlights - homogeneous shadow-free illumination
- ▼ flexible lightguides – directional illumination
- ▼ transmitted light illuminators homogenous transmitted light illumination

SYMBOLS USED ON THE DEVICE

CE Mark of conformity: Confirms conformity with EU directive.

Lights, suitable for direct attachment to normal inflammable mounting surfaces.

This symbol indicates the year of manufacture.

SAFETY INFORMATION

Two warning symbols are used throughout this document:

Caution: Risk of electrical shock!

Caution: Dangerous area. Warning: Refer to accompanying documentation!

This section contains safety information which must be observed strictly when using this device. It is in your own best interest to pay attention to all warnings and safety information on the unit itself and in this manual.

Target audience

People working with the device must read the sections which are relevant to their work. This applies especially to the chapter entitled "Safety information".

Duty of observation with respect to product

The operator must report all operational irregularities or changes to components of relevance to safety to the responsible supervisor or the manufacturer immediately.

Storage location of operating instructions

We recommend storing these operating instructions near the device to ensure quick access by operating staff when information is required.

Legal provisions

National and local safety and accident prevention regulations which are in force must be strictly observed in addition to the operational guidelines issued by the operating entity.

No technical modifications whatsoever are to be made to this device under any circumstances!

Refer to specifications in "Technical Data" section for binding operational limits.

Do not use this unit near water or in any space with extreme humidity. In order to prevent electric shock, do not expose this appliance to rain or high humidity.

Never spill liquid on the unit or insert foreign objects into the unit. This could result in electric shock or damage to the unit.

Do not place flammable materials on or near the unit at any time. Keep the unit away from any sources of heat. The device has not been approved for operation in areas subject to explosion hazards!

The device may only be operated using the mains voltage indicated.

Never open the device or any components unless instructed expressly to do so somewhere in these instructions.

Never look directly into the LEDs or into the end of the light guide as this may result in injury or damage of your eyes.

Cleaning

Before cleaning, disconnect the product from the electrical supply and let it cool off for at least 15 minutes. For cleaning, use only a cloth dampened with (soap) water or 5% ethanol or isopropylene alcohol, but never use any combustible or flammable liquids.

Allow the device to dry sufficiently before startup. If fluids accidentally enter the device, unplug the power supply cable immediately and take the unit to the service location for repair. Carefully remove any dust that may have collected at the vent openings.

The use of high concentrations or an excessive cleaning can lead to scratches or discoloration of the panel!

Replacement parts

Use only original spare parts. If this is not done, it can lead to personal injury and material damage. Refer to "Spare Parts" section for item numbers. If necessary contact your local dealer.

Responsibility

As the ultimate legal entity, the operating institution is responsible for ensuring the proper use of the device and for providing other operators with the necessary information and it specifies the competencies required to operate the device.

Disconnect the mains cable when the unit is not being used for an extended period of time. Only use original mains cables. Route the cable so that it cannot be pinched or severed. For reasons of safety only use the grounded 3-pin mains cable.

Repairs which are not described in this document must only be carried out by authorised workshops. The manufacturer is not liable for any damage resulting from a failure to comply with the above instructions!

Technological state

Photonic LED illumination devices constitute state-of-the-art technology and employ recognized safety standards.

This instrument was designed and built in accordance with the following regulations and standards:

▼ EEC/89/336: Electromagnetic compatibility (CE)

INSTALLATION AND CONNECTION

Remove the individual system components from their packaging and place them onto a sufficiently large horizontal surface.

The device may only be operated at the rated mains voltage. Only connect the unit to grounded sockets.

Keep air vents free!

The cold light source unit must only be operated from the power mains using the power supply provided!

Using the hollow pin plug, insert the power pack into the socket of the cold light lamp (6), providing the desktop power pack with a supply voltage. If necessary, use the ESD socket (5) for potential equalisation.


The connection from the ESD socket to the light guide (via the light guide connection) is of low impedance. In order to restrict the flow of current in the event of a malfunction, protecting the user and the equipment, it may be necessary to make provision for an additional protective resistor in the ESD connection cable.

Insert the requisite fibre optic product, e.g. ring light, gooseneck or transmitted light illuminator onto the light guide connection (4).

Never look at the lighting devices or into the end of the light guide as this may result in damage to your eyes.

Installation of the optional counter weight

When working with goosenecks lightguides an counter weight may be necessary. Loosen the 2 screws of the cover, insert the counter weight and fix it with the 2 screws again. Store the cover on a safe place for later user.

OPERATION AND CONTROLS

The LED cold light source can be switched on with (10).

Following a short self test, indicated by the two indicator LEDs (1 and 2) blinking, the cold light source is ready for operation. This state is indicated by the green LED (1) remaining illuminated.

Brightness control

The brightness can be varied by turning the button (3). Please note that increasing the brightness reduces service life. Therefore, only ever work with the brightness that you really actually need.

The current brightness setting can be read on a scale. This setting is stored when the device is switched off using the button and is available when the device is switched on the next time.

Temperature monitoring

The red LED (2) indicates if the operating temperature is too high as a result of environmental influences or a less than ideal placement of the device, resulting in a high ambient temperature. The stages are as follows:

The operating temperature approaches the permissible limit:

The red LED (2) on the front panel begins to pulsate.

If it is possible to lower the ambient temperature by taking suitable action, the LED goes off.

If this indicator is overlooked by the user, the LED starts to blink slowly and the brightness of the cold light source is reduced by approximately 30%. This enables the user to keep doing their work while the device is cooling down. If the specified operating temperature can be reached in this way, the LED control lamp (2) goes off and the LED cold light source returns to full power.

If both indicators are overlooked and/or cooling down is insufficient, the device is deactivated until the permissible operating temperature is reached. The red LED control lamp (2) is illuminated.

Fault message

If one or both LED control lamps (1 and 2) flash rapidly during operation the device is defective. In this situation, disconnect the device from the mains, let it cool down and restart it.

If the problem still persists and one or both LEDs start blinking rapidly again, remove the device from service and contact your service department or your dealer.

Never open the device or its components.


Service life

LED chips typically have a service life in excess of 100,000 hours. The actual service life depends to a great extent upon the ambient temperature and the operating current of the LEDs. A higher operating current produces not only more light but also more heat, which reduces the life span of the LEDs accordingly.

As a rule, the brightness of LEDs also decreases over the life cycle.

To increase the service life of the LEDs, the brightness should only ever be set to the actual level required. In this way, service life can extend to over 50,000 hours.

Overview of LED indicator lamps

LED status		Operational state of the device	
LED state/lighting behaviour		 GREEN POWER	 RED TEMP.
Inactive		Device is disconnected from power source and therefore inactive.	
Enabled		Device enabled	Device is overheating and switches off until it has cooled down to the specified operating temperature.
Rapid blinking		Fault message, malfunction	Fault message, malfunction
Slow blinking		Blinking during self test at device start-up.	
		Insufficient or excessive voltage. (only use power pack supplied)	The light output of the device is reduced by approximately 30% in order for it to cool down to the permissible operating temperature (early warning level 2).
Pulsating lamp		Device is in standby mode	Operating temperature is approaching the permissible limit (early warning level).


The development of this product is undergoing continual refinement. Although we make every effort to keep these operating instructions as current as possible, we reserve the right to modify device specifications at any time without prior notification. The operating instructions are prepared subject to possible errors and misprints.

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Leica Microsystems (Schweiz) AG