

Living up to Life

*Leica*  
MICROSYSTEMS

# Leica DMi8

The modular system: Stands,  
components, accessories




JANUARY 2015

# Leica DMI8

## Modular system

Leica DMI8 – Introduction .....	3
Dimensions and Technical Data .....	4
Stands and Optic carriers.....	6
Microscope Control-Elements and Microscope Function-Keys .....	7
Leica CTR Boxes .....	8
Tubes.....	9
Eyepieces .....	9
Focusing and framing graticules .....	10
Observation and Documentation Ports .....	11
Laser Safety .....	12
Focus and Objective Nosepiece .....	13
Objectives.....	14
Tube lens and Magnification changer.....	14
T-Houses and Infinity Port.....	15
T-House Accessories .....	15
Fluorescence Equipment .....	16
Light Sources, Lamp Housings, Supply Units for fluorescence.....	18
External Filter Wheel EFW.....	19
Transmitted Light Axis.....	20
LED Light Sources (for transmitted light) .....	20
Transmitted Light Axis Filters .....	21
Transmitted Light Filters.....	21
Condensers and Accessories.....	22
Integrated Phase Contrast (IPH) and Integrated Modulation Contrast (IMC) .....	26
Transmitted Light Polarization Contrast.....	28
Transmitted Light Differential Interference Contrast (DIC) .....	28
C-Mount Adapter .....	30
Stages and Specimen Holders .....	31
Accessories .....	34
Digital Image Documentation .....	35
Software .....	35
Micromanipulators .....	35
Anti-vibration.....	35
System overview Leica DMI8 Automated .....	36
System overview Leica DMI8 Coded/Automated .....	38



# Leica DMi8 – Introduction

The Leica DMi8 is the highly modular inverted research microscope. It is designed for all common microscope applications and techniques. All contrast methods such as brightfield, darkfield, phase contrast, differential interference contrast (DIC), fluorescence and modulation contrast are integral to the microscope and can be adapted or changed quickly and easily. Multiple illumination and imaging beam paths, as well as HCS optics, modular accessories and a comprehensive range of peripherals complement the Leica inverted research stands.

## Basic stand

The basic stand is the solid core of the microscope. It can be equipped with various focusing systems, objective turrets, stage mounts, and mounts for transmitted and incident light units. All Leica DMi8 are available with or without fluorescence. All Leica DMi8 can be upgraded with an infinity port.

To control or to see the status of the Leica DMi8 and its components, several control elements like touch screens or function keys are attachable.

Camera ports with up to 19mm FOV are standard on the Leica DMi8

The system diagrams (see page 36) provide a good overview of this modular system. Thanks to its modular design principle, you can modify and/or extend your system to suit your requirements at any time. Ergonomic considerations were given a high priority in all of the stand designs, such as the convenient positioning of important controls and the availability of ergomodules or ergonomically designed components.

The functional and rugged design of the stands ensure ease of use and image stability for a wide range of applications up to the highest magnifications.

Optics of the highest quality ensure brilliant images with high contrast and resolution for any conceivable application.

## Note:

A Leica DMi8 Microscope stand is defined as a combination of 11889xxx /11525xx article numbers such as:

- microscope body
- T-Houses (w/o Infinity Ports)
- optic carrier
- camera ports
- fluorescence module
- magnification changers
- transmitted light arms
- front modules
- etc.

11889xxx article numbers always come as an integral part of the aligned complete system. These articles are not possible to order as an individual component.



# Dimensions and Technical Data

## Compact Leica CTR electronics box

For indoor use only.  
Supply voltage: 100 – 240 VAC  
Frequency: 50 / 60 Hz  
Power consumption: max. 150 VA  
Fuses: 3.15 A, slow-blowing,  
Breaking capacity H, 250  
VAC  
Size: 5x20 mm  
Ambient temperature: 15° - 35°C  
Relative humidity: 90% up to 30°C,  
non-condensing  
Protection class: I  
Overvoltage category: II  
Pollution degree: 2

## Advanced/Advanced+ Leica CTR electronics box

For indoor use only.  
Supply voltage: 100 – 240 VAC  
Frequency: 50 / 60 Hz  
Power consumption: max. 290 VA  
Fuses: 6.3 A, slow-blowing,  
Breaking capacity H, 250  
VAC  
Size: 5x20 mm  
Ambient temperature: 15° - 35°C  
Relative humidity: 90% up to 30°C,  
non-condensing  
Protection class: I  
Overvoltage category: II  
Pollution degree: 2

## Leica DMi8 manual

For indoor use only.  
Supply voltage: 100 – 240 VAC  
Frequency: 50 / 60 Hz  
Power consumption: max. 55 VA  
Fuses: 1.6 A, slow-blowing,  
Breaking capacity H, 250  
VAC  
Size: 5x20 mm  
Ambient temperature: 15° - 35°C  
Relative humidity: 90% up to 30°C,  
non-condensing  
Protection class: I  
Overvoltage category: II  
Pollution degree: 2


## Leica DMi8 coded

For indoor use only.  
Supply voltage: 100 – 240 VAC  
Frequency: 50 / 60 Hz  
Power consumption: max. 55 VA  
Fuses: 1.6 A, slow-blowing,  
Breaking capacity H, 250  
VAC  
Size: 5x20 mm  
Ambient temperature: 15° - 35°C  
Relative humidity: 90% up to 30°C,  
non-condensing  
Protection class: I  
Overvoltage category: II  
Pollution degree: 2


## Leica DMi8 automated

For indoor use only.  
Supply voltage: 100 – 240 VAC (→ Leica CTR)  
Frequency: 50 / 60 Hz (→ Leica CTR)  
Power consumption: see Leica CTR  
Fuses: see Leica CTR  
Ambient temperature: 15° - 35°C  
Relative humidity: 90% up to 30°C,  
non-condensing  
Protection class: I (→ Leica CTR)  
Overvoltage category: II (→ Leica CTR)  
Pollution degree: 2 (→ Leica CTR)



**Leica DMI8 (with tiltable Transmitted Light arm)**



Weight: approx. 19 kg





**Leica CTR compact**



Weight: approx. 2 kg

**Leica CTR Advanced/Advanced+**



Weight: approx. 4 kg

# Stands and Optic carriers



Fig. 1: Stand Leica DMI8

## Stands

### Leica DMI8 stand motorized and automated functions for

motorized side ports (no downgrade to coded side ports)  
Upgradable from manual coded to fully automated.

- for research use only
- for regulated market (IVD, clinical use)

11889113  
11889013

### Leica DMI8 stand motorized and automated functions for

coded left side ports (no upgrade to motorized side ports)  
Upgradable from manual coded to fully automated.

- for research use only
- for regulated market (IVD, clinical use)

11889112  
11889012

### Leica DMI8 stand manual coded functions for

coded left side ports.

Including coded nosepiece and manual focus.

Including power supply and lamp-module.

Including 1x tube lens.

Upgradable with manual coded components.

- for research use only
- for regulated market (IVD, clinical use)

11889111  
11889011

### Leica DMI8 stand manual functions for

manual left side ports.

Including manual nosepiece and manual focus.

Including power supply and lamp-module.

Including 1x tube lens.

Including all control elements

Upgradable with manual components.

- for research use only
- for regulated market (IVD, clinical use)

11889110  
11889010

## Optic Carrier

Optic carrier

11888032

Optic carrier for IMC

11888033

Optic carrier with Bottom Port

11888034

Optic carrier for IMC with Bottom Port

11888035

When choosing an “Optic carrier for IMC” you need to select an additional Front Module (see page 26)

### C- Mount Adapter for Bottom Port (see page 30)

### Ergonomic height compensation plate

A height compensation plate was developed to raise the viewing height by 23 mm or to raise the side camera ports for oversize cameras or spinning disks, or to use a microscope with an inactive bottom port on workbenches without openings.

11525200




Fig. 2: Ergonomic height compensation plate

# Microscope Control-Elements and Microscope Function-Keys

## Control Elements, left hand side of the microscope

Control left: Light	11889066
Control left: Light and shutter	11889067
Control left: Light, shutter and Fluor/TL buttons	11889068
Control left: Light, shutter, Fluor/TL buttons and diaphragms	11889069

## Control Elements, right hand side of the microscope

Control right: empty	11889070
Control right: 4x function keys	11889071
Control left: 4x function keys an 3x focus keys	11889072

## Control Elements, front of the microscope

Control front: empty	11889064
Control front: 6x status LED	11889063
Control front: 12x status LED	11889061
Control front: 6x function keys for objectives	11889062
Control front: 12x function keys for objective & fluorescence-cubes	11889060

## Touch-Screen

On site microscope, high resolution 6" touch-screen, tiltable and intensity adjustable

### Leica STP4000

External, wired, high resolution 6" touch –screen, intensity adjustable without xyz-control

### Leica STP8000

External, wired, high resolution 6" touch –screen, intensity adjustable, with xyz-controls for focus and motor stages, with 11 programmable function keys

### Leica SmartMove

x/y/z-Ergo control panel for electronic focus and motor stage. With 4 programmable function keys.



Fig. 3: Control elements left hand side



Fig. 4: Control elements right hand side



Fig. 5: Touch-Screen



Fig. 6: Leica STP8000



Fig. 7: Leica SmartMove

# Leica CTR Boxes

Power consumption of Leica CTR Boxes (see page 4).



Fig. 8: Leica CTR compact

## Leica CTR compact

For all Leica DMI8 (11889012,11889013,11889112,11889119) without motorized stages, without high-speed and without objectives with motorized correction rings

11525206

## Leica CTR advanced

For all Leica DMI8 (11889012,11889013,11889112,11889119) with 1 Master-Module incl. 1x serial, 2x USB, 3 x I<sup>2</sup>C, upgradable with 6x CTR boards for motorized stages, for high-speed sequencer or for objectives with motorized correction rings

11525207



Fig. 9: Leica CTR advanced

## Leica CTR advanced without lamp module

For all non motorized Leica DMI8 (11889010,11889011, 11889110,11889111) with 1 Master-Module incl. 1x serial, 2x USB, 3 x I<sup>2</sup>C, upgradable with 6x CTR boards. For motorized stages, for high-speed sequencer or for objectives with motorized correction rings

11525208

## Leica CTR advanced +

For all Leica DMI8 (11889012,11889013,11889112,11889119) with 2 Master-Module incl. 1x serial, 4x USB, 8 x I<sup>2</sup>C, upgradable with 6x CTR boards for motorized stages, for high-speed sequencer or for objectives with motorized correction rings

11525209



Fig. 10: Leica CTR advanced+

## CTR board XY basic

CTR Board xy basic, 15 pins, to control 3-plate motor stages (rack and pinion)

11525210

## CTR board XY advanced

CTR Board xy advanced, 25 pins, to control Scanning stages

11525211

## CTR board MotCORR

CTR Board MotCORR to control 1 or 2 objectives with motorized correction rings

11525212

## CTR board high speed sequencer

CTR Board HS/Sequencer to control high speed experiments with sequencer. 4x TTL, 1x Analog, 2x Sys/Camera

11525213

## Cable adapter 15/25 pins

Using a 3-plate stage in combination with CTR Board xy advanced, 25 pins (11525211)

11505237



# Tubes

The following applies to all tubes: Field of view 25 mm, eyepiece diameter 30 mm, a interpupillary distance range 55–75 mm and a viewing angle of 45° (Fixed) or 30–45° (Ergo)

<b>Binocular fixed tube</b>	11889025
<b>Binocular ergonomic tube</b>	11889026
<b>Binocular ergonomic tube with Bertrand lens</b> (for observation of back focal plane)	11889031
<b>Trinocular ergonomic tube</b> Binocular observation tube with side camera port and variable light path, 100% visual / 0% camera, and 50% visual / 50% camera	11889030
<b>Trinocular ergonomic tube with Bertrand lens</b> Binocular observation tube with side camera port and variable light path, 100% visual / 0% camera, and 50% visual / 50% camera and Bertrand lens (for observation of back focal plane)	11889029
<b>Trinocular ergonomic tube</b> Binocular observation tube with side camera port and variable light path, 100% visual / 0% camera, and 0% visual / 100% camera	11889028
<b>Trinocular ergonomic tube with Bertrand lens</b> Binocular observation tube with side camera port and variable light path, 100% visual / 0% camera, and 0% visual 100% camera with Bertrand lens (for observation of back focal plane)	11880027

## C- Mount Adapter for Top Port (see page 30)

# Eyepieces

Eyepieces for eyeglass wearers are available with adjustable lenses (dioptric equalization from –6.8 to +4.2 or –6 to +5).

M eyepieces are designed to accommodate a variety of graticules. 10x eyepieces are standard. All eyepieces have removable or fold-down eyecups and can be used with or without eyeglasses.

## Eyepieces with FOV 20

- Eyepiece HC PLAN 10x/20 BR. 11507801
- Eyepiece HC PLAN 10x/20 BR.M 11507802

## Eyepiece with FOV 22

- Eyepiece HC PLAN S 10x/22 Br.M 11507807

## Eyepiece with FOV 25

- Eyepiece HC PLAN S 10x/25 Br.M 11507808

## Special eyepieces with high magnification

- Eyepiece HC PLAN 12.5x/16 BR.M 11506515
- Eyepiece 16x/14B, adjustable 10445301
- Distance ring for eyepieces 16x/14B 11506808




Fig. 11: Trinocular ergonomic tube with Bertrand lens




Fig. 12: Eyepiece HC PLAN 10x/20

# Focusing and framing graticules


- Graticule 10 mm = 100 parts, D = 26 mm 11 506 950
- Graticule 10 mm = 200 parts, D = 26 mm 11 506 951
- Crosshair graticule, D = 26 mm 11 506 953
- Crosshair graticule with graduation  
10 mm = 100 parts, D = 26 mm 11 506 952
- Graticule with grid 10 x 10 mm,  
0.1 mm graduation, D = 26 mm 11 506 954
- Graticule with grid 10 x 10 mm,  
1 mm graduation, D = 26 mm 11 506 955




Graticule 11 506 950




Graticule 11 506 951




Graticule 11 506 952



Graticule 11 506 953



Graticule 11 506 954



Graticule 11 506 955

# Observation and Documentation Ports

## Camera ports for Leica DMi8 (11889010,11889110)

Manual left camera side port 100/0	11889043
Manual left camera side port 80/20	11889042
Manual left camera side port 100/0 with interlock	11889087

## Camera ports for Leica DMi8 (11889011,11889111,11889012,1189112)

Coded left camera side port 100/0	11889045
Coded left camera side port 80/20	11889044
Coded left camera side port 100/0 with interlock	11889088

## Camera ports for Leica DMi8 (11889013,11889113)

### Motorized side port, left 11889046

You may select at least one and up to three different prisms from the following for this version:

Side port prism, 100% left	11888259
Side port prism, 80% left	11888262
Side port prism, 50% left	11888264
Side port prism, dichroic 620nm, left	11888260

### Motorized side port, left and interlock 11889077

Includes the prism 100% left.

You may select no or one additional prism from the following for this version:

Side port prism, 80% left	11888262
Side port prism, 50% left	11888264

### Motorized side port, right 11889047

You may select at least one and up to three different prisms from the following for this version:

Side port prism, 100% right	11888258
Side port prism, 80% right	11888261
Side port prism, 50% right	11888263

### Motorized side port, right and interlock 11889078

Includes the prism 100% right. You may select no or one additional prism from

Side port prism, 80% right	11888261
Side port prism, 50% right	11888263



Fig. 13: Camera Ports

<b>Motorized side port, right and left</b>	11889048
You may select at least one and up to three different prisms from the following for this version: You have to select at least one right and one left prism!	
Side port prism, 100% right	11888258
Side port prism, 100% left	11888259
Side port prism, 80% right	11888261
Side port prism, 80% left	11888262
Side port prism, 50% right	11888263
Side port prism, 50% left	11888264
Side port prism, dichroic 620nm, left	11888260
 <b>Motorized side port, right and left and interlock</b>	 11889079
Includes the prism 100% left Includes the prism 100% right No other choice for more prisms	
 <b>No side port</b>	
This version requires the port compensation module	11888256
 <b>C- Mount Adapter for side port (see page 30)</b>	

## Laser Safety

<b>Laser safety Kit for DMi8. All interlocks</b>	
assembled in the microscope at factory site	11889086
 <b>Laser safety Hood for DMi8 with interlock</b>	
to use in combination with a tiltable TL-arm and Condensers with free working distance not higher than 28 mm	11889065
 <b>UV protection Hood for DMi8 without interlock</b>	
to use in combination with a tiltable TL-arm and Condensers with free working distance not higher than 28 mm	11522087

# Focus and Objective Nosepiece

## Manual Focus

Manual focus system with 12 mm travel range.

Tactile Coarse / Fine manual Focus Drive with coarse and fine / knob on each side of the microscope

11889054

## Motorized Focus

Motorized focus system with 12 mm travel range.

Tactile Coarse / Fine motorized Focus Drive with 4 gears and 5 sensitivity levels (0.05  $\mu$ ; 0.1  $\mu$ ; 0.7  $\mu$ ; 1.5  $\mu$ ; 5.0  $\mu$ ). Electronic focus repositioning and electronic parfocality, with coarse / fine knob on each side of microscope

11889055

## Motorized Closed Loop Focus

Motorized closed loop focus system with encoded 12 mm travel range.

Reproducibility < 20 nm bidirectional. Tactile Coarse / Fine motorized Focus Drive with 4 gears and 5 sensitivity levels (0.05  $\mu$ ; 0.1  $\mu$ ; 0.7  $\mu$ ; 1.5  $\mu$ ; 5.0  $\mu$ )

Electronic focus repositioning and electronic parfocality, with coarse / fine knob on each side of microscope

11889056

## Motorized Closed Loop Focus with AFC

Adaptive Focus Control (AFC) actively keeps the focus position over time.

In combination with Closed Loop Focus system with encoded 12 mm travel range. Reproducibility < 20 nm bidirectional. Tactile Coarse / Fine motorized Focus Drive with 4 gears and 5 sensitivity levels (0.05  $\mu$ ; 0.1  $\mu$ ; 0.7  $\mu$ ; 1.5  $\mu$ ; 5.0  $\mu$ )

Electronic focus repositioning and electronic parfocality, with coarse / fine knob on each side of microscope.

11889073

## Motorized Focus with AFC

Adaptive Focus Control (AFC) actively keeps the focus position over time.

Tactile Coarse / Fine motorized Focus Drive with 4 gears and 5 sensitivity levels (0.05  $\mu$ ; 0.1  $\mu$ ; 0.7  $\mu$ ; 1.5  $\mu$ ; 5.0  $\mu$ )

Electronic focus repositioning and electronic parfocality, with coarse / fine knob on each side of microscope

11889074

## Manual Coded Nosepiece

Coded 6-fold objective nosepiece

with 25 mm threads and 45 mm parfocal distance

11889049

## Motorized Nosepiece

Motorized coded 6-fold objective nosepiece

with 25 mm threads and 45 mm parfocal distance

11889050

## Spill Protection

For motorized nosepiece the spill protection seals the revolver against leakage. All liquids will be collected in recipient.

11525215



Fig. 14: Objective Nosepiece




Fig. 15: Leica Adaptive Focus Control (AFC)



Fig. 16: Spill Protection

# Objectives

Based on the Leica principle of infinity distance correction of optics, the microscope objectives are infinity corrected for tubelens systems with 200 mm reference focal lengths. The calibration length is 45 mm. The objectives are divided into 4 correction classes:



Fig. 17: Objectives

Objective	Class	Field of view performance
Achromatic	HI PLAN	up to 20
Planachromatic	NPLAN	up to 22
Semi-apochromatic	HC/X PLFLUOTAR	up to 25
Apochromatic	HC/X PL APO	up to 25

When selecting the objectives, consider the intended use with regard to specimen covering, etc. For more detailed explanations, please refer to the appendix of the objective list.

<http://www.leica-microsystems.com/objectives>

## Tube lens and Magnification changer

**Fixed 1x tube lens** 11889023

**Motorized magnification changer**  
including 1x tube lens at position 1 11889024

**Magnification levels for motorized magnification changer**  
In addition to the 1x tube lens, one or two lenses must be installed on the disk of motorized magnification changer.

1.6x tube lens for motorized magnification changer 11888377  
2.0x tube lens for motorized magnification changer 11888376

**Manual Magnification Changer**  
for Eyepiece and Top Port. Find detailed information in the chapter Front-Modules (see page 26)

**Front Module for IMC/IPH with integrated manual 1.6 x Magnification changer (for manual stands, only)** 11889075

**Front Module for IMC/IPH with integrated coded 1.6 x Magnification changer** 11888076

# T-Houses and Infinity Port

Every Leica DMI8 has to be combined with a T-House. The T-Houses determine the contrasting techniques, coupling devices, and modalities of the microscope. The T-House always acts as rest for stages and transmitted light arms.

## Leica DMI8 T-House(0)

For Microscopes without fluorescence methods to hold all different transmitted light arms 11889036

## Leica DMI8 T-House(1) with one Infinity Port

Infinity Port for Microscopes with one external fluorescence module to the rear and to hold all different transmitted light arms 11889037

## Leica DMI8 T-House(2) with integrated manual Fluorescence Axis

with integrated manual fluorescence module including standard lamp-mount (e.g. for Hg lamp-house) to the rear, manual Shutter, Iris Field-Diaphragm, interface for optional manual Fluorescence Intensity Manager (FIM) and to hold all different transmitted light arms. 11889038

## Leica DMI8 T-House(3) with two Infinity Ports

Two infinity ports for Microscopes with 2 external in-coupling devices or fluorescence modules to the rear and left hand side and to hold all different transmitted light arms. 11889039

## Leica DMI8 T-House(4) as TIRF/GSD adapter

for microscopes with TIRF or GSD device and to hold all different transmitted light arms. 11889040



Fig. 18: T-House

# T-House Accessories

## Attenuator

Manual attenuator with 5 positions for T-House(2) with integrated manual Fluorescence Axis 11525205

## T-House Splitter 50/50

Manual 50/50 splitter for T-House(3) for use of both the infinity ports simultaneously 11525377

## T-House Splitter 100/0

Manual 100/0 splitter for T-House(3) for use of both the infinity ports consecutively and separately 11525388

# Fluorescence Equipment

Epi-fluorescence with the Leica DMi8 series features a completely new incident illumination system and comes with

- manual or motorized 6 fold filter turret
- filter turret with fast, magnetic filter exchange
- filter turret with access panels on both sides
- fluorescence axis: manual, coded or motorized, external or integrated
- manual or motorized diaphragms
- manual or motorized FIM (optional)
- manual or motorized diaphragms
- external or integrated light guide adapter

The optional fast internal filter wheel (IFW) supports switching between up to 3 different excitation wavelengths within 50 ms. This permits extremely fast switching between different emissions without moving filter blocks when using a dual or a triple-pass filter block.

The Excitation Manager is available together with the fast internal filter wheel (IFW). It supports the tuning of emissions over 15 steps, two excitation wavelengths, a variety of intensities and the mixing of various colors.

The filter systems (cubes) have magnets to guide them into the fluorescence turret without tools. All cubes are designed to guarantee zero pixel shift and completely suppress stray light to ensure a dark background for the fluorescent image.

Always recommended:

## **Fluor-Protection Shield**

for all the Leica DMi8 microscopes with Fluorescence Accessories for UV Protection. Mounted between stage and tube 11525114

## **Fluorescence Turrets**

### **Motorized 6-position fluorescence turret**

Operation via function-keys at the microscope, touch-screen or software depending on the configuration 11 889 022

### **Manual coded 6-position fluorescence turret**

Operated from both sides. The active filter position is indicated via LED at the front panel, touch Screen or software depending on the configuration 11 889 021

### **Manual 6-position fluorescence turret**

only for Leica DMi8 11889010, 11889110  
Manual 6-fold fluorescence turret. Operated from both sides. A color code indicates the active position 11889020




Fig. 19: Fluorescence Turret



## External Fluorescence Axis

The interface for external fluorescence axis are the Infinity Ports of T-House(1) and T-House(3)

### Motorized external fluorescence-axis

The module includes

- internal filter wheel (IFW)
- light fiber connection
- motorized diaphragms
- slot for Leica fast external filter wheels

11889053

### Motorized external fluorescence-axis

The module includes

- light fiber connection
- motorized diaphragms
- slot for Leica fast external filter wheels

11889052

### Manual external fluorescence-axis

The module includes

- light fiber connection
- access to field and aperture plane
- shutter and diaphragm
- interface for structured illumination (OptiGrid®)
- slot for fixed apertures and intensity attenuator

11889051

### Aperture Slider

Manual aperture slider for manual external fluorescence axis with 6 openings (apertures)

11525203

## Internal Fluorescence Axis

The T-House (2) comes with the integrated internal fluorescence axis.

### Leica DMi8 T-House(2) and integrated manual Fluorescence Axis

with integrated manual fluorescence module including standard lamp-mount (e.g. for Hg lamp-house) to the rear, manual shutter, iris field-diaphragm, interface for optional manual Fluorescence Intensity Manager (FIM) and to hold all different transmitted light arms.

11889038




Fig. 20: Motorized external fluorescence axis



Fig. 21: Manual external fluorescence axis




Fig. 22: Manual internal fluorescence axis

# Light Sources, Lamp Housings, Supply Units for fluorescence

## External Light Sources



Fig. 23: Leica EL6000

### Leica EL6000

- External light source 11504115
- Liquid light guide, 2 m 11504116
- Cable for shutter control Leica EL6000 I<sup>2</sup>C 11500336
- 1-inch fiber optics adapter for lamp mounts  
necessary when using T-House(2) 11504117

### Chroma Photofluor II

- External Light source 11504183
- Liquid light guide 11504184
- 1-inch fiber optics adapter for lamp mounts  
necessary when using T-House(2) 11504117



Fig. 24: Lamp housing 106z – Hg 100 W (6 lens)

### Lamp housing 106z – Hg 100 W (6-lens)

(In combination with T-House(2) only), with centerable lamp mount for Hg 100 W lamp, with 1.5 m power cable, without burner.

6-lens, achromatic 1-inch collector, with heat-absorbing filter

- left hand operation 11504106
- right hand operation 11504114

### Lamp housing 106z – Xe 75 W (6-lens)

(In combination with T-House(2) only) with centerable lamp mount for Xe 75 W lamp, with 1.5 m power cable, without burner. 6-lens, achromatic 1-inch collector.

With heat-absorbing filter, face protection, protective gloves 11504105

## Lamps and burners

- High-pressure mercury burner Hg 100 W/2 11500321
- Lamp HXP R120/45C-Vis for Leica EL6000 11504120
- High-pressure xenon burner Xe 75 W 11500139

## Supply units

### Supply unit Hg 100 W

With power supply cord, automatic switching to power supply voltage 90 V–250 V 50/60 Hz with operating hours display. 11500334

### Supply unit Xe 75 W

With power supply cord, automatic switching to power supply voltage 90 V–250 V 50/60 Hz with operating hours display. 11500335



Fig. 25: Leica SFL4000

### LED light sources (not available in the US)

- Leica SFL100, 365nm 11504196
- Leica SFL100, 470nm 11504138
- Leica SFL100, 530nm 11504195
- Leica SFL4000 11504139

# External Filter Wheel EFW

## Fast external Filter wheel slider

to insert in an external fluorescence axis of Leica DMi8.  
The slider comes without filter set.

11525361

## Cable connection

for a fast external filter wheel.  
For each EFW a connecting cable is necessary.

Cable EXT - I<sup>2</sup>C, 200 cm

11504132

Cable EXT - I<sup>2</sup>C, 50 cm

11525217

## External Holder for 2 EFW slider

External device to hold one or two external filter wheels (EFW)

11640266

## Lightguide coupler.

1-inch fiber optics adapter for lamp mounts necessary when using  
T-House(2) including light guide coupler for one Leica DMi8 EFW

11504121

## C-Mount - 0.7x HC for EFW. This c-mount can hold one

Leica DMi8 EFW slider for emission application

11541545

## C-Mount – 1,3x HC for EFW. This c-mount can hold one

Leica DMi8 EFW slider for emission application

11541547

## EFW Filter Sets

for the sliders

- Triple-Band Set: Filter cube Set D/F/TX-P 11525347
- Fura2-Set: Filter cube Set FURA 2 11525348
- CFP/YFP/mCherry Filter cube Set C/Y mC-P 11525349
- FRET-Sedat (2 Wheels): Filter cube set C/Y FRET 11525350
- Quad-Sedat (2 Wheels): Filter cube Set QUAD-S 11525353
- Neutral Density (ND) set  
– External Fluorescence Intensity Manager 11504129



Fig. 26: Leica External Filterwheel (EFW)

# Transmitted Light Axis

The transmitted light illumination unit essentially consists of an illumination source and condenser carrier. Excellent light utilization is ensured by the LED Lamp in replaceable Leica lamp housings. Optimal and homogeneous illumination is a prerequisite for all transmitted light contrast methods.



Fig. 27: Fixed basic transmitted light arm

## Fixed basic transmitted light arm

Basic fixed transmitted light arm including LED illumination, condenser base and aperture diaphragm. To hold S50/0.50 and S80/0.30 condenser lenses.

With:

- 70 cm cable (11889010,1189110,11889011,11889111) 11525105
- 185 cm cable (11889012,1189112,11889013,11889113) 11525116

## Tilting manual transmitted light axis

With integrated tilting mechanism for specimen clearance, integrated manual filter magazine for 2 replaceable filter positions (one preconfigured with manual shutter), Condenser Quick-Changer for all manual condensers, lamp housing adapter for LED, integrated duct for the lamp housing cable:

- without Field Diaphragm 11525106
- with manual Field Diaphragm 11525107
- with manual Field Diaphragm for Atomic Force App. 11525108



Fig. 28: Tilting transmitted light arm

## Tilting manual coded transmitted light axis

with integrated tilting mechanism for specimen clearance, integrated manual filter magazine for 2 replaceable filter positions (one preconfigured with manual shutter), Condenser Quick-Changer for all coded/motorized condensers, lamp housing adapter for LED, integrated duct for the lamp housing cable

- without Field Diaphragm 11525109
- with manual Field Diaphragm 11525110

## Tilting motorized transmitted light axis

with integrated tilting mechanism for specimen clearance, integrated manual filter magazine for 2 replaceable filter positions (one preconfigured with manual shutter), Condenser Quick-Changer for all coded/motorized condensers, lamp housing adapter for LED, integrated duct for the lamp housing cable

- with motorized Field Diaphragm 11525111



Fig. 29: LED lamp house

# LED Light Sources (for transmitted light)

## LED lamp house

- 70 cm cable 11525100
- 185 cm cable 11525102

## LED lamp house with TTL shutter

- 70 cm cable 11525101
- 185 cm cable 11525103
- 285 cm cable 11525104

# Transmitted Light Axis Filters

## Light filter Ø 40 mm, unframed

for tilting light axis.

2 transmitted light filters can be swung in on the tiltable transmitted light illumination arms. A broad selection of filters is available to optimize illumination for observation and documentation. All illumination arms come with a light stop pre-installed in one position. Users can replace it with one of the following filters, however.

- DLF, daylight filter 11521577
- Panchromatic green filter 11521582
- Neutral filter N 16 (6.3%) 11521579
- Neutral filter N 4 (25%) 11521580
- Neutral filter N 2 (50%) 11521581
- Green filter VG 9, narrow band filter 11521583

# Transmitted Light Filters

## Ø 32 mm in holder with handle

- Panchromatic green filter for black/white photography 11512077
- VG 9, green filter for contrast enhancement (B/W) 11563122
- IL 546 nm (Polarization microscopy, interferometry) 11563155
- Neutral filter N2 (50%) 11543092
- Neutral filter N4 (25%) 11543093
- Neutral filter N16 (6.3%) 11543184
- Polarizer 11505087
- Lambda plate 11513908
- Quarter lambda plate 11513570
- Polarizer with protective filter 11513711

## Ø 32 mm without holder

to be placed on the condenser.

- DLF, daylight filter 11504046
- ALF, artificial light filter 11504047
- Panchromatic green filter 11504011
- Green filter VG 9 11504004
- Neutral gray filter N16 6.3% transmission 11504005
- Neutral gray filter N4 25% transmission 11504006
- Neutral gray filter N2 50% transmission 11504007
- Diffusion filter N 11504012
- Interference green filter VSS 546 11504010
- Blue glass filter BG 20 11504009



Fig. 30: Filter with handle




Fig. 31: Transmitted light filter

# Condensers and Accessories

A variety of condensers are available for the Leica DMI8

Three different types of condenser systems are available.

## Type 1:

This condenser system features a 7-position condenser disk to accommodate light rings, IC prisms or slit diaphragms (3 positions with diameter 23 mm and 4 positions with diameter 32 mm) and mount for a polarizer.

Incl. height adjustment and field diaphragm adjustment, if suitable.

Incl. adjustment tools and tool-box.

Available as motorized and manual coded version.

## Type 2:

This condenser system features a 5-position condenser disk to accommodate light rings, IC prisms or slit diaphragms (2 positions with diameter 23 mm and 3 positions with diameter 43 mm), and mount for a polarizer.

Available manual version

## Type 3:

This condenser system features a 3-position condenser slider to accommodate light rings, or slit diaphragms condenser Slider for IMC.

Available manual version

Type 1 and 2 are used in combination with the tilting transmitted light axis.

Type 3 is used in combination with the fixed transmitted light axis.

All the condenser bases can be individually equipped with different condenser lenses.

The condensers

- can be used for Brightfield, Phase-Contrast, Polarization. There are exceptions for IMC, IPH and DIC
- Additional features
  - aperture diaphragm
  - mount for interchangeable condenser heads

## S70 Condenser System (Type 1)

The motorized S70/0.30 condenser features a free working distance of 70 mm, a numerical aperture of 0.30 and is particularly suitable for specimens in high-volume containers. Including motorized aperture diaphragm.

Possibility to mount motorized polarizer.

Contrasting Techniques: BF, PH, Pol, DIC

Magnifications: 4x up to 40x

### Motorized S70/0.30 with fixed condenser head

- Motorized condenser base incl. head S70/0.30 11525370
- Light ring set for phase contrast S70/0.30 11521506
- Motorized polarizer holder including polarizer 11522120



Fig. 32: Condenser S70/0.30

## S23 or S28 Condenser System (Type1)

The manual/coded or motorized S28/0.55 or S23/0.53 condensers feature a free working distance of 28 mm or 23 mm, a numerical aperture of 0.55 or 0.53 and are particularly suitable for thicker specimens (living cells) for highest resolution. This condenser is also suitable for dark field contrast at objective apertures of up to 0.40. Suitable inserts are available for this condenser for integrated modulation contrast / integrated phase contrast.

Including motorized or manual aperture diaphragm.

Possibility to mount motorized or manual polarizer

Contrasting Techniques: BF, PH, DF, Pol, DIC, IMC, IPH

Magnifications: 4x up to 100x.



Fig. 33: S28 condenser

### Manual S28/0.55 with fixed condenser head

- Coded manual condenser base with fixed mount 11525373
- Condenser head S28/0.55 11505234
- Light ring set for PH and DF contrast S23–S28 11521505
- Manual polarizer holder including polarizer 11522103

### Manual S28/0.55 with manually movable condenser head

- Coded manual condenser base with movable mount 11525374
- Condenser head S28/0.55 11505234
- Light ring set for PH and DF contrast S23–S28 11521505
- Manual polarizer holder including polarizer 11522103

### Motorized S28/0.55 with fixed condenser head

- Motorized condenser base with fixed mount 11525372
- Condenser head S28/0.55 11505234
- Light ring set for PH and DF contrast S23–S28 11521505
- Motorized polarizer holder including polarizer 11522120

### Motorized S28/0.55 with motorized condenser head

- Motorized condenser base with movable mount 11525371
- Condenser head S28/0.55 11505234
- Light ring set for PH and DF contrast S23–S28 11521505
- Motorized polarizer holder including polarizer 11522120

### Manual S23/0.53 with fixed condenser head

- Coded manual condenser base with fixed mount 11525373
- Condenser head S23/0.53 11521500
- Light ring set for PH and DF contrast S23–S28 11521505
- Manual polarizer holder including polarizer 11522103

### Manual S23/0.53 with manually movable condenser head

- Coded manual condenser base with movable mount 11525374
- Condenser head S23/0.53 11521500
- Light ring set for PH and DF contrast S23–S28 11521505
- Manual polarizer holder including polarizer 11522103

### Motorized S23/0.53 with fixed condenser head

- Motorized condenser base with fixed mount 11525372
- Condenser head S23/0.53 11521500
- Light ring set for PH and DF contrast S23–S28 11521505
- Motorized polarizer holder including polarizer 11522120



Fig. 34: S28 condenser head

### Motorized S23/0.53 with motorized condenser head

- Motorized condenser base with movable mount 11525371
- Condenser head S23/0.53 11521500
- Light ring set for PH and DF contrast S23–S28 11521505
- Motorized polarizer holder including polarizer 11522120



Fig. 35: Condenser base for S1–S28

### S1 Condenser System (Type1)

The two manual/coded or motorized S1 condensers feature a free working distance of 1 mm, a numerical aperture of 0.90 dry or 1.40 oil and are particularly suitable for highest resolution and magnifications of up to 100x.

This condenser is also suitable for dark field contrast at objective apertures of up to 0.70. Including motorized or manual aperture diaphragm.

Possibility to mount motorized or manual polarizer.

Contrasting Techniques: BF, PH, DF, Pol, DIC

Magnifications: 1.25 up to 100x

### Manual S1/0.90 or S1/1.40 with fixed condenser head

- Coded manual condenser base with fixed mount 11525373
  - Spacer ring for condenser head S1 11521502
  - Condenser head S1/0.90 dry 11505150
- or
- Condenser head S1/1.40 oil (no phase contrast) 11551004
  - Light ring set for PH and DF contrast S1 11521504
  - Manual polarizer holder including polarizer 11522103

### Manual S1/0.90 or S1/1.40 with movable condenser head

- Coded manual condenser base with movable mount 11525374
  - Spacer ring for condenser head S1 11521502
  - Condenser head S1/0.90 dry 11505150
- or
- Condenser head S1/1.40 oil (no phase contrast) 11551004
  - Light ring set for PH and DF contrast S1 11521504
  - Manual polarizer holder including polarizer 11522103



Fig. 36: Light rings

### Motorized S1/0.90 or S1/1.40 with fixed condenser head

- Motorized condenser base with fixed mount 11525372
  - Spacer ring for condenser head S1 11521502
  - Condenser head S1/0.90 dry 11505150
- or
- Condenser head S1/1.40 oil (no phase contrast) 11551004
  - Light ring set for PH and DF contrast S1 11521504
  - Motorized polarizer holder including polarizer 11522120

### Motorized S1/0.90 or S1/1.40 with motorized condenser head

- Motorized condenser base with movable mount 11525371
  - Spacer ring for condenser head S1 11521502
  - Condenser head S1/0.90 dry 11505150
- or
- Condenser head S1/1.40 oil (no phase contrast) 11551004
  - Light ring set for PH and DF contrast S1 11521504
  - Motorized polarizer holder including polarizer 11522120



## S70 Condenser System (Type 2)

The manual S70/0.30 condenser features a free working distance of 70 mm, a numerical aperture of 0.30 and is particularly suitable for specimens in high-volume containers. Including manual aperture diaphragm

Possibility to mount manual polarizer

Contrasting Techniques: BF, PH, Pol, DIC

Magnifications: 4x up to 40x

### Manual S70/0.30 with fixed condenser head

- Manual condenser base incl. head S70/0.30 11525089
- Light ring set for phase contrast S70/0.30 11522090
- Manual polarizer in round insert 11522094
- 32/43 mm adapter for Wollaston prisms 11522093
- Spare disc for manual S40 or S70 condenser 11522092

## S40 Condenser System (Type 2)

The manual S40/0.50 condenser features a free working distance of 40 mm, a numerical aperture of 0.50 and is particularly suitable for specimens in high-volume containers. Including manual aperture diaphragm

Possibility to mount manual polarizer

Contrasting Techniques: BF, PH, Pol, IMC, IPH

Magnifications: 4x up to 63x

### Manual S40/0.50 with fixed condenser head

- Manual condenser base incl. head S40/0.50 11525088
- Light ring set for phase contrast S70/0.30 11522091
- Manual polarizer in round insert 11522094
- Spare disc for manual S40 or S70 condenser 11522092

## S80 Condenser System (Type 3)

For fixed transmitted light arm. The manual S80/0.30 condenser features a free working distance of 80 mm, a numerical aperture of 0.30 and is particularly suitable for specimens in high-volume containers. A 4-position slider accommodates light rings. Including manual aperture diaphragm.

### S80/0.30 fixed condenser head

- Manual condenser lens S80/0.30 11521251
- 4 Position Slider for BF PH0,PH1 and PH2 11521255

## S40 Condenser System (Type 3)

For fixed transmitted light arm. The manual S40/0.45 condenser features a free working distance of 40 mm, a numerical aperture of 0.45 and is particularly suitable for specimens in high-volume containers. A 4-position slider accommodates light rings. Including manual aperture diaphragm.

### S40/0.45 fixed condenser head

- Manual condenser lens S40/0.45 11521252
- 4 Position Slider for BF PH0,PH1 and PH2 11521253



Fig. 37: Condenser S40/0.50



Fig. 38: Condenser heads and spacer ring



Fig. 39: Condensers S40/S80

# Integrated Phase Contrast (IPH) and Integrated Modulation Contrast (IMC)



Fig. 40: IMC Slider



Fig. 41: Condenser Modulators



Fig. 42: IPH Slider

Not all stand versions are ready for IMC. Prerequisites are:

- Optic Carrier with IMC optics
- Front Module
- S23, S28 or S40 Condenser in combination with tiltable TL arm
- S40 or S80 Condenser in combination with fixed TL arm (no IPH)

## Front-Modules

The Front Module is an interface in the microscope stand for the modulation contrast or phase contrast slider. The 2 front modules are identical, but one features an additional manual magnification changer.

The manual magnification changer affects the top camera port and is not available in combination with the motorized magnification changer.

- Front Module for IMC/IPH 11888384
- Front Module for IMC/IPH with integrated manual 1,6 x Magnification changer (for manual stands, only) 11889075
- Front Module for IMC/IPH with integrated coded 1,6 x Magnification changer 11889076

## Integrated phase contrast

(only in combination with tiltable TL-arm)

Phase rings must be placed in the phase contrast slider for the different back focal planes of the objectives. For the back focal planes of the objectives, please refer to the objective list, see <http://www.leica-microsystems.com/objectives>

- IPH Module for Phase rings 11522065
- Phase Contrast Ring B 11522080
- Phase Contrast Ring C 11522064
- Light ring (set) for S23/S28 11521505
- Light ring (set) for S40 (5 – 63x incl. 32x) 11522091

## Integrated modulation contrast

(in combination with fixed TL-arm)

Modulators have been placed in the modulation contrast slider for the different back focal planes of the objectives. For the back focal planes of the objectives, please refer to the objective list <http://www.leica-microsystems.com/objectives>.

The S40 condensers requires:

- IMC Condenser Slider S40/0.45 11525119
- IMC Condenser Slider S80/0.30 11525120

- Slider for modulation Contrast for objectives with back focal planes B, C or D 11522152

For IMC with in combination with fixed TL-arm no Polarizer is required

**Integrated modulation contrast**

(in combination with tiltable TL-arm)

Modulators have been placed in the modulation contrast slider for the different back focal planes of the objectives. For the back focal plane of the objectives, please refer to the objective list <http://www.leica-microsystems.com/objectives>.

A bright-field share can be added to the modulation contrast via the polarizer, reducing the relief contrast.

- Slider for modulation Contrast for objectives with back focal planes B and C 11522075

The S23 and S28 condensers require inserts:

- IMC S23/28 Condenser Modulators Complete Set 11522074
- IMC S23/28 Condenser Modulator for 5x (Ø 23 mm) 11522108
- IMC S23/28 Condenser Modulator for 10x (Ø 23 mm) 11522109
- IMC S23/28 Condenser Modulator for 20x (Ø 32 mm) 11522110
- IMC S23/28 Condenser Modulator for 32x (Ø 32 mm) 11522127
- IMC S23/28 Condenser Modulator for 40x (Ø 32 mm) 11522111
- IMC S23/28 Condenser Modulator for 63x (Ø 32 mm) 11522112

The S40 condensers require inserts:

- IMC S40 Condenser Modulators Complete Set 11522113
- IMC S40 Condenser Modulator for 5x (Ø 23 mm) 11522099
- IMC S40 Condenser Modulator for 10x (Ø 23 mm) 11522095
- IMC S40 Condenser Modulator for 20x (Ø 43 mm) 11522096
- IMC S40 Condenser Modulator for 32x (Ø 43 mm) 11522126
- IMC S40 Condenser Modulator for 40x (Ø 43 mm) 11522097
- IMC S40 Condenser Modulator for 63x (Ø 43 mm) 11522098

- Adapter 32/43 mm for Wollaston Prisms 11522093

A polarizer is required for IMC

**For motorized condenser (Type1)**

- Motorized polarizer holder including polarizer for S23/28 condenser 11522120

**For manual condenser (Type 1)**

- Manual polarizer holder including polarizer for S23/28 condenser 11522103

**For manual condenser (Type 2)**

- Manual polarizer in round insert for S40 condenser 11522094

**Focusing telescope**

A focusing telescope is required when adjusting phase contrast, modulation contrast or differential interference contrast in order to view the rear focal plane of the objective.

The following tubes:

- Binocular ergonomic tube with Bertrand lens
- Trinocular ergonomic tube with Bertrand lens are already equipped with a focusing telescope.

- Focusing telescope 11 505 070



Fig. 43: Focusing telescope

# Transmitted Light Polarization Contrast



Fig. 44: Polarizer and Analyzer

A revolving polarizer and an analyzer are required for transmitted light polarization contrast. Low-strain objectives (marked with a P in the objective list) enhance the quality of the polarization contrast.

A polarizer is required for Polarization Contrast

### For motorized condenser (Type1)

- Motorized polarizer holder including polarizer for S23/28 condenser 11522120

### For manual condenser (Type 1)

- Manual polarizer holder including polarizer for S23/28 condenser 11522103

### For manual condenser (Type 2)

- Manual polarizer in round insert for S40 condenser 11522094

### Analyzer ICT/P

The analyzer is located below the objective turret. To insert the analyzer, remove the cover for the unused objective DIC disk opening. (Component of stand 11-020-437-101-013)

- On slider 30 x 5 mm, fixed orientation 11522046

Alternatively for stands with fluorescence axis:

### Analyzer Block

This analyzer is built into an empty fluorescence filter block and can be rotated into the beam path using a motorized fluorescence disk.

The swing direction is fixed at 0° (east-west). 11525300

# Transmitted Light Differential Interference Contrast (DIC)

For transmitted light differential interference contrast a revolving polarizer, a fixed analyzer and a set of Wollaston prisms are required. The Leica DMI8 features both manual and motorized polarizer, analyzer, objective and condenser prism functions. The illumination-side IC prisms are inserted in the condenser disk and objective-side IC prisms in the objective prism disk.

The Leica DMI8 realizes DIC at working distances of up to 70 mm (S70 condenser). Objectives with magnifications from 5x to 100x can be used for DIC depending to the condenser.

Fast (automatic) switching between DIC, brightfield and phase contrast is possible at all times without the need for DIC prisms to remain in the beam path. For valid combinations of prisms, condensers and objectives, please refer to the objective list. <http://www.leica-microsystems.com/objectives> Coded objective turrets recognize the installed objectives. The motorized objective prism disk selects the correct objective prism and sets the bias. The analyzer is automatically positioned in the beam path by the fluorescence disk. In addition, the luminous intensity, aperture diaphragm – and in the case of a motorized transmitted light axis, the field diaphragm – are automatically set to the required values. The condenser head automatically swings in and out as necessary. The user can adjust and overwrite the pre-set values at any time. Note: The S40 Condenser is not designed for DIC.

**DIC Turret**

- Manual DIC objective system prism disk, 4 positions 11522123
  - Manual coded DIC objective system prism disk, 4 positions 11525117
  - Motorized coded DIC objective system prism, disk 4 positions 11525118
- A polarizer is required for Differential Interference Contrast

**For motorized condenser (Type1)**

- Motorized polarizer holder including polarizer for S23/28 condenser 11522120

**For manual condenser (Type 1)**

- Manual polarizer holder including polarizer for S23/28 condenser 11522103

**For manual condenser (Type 2)**

- Manual polarizer in round insert for S40 condenser 11522094

**Analyzer ICT/P**

The analyzer is located below the objective turret. To insert the analyzer, remove the cover for the unused objective DIC disk opening. (Component of stand 11-020-437-101-013)

- On slider 30 x 5 mm, fixed orientation 11522046

Alternatively for stands with fluorescence axis:

**Analyzer Cube**

This analyzer is built into an empty fluorescence filter cube and can be rotated into the beam path using a motorized fluorescence disk. The swing direction is fixed at 0° (east-west). 11525300

**Wollaston prisms**

**Objective prisms**

- IC objective prism B1 11555007
- IC objective prism C 11555009
- IC objective prism C1 11522038
- IC objective prism C2 11522039
- IC objective prism D 11555010
- IC objective prism D1 11555056
- IC objective prism E 11555046

**Condenser prisms (Ø 23 mm)**

- IC condenser prism K3 11521594
- IC condenser prism K11 11521545

**Condenser prisms (Ø 32 mm)**

- IC condenser prism K2 11555016
- IC condenser prism K3 11555017
- IC condenser prism K4 11555018
- IC condenser prism K5 11555019
- IC condenser prism K6 11521521
- IC condenser prism K7 11521522
- IC condenser prism K8 11521523
- IC condenser prism K9 11555030
- IC condenser prism K10 11521524
- IC condenser prism K11 11521529
- IC condenser prism K16 11522037
- IC condenser prism K17 11555091

For the manual S70 condensers turrets with 43 mm openings a 32/43 adapter is necessary for each of the prisms (max 3)

- Reducing diameter adapter 32/43** 11522093



Fig. 45: Motorized DIC objective prism disk




Fig. 46: Manual DIC objective prism disk



Fig. 47: DIC condenser prisms



Fig. 48: Objective prism D

# C-Mount Adapter

You can adapt analog and digital cameras to all ports with documentation output. The C and B-mount adapters are aligned to the dimensions of the holder thread. The various fixed and variable magnification factors allow adjustment of the rendering of the microscopic image on the camera chip. In order to display the largest possible portion of the field of view on the monitor, the magnification factor of the adapter must fit the chip size of the camera. If the magnification is too low, there will be a lack of uniformity to the illuminated area (shading) and/or vignetting.

	Recorded picture diagonal in mm with				Order No.
	1-inch camera	2/3-inch camera	1/2-inch camera	1/3-inch camera	
<b>Without zoom magnification, for 1-chip cameras only:</b>					
C-mount adapter 1x HC	16	11	8	6	11541510
C-mount adapter 0.7x HC	–	15.7	11.4	7.8	11541543
C-mount adapter 0.55x HC	–	–	14.5	10.9	11541544
C-mount adapter 0.35x HC	–	–	–	17.1	11541512
<b>With variable magnification level (Vario TV adapter) for 1–3 chip cameras:</b>					
C-mount 0.32–1.6x HC	–	–	19*–5	18–3.8	11541517
<b>Without variable magnification level, for 1-3 chip cameras:</b>					
C-mount adapter 1x	–	–	16	12	11543706
B-mount adapter 1x	–	–	16	12	11543702
C-mount adapter -0.7x for EFW					11541545
C-mount adapter 1x for sCMOS					11541546
C-mount adapter -1.3x for EFW					11541547

\* available beginning with Vario factor 0.42 x!



Fig. 49: TV adapter

# Stages and Specimen Holders

A wide range of specimen stages are available. The most popular are:

- Fixed stage (248x204 mm) normal, heatable / temperature controlled
- Fixed micromanipulation stage (248x112 mm) heatable / temperature controlled
- Manual or motorized 3-plate stages
- Scanning stage
- Linear Motor Stage

A complete description of the stage portfolio, with inserts frames and accessories can be found in the **Live on stage Brochure**. Here a list of components:

## Fixed stage (248 x 204 mm) 11522078

88 mm Inserts with different openings for fixed stages, slim 3-plates stages and 160 x 110 mm plates.

- Insert with 5 mm opening 11522083
- Insert with 10 mm opening \*) 11522084
- Insert with 20 mm opening 11522085
- Insert with 40 mm opening 11522086
- Glass Heating Insert (Tokaihit MATS) 11532456

\*) Insert with 10 mm opening is already included in stages

11522078 • 11522015 • 11522020 • 11522069

## Fixed Heating Stage (248 x 212 mm) 11533025

## Fixed Cooling Stage (248 x 212 mm) 11522013

## Object guide for fixed regular stages

(for 11522078, 11522012 and 11522013)

Attachable object guides for all regular fixed stages measuring 248 mm x >200 mm to accommodate a variety of application inserts. 11522 014

## Inserts

for attachable object guide 11 522 014

- Holder for tissue culture plates (136 x 92 mm) 11520584
- Holder for Terasaki Plates, 60 or 72 wells 11520585
- Holder for flasks, bottles or plankton chambers 11520586
- Holder for flasks, bottles or plankton chambers 11520587
- Holder-Titer Trays, 96 or 120 wells 11520589
- Holder for Petri Dish Ø 88 mm / height 16 mm 11520590
- Holder for slides 76 mm x 26 mm 11520593
- Holding frame for plankton chamber with a basic area (102.5 x 28–50.5 mm) 11520595
- Universal Holding frame M for Petri dishes (24–68 mm) or glass slides 11533041
- for Uthermol™ Counting chambers 11532494
- Universal Holding frame MX for large Petri dishes (87–92 mm) and multiwells 11520689
- Universal Holding frame M-Duo for 1 or 2 Petri dishes (24–56mm) and/or 1 glass slides 11531798
- Heatable Universal Holding frame MH 2000\* for various dishes (24–68 mm) or microscopy-slides 11533045
- Heatable Universal Holding frame MH-L 2000\* for microscopy-slide sized vessels 11533046
- Heatable Universal Holding frame MH-R 2000\* for round dishes (24–68 mm) 11533047
- Cooling Insert X Universal Mounting frame KX 11532510




Fig. 50: Fixed stage



Fig. 51: Fixed heating stage



Fig. 52: Fixed cooling stage



Fig. 53: Universal holding frame M




Fig. 54: Various inserts for attachable mechanical stage




Fig. 55: Fixed micromanipulation stage




Fig. 56: Heating micromanipulation stage

<b>Slim Fixed Stage (248 x 112 mm)</b>	11522015
<b>Slim Fixed Heating Stage (248 x 112 mm)</b>	11533026
<b>Slim Fixed Cooling Stage (248 x 112 mm)</b>	11522017

<b>Object guide for slim fixed stages</b>	11522018
---	----------

<b>Inserts for attachable object guide</b>	11522018
• Holding frame for Petri dishes with Ø 30 mm	11522042
• Holding frame for Petri dishes with Ø 50 mm	11522043
• Holding frame for glass slides 76 mm x 26 mm	11522044

<b>3-plate-stage with rack and pinion</b>	
• Manual	11522076
• Motorized	11522068

<b>Slim 3-plate-stage with rack and pinion</b>	
• Manual	11522077
• Motorized	11522069

<b>Rotating Gliding Stage</b>	11533156
<b>Scanning Stage 127 x 83</b>	11522100
<b>SCAN Plus IM 130 x 85</b>	11525407
<b>LMT200 ITK Linear Motor Stage 120 x 80</b>	11525376

**Inserts**

for the manual 3-Plate stage 11522076,  
 for the motorized 3-Plate stage 11522068,  
 for the Scanning stages 11522100, 11522129 11532536,  
 for Linear Motor stage 11525376

• Holder for slides 3" x 1" (76 x 26 mm)	11531433
• Glass stage plate with 20 mm round opening	11522045
• Holder for Micro-Titer Trays 96 or 120 wells (tray size: 127 x 85 mm)	11531434
• Holder for Terasaki Trays 60 or 72 wells (tray size: 56 x 82 mm)	11531435
• Holder for Hamax Trays (tray size: 93 x 66 mm)	11531436
• Holder for Petri Dish Ø 36 mm	11531437
• Holder for Petri Dish Ø 54 mm	11531438
• Holder for Petri Dish Ø 65 mm	11531439
• Holder for Petri Dish Ø 88.5 mm	11531440
• Adjustable Universal Holder for Petri dishes (20-68 mm) or glass slides	11531441
• Universal Mounting frame K for Petri dishes (24–68 mm) or glass slides (with depressions at sides for micromanipulation)	11600234
• Universal insert frame KX for Petri dishes 87–92 mm or multiwells	11532338
• Universal Mounting frame K-Duo for 1 or 2 Petri dishes and/or 1 glass slide	11532514
• Mounting frame for 88 mm round inserts	11522063
• Metal Plate for 88 mm inserts lowered by 4 mm	11600237
• Heatable Universal Holding frame KH 2000 for	



- various dishes (24-68 mm) or microscopy-slides 11533048
  - Heatable Universal Holding frame KH-L 2000 for  
microscopy-slide sized vessels (47 x 11 mm opening) 11533049
  - Heatable Universal Holding frame KH-R 2000  
for round dishes (24–68 mm) 11533050
  - Heating Insert P 2000 for round dishes (24-68 mm) Compatible  
with incubator S-2. Round opening (32 mm).  
Cover with glass element for DIC 11533027
  - Heating Insert P 2000 for Lab-Tek™-type.  
Compatible with incubator S-2. Rectangular opening (46 x 21 mm).  
Cover with glass element for DIC 11533080
- Both Heating Inserts P are available as non heated versions
- Insert N for Lab-Tek™ type 11533037
  - Cooling Insert X Universal Mounting frame KX  
Especially for cooling down to 0°C. Round opening (8 mm) and  
especially suited for micromanipulation 11532510

### Accessories for scanning stages

- Cable for scanning stage to XY-advanced board 11525218

### Stage micrometer

- Transmitted light 2 mm = 200T, glass carrier with scale  
1 scale interval = 10 µm 11513106
- Incident light 10 mm = 100T for overview objectives (e.g. 1.25) 11519963

### Control Units for Heating Stages and Incubators CO<sub>2</sub> Incubation

See the **Live on Stage Brochure** for complete description.

# Accessories



Fig. 57: Screw cap

## Immersion oil, 10 ml

Free of natural fluorescence as per ISO 8036/1, refraction index  $n_e^{23} = 1.5180 \pm 0.005$ , dispersion  $v_e^{23} = 44 \pm 2$

11513859

## Immersion oil, 20 ml

as per ISO 8036/1, refraction index  $n_e^{23} = 1.5180 \pm 0.005$ , dispersion  $v_e^{23} = 44 \pm 2$

11513860

## Immersion oil, 250 ml

as per ISO 8036/1, refraction index  $n_e^{23} = 1.5180 \pm 0.005$ , dispersion  $v_e^{23} = 44 \pm 2$

11513861



Fig. 58: DIC cover

## Stage micrometer

Transmitted light 2 mm = 200 parts

11513106

## Stage micrometer

Incident light 1 mm = 100 parts

11563011

## Focusing telescope dia. 30.0 mm

11505070

## Hg 100 W/2 burner

11500321

## Xe high-pressure burner

11500139



Fig. 59: Analyzer opening cover

## Screw cap for empty objective positions

Component of stand

11 020-422-570-000

## Cover for unused objective DIC disk opening

Component of stand

11 090-144-020-058

## Dust and light protection cover for polarizer opening

Component of stand

11 020-437-101-013

## Dust and light protection cover for camera port openings

Component of stand

11 020-387-556-009



Fig. 60: Camera port cover

# Digital Image Documentation

## Leica digital camera system DC

Monochrome and color digital cameras for all applications  
(see [www.leica-microsystems.com](http://www.leica-microsystems.com))

## Software

### Leica LAS X

(see [www.leica-microsystems.com](http://www.leica-microsystems.com))

## Micromanipulators

<b>Micromanipulator</b> (for right side of microscope)	11520137
<b>Micromanipulator</b> (for left side of microscope)	11520138
<b>Large baseplate for Leica DMI8</b> to assemble the microscope and manipulators	11525214
• Extension for Leica DMI8	11521630
• Single instrument holder	11520142
• Dual instrument holder	11520143
• 3 instrument sleeves	11520145

While additional micromanipulation accessories such as microtools, pullers, grinders, microforges and anti-vibration stages are not a part of the Leica program, they can be supplied by your Leica distributor on request. Information can be provided on the adaptation of non-Leica micromanipulation tools for in-vitro fertilization, microinjection, microdissection or patch clamping.

## Anti-vibration

Large antivibration table for Leica DMI8	11525405
System Desk	11640255



Fig. 61: Leica DFC Microscope Camera




Fig. 62: Leica DMI8 with mechanical micromanipulators

# System overview


## Leica DMi8

### Automated


#### Eyepieces



#### Camera Adapter




#### 2a Condensers




#### Bino FixTube

#### Bino ErgoTubes


#### Trino ErgoTubes




#### 1 T-House




#### Mot. Side Ports



#### Coded Side Ports



2 Transmitted-Light Arms




Lasersafety

Lasersafety Kit 11889086  
 Lasersafety Hood with interlock 11889065

**Laser Safety musts:**  
 a) Sideport 0/100 % only  
 b) Condensers S1, S23 or S28 only  
 c) Tilttable TL-Arm

DIC systems



Optic Carrier

**without Bottom Port**  
 w/o IMC 11889032  
 prep. for IMC 11889033

**with Bottom Port**  
 w/o IMC 11889034  
 prep. for IMC 11889035

Front Module 11888384  
 Front Module + coded 1.6x Mag.Changer 11889076

Ergonomic Base Plate

Height Compensation Plate 23mm 11525200

Fluo Turret

Coded 11889021  
 Motorized 11889022

Filter cubes (see www.leica-microsystems.com)

11525300  
 Analyser

6 Focus

Manual Focus 11889054  
 Motorized Focus 11889055  
 Closed Loop Focus 11889056  
 Motorized Focus + AFC 11889074  
 Closed Loop Focus + AFC 11889073

**Coded Nosepiece** 11889049  
**Coded Nosepiece** 11889049  
**Motorized Nosepiece** 11889050

Tube Lens/Mag.Changer

Motorized Mag. Changer 11889024

Tube Lens 1.6x 11888377  
 Tube Lens 2.0x 11888376  
 Fixed Tube Lens 1x 11889023

Fluo Axis

Man. external 11889051  
 Mot. external 11889052  
 Mot. external with IFW 11889053  
 Cage System 30 mm 11525362

CTR Boxes

CTR Box compact 11525206

CTR Box advanced 11525207

CTR Box advanced + additional Sequencer possible 11525209

CTR Box advanced without power supply for manual stands with mot. stage 11525208

Cable Adapter 15/25 Pins 11505237  
 XY Basic Board 11525210  
 XY Advanced Board 11525211  
 MotCorr Board 11525212  
 Sequencer 11525213

Fluorescence lamp housings

Power supply EBQ 100-04-L 11500334

Lamp housing 106z L 6 Lenses Hg 100 W 1 Inch collector 11504106

Gel fiber-optics, 2 m 11504116

Leica EL6000 11504115

1-inch fiber-optics adapter 11504136

11504196  
 11504138  
 11504195  
 11504139

Leica SFL100, 365 nm (not available in the US)  
 Leica SFL100, 470 nm (not available in the US)  
 Leica SFL100, 530 nm (not available in the US)  
 Leica SFL4000 (not available in the US)

Control Elements

Leica SmartMove 11525115

Leica STP8000 11525113

Leica STP4000 11525220


TouchScreen 11889059

# System overview


## Leica DMi8

### Coded/Manual


#### 2a Condensers



#### Camera Adapter




#### Eyepieces





#### Bino FixTube

#### Bino ErgoTubes


#### Trino ErgoTubes




#### Lasersafety



#### 1 T-House



#### Optic Carrier



#### Optic Carrier




#### Coded Fluo Turret




#### Manual Fluo Turret




#### Coded Side Ports




#### Ergonomic Base Plate




#### Manual Side Ports




2 Transmitted-Light Arms




LED Lamp House



DIC systems



Fluorescence lamp housings



The statement by Ernst Leitz in 1907, “With the User, For the User,” describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Leica Microsystems operates globally in three divisions, where we rank with the market leaders.

#### LIFE SCIENCE DIVISION

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems’ customers at the leading edge of science.

#### INDUSTRY DIVISION

The Leica Microsystems Industry Division’s focus is to support customers’ pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

#### MEDICAL DIVISION

The Leica Microsystems Medical Division’s focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

Leica Microsystems – an international company with a strong network of worldwide customer services:

Active worldwide		Tel.	Fax
Australia · North Ryde	+61	2 8870 3500	2 9878 1055
Austria · Vienna	+43	1 486 80 50 0	1 486 80 50 30
Belgium · Diegem	+32	2 790 98 50	2 790 98 68
Brazil · São Paulo	+55	11 2764-2411	11 2764-2400
Canada · Concord/Ontario	+1	800 248 0123	847 405 0164
Denmark · Ballerup	+45	4454 0101	4454 0111
France · Nanterre Cedex	+33	811 000 664	1 56 05 23 23
Germany · Wetzlar	+49	64 41 29 40 00	64 41 29 41 55
India · Mumbai	+91	226 1880 200	226 1880 333
Italy · Milan	+39	02 574 861	02 574 03392
Japan · Tokyo	+81	3 6758 5630	3 5155 4333
Korea · Seoul	+82	2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31	70 4132 100	70 4132 109
People’s Rep. of China · Hong Kong · Shanghai	+852 +86	2564 6699 21 6039 6000	2564 4163 21 6387 6698
Portugal · Lisbon	+351	21 388 9112	21 385 4668
Singapore	+65	6550 5999	6564 5955
Spain · Barcelona	+34	93 494 95 30	93 494 95 32
Sweden · Bromma	+46	8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41	71 726 34 34	71 726 34 44
Turkey · Istanbul	+90	216 504 0100	216 504 0110
United Kingdom · Milton Keynes	+44	800 298 2344	1908 577640
USA · Buffalo Grove/Illinois	+1	800 248 0123	847 405 0164