Living up to Life

Leica
MICROSYSTEMS

Leica DMi8 A Leica DMi8 C Leica DMi8 M

Instructions for Use • Gebrauchsanweisung • Mode d'emploi

Leica Microsystems CMS GmbH • Manual 11934096, Revision 1.0 from 2015-01-28

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Leica Microsystems CMS GmbH • Manual 11934096, Revision 1.0 from 2015-01-28

Copyright

All rights to this documentation are held by Leica Microsystems CMS GmbH. Reproduction of text or illustrations (in whole or in part) by print, photocopy, microfilm or other method (including electronic systems) is not allowed without express written permission from Leica Microsystems CMS GmbH.

The instructions contained in the following documentation reflect state-of-the-art technology. We have compiled the texts and illustrations as accurately as possible. Still, we are always grateful for comments and suggestions regarding potential mistakes within this documentation.

The information included in this manual may be changed without prior notice.

Revision 1.0, published on 2015-01-28 by:

Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
D-35578 Wetzlar (Germany)
<http://www.leica-microsystems.com>
Responsible for contents:
Marketing CMS

Contents

1	Important Notes about this Manual	4
1.1	Text Symbols and Their Meaning	4
2.	Intended Use of the Microscopes	6
3.	Safety Notes	8
3.1	General safety notes	8
3.2	Electrical safety	8
3.3	Laser safety	11
3.3.1	Safety switch (interlock) on the microscope	11
3.3.2	Screwed-on Filter Turret Cover	12
3.3.3	Laser coupling by the user	12
3.3.3.1	Safety hazard due to folded back transmitted light illuminator arm	13
3.3.3.2	Safety hazard at eyepiece	13
3.3.3.3	Safety hazard in the condenser area during operation	13
3.3.3.4	Safety hazard above the condenser	14
3.3.3.5	Safety hazard in the specimen area	14
3.3.3.6	Safety hazard due to openings at microscope	15
3.3.3.7	Note regarding the safety hazards in Section 3.3.3.1 and Section 3.3.3.2	16
3.3.4.	Exception for Class 1 Lasers	17
3.4	Storage and transport	18
3.5	Ambient conditions	18
3.6	Notes on handling light sources	19
3.7	Notes on handling the touch screen	20
3.8	Notes on handling immersion oil	20
3.9	Notes on handling acids and bases	20
3.10	Disposal	20
4.	Leica DMI8 series instrument overview	21
4.1	Microscope variants (examples)	21
4.2	Dimensions (specified in mm)	22
4.3	Specifications ⁺	25
4.4	Controls	29
4.4.1	Control panel fields	30
4.4.2	Additional Control Elements	32
4.4.3	The touch screen	36
4.4.4	The Leica SmartMove control element	38
4.4.5	Leica STP4000/8000 external control panel	39

5.	Operation	40	7.	Troubleshooting	86
5.1	Switching on the microscope	41	8.	Care of the Microscope.....	89
5.2	Illumination	43	8.1	Dust cover	89
5.2.1	LED lamp housing	43	8.2	Cleaning.....	89
5.2.2	Leica EL6000 (fluorescence)	43	8.3	Handling acids and bases	90
5.3	Objectives	44	9.	Essential Consumable and Spare Parts	91
5.4	Tube.....	47	10.	Appendix.....	92
5.5	Eyepieces.....	48	10.1	Assigning the variable function keys	92
5.6	Stages and object displacement.....	49	10.2	Abbreviation glossary	94
5.6.1	Object displacement	49	10.3	Overview of touch screen pictograms.....	95
5.6.2	Adjusting the Step Sizes on Motorized Stages.....	50	10.4	Additional information.....	96
5.6.3	Saving and moving to stage positions.....	50	11.	Index	97
5.7	Focusing	52	12.	EU Declaration of Conformity.....	99
5.7.1	Focusing the image	52			
5.7.2	Focus thresholds.....	53			
5.7.3	Adjusting the step sizes.....	54			
5.8	Selecting the ports	55			
5.9	Magnification Changer	56			
5.10	Aperture diaphragm and field diaphragm	57			
5.10.1	Manual diaphragms for incident light	57			
5.10.2	Manual diaphragms for transmitted light	58			
5.10.3	Motorized diaphragms.....	59			
5.11	Incident light reflector disk	60			
5.12	Filter	62			
5.12.1	Filter slider for incident light	62			
5.12.2	Transmitted light filter	62			
5.13	The memory function on the touch screen	63			
6.	Operation Contrast Methods	64			
6.1	Bright field (IL)	65			
6.2	Dark field (IL)	66			
6.3	Polarization (IL)	67			
6.4	Differential interference contrast (ICR)	68			
6.5	Oblique light (IL)	70			
6.6	MACRO overview magnification	72			
6.7	Fluorescence	73			
6.8	Bright field (TL).....	75			
6.9	Phase Contrast (TL).....	77			
6.10	Dark field (TL)	79			
6.11	Polarization (TL)	81			
6.12	Differential interference contrast – ICT (TL)	83			
6.13	Combination method	85			

1 Important Notes about this Manual

Caution!

This User Manual is an essential component of the product; it must be read carefully before the product is assembled, put into operation or used, and must be kept for future reference.

This User Manual includes important instructions and information related to the operating safety and maintenance of the instrument.

1.1 Text Symbols and Their Meaning

(1.2)

Numbers in parentheses refer to illustrations; for example, (1.2) refers to Fig. 1, item 2.

→ p. 20

Numbers with an arrow, e.g. → p. 20, refer to a specific page of this manual.

WARNING indicates a potentially hazardous situation with a medium degree of risk that, if not avoided, could result in death or serious injury.

CAUTION indicates a potentially hazardous situation with a low degree of risk that, if not avoided, may result in minor or moderate injury.

Caution!

In this manual, additional safety notes are indicated with the triangle symbol shown here, and have a gray background.

!

Caution! The instrument and accessories can be damaged when operated incorrectly.

Warning of hazardous electrical voltage! Danger of electric shock!

Danger due to hot surface

Warning of optical radiation! Never look directly into the light beam! Wear safety goggles!

Warning of electromagnetic field

Warning of UV radiation! Wear safety goggles!

Warning of hazardous laser radiation – take appropriate safety precautions!

Warning of hand injuries

Instructions for disposing of the instrument, its accessories and consumables.

Connection to ground

Explanatory note

*

Item not contained in all equipment configurations

Manufacturing date, e.g. 11 / 2011 for November 2011

China RoHS 50 years EFUP
(Environmentally friendly use period)

2. Intended Use of the Microscopes

The Leica DMI8 M, DMI8 C and DMI8 A microscopes are inverted light microscopes and intended for routine and research applications.

The Leica DMI8 series can be used universally. All contrast methods, such as bright field, dark field, polarization, phase contrast, DIC (ICR/ICT), oblique light or fluorescence, are integral components of the microscope and can be adapted or exchanged quickly and easily. Variable illumination and imaging beam paths as well as HCSoptics, modular accessories and a wider range of available peripherals supplement the inverted Leica research stand from Leica Microsystems CMS GmbH.

Leica DMI8 microscopes:

Leica DMI8 M (manual)	11889092
Leica DMI8 C (coded)	11889117
Leica DMI8 A (automated)	
with coded side port	11889121
Leica DMI8 A (automated)	
with motorized side port	11889122

Directives of the European Community (EU Directives)

The requirements in EU Directives 2006/95/EC and 2014/35/EU regarding the safety of electrical equipment, and 2004/108/EC and 2014/30/EU regarding electromagnetic compatibility have been met for the microscopes described here for use in industrial environments.

The automated variant DMI8 A meets the requirements of the Machinery Directive 2006/42/EC.

Reasonably foreseeable misuse

The following are prohibited:

- To use the microscope for any purpose not in accordance with the Declaration of Conformity (e.g. the use of medical products in accordance with *EU Directive 93/42/EEC*)
- Making use of auxiliary equipment to use the area between the microscope stage and objectives as a clamp or holder (working like a vice)
- Operating the microscope in an inclined position
- Cleaning the microscope in a way other than specified in the manual
- Modifying the built-in safety circuits in the instrument
- Allowing unauthorized personnel to open the instrument
- Using cables that Leica Microsystems CMS GmbH has not provided or permitted
- Using the swiveling touch screen of the microscope as a transport handle
- Using combinations with non-Leica components that go beyond the scope of the manual

Caution!

The manufacturer assumes no liability for damage caused by, or any risks arising from, using the microscope for other purposes than those for which they are intended or not using them within the specifications of Leica Microsystems CMS GmbH.

In such cases, the Declaration of Conformity shall be invalid.

WARNING

Notes on handling laser devices

In the standard design without additional laser safety precautions, the microscopes are not suitable for coupling laser radiation (such as to the camera ports), since this radiation poses a hazard for the user (eye injury in particular).

For use of the microscope with lasers, Leica Microsystems CMS GmbH offers special microscope variants with additional safety features. Laser couplings require corresponding safety devices that have to be inspected and installed by trained personnel. For further information, please contact your authorized Leica Microsystems CMS GmbH representative.

3. Safety Notes

3.1 General safety notes

These instruments of protection class 1 are built and inspected in accordance with
EN 61010-1/IEC 61010-1, Safety requirements for electrical equipment for measurement, control and laboratory use.

They also meet the requirements for *EN / IEC 61326-1, Electrical equipment for measurement, control and laboratory use – EMC requirements.*

The products were also tested in accordance with *EN 62471 / IEC 62471 Photobiological safety of lamps and lamp systems.*

Caution!

In order to maintain this condition and to ensure safe operation, the user must follow the instructions and warnings contained in this user manual.

Caution!

The instruments and accessories described in this user manual have been tested for safety and potential hazards.

The responsible Leica Microsystems CMS GmbH affiliate or the main plant in Wetzlar must be consulted whenever the instrument is altered, modified or used in conjunction with non-Leica components that go beyond the scope of this manual.

Unauthorized alterations to the instrument or noncompliant use shall void all rights to any warranty claims!

3.2 Electrical safety

General specifications

Leica CTR compact electronics box

For indoor use only.

Supply voltage:	100 – 240 VAC
Frequency:	50 / 60 Hz
Power consumption:	max. 150 VA
Fuses:	3.15 A, slow-blowing, Breaking capacity H, 250 VAC, Size: 5x20 mm
Ambient temperature:	15° - 35°C
Relative humidity:	90% up to 30°C, non-condensing
Protection class:	I
Overvoltage category:	II
Pollution degree:	2

Leica CTR advanced / advanced+ electronics box

For indoor use only.

Supply voltage:	100 – 240 VAC
Frequency:	50 / 60 Hz
Power consumption:	max. 290 VA
Fuses:	6.3 A, slow-blowing, Breaking capacity H, 250 VAC Size: 5x20 mm
Ambient temperature:	15° - 35°C
Relative humidity:	90% up to 30°C, non-condensing
Protection class:	I
Overvoltage category:	II
Pollution degree:	2

Leica DMI8 M (manual)

For indoor use only.

Supply voltage:	100 – 240 VAC
Frequency:	50 / 60 Hz
Power consumption:	max. 55 VA
Fuses:	1.6 A, slow-blowing, Breaking capacity H, 250 VAC Size: 5x20 mm
Ambient temperature:	15° - 35°C
Relative humidity:	90% up to 30°C, non-condensing
Protection class:	I
Overvoltage category:	II
Pollution degree:	2

Leica DMI8 C (coded)

For indoor use only.

Supply voltage:	100 – 240 VAC
Frequency:	50 / 60 Hz
Power consumption:	max. 55 VA
Fuses:	1.6 A, slow-blowing, Breaking capacity H, 250 VAC Size: 5x20 mm
Ambient temperature:	15° - 35°C
Relative humidity:	90% up to 30°C, non-condensing
Protection class:	I
Overvoltage category:	II
Pollution degree:	2

Leica DMI8 A (automated)

For indoor use only.

Supply voltage:	100 – 240 VAC (→ Leica CTR)
Frequency:	50 / 60 Hz (→ Leica CTR)
Power consumption:	see Leica CTR
Fuses:	see Leica CTR
Ambient temperature:	15° - 35°C
Relative humidity:	90% up to 30°C, non-condensing
Protection class:	I (→ Leica CTR)
Overvoltage category:	II (→ Leica CTR)
Pollution degree:	2 (→ Leica CTR)

Leica EL6000*

For indoor use only.

Supply voltage:	100 – 240 VAC
Frequency:	50 / 60 Hz
Power consumption:	max. 200 VA
Fuses:	2.5 A, slow-blowing, Breaking capacity H, 250 VAC Size: 5x20 mm
Ambient temperature:	0° - 40°C
Relative humidity:	90% up to 30°C, non-condensing
Protection class:	I
Overvoltage category:	II
Pollution degree:	2
(* See provided User Manual)	

The power plug may only be plugged into an outlet equipped with a grounding contact.

Do not interfere with the grounding function by using an extension cord without a ground wire. Any interruption of the ground wire inside or outside of the instrument, or release of the protective conductor terminal, can cause the instrument to become hazardous. Intentional ground interruption is not permitted!

Caution!

By connecting it to the grounding connection (ground screw at the back of the electronics box) ancillary equipment with its own and/or extra power supply may be brought to the same ground wire potential. For connections without a ground protection conductor, Leica Microsystems CMS GmbH Service must be consulted.

Caution!

Never use any fuses as replacements other than those of the types and the current ratings listed here. Using other fuses or bridging the fuse holder is not permitted. The use of other fuses poses a fire hazard.

Caution!

The microscope's electrical accessory components are not protected against water. Water can cause electric shock.

Caution!

Protect the microscope from excessive temperature fluctuations. Such fluctuations can lead to the accumulation of condensation, which can damage the electrical and optical components.
Operating temperature: 15-35°C.

Before exchanging the fuses or lamps, be absolutely certain to switch off the main power switch and remove the power cable.

Caution!

Do not plug in or unplug data lines and control circuits unless the instrument is switched off; otherwise, the instrument may be damaged.

WARNING

By definition, the main circuit breaker of this instrument is the connection between the power cable and device port. The user must ensure unobstructed access to the main circuit breaker at all times.

WARNING

Only use power cables, or alternative cables with VDE/HAR codes that minimally fulfill the following conditions: 3x0.75 mm² and 10A/250V.

3.3 Laser safety

3.3.1 Safety switch (interlock) on the microscope

Two safety switches are integrated into the microscope with the laser safety kit option. These must be incorporated into the laser safety concept so that the light path of the laser beam is interrupted when the interlock switch is opened.

Interlock switch (1) and (2)


Interlock switch 1:

An interlock switch is integrated on the transmitted light illuminator arm and activated when the transmitted light illuminator arm is folded back (interlock switch opens).

Interlock switch 2:

An additional interlock switch is also installed in the microscope, which is activated (interlock switch opens) when the beam path to the eyepiece is released.

The interlock switches integrated into the microscope lead to two equivalent interlock connectors on the back of the microscope. The interlock connector of the microscope must be connected to your laser beam interrupter.

Interlock connectors on the rear side of the stand

Technical data for the interlock connector:

Plug:

Manufacturer: Binder
Type: 09-0097-20-05
(5-pin flange plug,
711 series subminiature round plug connector)

Compatible coupling:

Manufacturer: Binder
Type: 99-0096-102-05
(5-pin cable socket,
711 series subminiature round plug connector)

Type: 79-1414-12-05
(5-pin cable socket with 2 m cable,
702 series subminiature round plug connector)

Type: 99-0414-00-05
(5-pin cable socket,
712 series subminiature round plug connector)

Pin assignment:

Switching contact between pin 1 and pin 5.

Maximum permitted electrical characteristic values of the integrated interlock switch:

$U_{max} = 40 \text{ V DC}$
 $I_{max} = 0.2 \text{ A}$

3.3.2 Screwed-on Filter Turret Cover

As an additional safety precaution on laser microscope systems, the incident light reflector disk cover is secured using 2

DIN 912 M4x6 lens head Allen screws
(Order No. 11 703 121 630 000)

Store the screws in a safe place. Never operate the microscope when the cover is not screwed into place!

Screws for securing the reflector turret cover

3.3.3 Laser coupling by the user

Due to the danger posed to the user by laser radiation, the standard version of Leica DMI8 microscopes is not suitable for coupling laser radiation¹⁾ without additional laser safety measures.

If you as the user modify a microscope by coupling a laser beam into it, you assume responsibility as the manufacturer of the complete system.

For laser safety reasons, the microscope has the basic requirement of being equipped with a transmitted light illuminator arm.

We do not make provisions for you opening or working on the instrument, therefore you are fully responsible for your modifications complying with the basic safety requirements and laser safety precautions set forth by the European Union directives, and with the legal and technical directives and standards of your country.

If you modify your microscope by coupling a laser beam independently, when determining and implementing laser safety measures also take note of the microscope's hazard areas described in the following, as these areas pose the danger of serious and irreversible eye damage and skin injuries!

Danger of serious and irreversible eye damage from laser radiation!

By all means avoid letting direct or indirect laser radiation strike the eyes or skin! The laser light can cause serious eye damage and skin damage!

If you want to couple another laser beam into an existing laser microscope system from Leica Microsystems, you must first also check whether the already integrated shutter (laser beam interrupter) is suitable for the future higher laser power. Furthermore, you have to verify whether the existing laser warning messages on the laser microscope system and the existing laser safety precautions will still be sufficient after your modifications.

1) Laser radiation of Class 1 lasers is exempt. Further details regarding this topic in Section 3.3.4

3.3.3.1 Safety hazard due to folded back transmitted light illuminator arm

If you fold back the transmitted light illuminator arm, the laser radiation escaping from the lenses is radiated directly into the open space.

Important safety note

As a safety measure, we recommend the microscope models with integrated interlock connectors as mentioned in Section 3.3.3.7.

Folded back transmitted light illuminator arm; the light output out of the objective is not shadowed by the condenser base.

3.3.3.2 Safety hazard at eyepiece

The laser radiation that you coupled into the microscope can disperse from the eyepiece and directly hit the user's eye!

Important safety note

As a safety measure, we recommend the microscope models with integrated interlock connectors as mentioned in Section 3.3.3.7.

Eyepieces

3.3.3.3 Safety hazard in the condenser area during operation


If you couple a high power laser beam into the microscope, the laser radiation, broadened by the lens, that passes the condenser laterally, can exceed the permitted limit values for eyes and skin.

Increased laser radiation can beam laterally past the condenser

3.3.3.4 Safety hazard above the condenser


The laser beam path between the condenser base and the illuminated field diaphragm is without obstruction.

Laser beam path between condenser base and illuminated field diaphragm without obstruction

3.3.3.5 Safety hazard in the specimen area

The laser radiation is accessible without obstruction in the specimen area of the microscope!

Specimen area

3.3.3.6 Safety hazard due to openings at microscope

Make sure, as part of your risk assessment, that all openings at the microscope are tightly sealed and can only be opened with tools to prevent laser radiation from escaping!

The openings to consider, for example: top port next to the eyepieces, left and right side port, connection for transmitted light lamp housing at the illuminator arm, lamp connection for fluorescence/connection for incident light lamp housing, cover for fluorescent turret disk.

The material that is used to seal off the openings has to be suitable for the intensity of the coupled laser beam.

Furthermore, the terminals that are not occupied by objectives have to be closed with caps at the objective turret.

Openings at the stand that have to be tightly sealed

- a Cover for fluorescent turret disk
- b Transmitted Light lamp housing
- c Incident light lamp housing
- d Side port right
- e Left side port
- f Objective turret
- g Top port
- h Objective turret openings

3.3.3.7 Note regarding the safety hazards in Section 3.3.3.1 and Section 3.3.3.2

Important safety note

The overall responsibility for laser safety of the laser microscope system assembled by you remains with you as the responsible body, even if you use the microscopes from Leica Microsystems CMS GmbH listed below in conjunction. Within the context of your system composition and modifications, you have to verify whether the component selected by you from this chapter is effective and suitable for your laser protection concept.

We recommend that for your laser applications you use a microscope variant with the following selection options because the microscope is already equipped with interlock switches at the factory:

Leica No. 11889086
Laser Safety Kit

and also one of the following options:

Leica No. 11889077
Left motorized side port, 100% and interlock

Leica No. 11889078
Right motorized side port, 100% and interlock

Leica No. 11889079
Right/left motorized side port, 100% and interlock

Leica No. 11889087
Left side port 100/0, manual with interlock

Leica No. 11889088
Left side port 100/0, coded with interlock

Note:

The fixed transmitted light illuminator arm (order No. 11525105 (with short cable) or order No. 11525116 (with long cable) is not suitable for laser applications and therefore is not provided for these applications.

a) Safety function for the safety hazard mentioned in Section 3.3.3.1.

The special microscope variants have a built-in interlock switch on the transmitted light illuminator arm that is activated when folding back the transmitted light illuminator arm (interlock switch opens).

The interlock connector of the microscope must be connected to your laser beam interrupter, which interrupts the laser beam coupled into the microscope when opening the interlock switch or folding back the illumination arm.

b) Safety function for the safety hazard mentioned in Section 3.3.3.2.

An additional interlock switch is also installed in the microscope, which is activated (interlock switch opens) when the beam path to the eyepiece is released.

The interlock switch of the microscope must be connected to your laser beam interrupter, which interrupts the laser beam coupled into the microscope when opening the interlock switch or releasing the beam path towards the eyepiece.

Both interlock switches a) and b) are routed to a shared interlock connector on the microscope; you have to connect your laser beam interrupter to this connector.

When doing so, note the electrical characteristic values of the integrated interlock switches in the microscope. You must not exceed these values by connecting a laser beam interrupter, as otherwise the interlock switches will be damaged and no longer function safely and reliably.

Maximum permitted electrical characteristic values of the integrated interlock switch:

$U_{\max} = 40 \text{ V DC}$

$I_{\max} = 0.2 \text{ A}$

The integrated interlock switches in the microscope and their wiring with a laser beam interrupter/laser ensure that the laser radiation you coupled into the microscope cannot get into the user's eyes through the eyepieces.

Important safety note

In addition, observe the safety hazards described in Section 3.3.3.3 to 3.3.3.6!

3.3.4. Exception for Class 1 Lasers

Coupling a laser of Class 1 in accordance with IEC/EN 60825-1 into a microscope does not result in any danger to the user.

The listed safety hazards apply only for lasers of Class 1M, 2, 2M, 3R, 3B and 4.

Note that the rules and regulations of the United States FDA 21CFR1040.10 Class I (CFR – Code of Federal Regulation Title 21, revised as of April 1, 2011) do not distinguish between Class 1 and 1M. Thus the statement that there are no safety hazards on the microscope do not apply for the laser classification in accordance with FDA Class 1.

3.4 Storage and transport

Transport

Caution!

Transport and storage at -20 °C – +85 °C and max. 90% humidity.

For shipping or transporting the microscope and its accessory components, the original packaging should be used.

As a precaution to prevent damage from vibrations, the following components should be disassembled and packaged separately:

- Unscrew the objectives.
- Secure the objective turret, e.g. with adhesive tape.
- Remove the eyepieces.
- Remove the condenser.
- Remove the specimen stage.
- If necessary, remove the transmitted light illuminator arm.
- Remove the lamp housings.
- Secure the magnetic reflector turret cover, e.g. with adhesive tape.
- Remove all moving or loose parts
- Protect the microscope from dust and dirt using a plastic bag.

Weight

At a full load, the microscope weighs more than 18 kg. The user must take corresponding precautions for transporting the instrument.

Caution!

For transport, always remove all components listed under "Transport!".

3.5 Ambient conditions

Caution!

This microscope may not be used at elevations over 2000 m above sea level.

Caution!

Do not use this instrument near sources of high electromagnetic radiation (for example, unshielded, intentionally operated ultra-high frequency sources), because these can disrupt proper operation.

We recommend assessing the electromagnetic environment before operation of this instrument and then giving corresponding instructions.

Installation location

Work with the microscope should be performed in a dust-free room, which is free of oil vapors, other chemical vapors and extreme humidity. At the workstation, large temperature fluctuations, direct sunlight and vibrations should be avoided. These conditions can distort measurements and longer microscopic scans.

Permitted ambient conditions:

Temperature	15-35°C
Relative humidity	maximum 90% up to 30°C non-condensing

Microscopes in warm and warm-damp climatic zones require special care in order to prevent fungus contamination. For additional instructions, see the chapter on "Care of the Microscope" → p. 89.

Caution!

The microscope must be set up so that the power voltage input is not obstructed and the instrument can quickly be removed from the power supply if necessary.

Caution!

Electrical components must be placed at least 10 cm away from the wall and from flammable substances.

3.6 Notes on handling light sources

CAUTION

Light sources pose a potential irradiation risk (glare, UV radiation, IR radiation). Therefore, lamps may only be operated in closed housing and in an installed condition. Never look directly into the beam path (blinding hazard).

CAUTION

Do not switch on the light source until the light guide is connected to the microscope.
Uncontrolled light output from the light guide poses a blinding hazard!
There is also a burn hazard at the light guide output!

CAUTION

Make sure to follow the safety notes in the operating manual for the light source.

CAUTION

Switch on the light source only if it is firmly connected to the microscope.
Uncontrolled light output poses a blinding hazard!

CAUTION

Do not look into the eyepieces when toggling the contrast method!
During the toggling procedure, the maximum radiant power of the light source may be present on the eyepieces for a short time and can temporarily blind the user!

Lamps and lamp housing may be hot!
They must be placed at least 10 cm away from the wall
and away from flammable substances.
Supply and data lines in particular may not come into
contact with lamp housings!

3.7 Notes on handling the touch screen

Caution!

Touch the touch screen using your finger only. Never
use a pen or other hard, sharp or pointed objects.

3.8 Notes on handling immersion oil

Caution!

When using immersion oil, take care to avoid skin con-
tact! Ask the supplier for a safety data sheet!

3.9 Notes on handling acids and bases

For examinations using acids or other aggressive chemi-
cals, particular caution must be taken.

Caution!

Be absolutely certain to avoid direct contact with these
chemicals.

3.10 Disposal

Once the product has reached the end of its service life,
please contact Leica Microsystems GmbH Service or Sales
about disposal.

Please observe and comply with the national and federal
laws and regulations that are equivalent to EC directives
such as WEEE.

Note!

Like all electronic devices, the microscope, its acces-
sory components and consumables must never be dis-
posed of with general household waste.

4. Leica DMi8 series instrument overview

4.1 Microscope variants (examples)

Leica DMi8 A with at least one motorized component


Leica DMi8 C with coded components

Leica DMi8 M with manual components


4.2 Dimensions (specified in mm)

Leica CTR compact

Weight: approx. 2 kg



Leica CTR advanced/advanced+

Weight: approx. 4 kg



4. OVERVIEW OF THE INSTRUMENT

Leica DMi8 M/C

Weight: approx. 15 kg

Leica DMI8 A (with tiltable transmitted light illuminator arm)

Weight: approx. 19 kg

Height compensation plate*

In order to enlarge the viewing height by 23 mm, raise the side camera ports for very large cameras or spinning disks, or to use a microscope with an inactive bottom port even without a hole in the work table, a height compensation plate (order No. 11525200) was developed.

4.3 Specifications⁺⁾

Contrast Methods	- Incident light (IL): BF, DF, POL, DIC (ICR), oblique light, Fluo, Macro Additionally: MACRO overview magnification (motorized or manual) - Transmitted light (TL): BF, DF, POL, ICT - Combination (IL/TL): BF/BF
Incident light axis	<u>Manual:</u> - Lamp housing receptacle - Centerable iris diaphragms with aperture and field diaphragm function - Holder for polarizer slider - Holder for light filter <u>Motorized:</u> - Lamp housing receptacle - Motorized pinhole diaphragm disks with aperture and field diaphragm function - Holder for polarizer slider - Holder for light filter
Transmitted light axis	<u>Tiltable transmitted light illuminator arm</u> - Manual and coded - With integrated mechanical tilting mechanism with sufficient room for specimens and micro-manipulators - For condensers $\leq$ S28 with an integrated manual or motorized illuminated field diaphragm - Filter magazine for 2 replaceable filters - With condenser quick change - Manual shutter - Lamp housing receptacle for replaceable LEDs - With integrated cable duct - Condenser stop (to prevent from being pushed down accidentally) <u>Fixed transmitted light illuminator arm</u> - With integrated LED illumination - Manual adjustment of brightness and aperture diaphragm - Automatic brightness adjustment when changing the contrast method - Automatic shutoff (adjustable) - Integrated condenser holder - Receptacle for transmitted light filter
Tube	- Ergonomic with or without photo output on the left side 2 switching positions: 100% VIS and 50%VIS/50%CAM or 2 switching positions: 100% VIS and 0%VIS/100%CAM - Interpupillary distance control - Height and angle adjustment (30° - 45°)

⁺⁾ Depending on configuration

Magnification Changer	<u>Manual</u> • 2 switch positions Magnification factors of 2x or 1.5x • Works on the tube port and eyepieces <u>Motorized</u> • 3 or 2 switch positions Magnification factors of 1x; 1.5x; 2x or 1x; 1.5x or 1x; 2x • Works on all camera ports and eyepieces
Objective turret	<u>Manual</u> • 6-position for objectives with M25 threads and a 45mm parfocalizing distance <u>Coded</u> • 6-position for objectives with M32 threads and a 45mm parfocalizing distance • 6-position for objectives with M25 threads and a 45mm parfocalizing distance <u>Motorized</u> • 6-position for objectives with M32 threads and a 45mm parfocalizing distance • 6-position for objectives with M25 threads and a 45mm parfocalizing distance All turrets can be equipped with a manual DIC prism slide holder or with a motorized or manual/coded Wollaston prism turntable.
Incident light reflector disk	<u>Manual:</u> • 6-position <u>Coded:</u> • 2-position • 6-position <u>Motorized:</u> • 6-position
Stages	<u>Fixed stage:</u> • Stage plate coated in ceramic (248 mm x 204 mm) <u>Manual 3-plate cross-stage</u> • Travel range: 100 mm x 50 mm • Optional different inserts for various applications, Size of the inserts: 160 mm x 110 mm (compatible with manual stage and motorized stage) and round insert disks (diameter 110 mm) with a center opening of 5, 10, 20 and 40 mm <u>Motorized 3-plate cross-stage:</u> • Travel range: 100 mm x 45 mm • Optional different inserts for various applications, Size of the inserts: 160 mm x 110 mm (compatible with manual stage and motorized stage) and round insert disks (diameter 110 mm) with a center opening of 5, 10, 20 and 40 mm

Stages	<u>Scanning stage (motors at bottom)</u> • 2 mm screw pitch • Travel range: 100 mm x 45 mm • Optional different inserts for various applications, Size of the inserts: 160 mm x 110 mm (compatible with manual stage and motorized stage) and round insert disks (diameter 110 mm) with a center opening of 5, 10, 20 and 40 mm
Condensers	<u>For tiltable transmitted light illuminator arms:</u> • Motorized and coded • Motorized aperture diaphragm • Contrast method: BF, DF, PH, ICT, Pol • Automatic method toggling • Condenser disk with 7 positions for contrast methods • Condenser housing for S1–S28 • Condenser heads: S1/1.4 oil, S1/0.9 dry, S23/0.53, S28/0.55 • Condenser heads can be swung out • All condensers for magnifications of 1.25x to 100x • With or without motorized or manual polarizer (optional) • With motorized or coded Wollaston prism disk (optional) <u>For the fixed transmitted light illuminator arm</u> • S40/0.50 condenser • S80/0.30 condenser
Z-drive	<u>Manual</u> • 12 mm travel path (2 mm below, 10 mm above the stage) • Handwheels with 3 increments for coarse, medium and fine focusing and an adjustable upper focus stop or • Handwheels with 2 increments for coarse and fine focusing <u>Motorized</u> • 12 mm travel path (2 mm below, 10 mm above the stage) • Maximum travel speed: 5 mm/s • 5 focus increments: 0.05 μm; 0.1 μm; 0.7 μm; 1.5 μm; 5.0 μm • Electronic focus repositioning • Automatic lowering before objective change • Electronic parfocality • Handwheels, each with 2 increments for coarse and fine focusing (4 increments)
Observation ports	<u>Manual</u> • Left side port (80% or 100% transmission) • Optional: Top port with 2 switch positions • 100% to the eyepieces • 50% to the eyepieces / 50% to the camera port <u>Coded:</u> • Left side port (80% or 100% transmission) • Optional: Top port with 2 switch positions • 100% to the eyepieces • 50% to the eyepieces / 50% to the camera port

Observation ports	<u>Motorized</u> • Left side ports (100%, 80% or 50% transmission) • Optional: Right side ports (100%, 80% or 50% transmission) • Optional: Top port with 2 switch positions • 100% to the eyepieces • 50% to the eyepieces / 50% to the camera port
Controls	<u>Manual:</u> • Rotary knob for setting the brightness • Switch rods for operating the aperture diaphragm, field diaphragm and oblique illumination <u>Coded:</u> • Rotary knob for setting the brightness • Button for toggling between transmitted light/incident light • Switch rods for operating the aperture diaphragm, field diaphragm and oblique illumination • Switch rods for operating the port and manual magnification changer <u>Motorized:</u> • Touch screen (only for Leica DMI8 A) • Button for toggling between TL/IL • 2 buttons for motorized aperture diaphragm adjustment • 2 buttons for motorized illuminated field diaphragm adjustment • 1 button for the oblique illumination • 4 variable function keys • 3 operating buttons for focus thresholds (only for electronic focusing) • 2 handwheels for focusing • Rotary knob for setting the brightness • 6 buttons for reflector cubes • 6 buttons for objectives • Leica SmartMove: ergonomic control element for monitoring x,y,z and 4 additional variable function keys • Leica STP4000: operation via external touch screen • Leica STP8000: operation via touch screen using control element for monitoring x,y,z and additional variable function keys • Switch rods for operating the aperture diaphragm, field diaphragm and oblique illumination • Switch rods for operating the port and manual magnification changer
Electronics box	Separate unit for controlling motorized elements of the microscope <u>Leica CTR compact:</u> • for motorized stages <u>Leica CTR advanced:</u> • for scanning stages

Interfaces	• Service port for saving/loading microscope settings via the touch screen <u>Leica DMI8 C:</u> • 1x I²C for external/internal peripheral devices <u>Leica DMI8 A:</u> • 6 x I²C for external/internal peripheral devices • CTR boxes

4.4 Controls

Depending on the configuration of your microscope, various controls are available to you. Accordingly, the control panel fields on the left and right sides of the stand contain a different number of buttons. However, the functions themselves are assigned to specific buttons.

The controls are displayed below using typical examples for the manual, coded and fully automated microscope variant.


The touch screen is available as a control element on the Leica DMI8 with at least one motorized component.

The external controls are the Leica SmartMove control and the Leica STP4000 and Leica STP8000 control panels.

4.4.1 Control panel fields

Control panel field on the left microscope side


- 1 Motorized aperture diaphragm adjustment
- 2 Motorized illuminated field diaphragm adjustment
- 3 Motorized toggling between the TL/IL (Fluor) illumination axes
- 4 Adjusting the brightness
An LED indicates the illumination method which is currently affected by the brightness adjustment.
- 5 Adjusting the oblique illumination
One LED shows the status:
LED off = oblique illumination switched on
LED on = oblique illumination switched off
- 6 Focusing

Maximum system design of the left control panel field
The number of buttons may vary depending on the configuration.

Control panel field on the right microscope side

- 1 Variable function keys (freely configurable)
The function keys are preassigned at the factory (see "Identification Sheet" and p. 92).
- 2 Focus threshold (the LED illuminates when the focus threshold is set)
- 3 Setting the focus position and focus threshold
- 4 Focus position (the LED illuminates when the focus position is set)
- 5 Focusing

Maximum system design of the right control panel field
The number of buttons may vary depending on the configuration.


Control panel fields on the front side

Swiveling touch screen for adjusting all automated components on the Leica DMI8 A
Adjustable angle for better readability.
Refer also to Chapter "The Touch Screen" p. 36.


Coded objective turret

- 1 6 LEDs for indicating the current objective
- 2 Label fields ⁺⁾


Coded objective turret and coded incident light reflector disk

- 1 6 LEDs for indicating the current reflector cube
- 2 Label fields ⁺⁾
- 3 6 LEDs for indicating the current objective

Motorized objective turret

- 1 6 objective changer buttons with LEDs for indicating the current objective
- 2 Label fields ⁺⁾

Motorized objective turret and motorized incident light reflector disk

- 1 6 reflector changer buttons with LEDs for indicating the current reflector cube
- 2 Label fields ⁺⁾
- 3 6 objective changer buttons with LEDs for indicating the current objective


⁺⁾ Stickers for labeling are included in the delivery package

4.4.2 Additional Control Elements


Right side of Leica DMi8 with manual incident light axis

- 1 Manual Magnification Changer
- 2 Stage movement
- 3 Manual reflector change
- 4 Illuminated field diaphragm adjustment
- 5 Holder for filter slider
- 6 Centering of aperture diaphragm
- 7 Centering illuminated field diaphragm
- 8 Identification label
- 9 On/off switch (only for Leica DMi8 M)
- 10 Polarizer
- 11 Toggling the left side port
- 12 Analyzer
- 13 Changing objectives
- 14 Focusing

Left side of Leica DMi8 with manual incident light axis

- 1 Aperture diaphragm adjustment
- 2 Slider for oblique illumination
- 3 Cover for left side port
- 4 Focusing
- 5 Brightness adjustment
- 6 Adjusting the oblique illumination
- 7 Manual reflector change
- 8 Top port toggle

Right side of Leica DMi8 A

- 1 Tube setting
- 2 Polarizer
- 3 Centering of aperture diaphragm
- 4 Centering illuminated field diaphragm
- 5 Analyzer
- 6 Holder for color filter
- 7 Cover for right side port
- 8 Focusing
- 9 Focus thresholds
- 10 Variable function keys

DIC setting

- 1 DIC fine adjustment
- 2 Selecting the objective prisms
- 3 Identifying letter to match objective

Manual alternative: objective prism slider


Manual DIC objective prism disk

Motorized DIC objective prism disk

Left side of the Leica DMi8 with a tiltable transmitted light illuminator arm

- 1 Illuminated field diaphragm adjustment (motorized or manual)
- 2 Condenser height adjuster
- 3 Swinging in the polarizer
- 4 Aperture Diaphragm Adjustment
- 5 Condenser setting
- 6 Condenser height stop


Right side of the Leica DMi8 with a tiltable transmitted light illuminator arm

- 7 Swinging the transmitted light filter in and out
The front position is assigned at the factory using a shutter (dark stop), which is interchangeable with a filter.

Left side of the Leica DMi8 with a fixed transmitted light illuminator arm

- 1 Filter receptacle
- 2 Aperture Diaphragm Adjustment
- 3 Stop lever for condenser height adjustment
- 4 Light ring slide

4.4.3 The touch screen

The touch screen adjusts all motorized components of the Leica DMI8 A automatically via buttons. After the instrument is switched on, the display shows the current microscope status. The display and functions that can be operated using the touch screen depend on the features of the individual microscope.

The meanings of the most important pictograms are listed in the Appendix → p. 95.

The display shows various levels of operation:

- Level 1: Navigation bar
- Level 2: Menu bar
- Level 3: Control panel field

When you select a pictogram in the navigation bar (on the left-hand border of the display), the corresponding menu bar (on the top border of the display) is shown. The display on the control panel changes accordingly. For some functions, an additional information field appears above the operating keys.

The navigation bar allows you to quickly move between navigation items:


Basic microscope settings


Contrast methods

Magnification

Stage and focus controls


Optional

Touch screen settings

Initialization

Touch screen after initialization

**Note:**

Before working with the touch screen, we recommend calibrating it first.

Calibration

- Press the  button.
- Press the  button.
- Follow the instructions and press all 4 corners of the working area using your finger. Press and hold the calibration cross until the next cross is displayed. This adjusts the sensitivity of the touch screen and the pressure surface to the user.

Other touch screen settings:**Language selection**


The current language for the display texts is shown, and you can change it by selecting the corresponding flag.

Saving the microscope settings

In Microscope Settings, press **Save** to write the current microscope settings to a USB stick, inserted into the left side of the microscope (Service port). Pressing **Load** restores all settings.

Display brightness

To prevent unwanted interference from light, you can switch off the display using the **Switch off** key. To switch it back on, simply touch any location on the display. You can also change the brightness incrementally using the **Brightness** key.

**Version number**

This button shows you the version number of the touch screen software and that of the master firmware of the connected microscope.

Please have these available for any inquiries.

**Note:**

The meaning of the other menu items of the navigation bar is explained in the respective chapters.

4.4.4 The Leica SmartMove control element

Rotary knobs on the SmartMove

The rotary knobs 1 and 2 move the specimen stage in the x and y-direction.

The image is focused using the rotary knob 3.

You can adjust the height of the rotary knobs by turning the wheel 4 to a working height that is convenient for you.

Variable function keys on the SmartMove

The variable function keys are assigned functions at the factory that are appropriate to the features of your microscope. They are labeled accordingly. For details on key assignments, please refer to the included identification sheet. For information on the abbreviations used, please refer to the list → p. 94.

Note:

The key assignment can be changed using the LAS X software.

Note:


To lock the controls, press and hold the two front keys (5) simultaneously for 5 seconds. Repeat this process to unlock the controls. The status LED (6) indicates whether the Leica SmartMove is in normal mode or locked mode.

Red LED = Entries blocked

Green LED = Normal mode

Leica SmartMove control

- 1 Movement in x-direction
- 2 Movement in y-direction
- 3 Focus setting
- 4 Individual setting of the knob height position
- 5 Variable function keys (preset at the factory)
- 6 Status LED

4.4.5 Leica STP4000/8000 external control panel


The Leica STP4000 enables remote control of all automated microscope functions using an external touch screen. The Leica STP8000 control panel also enables monitoring in X, Y and Z and has user-definable function keys. Separate instructions are provided for operation of the Leica STP4000/8000.

Leica STP4000


Controls of the Leica STP8000

- 1 Touch screen
- 2 Information key
- 3, 4 Variable function keys, user-programmable
- 5 Fine focus adjustment
- 6 Coarse focus adjustment
- 7 Movement in y-direction
- 8 Movement in x-direction

Touch screen for STP4000/8000

5. Operation

Caution!

The microscope may only be set up, installed and put into operation by authorized Leica personnel!

Caution!

The microscope may only be operated by personnel who have been instructed by an authorized Leica employee!

Note:

If at all possible, avoid touching optical surfaces such as the lens surfaces on the objective. If fingerprints do appear on the glass surfaces, remove them with a soft leather or linen cloth. Even small traces of finger perspiration can damage the surfaces in a short time. For additional instructions, see the chapter entitled "Care of the Instrument" → p. 89.

Note:

Operating the motorized components using the Leica STP4000/8000 control panel is described in separate Instructions for Use.

Note:

For information about using the LAS software, refer to the corresponding online help.

5.1 Switching on the microscope

Leica DMi8 M

- Switch on the microscope at the on/off switch on the right side of the stand. The indicator lamp is illuminated in green during operation.

Leica DMi8 C and Leica DMi8 A

- Switch on the electronics box on the on/off switch. The indicator lamp is illuminated in green during operation.

Note:

If you have connected a PC, switch on the **electronics box first**, and **then the computer**.

All motorized microscope components will then run through an initialization phase.

Note:

At initialization, all of the user's most recent settings are restored.


Caution:

The focus position and lower threshold are stored when switching off the microscope.

After the initialization is completed, the touch screen shows the status screen with the current microscope settings.

For fluorescence:

- Switch on the Leica EL6000 compact light source on the on/off switch. The indicator lamp is illuminated during operation.

Do not switch on the light source until the light guide is connected to the microscope.
Uncontrolled light output from the light guide poses a blinding hazard!
There is also a burn hazard at the light guide output!

Caution!

When using other light sources, observe the respective instructions for the light source.

5.2 Illumination


Caution!

When using other light sources than those described here, observe the respective instructions for the light source.

5.2.1 LED lamp housing

- If necessary, switch to the incident light axis or transmitted light axis.
- Set the light intensity using the knob on the left side of the stand.


Adjusting the brightness using the touch screen:

- If necessary, switch to the incident light axis with the **IL** button or to the transmitted light axis with the **TL** button.
- Adjust the intensity using **+** and **-**.
The setting can be made in coarse and fine steps in a range from 0 to 255. Pressing the middle key, which also displays the current value, switches between coarse and fine adjustment.

The light source can be shut off using the  button (shutter function).

- The intensity is individually adjusted and stored for each objective and contrast method.

5.2.2 Leica EL6000 (fluorescence)

- If necessary, switch to the incident light axis.
- Set the light intensity on the Leica EL6000 compact light source.
(See separate Instructions for Use).

CAUTION

Do not switch on the light source until the light guide is connected to the microscope.

Uncontrolled light output from the light guide poses a blinding hazard!

There is also a burn hazard at the light guide output!

5.3 Objectives

! Caution:

If necessary, first lower the focus to prevent damaging the objective front lenses!

5.3.1 Changing objectives


Leica DMi8 M with manual objective turret

- Configure an objective by rotating the objective turret on the knurled disk.
During the manual objective change, be sure that the turret locks into place.

Leica DMi8 C/A with manually coded objective turret

- Configure an objective by rotating the objective turret on the knurled disk.
During the manual objective change, be sure that the turret locks into place.
The LED assigned to the selected objective on the front side of the microscope is illuminated.

a Keypad for coded objective turret

b Keypad for coded objective turret and coded reflector turret

Leica DMi8 A with motorized objective turret and front operating panel with function keys

- Configure an objective by pressing the corresponding objective changer button on the keypad of the front side of the microscope.
The LED of the objective changer button is illuminated.

a Keypad for motorized objective turret

b Keypad for motorized objective turret and motorized reflector turret

Leica DMI8 A

The objectives can be swung into the beam path by motorized action using variable function keys on the stand or Leica SmartMove or using the touch screen.

Note:

The positions of the objectives in the objective turret have been specified at the factory and must be observed when installing the objectives.
(See Installation Manual).

When selecting an objective, the microscope automatically selects:

- The most recently configured position of the illuminated field diaphragm
- The most recently configured position of the aperture diaphragm
- The most recently configured light intensity in the respective contrast method (with DIC fine adjustment where applicable)
- The most recently configured focus (including parfocality) for microscopes with motorized focus

Adjusting the objectives using the touch screen:

- Use the  tab to switch to the Magnification menu.


The information field on the top edge of the screen shows the current total magnification at the specified observation port.


Pressing the  key allows you to switch between details for the various observation outputs.

The objective keys are labeled with the corresponding objective magnification. The objectives permitted for the selected contrast method are marked with a white triangle. Immersion objectives are also marked with drop symbol.

- To switch between objectives, press the corresponding key, which is then color-coded.

Function keys on the right side of microscope

Immersion objectives

- For **immersion objectives** use the appropriate immersion medium.

OIL: Use optical immersion oil in accordance with DIN/ISO standards

only

Cleaning → p. 89

W: Water immersion

IMM: Universal objective for water, glycerin,
Oil immersion

Caution!

Follow safety data sheet for immersion oil!

- Select an immersion objective.
The objective turret is lowered to the lower threshold. This makes it possible to apply the immersion fluid when switching from the dry objective to the immersion objective. Conversely, the immersion fluid can be removed. The current objective remains in the beam path.
- Press the key for the desired immersion objective, which is now approached.

Note:

For lockable immersion objectives, press the front part upwards as far as it will go (approx. 2 mm). Then, after a gentle turning motion to the right, the objective is locked.

For objectives with a corrective mount, adapt the objective to the thickness of the cover slip by turning the knurled knob.

Note:

If objectives are retrofitted, they have to be taught in using the LAS software. Afterwards, the parfocality should also be taught in again.

See the instructions for the LAS software.

5.4 Tube

Note:

Close any unused tube openings, as otherwise stray light can interfere with observation.

Adjusting the interpupillary distance

- Adjust the viewing distance of the eyepiece tubes so that a congruent total image is seen.

Adjusting the viewing angle


- For the ergo tubes, the viewing angle can be adjusted from 30° to 45° by tilting the binocular eyepiece.

Beam splitting in photo tubes

The light division is set manually by pulling out a control bar.

	Observation(VIS)	Photo
☐ 100 %	100 %	0 %
☐ 0 %	0 %	100 %
Alternatively:	50 %	50 %

5.5 Eyepieces

Note:

The eyepiece's aperture protector must be removed or folded back during microscopy while wearing eyeglasses. We recommend removing bifocals and spectacles with progressive-addition lenses when using the microscope.

- For the adjustable tubes with documentation output, choose the 100% VIS position.

Eyepieces with inlaid graticule

- Focus the graticule by adjusting the eyelens in the eyepiece.
- Focus on the object through this eyepiece.
- Then, close that eye and focus on the object by adjusting the second ocular only.

Correction for vision problems

- With your right eye, look through the right eyepiece and bring the specimen into sharp focus.
- Then, with your left eye, view the same position of the specimen and rotate the left eyepiece tube until this position is brought into sharp focus. While doing so, do not change the adjustment of the z-position!

Note:

We recommend using the LAS software to teach in eyepieces that are not included in the scope of delivery or are retrofitted. This ensures that the total magnification specified on the touch screen is correct. See the instructions for the LAS software.

5.6 Stages and object displacement

A wide variety of stages is available in variants with fixed stage, manual and motorized 3-plate cross-stage and scanning stage with various mechanical stages and inserts.

5.6.1 Object displacement

Fixed stage

The object displacement is carried out manually by moving the specimen on the specimen stage or using the mechanical stage.

Manual 3-plate cross-stage

The object displacement in the x- and y-direction takes place using a coaxial drive.

Motorized 3-plate cross-stage

Object displacement in the x- and y-direction can be performed by the Leica SmartMove control element, the Leica STP8000 main switch board or the LAS software.

Scanning stage/linear motorized stage

Also read the separately provided instruction manual.



Note:

When using other external stages, observe the respective instructions for the external stages.

In spite of the low speed, there is a risk of pinching and crushing injuries when moving the scanning stage. Never put objects or body parts in the travel path!


Stage movement using Leica SmartMove

5.6.2 Adjusting the Step Sizes on Motorized Stages

- On the touch screen, use the  key for switching to stage and focus controls. The X/Y coordinates of the current stage position are displayed.


- Select the **Stage** menu using the  key.
- You can change the travel speed of the stage by switching between the **Coarse** and **Fine** step sizes. If the **Coarse** value is selected, the travel speed is the same for all objectives. The **Fine** value is adapted to the respective objective.

5.6.3 Saving and moving to stage positions

In the LAS software, the programmed working area of the stage is marked by a dashed line frame.

The **Limits**  key resets the working area to the maximum travel range of the stage.

Up to 5 stage positions and one load position (L) can be stored.

- Selecting stored positions using the + and – keys

- Saving the current stage position under this number using the key

If a stage position was previously saved in this location, it will be overwritten.

- Moving to the current stored position using the  key
When the load position is approached, the objective turret is lowered to the lower threshold.

Note:

If a stored position is outside of the stage's working area, for example, as the result of subsequently changing the working area using the LAS software, you must restore the stage's maximum positioning range via the **Limits** button

to reach the position.

The stored positions will be saved until the microscope is switched off and are deleted after a restart.

5.7 Focusing

Note:

Alternatively, handwheels with two increments for coarse and fine focusing or with three increments for coarse, medium and fine focusing are available for the Leica DMI8 M, as is an adjustable focus stop.

5.7.1 Focusing the image

Leica DMI8 manual focus

- Focus the image by turning the focus dials on the left and right sides of the stand.

Speed switch (optional)

Fine focusing can be adjusted in two speeds. We distinguish between the fine and super-fine switching stages.

- To switch speeds, press the left focus knob to the right or the right focus knob to the left.


Leica DMI8 A with motorized focus

- Focus the image using the rear rotary knob on the Leica SmartMove control element.
- Alternatively, the focus dials on the left and right sides of the stand can also be used.

Note:

For coarse and fine focusing with the handwheel, 2 focus stages each can be selected using the LAS software.

5.7.2 Focus thresholds

Leica DMi8 A with motorized focus

Note:

The parfocality is taught in at the factory. As a result of screwing in the objectives during assembly, it might be necessary to teach in the parfocality again. We recommend checking the parfocality before setting the thresholds and, if necessary, teaching them in again using the LAS software.

Note:

The focus position and stops are stored by the microscope and retained from one session to the next when power is switched off.

Setting thresholds

Lower focus threshold:

- Press the **SET** key and the key

When the threshold is set, the LED is illuminated.

Pressing the same key combination again deletes the threshold.

The **lower threshold** is the same for all objectives and cannot be exceeded.

Focus position:

- Press the **SET** key and the key


When the focus position is set, the LED is illuminated.

Pressing the same key combination again deletes the focus position.

The **focus position** cannot be exceeded.

It should be defined for the dry objective with the highest magnification level. For all other objectives, it is set automatically, taking into account the parfocality adjustment and the working distance.

Right side of Leica DMi8 with motorized focus

- 1 Focusing
- 2 Focus stop
(the LED illuminates when focus threshold is set)
- 3 Setting the focus position and focus threshold
- 4 Focus position
(the LED illuminates when the focus position is set)

Moving to the thresholds

Lower focus threshold:

- Press and hold the key

Focus position:

Press and hold the key

Note:

When moving to the thresholds, the keys must be held until the position has been reached.



Setting the focus via the touch screen

- On the touch screen, use the  key for switching to stage and focus controls.
- Select the **Focus Drive** menu using the  key. The current Z position is shown in the information field.
- When the focus position is set, it can be approached directly using the  key.

When the lower threshold is set, use the  key to reach the lower threshold.

The Z drive continues only as long as the button is held depressed. The process can be interrupted at any time to prevent a collision of the objective and specimen in case of an incorrectly set lower stop.

- To set or clear the focus position or the lower threshold, select the corresponding position under **Set/Clear**, then press the **Set** or **Clear** key.

5.7.3 Adjusting the step sizes

Leica DMI8 A with motorized focus

- Use **Coarse** and **Fine** to switch between coarse and fine focusing at the stand handwheel or the focus knob of the Leica SmartMove as applicable.

The **Fine** value is adapted to the respective objective. The values are predefined logically. You can change the allocation using the LAS software.

If the **Coarse** value is selected, the travel speed is the same for all objectives. **Coarse** corresponds to the maximum speed.

Note:

The allocation of a certain step size to an objective applies not only for the z-drive, but also for defining the stage step size, which is allocated to this objective after selecting **Fine**.


5.8 Selecting the ports

- The top port can be operated manually using the switch rod (see Tube Operation).

Leica DMi8 M and Leica DMi8 C

The manual switch rod on the right side of the stand enables or disables the left side port.



	VIS	LEFT
	100 %	0 %
	20 %	80 %
Alternatively:	0 %	100 %

Leica DMi8 A

Adjusting the ports using the touch screen:

- Use the  button to switch to the Magnification menu.
- Select the left side port using the corresponding button or point 100% of the light at the eyepieces. The middle button distributes 50% of the light to the left side port and 50% to the eyepieces. Other options are possible, depending on the configuration.

Note:

Side ports not in use must be covered with cover caps.

5.9 Magnification Changer

Leica DMi8 M and Leica DMi8 C

As an option, a mechanical Magnification Changer can be used.

Magnification factors: 1x ; 1.5x or 1x ; 2x

- Using the slide bar, toggle between 1x and the magnification factor.

Note:

The mechanical Magnification Changer works on the eyepiece and the top port.

Leica DMi8 A

As an option, a motorized Magnification Changer with 2 or 3 switch positions can be used.

Magnification factors:

For 2 switch positions: 1x ; 1.5x or 1x ; 2x

For 3 switch positions: 1x ; 1.5x ; 2x

Setting the magnification changer via the touch screen:

- Use the  button to switch to the Magnification menu.
- Select the magnification factor using the corresponding key.


Note:

The selected factor is included in the calculation of the total magnification.

Note:

The motorized Magnification Changer works on all ports.

5.10 Aperture diaphragm and field diaphragm

5.10.1 Manual diaphragms for incident light

Leica DMi8 with manual aperture diaphragm

- The manual aperture diaphragm is adjusted using a slide bar on the left side of the stand.


Leica DMi8 with manual illuminated field diaphragm

- The manual illuminated field diaphragm is adjusted using a slide bar on the right side of the stand.

Note:

The aperture diaphragm and illuminated field diaphragm are adjusted when setting up the instrument (see the installation instructions).

Adjusting the aperture diaphragm

Adjusting the illuminated field diaphragm

5.10.2 Manual diaphragms for transmitted light

Leica DMi8 with manual aperture diaphragm

For a tiltable transmitted light illuminator arm:

- The manual aperture diaphragm is adjusted on the condenser.

For a fixed transmitted light illuminator arm:

- The manual aperture diaphragm is adjusted on the illuminator arm.

Leica DMi8 with manual illuminated field diaphragm

For a tiltable transmitted light illuminator arm:


- The manual illuminated field diaphragm is adjusted on the illuminator arm.


Note:

The illuminated field diaphragm is available only in conjunction with $\leq$ S28 condensers.

Leica DMi8 with a fixed transmitted light illuminator arm

Leica DMi8 with tiltable transmitted light illuminator arm

5.10.3 Motorized diaphragms

Both diaphragms have been set to suitable values for the current objective and contrast method at the factory.

- Using the **Aperture** (aperture diaphragm) or **Field** (field diaphragm) keys, the motorized diaphragms can be changed at any time. The display on the touch screen changes accordingly.

Then, the function keys are assigned to the currently active axis for incident light (IL) or transmitted light (TL).

Note:

The old values will be overwritten by the current ones!


Note:

When using **PH** or **DF**, the aperture diaphragm is open all the way and cannot be closed.

Adjusting the diaphragms using the touch screen

- On the touch screen, use the  button to switch to the basic microscope settings.
- Select the **Status/Illumination** menu using the  key.
- Adjust the aperture and field diaphragms using the **+** and **–** keys.

5.11 Incident light reflector disk

Note:

Observe the Installation Manual when equipping the reflector turret or replacing individual reflector cubes.

Leica DMi8 with manual incident light reflector disk

- The reflector cubes are manually swung into the beam path by turning the turret disk.
- A color mark on the edge of the disk indicates the current reflector cube. If the shipped configuration is changed by inserting another reflector cube, note the color of the position where the filter cube was inserted so you can allocate it later.

Changing reflector cubes:

- Remove the cover from the turret disk. The cover is secured in place by a magnetic holder and can be detached from the front.
- Pull the reflector cube forward to remove it and slide the new reflector cube straight into the receptacle.
- Put the cover back in place.

Note:


If your microscope is equipped with a 2x reflector disk, to remove or insert the reflector cube you have to turn the reflector disk slightly to the right or left. You cannot remove the reflector cube when it is in the engaged position, i.e. directly from the front.

Changing reflector cubes for laser applications

WARNING

For use of the microscope with lasers, Leica Microsystems CMS GmbH offers special microscope variants with additional safety features. Laser couplings require corresponding safety devices that have to be inspected and installed by trained personnel.

Empty position for mounting a reflector cube. The reflector cube is pushed in straight from the front until it clicks into place.

- The cover of the incident light reflector disk is secured using 2 DIN 912 M4x6 lens head Allen screws (Order No. 11 703 121 630 000)
Detach both screws.

Store the screws in a safe place. Never operate the microscope when the cover is not screwed into place!

- Change the reflector cubes as described for the Leica DMi8 with manual incident light reflector disk.
- Close and refasten the cover.

Leica DMi8C with manually coded incident light reflector disk

- The reflector cubes are manually swung into the beam path by turning the turret disk.
The LED allocated to the reflector cube and located on the front side of the microscope is illuminated.

Leica DMi8 A with incident light reflector disk and front operating panel with function keys

- Set a reflector cube by pressing the corresponding reflector cube changer button on the keypad on the front side of the microscope.
The LED of the reflector cube changer button illuminates.

Note:

For the Leica DMi8 A with touch screen, selection of a reflector cube can also be moved to one of the freely assignable function keys on the right side of the stand.

Changing reflector cubes → p. 60.

For coded reflector turret disk

For motorized reflector turret disk

5.12 Filter

5.12.1 Filter slider for incident light

- Insert the filter slider into the receptacle until it locks into place.

5.12.2 Transmitted light filter

For a tiltable transmitted light illuminator arm:

- The transmitted light filter and the dark stop are swung in and out on the right side of the transmitted light illuminator arm.

Note:

The front position is assigned at the factory using a shutter (dark stop), which is interchangeable with a filter.

For a fixed transmitted light illuminator arm:

- Insert the filter into the filter receptacle on the transmitted light illuminator arm.

Leica DMi8 with tiltable transmitted light illuminator arm:
Transmitted Light filter

Leica DMi8 with a fixed transmitted light illuminator arm:
Transmitted Light filter

5.13 The memory function on the touch screen

Up to 6 user-specific settings can be saved via the touch screen. In doing so, combinations of contrast method and objective are stored. One exception is fluorescence mode. In this mode, the combination of fluorescence filter and objective is stored.

- On the touch screen, use the  button to switch to the basic microscope settings.

- Select the **Memory Function** menu using the  key.

Teaching in a combination


- Set the desired combination.
- Press the key to be allocated.
- Then, press the **Set** key.

Deleting a saved combination

- Press the corresponding key.
- Then, press the **Clear** key.

Setting a saved combination

- Press the corresponding key.
- Finally, press **GoTo**.

6. Operation Contrast Methods

Do not look into the eyepieces when toggling the contrast method!

During the toggling procedure, the maximum radiant power of the light source may be present on the eyepieces for a short time and can temporarily blind the user!

Leica DMi8 M and Leica DMi8 C

The contrast method on the manual or coded Leica DMi8 can be operated using the manual focus, manual objective turret, or rotary knobs and slide bars on the stand.

- If necessary, first switch to the corresponding axis.



Leica DMi8 A

All contrast methods on the automated Leica DMi8 can be selected and operated via the touch screen, variable function keys or LAS software. The only exceptions are methods involving components that have to be operated manually (e.g. systems with manual analyzers). The following describes operation of the automated Leica DMi8 via the touch screen. Refer to the separate instructions for operation via the software.

- If necessary, switch to the corresponding axis.

- On the touch screen, use the  key to configure the contrast methods.

- Select the desired contrast method.


Note:

If all positions of the reflector turret disk have been assigned, the "DAPI" filter cube can be replaced with the "DAPI-TL" filter cube using LAS software. Transmitted Light contrast methods are possible with this filter cube.

6.1 Bright field (IL)


Leica DMI8 with manual incident light methods

- If necessary, insert the BF (= bright field) or Smith reflector on the reflector turret.
- Swing all remaining optical components (such as the filter slider, analyzer, polarizer or IC prisms) out of the beam path.
- Insert a suitable specimen.
- Select your objective.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.
- Insert the aperture diaphragm and illuminated field diaphragm.

Leica DMI8 with motorized incident light methods


- Select the contrast method **BF** (bright field).
- Insert a suitable specimen.
- Select a suitable objective in the magnification menu.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Set the brightness with the intensity rotary knob.

6.2 Dark field (IL)


Leica DMi8 with manual incident light methods

- Use the reflector turret to set the dark field position.
- Select a dark field objective.
- Fully open the aperture.
- Swing all remaining optical components (such as the filter slider, analyzer, polarizer or IC prisms) out of the beam path.
- Insert a dark field specimen.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.

Leica DMi8 with motorized incident light methods

- Select the contrast method **DF** (dark field).
- Insert a dark field specimen.
- Select a suitable objective in the magnification menu.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.


Note:

When the dark field method is selected, the aperture diaphragm is opened fully and cannot be adjusted.

6.3 Polarization (IL)


Leica DMI8 with manual polarizer and analyzer

- Slide the polarizer into the beam path.
- Insert the analyzer on the right side of the stand up to the click stop.
- Bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Select a suitable objective.
- Insert a specimen.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Set the brightness with the intensity rotary knob and approach a suitable objective.

Leica DMI8 with motorized incident light methods



- Select the **POL** (polarization) contrast method.
- The analyzer cube is automatically brought into the beam path.
- Slide the polarizer into the beam path.
- Bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Select a suitable objective in the magnification menu.
- Insert a specimen.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Set the brightness with the intensity rotary knob and approach a suitable objective.

6.4 Differential interference contrast (ICR)

Leica DMI8 with manual incident light methods

- Use the reflector turret to set the "ICR" reflector cube (contains polarizer and analyzer).
Alternatively:
Manually slide the ICR polarizer and the IC/P analyzer on the stand into the beam path.
- Select a suitable objective.
- Select the compatible objective prism (using the objective prism slider or objective turret.)
- Insert a specimen.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.
- For fine adjustment, use the adjusting screw below the objective turret or the knurled screw on the objective prism slider.

Manual DIC objective prism disk

Leica DMI8 with motorized incident light methods

- Select the **DIC** (Differential Interference Contrast) contrast method.
- The "ICR" reflector cube (contains polarizer and analyzer) and the compatible condenser prism are automatically brought into the beam path. The corresponding objective prism is also approached automatically.
- Place a DIC specimen on the stage.
- Select a suitable objective in the magnification menu.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.
- The adjustment is carried out via the + and – buttons on the touch screen.


Manual alternative:


- Slide the polarizer into the beam path manually.
- Then, also manually insert the analyzer on the right side of the stand until it locks into place.

Guide the objective and condenser prisms manually until a correct combination is displayed on the touch screen.

- For fine adjustment use the knurled knob below the objective turret.

Motorized DIC objective prism disk

6.5 Oblique light (IL)

Leica DMI8 with manual incident light methods

- First, configure the bright field basic setting as described in Item 6.1.

- Enable oblique light mode by pressing the

button. The LED is off.

- The desired effect is set using a slider.

- To disable oblique light mode, press the

button again. The LED illuminates.


Leica DMI8 with motorized incident light methods

- Select a suitable incident light contrast method, for example, **IL-BF**.

- On the touch screen, use the  button to switch to the basic microscope settings.


- Select the **Status/Illumination** menu using the  key.



- Enable oblique light mode by pressing the **Oblique** button. When oblique light is switched on, the button is red.

- The desired effect is set in the range from -120 to +120 using the + and – buttons.

- To disable oblique light mode, press the **Oblique** button again.

6.6 MACRO overview magnification

MACRO overview magnification enables you to illuminate a field with a diameter of approximately 40 mm. This is done using a special objective (MACRO objective) and the MACRO filter cube.

Note:

MACRO overview magnification can be used only in conjunction with the bright field incident light contrast method.

Note:

The MACRO objective is focused by trained Leica personnel when the instrument is installed. Due to the high depth of field of the MACRO objective, the image of the specimen is prefocused. Refocusing with the focus dial is not necessary.

Note:


The MACRO objective is parfocal with the Leica objectives.

Leica DMI8 with manual incident light methods

- Swing in the MACRO objective.
- Use the reflector turret to set the MACRO filter cube.
- Insert a specimen.
- Configure the brightness on the intensity rotary knob.

Leica DMI8 with motorized incident light methods

- Press the **MACRO** button.
The MACRO objective and MACRO filter cube are approached automatically.
- Insert a specimen.
- Configure the brightness on the intensity rotary knob.

6.7 Fluorescence


Note:

To replace filter cubes, refer to Chapter 5.11, Reflector turret disk →. p. 60.

Leica DMi8 with manual incident light reflector disk

- Where appropriate, pull the filter slider and analyzer slider out completely to open the beam path.
- Set the light intensity.
- The filter cubes are swiveled into the beam path manually by turning the incident light reflector disk.
- A color mark on the edge of the disk indicates the current filter cube. If the shipped configuration is changed by inserting another filter cube, note the color of the position where the filter cube was inserted so you can allocate it later.

Leica DMI8 with manually coded incident light reflector disk

- The filter cubes are swiveled into the beam path manually by turning the incident light reflector disk. The LED allocated to the filter cube and located on the front side of the microscope is illuminated.

Leica DMI8 manually coded with motorized incident light reflector disk

- Set a filter cube by pressing the corresponding reflector changer button on the keypad on the front side of the microscope. The LED of the reflector changer button illuminates.

For coded incident light reflector disk


For motorized incident light reflector disk

Leica DMI8 with motorized incident light reflector disk

Operating the contrast methods using the touch screen:

- On the touch screen, use the  key to configure the contrast methods.
- For incident light, select the corresponding contrast method.
- Where appropriate, pull the filter slider and analyzer slider out completely to open the beam path.
- Insert a specimen.
- Select a suitable objective in the magnification menu.

Note:

For the Leica DMI8 A with touch screen, selection of a filter cube can be moved to one of the freely assignable function keys on the right side of the stand.

- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Set the light intensity.

6.8 Bright field (TL)

Leica DMi8 with manual transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- If necessary, use the reflector turret to set the TL bright-field position or an empty position, or select the "DAPI-TL" filter cube.

For a tiltable transmitted light illuminator arm:

- Select the bright-field position on the condenser **BF**.

For a fixed transmitted light illuminator arm:

- For correct height adjustment of the S80/0.30 and S40/0.45 condensers, the stand has marks (1). These marks pertain to a liquid height of 15 mm. For stands with double marks, the upper and lower mark represent a range with different liquid heights. Move the condenser holder while observing the specimen to attain an optimal image.


Note:

For the stands with 3-plate mechanical stages, position the illuminator arm 25 mm below the marks, as the illumination axis is installed 25 mm higher on an adapter than on the fixed stage.

- Press the stop lever (2) and adjust the transmitted light illumination arm until the upper edge of the carrier is even with the corresponding condenser height mark.
- Swing all remaining optical components, such as the analyzer, polarizer or IC prisms, from the beam path.
- Insert a transmitted light specimen.
- Select your objective.
- Focus the image using the focus dials.

Condenser setting for tiltable transmitted light illuminator arm - configuring **BF** position

Condenser setting for fixed transmitted light illuminator arm

- 1 Marks
- 2 Stop lever for condenser adjustment
- 3 Aperture diaphragm

- Configure the brightness on the intensity rotary knob.
- For an optimal field diaphragm setting, check the Koehler illumination.

Leica DMI8 with motorized transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- Select the contrast method **BF** (bright field).

The bright-field position is approached for the motorized condenser.

An empty position or the filter cube "DAPI-TL" on the fluorescent turret disk is approached automatically.

- Insert a transmitted light specimen.
- Select a suitable objective in the magnification menu.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.


Note:

You can find additional information regarding the correct settings for Koehler illumination in Leica Science Lab.

Link:

<http://www.leica-microsystems.com/science-lab/>

6.9 Phase Contrast (TL)

Leica DMi8 with manual transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- If necessary, use the reflector turret to set the TL bright-field position or an empty position, or select the "DAPI-TL" filter cube.
- Select a phase contrast objective. Objectives that are suitable for phase contrast are engraved with PH.

For a tiltable transmitted light illuminator arm:

- Select the corresponding light ring on the condenser.

For a fixed transmitted light illuminator arm:

- Adjust the condenser height.
- Guide the phase contrast light ring slide with the marking "TOP LEFT" on the top left laterally into the receptacle until it clicks into place.
(Caution: there are various slide bars for the condensers!)
- Use the light ring of the light ring slide that corresponds to the magnification (5, 10/20 or 40).

Note:

The phase contrast light ring slide is coded. As soon as the slider is taken out of the bright-field position (available opening) and engaged in a phase contrast position (any light ring), the brightness increases. Conversely, switching from phase contrast to bright field decreases the brightness.


- 1 Slider for light rings (S80/0.30)
- 2 Slider for light rings (S40/0.45)
- 3 Click stop
- 4 Light rings
- 5 Bright-field position
- 6 "TOP LEFT" marking

Condenser setting for tiltable transmitted light illuminator arm - configuring light ring, e.g. **PH1**

Phase contrast for fixed transmitted light illuminator arm

Light ring slide for fixed transmitted light illuminator arm

- Completely open the aperture on the condenser.
- Swing all remaining optical components, such as the analyzer, polarizer or IC prisms, from the beam path.
- Insert a phase contrast specimen.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.

Leica DMI8 with motorized transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- Select the **PH** (phase contrast) contrast method.


The correct light ring on the motorized condenser is swung inwards.

- Insert a transmitted light specimen.
- Select a suitable objective in the magnification menu. Objectives that are suitable for phase contrast are engraved with PH.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.

Note:

When selecting the phase contrast method, the aperture diaphragm is opened fully and can not be adjusted.

6.10 Dark field (TL)

Leica DMI8 with manual transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- If necessary, use the reflector turret to set the TL bright-field position or an empty position, or select the "DAPI-TL" filter cube.
- Select a dark field objective.
- Select the corresponding dark field stop (light ring) on the condenser.
- Completely open the aperture on the condenser.
- Swing all remaining optical components, such as the analyzer, polarizer or IC prisms, from the beam path.
- Insert a dark field specimen.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.


Condenser setting for tiltable transmitted light illuminator arm - configuring dark field stop (light ring)

Leica DMI8 with motorized transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- Select the contrast method **DF** (dark field).

The dark field ring (dark stop) is swung in for the motorized condenser.

- Insert a transmitted light specimen.
- Select a suitable objective in the magnification menu.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.

Note:

The maximum applicable objective aperture for the dark field is for the condenser S1 0.70 and for the condenser S23/S28 0.40.

Note:

When the dark field method is selected, the aperture diaphragm is opened fully and cannot be adjusted.


6.11 Polarization (TL)

Leica DMI8 with manual transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- If necessary, use the reflector turret to set the TL bright-field position or an empty position, or select the "DAPI-TL" filter cube.
- Select the bright-field position on the condenser.
- Swing all IC prisms out of the beam path.
- Select your objective.
- Swing the polarizer on the condenser into the beam path.
- Insert the analyzer on the right side of the stand up to the click stop.
- Bring the polarizer and analyzer into cross position until they reach maximum darkness.
- Insert a specimen.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.

Condenser setting for tiltable transmitted light illuminator arm - configuring **BF** position

Leica DMI8 with motorized transmitted light method

- If necessary, use the filter slider to set the TL bright-field position.
- Select the **POL** (polarization) contrast method.

Motor procedure:

- After the **POL** contrast method has been selected, the polarizer is swung in (if the microscope is equipped with these components). The analyzer cube is also automatically brought into the beam path.


Combined procedure:

- You also have the option of combining purely mechanical and motorized components, i.e. you can combine a mechanical analyzer and a motorized polarizer.

Manual procedure:

For a manual swing polarizer and a mechanical analyzer, proceed with the manual transmitted light method as described for the Leica DMI8.

- Select a suitable objective in the magnification menu.
- Insert a specimen.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.


6.12 Differential interference contrast – ICT (TL)

Leica DMI8 with manual transmitted light method


- If necessary, use the filter slider to set the TL bright-field position.
- If necessary, use the reflector turret to set the TL bright-field position or an empty position, or select the "DAPI-TL" filter cube.
- Select an objective.
- Select the corresponding Wollaston prism condenser on the condenser.
- Select the corresponding Wollaston prism objective on the objective turret.
- Swing the polarizer on the condenser into the beam path.
- Insert the analyzer on the right side of the stand up to the click stop.
- Insert a specimen.
- Focus the image using the focus dials.
- Configure the brightness on the intensity rotary knob.
- For fine adjustment, use the adjusting screw below the objective turret.

Condenser setting for tiltable transmitted light illuminator arm - configuring Wollaston prism

Objective prism slider

Manual DIC objective prism disk

Leica DMI8 with motorized transmitted light method



- If necessary, use the filter slider to set the TL bright-field position.
- Select the **DIC** (Differential Interference Contrast) contrast method.
- The polarizer located in the condenser and the fitting condenser prisms are automatically brought into the beam path. The corresponding objective prism and the analyzer cube are also approached automatically.
- Place a DIC specimen on the stage.
- Select a suitable objective in the magnification menu.
- Focus the image using the rotary knob on the Leica SmartMove or using the focus dial.
- Configure the brightness on the intensity rotary knob.
- The adjustment is carried out via the + and – buttons on the touch screen.

**Manual alternative:**

- Manually swing the polarizer on the condenser into the beam path.
- Then, also manually insert the analyzer on the right side of the stand until it locks into place.

Guide the objective and condenser prisms manually until a correct combination is displayed on the touch screen.

- For fine adjustment use the knurled knob below the objective turret.



6.13 Combination method

Leica DMI8 with motorized incident and transmitted light methods

- If necessary, use the filter slider to set the TL bright-field position.

- On the touch screen, use the  key to configure the contrast methods.
- Insert a specimen.
- Select a suitable objective in the magnification menu.
- For the combination, select the contrast method **BF-BF**.
- The illumination settings for the incident and transmitted light axes can be adjusted separately.
To do so, toggle between the illuminators using the **TL/Fluor** button.

7. Troubleshooting

Problem	Cause/remedy
Stand	
The microscope does not respond.	▶ Ensure that the AC outlet has power. ▶ Make sure that the electronics box or the microscope is connected to the power supply. ▶ Check the cable connections. ▶ Inform Service and have the supply unit fuse checked.
Illumination	
The image is completely dark.	▶ Check the connection to the lamp housings on the microscope. ▶ Make sure the lamps are not defective and are connected to the power supply. ▶ Select the corresponding illumination axis. ▶ Switch over to the active port. Check the position of the beam splitter (camera/tube). ▶ Inform the service department and have them check whether a fuse is defective.
Bright field	
The specimen cannot be brought into focus.	▶ Use the correct immersion medium. ▶ Make sure that you are using an objective with cover slip correction. ▶ Adjust the correction ring if present on the objective.
The image shows too little contrast.	▶ Close the aperture diaphragm. Aperture diaphragm open = low contrast, low depth of field, high resolution Aperture diaphragm closed = high contrast, high depth of field, low resolution
The image is not uniformly illuminated.	▶ Remove all unneeded filters from the light path. ▶ Check whether oblique light is switched on.

Problem	Cause/remedy
Dark field	
No definite DF contrast is possible.	‣ Be sure that a DF objective is being used. ‣ The light intensity is too low. ‣ Open the aperture diaphragm. For DF in the transmitted light: ‣ The objective aperture setting is too high: Maximum 0.7 for condenser S1 Maximal 0.4 for condenser S23/28 ‣ Check the condenser centering.
The image is not uniformly illuminated.	‣ Open the aperture diaphragm. ‣ Check whether oblique light is switched on. ‣ Remove all unneeded filters from the light path.
Phase contrast (TL)	
Phase contrast cannot be adjusted.	‣ The specimen is too thick, too thin or too brightly stained. ‣ Refractive indexes of the mounting medium and specimen are identical, so that there is no phase jump. ‣ The cover glass is not placed uniformly. ‣ Check that the correct light ring is positioned. ‣ Check that the correct light slider is positioned. ‣ The S40/0.45 condenser has been switched with the S80/0.30 condenser. ‣ Open the aperture diaphragm completely. ‣ Check the centering of the light rings.
Polarization	
Polarization contrast cannot be adjusted.	‣ Bring the polarizer and analyzer into cross position until they reach maximum darkness (for transmitted light without specimen). ‣ Check whether there may be ICR/ICT prisms in the beam path.

Problem	Cause/remedy
Transmitted-light interference contrast	
No interference contrast can be set (IL).	▶ Bring the polarizer and analyzer into cross position until they reach maximum darkness. ▶ Check whether the IC prism compatible with the objective being used has been selected.
No interference contrast can be set (TL).	▶ The specimen is too thick or too thin. ▶ The mounting medium or specimen are of birefringent material. Rotate the specimen. ▶ The difference in the refractive indices of the specimen and the mounting medium is too small. ▶ The cover slip is too thick. ▶ Check the Koehler illumination. ▶ Bring the polarizer and analyzer into cross position until they reach maximum darkness (without specimen). ▶ Check whether the compatible condenser prism and the corresponding objective prism are configured. ▶ Check the correct placement of the IC condenser prisms.
Fluorescence	
The image is completely dark (no fluorescence).	▶ Select the incident light axis (IL). ▶ Check the selected fluorescence filter – fluorochrome combination. ▶ Check the mirror housing to see if the correct path has been selected.
The fluorescence is too weak.	▶ Unsuitable specimen (stored incorrectly, too old, bleached). ▶ Unspecific filter combination. ▶ Objective with too low of a numerical aperture. ▶ Eyepiece magnification too high. ▶ Room where microscope is located too bright. ▶ Incorrect beam split on the trinocular tube. ▶ Secondary light due to reflection at condenser.

8. Care of the Microscope

Caution!

Unplug the power supply before performing cleaning and maintenance work!
Protect electrical components from moisture!

Microscopes in warm and warm-damp climatic zones require special care in order to prevent the build up of fungus. The microscope should be cleaned after each use, and the microscope optics should be kept painstakingly clean.

8.1 Dust cover

Note:

To protect against dust, cover the microscope and accessories with the dust cover after each use.

Let the microscope and lamp housings cool down before covering the stand with a dust cover. The dust cover is not heat-resistant. In addition, condensation may occur.

8.2 Cleaning

Caution:

Residual fiber and dust can create unwanted background fluorescence during fluorescence microscopy.

Cleaning coated parts

Dust and loose dirt particles can be removed with a soft brush or lint-free cotton cloth.

Clinging dirt can be cleaned as necessary with a low-concentrated soap solution, petroleum ether or ethyl alcohol. For cleaning coated parts, use a linen or leather cloth that is moistened with one of these substances.

Caution:

Thinners containing acetone, xylene or nitrogen can harm the microscope and thus must not be used.

Test cleaning solutions of unknown composition on a less visible area of the unit first. Be sure that coated or plastic surfaces do not become matted or etched.

Cleaning the specimen stage

Remove light-colored spots on the stage by rubbing with paraffin oil or acid-free Vaseline.

Cleaning glass surfaces and objectives

The cleaning of glass surfaces and objectives in particular should be carried out exclusively as outlined in the brochure "Cleaning of Microscope Optics". The information can be downloaded at:

<http://www.leica-microsystems.com/products/light-microscopes/industrial-materials/inverted-microscopes>

Select the microscope type and go to the "Download" page.

In case of questions, please contact our Technical Service.

Removing immersion oil

Caution!

Follow safety notes for immersion oil!

First, wipe off the immersion oil with a clean cotton cloth, and then re-wipe the surface several times with ethyl alcohol.

8.3 Handling acids and bases

For examinations using acids or other aggressive chemicals, particular caution must be taken.

Caution:

Never allow the optics and mechanical parts to come into direct contact with these chemicals.

9. Essential Consumable and Spare Parts

Order No. Material No.	Name	Used for
<u>Screw cap for unused objective receptacles</u>		
020-422-570-000	Screw cap M 25	Objective turret
11 020-422-580-028	Screw cap M 32	Objective turret
<u>Cover for unoccupied objective-DIC-disk opening</u>		
11 090-144-020-088	DIC cover	Microscope stand
<u>Dust and light stop for analyzer opening</u>		
11 020-530-100-013	Analyzer opening cover	Microscope stand
11 020-437-101-013	Polarizer opening cover	Microscope stand
<u>Dust and light stop for camera port openings</u>		
11 020-387-556-009	Port opening cover	Microscope stand
<u>Replacement eyecup (antiglare protection) for HC PLAN eyepiece</u>		
021-500-017-005	HC PLAN eyecup	10x/25 eyepiece
021-264-520-018	HC PLAN eyecup	10x/22 eyepiece
021-264-520-018	HC PLAN eyecup	10x/20 eyepiece
<u>Immersion oil in accordance with DIN/ISO standards, fluorescence-free</u>		
11 513 859	Type F, ISO 8036, Very low self-fluorescence. Expressly recommended for fluorescence applications and all APO objectives, 10 ml	OIL and IMM objectives and oil condenser heads
11 513 860	Type N, ISO 8036, Low self-fluorescence, 20 ml	
11 513 861	Type N, ISO 8036, Low self-fluorescence, 250 ml	
<u>Screws for the fluorescence turret cover</u>		
11 703 121 630 000	DIN 912 Allen screw, M4x6 lens head	Cover of the fluo turret disk for laser microscope systems
<u>Fuses</u>		
11 362 150 010 160	Fine-wire fuse T1.6A, 250 VAC, Breaking capacity H	Leica DMi8 M / Leica DMi8 C


Additional spare parts such as cover caps, stickers, screws, dust covers etc., on request.

10. Appendix

10.1 Assigning the variable function keys

The variable function keys can always be defined by the user via the LAS software.

Depending on the configuration, useful functions are already assigned to the function key upon delivery.

Preassigning the function keys:

Industry applications	
Microscope with Magnification Changer:	
Button 1	TL
Button 2	Photo
Button 3	Toggle Magnification Changer
Button 4	Toggle ports

Industry applications	
Microscope without Magnification Changer:	
Button 1	Switch through objectives to the right
Button 2	Photo
Button 3	Switch through objectives to the left
Button 4	Toggle ports

Possible assignments of the variable function keys to the stand and SmartMove:

Function button	Function
BF (IL)	Bright field incident light
DF (IL)	Dark field incident light
POL (IL)	Polarization incident light
ICR	Interference contrast incident light
OBL	Oblique light
CHANGE IL ●	Switches through all incident light contrast methods
BF (TL)	Bright field transmitted light
PH (TL)	Phase contrast transmitted light
ICT (TL)	Interference contrast transmitted light
DF (TL)	Dark field transmitted light
POL (TL)	Polarization transmitted light
CHANGE TL ●	Switches through all transmitted light contrast methods
INT ↑	Increases the brightness (Transmitted Light)
INT ↓	Reduces the brightness (Transmitted Light)
AP ↑	Opens aperture diaphragm (Transmitted Light)
AP ↓	Closes the aperture diaphragm (Transmitted Light)
FD ↑	Opens the field diaphragm (Transmitted Light)
FD ↓	Closes the field diaphragm (Transmitted Light)
CHG FLT	Switch TL filter
FLUO	Fluorescence (last filter cube)
CUBE 1-6	Chooses the fluorescence cube at position 1-6
CUBE ↻	Switches through the fluorescence cube clockwise (1 → 6)
CUBE ↺	Switches through the fluorescence cube counterclockwise (6 → 1)
CHG FLUO	Switches through all fluorescence cubes
INT FLUO ↑	Increases the brightness (fluorescence)
INT FLUO ↓	Reduces the brightness (fluorescence)
FD FLUO ↑	Opens the field diaphragm (fluorescence)
FD FLUO ↓	Closes the field diaphragm (fluorescence)
CHG FW	Toggles the filter wheel functions
IFW	Activates the internal filter wheel
COMBI ●	Combination method (PH fluorescence or DIC fluorescence)
CHG CB ●	Switches through all combination methods
CHG OBJ ↻	Switches through objectives clockwise
CHG OBJ ↺	Switches through objectives counterclockwise
DRY/IMM	Switch dry/immersion
CHG ↻ Mag.-Chg.	Switches through the Magnification Changer clockwise
XY PRECISE	Activate stage precisely
1x Tube lens	Activate lens 1x
Z COARSE	Enable coarse focus
XY PRECISE	Activate stage precisely
XY FAST	Activate stage quickly
Acquire Image/	
Photo	Capture an image with the camera


Additional abbreviations for objective names, filter name, magnifications.

10.2 Abbreviation glossary


AP	Aperture diaphragm
BF	Bright field
COMBI	Combination methods
CUBE	Fluorescence cube
DF	Dark field incident / transmitted light
DIC	Differential interference contrast
FD	Field diaphragm
FLUO/Fluor	Fluorescence axis (incident light)
ICR	Interference contrast, incident light
ICT	Interference contrast transmitted light
IFW	Internal filter wheel
IL	Incident light
INT	Intensity
OBL	Oblique light
PH	Phase Contrast
POL	Polarization, incident / Transmitted Light
TL	Transmitted Light
SET	Set focus position or lower threshold
$\overline{\uparrow}$ Z	Focus position
$\underline{\downarrow}$ Z	Lower threshold

10.3 Overview of touch screen pictograms


Navigation bar


Menus


Z-drive


Objectives


Illumination


Observation ports

Control

Various

10.4 Additional information

General microscopy

Refer to the Leica Science Lab

<http://www.leica-microsystems.com/science-lab/>

to find additional useful information and instructions on the microscope methods.

Under the category "Basics in Microscopy", you can find, for instance:

- The instructions for setting the correct Koehler illumination
- The procedure for contrast methods such as
 - Phase contrast
 - Modulation contrast
 - Differential interference contrast
- Fluorescence techniques
- The instructions for cleaning the microscope optics

and much more.

Objectives

For information on the objectives such as the objective classes and the marking of objectives, refer to

<http://www.leica-microsystems.com/de/produkte/objektive/>

Accessories

The "Live on Stage" brochure provides an overview of the accessories for live cell microscopy including the various specimen stages, specimen holders and inserts.

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3-plate cross-stage 49

A

Abbreviation Glossary 94
Accessories 96
Acids and bases 20, 90
Activating the oblique light 30
Adjusting the oblique light 33
Ambient conditions 18, 19
Analyzer 32, 67, 81, 83
Analyzer cube 84
Analyzer slot 33
Aperture diaphragm 59
Aperture diaphragm adjustment (IL) 33
Aperture diaphragm adjustment (TL) 35
Aperture diaphragm (IL) 57
Aperture Diaphragm Setting 30
Aperture diaphragm (TL) 58

B

Basic microscope settings 36
Beam splitting in photo tubes 47
Bright field (IL) 65
Bright field (TL) 75
Brightness Setting 30, 33, 43

C

Centering of aperture diaphragm (IL) 32, 33
Centering of illuminated field diaphragm (IL) 32, 33
Changing objectives 32, 44
Changing reflector cubes 60
Cleaning 89
Color filter 33
Combination method 85
Condenser 27
Condenser height adjuster 35
Condenser height stop 35
Condenser setting 35, 75, 77, 79
Contrast methods 25, 36, 64
Contrast Methods 19

Controls 28, 29
Correction for vision problems 48
Cover of the incident light reflector disk 12, 61

D

„DAPI“, „DAPI-TL“ filter cubes 64
Dark field (IL) 66
Dark field (TL) 79
Dark stop 79
Declaration of Conformity 99
DIC fine adjustment 34
Differential interference contrast (IL) 68
Differential interference contrast (TL) 83
Dimensions 22
Display brightness 37
Disposal 20
Dust protector 89

E

EC Declaration of Conformity 99
Electronics box 28
EU Directives 6
Eyepieces 48

F

Field diaphragm 59
Filter cube 73, 74
Filter holder (TL) 35
Filter slider 32, 62, 73
Fixed stage 49
Fluorescence 73
Focus dials 52
Focus drive 54
Focusing 30, 32, 33, 52
Focus position 30, 53, 54
Focus stages 52
Focus threshold 30, 33, 53
Fuses 91

G

Graticule 48

H

Height compensation plate 24
Holder for color filter 33

I

IC/P analyzer 68
ICR 68
ICR polarizer 68
ICT 83
Identification sheet 30
Illuminated field diaphragm adjustment 30
Illuminated field diaphragm adjustment (IL) 32
Illuminated field diaphragm adjustment (TL) 35
Illuminated field diaphragm (IL) 57
Illuminated field diaphragm (TL) 58
Illumination 43
Immersion objectives 45, 46
Immersion objectives, lockable 46
Immersion oil 20, 90, 91
Incident light axis 25
Incident light reflector disk 26, 31, 60, 73, 74
Installation location 19
Interfaces 29
Interlock 11
Interlock connector 16
Interlock switch 16
Interpupillary distance 47

K

Koehler illumination 96

L

Language selection 37
Laser class 1 17
Laser coupling 12, 60
Laser radiation 7, 13
Laser Radiation 14, 15
Laser safety 11
Laser Safety Kit 16
Laser safety precautions 7
LED lamp housing 43
Leica CTR electronics box 22, 28, 29
Leica EL6000 42, 43
Leica Science Lab 96
Leica SmartMove 29, 38, 49, 52
Leica STP4000/8000 29, 39
Light guide 42, 43
Light ring 77
Light ring slide 35, 77
Light sources 19, 42, 43
Linear motorized stage 49

M

MACRO filter cube 72
MACRO objective 72
MACRO overview magnification 72
Magnification 36, 45
Magnification Changer 26, 32, 56
Memory function on the touch screen 63

O

Object displacement 49
Objective changer buttons 31
Objective Prism 34, 68, 69, 84
Objective prism slider 68
Objectives 26, 44, 45, 96
Objectives with corrective mount 46
Objective turret 26, 31, 45
Oblique button 71
Oblique light (IL) 70
Observation ports 27, 28, 45
On/off switch 32, 41

P

Parfocality 46
Phase contrast (TL) 77
Pictograms 95
Polarization (IL) 67
Polarization (TL) 81
Polarizer (IL) 32, 67
Polarizer (TL) 35, 81, 83
Positioning range of the stage 50
Preassigning the function keys 92

R

Reflector change 32, 33
Reflector changer buttons 31, 74
Reflector cube, ICR 68

S

Safety Notes 8
Safety switch 11
Scanning stage 49
Selecting the ports 55
Service port 37
Side-Port 32, 33, 55
Software, LAS 38, 40, 46, 48, 49, 50, 51, 52, 54, 64, 92
Spare parts 91
Specifications 8
Speed switch (focus) 52
Stage and focus controls 36
Stage movement 32
Stage positions 50
Stages 26, 27, 49
Step sizes 50, 54
Storage 18
Switching on 41
Symbols 4

T

Technical data for the interlock connector 11
Threshold 54
Toggling between the illumination axes 30
Top-Port 33, 55
Total magnification 45, 48, 56
Touch Screen 20, 29, 30, 36
Touch screen calibration 37
Transmitted light axis 25
Transmitted light filter 62
Transmitted Light filter 35
Transport 18
Travel speed of the stage 50
Tube 25, 47
Tube setting 33
Type plate 32

V

Variable function keys 30, 33, 38, 45, 92, 93
Version number 37
Viewing angle 47

W

Weight 18, 22, 23, 24
Wollaston prism 83

Z

Z-drive 27

12. EU Declaration of Conformity

To download the EC Declaration of Conformity, use this link

<http://www.leica-microsystems.com/products/light-microscopes/industrial-materials/inverted-microscopes>

Select the microscope type and go to the "Download" page.

People's Republic of China
- Administrative Measures on the Control of Pollution Caused by Electronic Information Products -

部件名称 Name of the part	有毒有害物质或元素 Hazardous substances					
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (Cr ⁶⁺)	多溴联苯 (PBB)	多溴二苯醚 (PBDE)
印刷电路板 printed circuit boards	X	0	0	0	0	0
电子元件 electronic components	X	0	0	0	0	0
机械部件 mechanical parts	X	0	0	X	0	0
光学元件 optical components	X	0	X	0	0	0
电缆 cables	0	0	0	0	X	X
光源 light sources	0	X	0	0	0	0

- o : 表示该有毒有害物质在该部件中的含量均在SJ/T 11363-2006 标准规定的限量要求以下。
Indicates that the concentration of the hazardous substance in all materials in the parts is below the relevant threshold of the SJ/T 11363-2006 standard.
- x : 表示该有毒有害物质至少在该部件的某一材料中的含量超出SJ/T 11363-2006 标准规定的限量要求。
Indicates that the concentration of the hazardous substance of at least one of all materials in the parts is above the relevant threshold of the SJ/T 11363-2006 standard.

Note: The actual product may or may not include in all the part types listed above

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