


System diagram Leica DM500, DM750



Plan FOV 18 objectives

13 613 240 - Plan 4x/0.10 NA, 26.2mm *
13 613 241 - Plan 10x/0.22 NA, 7.8mm *
13 613 242 - Plan 40x/0.65 NA, 0.31mm *
13 613 243 - Plan 100x/1.25 NA, 0.10mm *, Oil
* working distance

HI Plan FOV 20 objectives

11 506 226 - HI Plan 4x/0.10 NA, 18.0mm *
11 506 228 - HI Plan 10x/0.25 NA, 12.0mm *
11 506 276 - HI Plan 20x/0.40 NA, 0.92 mm *
11 506 236 - HI Plan 40x/0.65 NA, 0.36mm *
11 506 237 - HI Plan 63x/0.75 NA, 0.31mm *
11 506 238 - HI Plan 100x/1.25 NA, 0.10mm *, Oil
11 506 322 - HI Plan 100x/0.80 NA, 0.3mm* dry objective
* working distance

HI Plan Phase FOV 20 objectives

11 506 230 - HI Plan Phase 10x/0.25 NA, 12.0mm *.
11 506 278 - HI Plan Phase 20x/0.40 NA, 0.92mm *.
11 506 240 - HI Plan Phase 40x/0.65 NA, 0.36mm *.
11 506 241 - HI Plan Phase 100x/1.25 NA, 0.10mm * Oil
* working distance

Replacements

13 613 563 Eyeguard pair
13 WFAG30001 Fuses for DM500, DM750
13 613 571 Thumbscrew pair for DM750 condenser mount
13 613 572 Set screw pair for DM750 condenser mount
13 613 584 Dust cover DM500, DM750
13 613 210-100 DM500 User documents
13 613 510-100 DM750 User documents
13 613 705-100 ICC50 User documents
13 613 573 Add on lens for condensers with DM750 (included with stand)
11 505 249 Rechargeable battery set

Analyzer Module Sliders

11 555 045 Fixed analyzer slider
11 555 079 180° Rotatable analyzer slider
11 555 080 360° Rotatable analyzer slider

Others

13 613 262 Captive thumbscrew for EZ tube
13 613 261 Rubber grips for X/Y controls
13 614 800 Immersion oil
11 513 106 Stage micrometer
11 505 091 Diffuser slider for using 2.5x objective with Abbe condenser
11 505 507 Auxiliary lens for using 2.5x objective with 11501183 flip top condenser
13 613 566 32mm filter holder
13 596 016 Hard carrying case
10 450 245 Soft carrying case
13 613 562 Simple pol kit
11 505 250 Solar panel

International Power Cords

13 613 900 U.S.
13 613 901 Europe Continental
13 613 902 U.K.
13 613 903 Switzerland
13 613 904 Denmark
13 613 905 Italy
13 613 906 Australia
13 613 907 China
13 613 908 Japan
13 613 909 Israel
13 613 910 South Africa
13 613 911 India
13 613 912 Argentina